Dimension® clinical chemistry system

Flex® reagent cartridge

SAL

See shaded sections: Updated information from 2017-12 version.


Issue Date 2019-04-22

Salicylate

Intended Use: The SAL method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to measure salicylates, a class of analgesic, antipyretic and anti-inflammatory drugs (including aspirin) in human serum. Salicylate test results may be used in the diagnosis and treatment of salicylate overdose and for monitoring salicylate levels during therapy.

Summary: Salicylates (aspirin, acetylsalicylic acid) have analgesic, antipyretic, and anti-inflammatory properties and have been used for centuries to relieve pain. Salicylate overdose may cause intoxication. Measurement of salicylate concentration is important for assessment of the severity of intoxication.¹

Principles of Procedure: The salicylate method is a modification of the Trinder colorimetric technique.² Under acidic conditions, ferric nitrate forms a complex with salicylate. The amount of salicylate-ferric complex formed is proportional to the salicylate concentration and is measured using a two filter (510, 700 nm) endpoint with sample-reagent blanking technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1, 2	Liquid	Surfactant	0.125 mg
3, 4, 5	Liquid	Nitric acid	7.15 µmol
6, 7, 8	Liquid	Ferric Nitrate	3.3 µmol

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

Risk and Safety:

H290, H319, H317
P272, P280, P305 + P351 + P338, P501

Warning!

May be corrosive to metals. Causes serious eye irritation. May cause an allergic skin reaction.

Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Nitric acid; Iron (III) nitrate nonahydrate; 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 5 days for wells 1 – 8

Specimen Collection and Handling: Normal procedures for collecting serum may be used for samples to be analyzed by this method.³

Lithium heparin and EDTA anticoagulants do not interfere with the SAL method at concentrations normally found in blood collection tubes.

Follow the instructions provided with your specimen collection device for use and processing.⁴

Complete clot formation should take place before centrifugation.⁵

Specimens should be free of particulate matter.

Specimens are stable for 7 days at 20 – 25 °C, 2 weeks at 2 – 8 °C and 6 months frozen at -20 °C or colder.⁶

Procedure

Materials Provided

SAL Flex® reagent cartridge, Cat. No. DF20

Materials Required But Not Provided

SAL Calibrator, Cat. No. DC38
Quality Control Materials

Test Steps

Sampling,^c reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume. Precise container filling is not required.

Test Conditions

Sample Volume	15 µL
Reagent 1 Volume	50 µL
Reagent 2 Volume	50 µL
Reagent 3 Volume	50 µL
Temperature	37 °C
Wavelength	510 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	0.2 – 100 mg/dL [0.01 – 7.24 mmol/L] ^d
Calibration Material	Salicylate Calibrator, Cat. No. DC38
Calibration Scheme	3 levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.0724) = [mmol/L]
Typical Calibration Levels	0.0, 20.0, 100.0 mg/dL [0, 1.45, 7.24 mmol/L]
Calibration Frequency	Every 3 months for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -1.330 C ₁ 0.666

d. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known salicylate concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of SAL in mg/dL [mmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.2 – 100.0 mg/dL [0.01 – 7.24 mmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 100.0 mg/dL [7.24 mmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water to obtain result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

Results less than 2.8 mg/dL [0.20 mmol/L] should be reported as "less than 2.8 mg/dL [0.20 mmol/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide for troubleshooting.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
20 mg/dL [1.45 mmol/L]	>1.0 mg/dL [0.07 mmol/L]
100 mg/dL [7.24 mmol/L]	>2.0 mg/dL [0.14 mmol/L]

Interfering Substances

Sodium azide at a concentration of 0.05% increases salicylate results by less than or equal to 65 mg/dL [4.71 mmol/L].

Lipemia (Intralipid®) of 600 mg/dL [6.78 mmol/L] and above at a SAL concentration of 33 mg/dL [2.36 mmol/L] tripped a test report message; therefore the magnitude of the interference could not be determined.

Bilirubin (unconjugated) of 20 mg/dL [342 µmol/L] at a SAL concentration of 33 mg/dL [2.39 mmol/L] decreased SAL results by 22%.

At salicylates concentration of 20 and 40 mg/dL [1.45 and 2.90 mmol/L], salicylic acid at 100 mg/dL increases the SAL result by 55%.

Therapeutic doses of ceftriaxone (Rocephin®) may produce falsely elevated results with this assay. A falsely elevated result should be interpreted with caution and confirmed by another method.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Rocephin® is a registered trademark of Hoffmann-La Roche (Roche) Pharmaceuticals, Nutley, NJ.

Therapeutic Range: 2.8 – 20.0 mg/dL [0.20 – 1.45 mmol/L]⁷

SAL concentrations less than 2.8 mg/dL [0.20 mmol/L] should be considered negative (serum blank).
Therapeutic salicylate concentrations vary significantly, depending on the individual. A range of 2.8 – 20 mg/dL [0.20 – 1.45 mmol/L] includes effective serum concentrations for many patients; however, some individuals are best treated at concentrations outside this range. Concentrations greater than 30 mg/dL [2.17 mmol/L] are often associated with toxic symptoms. Concentrations greater than 60 mg/dL [4.34 mmol/L] are often lethal.²

Specific Performance Characteristics^a

Material	Precision ^f		
	Mean	Standard Deviation (% CV)	
	mg/dL [mmol/L]	Within-run	Between-day
Serum Pool			
Level 1	25.9 [1.9]	0.22 [0.02] (0.84)	0.26 [0.02] (1.00)
Level 2	46.6 [3.4]	0.22 [0.02] (0.48)	0.37 [0.03] (0.79)
Level 3	88.6 [6.4]	0.33 [0.03] (0.38)	0.76 [0.06] (0.86)
Dade® IAC, TRI	9.8 [0.7]	0.17 [0.01] (1.79)	0.22 [0.02] (2.29)

e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
f. Specimens at each level were analyzed in duplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

Dade® is a registered trademark of Siemens Healthcare Diagnostics, Deerfield, IL 60015.

Method Comparison

Regression Statistics^g

Comparative Method	Slope	Intercept	Correlation Coefficient	n ^h
SAL Method on the on the aca® discrete clinical analyzer				
Site 1	1.03	-0.33	0.999	100
Site 2	1.03	-0.28	0.999	100
TDx® Salicylate				
Site 1	1.01	-0.02	0.997	100
Site 2	1.02	1.35	0.997	99

g. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.
h. Clinical specimens from patients undergoing salicylate therapy.

TDx® is a registered trademark of Abbott Diagnostics, Abbott Park, IL 60064.

Specificity

HIL Interference

The SAL method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	Salicylate Concentration mg/dL [mmol/L]	Bias ⁱ %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	32.9 [2.38]	<10
Bilirubin (unconjugated)	5 mg/dL [86 µmol/L]	33.0 mg/dL [2.39 mmol/L]	<10
Lipemia (Intralipid®)	200 mg/dL [6.78 mmol/L]	32.6 mg/dL [2.36 mmol/L]	<10

i. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the SAL method when present in serum in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 10% at a salicylate concentration of 20 mg/dL [1.45 mmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Ascorbic acid	3 mg/dL	0.17 mmol/L
Caffeine	10 mg/dL	515 µmol/L
Cimetidine	10 mg/dL	0.4 mmol/L
Ethanol	350 mg/dL	76 mmol/L
Ibuprofen	40 mg/dL	1.9 mmol/L
Nicotine	2 mg/dL	0.1 mmol/L
Propoxyphene	0.4 mg/dL	11.8 µmol/L

Analytical Sensitivity: 0.2 mg/dL [0.01 mmol/L]

The analytical sensitivity represents the lowest concentration of salicylate that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Tietz NW. Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Co., 1994:1187-1189.
2. Trinder P. Rapid determination of Salicylate in Biological Fluids, Biochem. J., 1954;57:301-303.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Banfi G, Bauer K, Brand W et al. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO/DIL/LAB/99.1 Rev 2, 2003, pp 1 - 62.
7. Tietz NW. Clinical Guide to Laboratory Tests. Philadelphia: W.B.Saunders Co., 1990:746.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbrikante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare