Dimension Vista® System

Flex® reagent cartridge

SAL

See shaded sections: Updated information from 2017-12 version.

Issue Date 2019-05-30

Salicylate

Intended Use: The SAL method is an *in vitro* diagnostic test for the quantitative measurement of salicylate in human serum on the Dimension Vista® System. Salicylate test results may be used in the diagnosis and treatment of salicylate overdose and for monitoring salicylate levels during therapy.

Summary: Salicylates (aspirin, acetylsalicylic acid) have analgesic, antipyretic, and anti-inflammatory properties and have been used for centuries to relieve pain. Salicylate overdose may cause intoxication. Measurement of salicylate concentration is important for assessment of the severity of intoxication.¹

Principles of Procedure: The salicylate method is a modification of the Trinder colorimetric technique.² Under acidic conditions, ferric nitrate forms a colored complex with salicylate. The amount of salicylate-ferric complex formed is proportional to the salicylate concentration and is measured using a two filter (510, 700 nm) endpoint with sample-reagent blanking technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1–4	Liquid	Surfactant	0.125 mg
5–8	Liquid	Nitric acid	7.15 µmol
9–12	Liquid	Ferric Nitrate	3.3 µmol

- a. Wells are numbered consecutively from the wide end of the cartridge.
b. Nominal value per well in a cartridge.

Risk and Safety:

H290, H319, H317
P280, P305 + P351 + P338, P501

Warning!
May be corrosive to metals. Causes serious eye irritation. May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture w/ 2-methyl-3(2h)-isothiazolone; Nitric acid

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2–8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 7 days for wells 1–12

Specimen Collection and Handling: Recommended specimen type: serum.

Serum can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.³

Follow the instructions provided with your specimen collection device for use and processing.⁴

Complete clot formation should take place before centrifugation. Serum should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁵ Specimens should be free of particulate matter.

Specimens are stable for 7 days at 20–25 °C, 2 weeks at 2–8 °C and 6 months frozen at -20 °C or colder.⁶

Procedure

Materials Provided

SAL Flex® reagent cartridge, Cat. No. K5020

Materials Required But Not Provided

CHEM 2 CAL, Cat. No. KC120

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

Test Conditions

Sample Volume (delivered to the cuvette)	6.4 µL
Reagent 1 Volume	21.4 µL
Reagent 2 Volume	21.4 µL
Reagent 3 Volume	21.4 µL
Temperature	37.0 °C
Reaction Time	0.5 minutes
Wavelength	510 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Calibration Material	CHEM 2 CAL, Cat. No. KC120
Calibration Scheme	2 levels, n=5
Units	mg/dL [mmol/L] ^c (mg/dL x 0.0724) = [mmol/L]
Typical Calibration Levels	Level 1 (Calibrator A): 0 mg/dL [0 mmol/L] Level 2 (Calibrator B): 110.0 mg/dL [8.0 mmol/L]
Calibration Frequency	Every 90 days for any one lot Calibration interval may be extended based on acceptable verification of calibration.
A new calibration is required:	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations

c. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze two levels of a Quality Control (QC) material with known salicylate concentration. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the activity of salicylate in mg/dL [mmol/L] using the calculation scheme described in your Dimension Vista® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 1.7–100.0 mg/dL [0.12–7.24 mmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 100.0 mg/dL [7.24 mmol/L] should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water to obtain results within reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 2.13 µL (dilution factor = 3) for serum. Refer to your Dimension Vista® Operator's Guide.

- Samples with results less than 1.7 mg/dL [0.12 mmol/L] will be reported as "less than 1.7 mg/dL [0.12 mmol/L]" by the instrument.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in salicylate results. Refer to your Dimension Vista® Operator’s Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported.

Interfering Substances

The SAL method was evaluated for interference according to CLSI/NCCLS EP7-A2.⁷ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

- Bilirubin (conjugated) at 10 mg/dL [0.006 mmol/L] depresses SAL results by 47% at 21 mg/dL of SAL.
- Bilirubin (unconjugated) at 10 mg/dL [0.006 mmol/L] depresses SAL results by 20% at 21 mg/dL of SAL.
- Turbidity (as Intralipid®) at 1000 mg/dL [11.3 mmol/L] increases SAL results by 19% at SAL concentration of 20 mg/dL.

‡ At salicylate concentration of 20 and 40 mg/dL [1.45 and 2.90 mmol/L], salicyluric acid at 100 mg/dL increases the SAL result by 55%.

‡ Sodium azide at a concentration of 0.05% increases salicylate results by less than or equal to 41 mg/dL [2.97 mmol/L].

Therapeutic Range: 2.8–20.0 mg/dL [0.20–1.45 mmol/L]⁸

SAL concentrations less than 2.8 mg/dL [0.20 mmol/L] should be considered negative (serum blank). Therapeutic salicylate concentrations vary significantly, depending on the individual. A range of 2.8–20.0 mg/dL [0.20–1.45 mmol/L] includes effective serum concentrations for many patients; however, some individuals are best treated at concentrations outside this range. Concentrations greater than 30 mg/dL [2.17 mmol/L] are often associated with toxic symptoms. Concentrations greater than 60 mg/dL [4.34 mmol/L] are often lethal.

Each laboratory should establish its own expected values for salicylate as performed on the Dimension Vista® System.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n=5 replicates at the following salicylate concentration are:

SAL concentration	Acceptable SD Maximum
6 mg/dL [0.4 mmol/L]	0.4 mg/dL [0.03 mmol/L]
82 mg/dL [5.9 mmol/L]	2.9 mg/dL [0.21 mmol/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Material	Precision ^{9, d}		
	Mean	Standard Deviation (% CV)	
	mg/dL [mmol/L]	Repeatability	Within-Lab
Multiqua® Control			
Level 1	6.2 [0.45]	0.1 [0.01] (2)	0.2 [0.01] (3)
Level 2	47.1 [3.41]	0.4 [0.03] (1)	1.2 [0.08] (3)
Level 3	82.0 [5.94]	0.7 [0.05] (1)	0.9 [0.06] (1)

- d. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.
Multiqua® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Method Comparison ¹⁰				
Comparative Method	Slope	Regression Statistics ^e		n
		Intercept	Correlation Coefficient	
Dimension® RxL System	1.01	mg/dL [mmol/L]	-0.4 [-0.03]	1.000
				47 ^f

- e. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.
f. The range of 47 salicylate values in the correlation study was 1.8–70 mg/dL [0.13–5.1 mmol/L].

Specificity

HIL Interference
The SAL method was evaluated for interference according to CLSI/NCCLS EP7-A2.⁷ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Salicylate mg/dL [mmol/L]	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer) 1000 mg/dL [0.62 mmol/L]	20 [1.45]	<10
Bilirubin (unconjugated)	2 mg/dL [34.2 µmol/L] 5 mg/dL [85.5 µmol/L]	21 [1.52]	-11 -18
Bilirubin conjugated)	2 mg/dL [34.2 µmol/L] 5 mg/dL [85.5 µmol/L]	21 [1.52]	<10 -19
Lipemia (Intralipid®)	800 mg/dL [9.0 mmol/L] 1000 mg/dL [11.3 mmol/L]	20 [1.45]	<10 19

* Analyte results should not be corrected based on this bias.

Non Interfering Substances

The following substances do not interfere with the SAL method when present in serum in the amounts indicated. Inaccuracies (biases) due to these substances are less than 10% at SAL concentration of 20 mg/dL [1.45 mmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	5 ng/mL	6.15 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	3000 mg/dL	33.9 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 1.7 mg/dL [0.12 mmol/L]

The analytical sensitivity represents the lowest concentration of salicylate that can be distinguished from zero. This sensitivity is defined as the mean value (n=20) plus two standard deviations of the level 1 (0 mg/dL [0 mmol/L]) CHEM 2 CAL.

‡ The Dimension Vista® SAL method (REF K5020) and the Dimension® SAL (REF DF20) method utilize the same reagents under equivalent reaction conditions. Interfering substances were tested using the Dimension® SAL (REF DF20) method and the results are representative of both methods.

Symbols Key: See Adjacent Panel.

Bibliography: See Adjacent Panel.

Dimension Vista®, Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

1. Tietz NW. Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Co., 1994:1187-1189.
2. Trinder P. Rapid determination of Salicylate in Biological Fluids, Biochem. J., 1954; 57:301-303.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection—Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Banfi G, Bauer K, Brand W et al. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO/DIL/LAB/99.1 Rev 2, 2003, pp 1-62.
7. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
8. Tietz NW. Clinical Guide to Laboratory Tests. Philadelphia: W.B. Saunders Co., 1990:746.
9. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
10. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFSS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare