Dimension Vista® System

Flex® reagent cartridge

ALTI

See shaded sections: Updated information from 2017-08 version.


Issue Date 2019-06-14

Alanine Aminotransferase

Intended Use: The alanine aminotransferase (ALT) method is an *in vitro* diagnostic test for the quantitative measurement of alanine aminotransferase activity in human serum or plasma on the Dimension Vista® System. Measurements of alanine aminotransferase are used in the diagnosis and treatment of certain liver diseases and heart diseases.

Summary: The Dimension Vista® ALTI method is an adaptation of the recommended alanine aminotransferase procedure of the IFCC as described by Bergmeyer.¹ The procedure is based on the principles outlined by Wroblewski and LaDue² but is modified to contain pyridoxal-5-phosphate (P5P) as an activator and to replace phosphate buffer with tris (hydroxymethyl) aminomethane. Significant elevations of alanine aminotransferase are found in diseases of the liver, such as hepatitis, necrosis, jaundice and cirrhosis. Alanine aminotransferase levels can be elevated even before clinical jaundice appears.³

Principles of Procedure: Alanine aminotransferase catalyzes the transamination of L-alanine to α-ketoglutarate (α-KG), forming L-glutamate and pyruvate. The pyruvate formed is reduced to lactate by lactate dehydrogenase (LDH) with simultaneous oxidation of reduced nicotinamide-adenine dinucleotide (NADH). The change in absorbance is directly proportional to the alanine aminotransferase activity and is measured using a bichromatic (340, 700 nm) rate technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 2	Liquid	Tris buffer	1200 mmol/L	Porcine muscle
3 – 6	Tablet ^c	LDH	7712 U/L	
(Reagent 1)		NADH	0.57 mmol/L	
		P5P	0.39 mmol/L	
7 – 10	Tablet ^c	α-KG	61 mmol/L	
(Reagent 2)				
11 – 12	Liquid	Alanine	1170 mmol/L	

- a. Wells are numbered consecutively from the wide end of the cartridge.
 b. Nominal value in hydrated cartridge.
 c. Tablet contains excipients.

Risk and Safety:

H317
 P280, P272, P302 + P352, P333 + P313, P501
Warning!
 May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 2-Chloracetamide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the Dimension Vista® System.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 10 days for wells 1 – 2, 11 – 12
 5 days for wells 3 – 10

Specimen Collection and Handling:

Recommended specimen types: serum and plasma (lithium heparin).

Serum and plasma specimens can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁴

Follow the instructions provided with your specimen collection device for use and processing.⁵

For Serum: Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁶

Separated samples are stable for 7 days refrigerated at 2 – 8 °C. For longer storage, specimens may be frozen for 1 month at -20 °C or colder.⁷ Avoid repeated freezing and thawing. Thawed frozen specimens which are turbid must be clarified by centrifugation prior to testing.

The purpose of specimen storage information is to provide guidance to the customer; however, customers may validate their own procedures for storing patient samples.

Procedure

Materials Provided

ALTI Flex® reagent cartridge, Cat. No. K2143

Materials Required But Not Provided

ENZ 2 CAL, Cat. No. KC321

Enzyme Diluent, Cat. No. 790035901

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

Test Conditions

Reaction Cuvette	
Sample Volume (delivered to the cuvette)	14.6 µL
Reagent 1 Volume	58.3 µL
Reagent 2 Volume	49.2 µL
Temperature	37.0 °C
Reaction Time	5.6 minutes
Wavelength	340 and 700 nm
Type of Measurement	Bichromatic Rate

Calibration

Calibration Material	ENZ 2 CAL, Cat. No. KC321
Calibration Scheme	2 levels, n=5
Units	U/L [µkat/L] ^d (U/L x 0.0167 = µkat/L)
Typical Calibration Levels	Level 1 (Calibrator A): 0 U/L [0 µkat/L] Level 2 (Calibrator B): 1050 U/L [17.54 µkat/L]
Calibration Frequency	Every 90 days for any one lot Calibration interval may be extended based on acceptable verification of calibration.
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations

d. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known alanine aminotransferase activity. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the activity of alanine aminotransferase in U/L [µkat/L] using the calculation scheme described in your Dimension Vista® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 6 – 1000 U/L [0.10 – 16.70 µkat/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 1000 U/L [16.70 µkat/L] are reported as "Above Assay Range" and should be repeated on dilution.

Autodilution (AD): The autodilution feature options for the Vista® ALTI method are:

Sample Type	Sample volume (µL)	Dilution Factor
Serum/plasma	4.16	3.5
	10	10
	10	20

Refer to your Dimension Vista® Operator's Guide for more information on the autodilution and special dilution features.

Manual Dilution: Dilute with Enzyme Diluent, Cat. No. 790035901, to obtain results within reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

- Samples with results less than 6 U/L [0.10 µkat/L] will be reported as "less than 6 U/L [0.10 µkat/L]" by the instrument.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in alanine aminotransferase results. Refer to your Dimension Vista® Operator's Guide for the meaning of report flags and comments. Any result containing flags and/or comments should be addressed according to your laboratory's procedure manual.

Venipuncture should occur prior to sulfasalazine and/or sulfapyridine administration due to the potential for falsely depressed results.

Interfering substances

Bilirubin (conjugated) at 40 mg/dL [684 µmol/L] decreases ALTI results at an activity of 70 U/L [1.17 µkat/L] by -12%.

Bilirubin (conjugated) at 60 mg/dL [1026 µmol/L] decreases ALTI results at an activity of 144 U/L [2.40 µkat/L] by -13%.

Triglycerides above 400 mg/dL [4.52 mmol/L] tripped a test report message; therefore the magnitude of the interference could not be determined.

Lipemia (Intralipid®) of 600 mg/dL [6.78 mmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values^a

Females: 13 – 56 U/L [0.22 – 0.93 µkat/L], n= 144

Males: 16 – 61 U/L [0.27 – 1.02 µkat/L], n= 125

The reference intervals were calculated non-parametrically and represent the central 95% of results determined from a population of healthy adults with normal total bilirubin, negative for Hepatitis C, and liver enzymes <1.5x normal limits.

Each laboratory should establish its own expected values for ALTI as performed on the Dimension Vista® System.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n=5 replicates at the following ALTI activity are:

ALTI Activity	Acceptable SD Maximum
24 U/L [0.40 µkat/L]	2.3 U/L [0.04 µkat/L]
100 U/L [1.67 µkat/L]	4.8 U/L [0.08 µkat/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Precision^{a, e}

Material	Mean U/L [µkat/L]	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Serum Pool 1	38 [0.64]	0.7 [0.01] (1.8)	1.9 [0.03] (4.9)
Serum Pool 2	143 [2.39]	1.4 [0.02] (1.0)	4.1 [0.07] (2.8)
Plasma Pool 1	139 [2.32]	1.1 [0.02] (0.8)	4.8 [0.08] (3.5)
Bio-Rad Multiquel® Control			
Level 1	24 [0.40]	0.6 [0.01] (2.6)	1.1 [0.02] (4.8)
Level 2	77 [1.28]	1.0 [0.02] (1.4)	1.4 [0.02] (1.9)
Level 3	175 [2.93]	1.7 [0.03] (1.0)	3.4 [0.06] (1.9)
MAS® chemTRAK® H Control			
Level 1	24 [0.40]	0.5 [0.01] (2.3)	0.7 [0.01] (2.9)
Level 2	100 [1.67]	1.1 [0.02] (1.1)	2.9 [0.05] (2.9)
Level 3	205 [3.43]	1.8 [0.03] (0.9)	2.1 [0.03] (1.0)

e. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Multiquel® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

MAS® and chemTRAK® are registered trademarks of Medical Analysis Systems Inc., Camarillo, CA 93012-8058, USA.

Method Comparison¹⁰
Regression Statistics^f

Comparative Method	Slope	Intercept U/L [µkat/L]	Correlation Coefficient	n
Siemens ADVIA® 1650 ALT	1.02	1.20 [0.02]	0.9996	118 ^g
Dimension® ALTI	1.00	-1.26 [-0.02]	0.9996	118 ^g

f. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.

g. The range of 118 values in the correlation study was 9 – 867 U/L [0.15 – 14.47 µkat/L].

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The ALTI method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	ALTI Activity U/L [µkat/L]	Bias* %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L]	46 [0.77]	<10
	1000 mg/dL [0.62 mmol/L]	132 [2.20]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	68 [1.14]	<10
	80 mg/dL [1368 µmol/L]	138 [2.30]	<10
Bilirubin (conjugated)	30 mg/dL [513 µmol/L]	70 [1.17]	<10
	40 mg/dL [684 µmol/L]	144 [2.40]	<10
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	47 [0.78]	<10
	200 mg/dL [2.26 mmol/L]	135 [2.25]	<10
	600 mg/dL [6.78 mmol/L]	47 [0.78]	--- ^h
	600 mg/dL [6.78 mmol/L]	135 [2.25]	--- ^h

* Analyte results should not be corrected based on this bias.

h. The interference testing at this level and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Non-Interfering Substances

The following substances do not interfere with the ALTI method when present in serum and plasma at the activity indicated. Inaccuracies (biases) due to these substances are less than 10% at ALTI activity of 55 U/L [0.92 µkat/L] and 165 U/L [2.8 µkat/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Amikacin	8 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1.0 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	50 mg/dL	3.62 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	400 mg/dL	4.52 mmol/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vancomycin	10 mg/dL	69 µmol/L

Limit of Detection and Limit of Blank^{12, i}

The Limit of Detection (LoD) for ALTI is 6 U/L [0.10 µkat/L], determined consistent with CLSI guideline EP17-A¹² and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 240 determinations, with 5 blank and 5 low level samples. The Limit of Blank (LoB) is 2.3 U/L [0.04 µkat/L].

i. LoD is the lowest concentration of analyte that can be detected reliably. LoB is the highest concentration that is likely to be observed for a blank sample.

Symbols Key: See Adjacent Panel

Bibliography: See Adjacent Panel

Dimension®, Vista®, ADVIA®, and Flex® are trademarks of Siemens Healthcare Diagnostics.

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©2011 Siemens Healthcare Diagnostics

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Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

1. Bergmeyer HU, Scheibe P, Wahlefeld AW. Optimization of methods for aspartate aminotransferase and alanine aminotransferase, Clin Chem, 1978; 24/1:58-73.
2. Wroblewski F, LaDue JS. Serum glutamic pyruvic transaminase in cardiac and hepatic disease, Proc Soc Exp Biol Med 1956; 91:569.
3. Tietz NW, Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Burtis CA, Ashwood ER, Bruns DE, Elsevier Saunders, 2006; pp 604-607.
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
7. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition, AACC Press, Washington, DC 1997; pp.3-10 to 3-12.
8. Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline - Third Edition*. CLSI document EP28-A3c [ISBN 1-56238-682-4]. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2008.
9. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative-Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
10. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
11. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
12. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EHGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

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All rights reserved.