Dimension Vista® System

Flex® reagent cartridge

UCFP

See shaded sections: Updated information from 2015-05 version.


Issue Date 2019-05-30

Urinary/Cerebrospinal Fluid Protein

Intended Use: The UCFP method is an *in vitro* diagnostic test for the quantitative measurement of protein in urine and cerebrospinal fluid on the Dimension Vista® System.

Summary: The UCFP method allows direct quantitation of proteins in urine and cerebrospinal fluid specimens. The UCFP method is an adaptation of pyrogallol red-molybdate method by Y. Fujita, I. Mori and S. Kitano.¹ Measurement of the protein content in urine is used in diagnosis and treatment of kidney diseases. Measurement of the protein content in cerebrospinal fluid is used in the diagnosis and treatment of central nervous system diseases.

Principles of Procedure: In the reaction sequence, pyrogallol red combines with sodium molybdate to form a red complex with maximum absorbance at 470 nm. The protein in the sample reacts with this complex in acid solution to form a bluish-purple colored complex, which absorbs at 600 nm. The absorbance at 600 nm is directly proportional to the concentration of protein in the sample. The analyte concentration is determined by calculation using a logit curve fit on a previously stored calibration curve.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 6	Liquid	Pyrogallol red in methanol	0.20 mM
		Sodium molybdate	0.35 mM
		Stabilizers and Surfactants	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

Risk and Safety:

H302, H371
P270, P264, P309 + P311, P301 + P312, P501

Warning!

Harmful if swallowed. May cause damage to organs.

Do not eat, drink or smoke when using this product. Wash hands thoroughly after handling. IF exposed or if you feel unwell: Call a POISON CENTER or doctor/physician. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Methanol

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 7 days for wells 1 – 6

Specimen Collection and Handling: Recommended specimen types: urine and cerebrospinal fluid. Normal procedures for collecting urine and cerebrospinal fluid may be used for samples to be analyzed by this method.

Random urine specimens may be used but timed 24-hr specimens are preferred. No preservative is required during 24-hr collection, but thereafter urine aliquots should be stored at 2 – 4 °C or frozen. Urine specimens must be free of any particulate matter before analysis.²

Specimens stored at 4 °C with no additives are stable for at least three days.

Cerebrospinal fluid specimens should be collected with care to avoid contamination with plasma proteins. Blood present in the cerebrospinal fluid invalidates the protein values since it reflects contamination with plasma proteins. Remaining cerebrospinal fluid specimens should be refrigerated or frozen.²

Procedure

UCFP Flex® reagent cartridge, Cat. No. K3026

Materials Provided

Materials Required But Not Provided

UCFP CAL, Cat. No. KC260
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

Test Conditions

Sample Volume (delivered to the cuvette)	3.68 µL
Reagent Volume	128.68 µL
Temperature	37.0 °C
Reaction time	2 minutes
Wavelength	600 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Calibration Material	UCFP CAL, Cat. No. KC260
Calibration Scheme	5 levels, n=5
Units	mg/dL [mg/L] ^c (mg/dL x 10) = [mg/L]
Typical Calibration Levels	Level 1 (System Water): 0 mg/dL [0.0 mg/L] Intermediate level/levels automatically prepared by the instrument: 25, 50, 125 mg/dL [250, 500, 1250 mg/L] Level 5 (Calibrator A): 250 mg/dL [2500 mg/L]
Calibration Frequency	Every 60 days for any one lot Calibration interval may be extended based on acceptable verification of calibration.
A new calibration is required:	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations

c. Système International d'Unités [SI units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known UCFP concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of Urinary/Cerebrospinal Fluid Protein in mg/dL [mg/L] using the calculation scheme described in your Dimension Vista® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5–250 mg/dL [50–2500 mg/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 250 mg/dL [2500 mg/L] should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water to obtain results within reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 10.0 µL (dilution factor = 10) for CSF/urine. Refer to your Dimension Vista® Operator's Guide.

- Samples with results less than 5 mg/dL [50 mg/L] will be reported as "less than 5 mg/dL [50 mg/L]" by the instrument.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in UCFP results. Refer to your Dimension Vista® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

Interfering Substances

Hemolyzed samples should be avoided since hemolysis increases UCFP results. Hemoglobin at 10 mg/dL [0.006 mmol/L] increases UCFP result by 22% at UCFP level of 48 mg/dL [480 mg/L].

‡ Samples containing amikacin, gentamicin, kanamycin, and tobramycin should be avoided since these substances falsely increase UCFP results.

‡ Neomycin sulfate at 15 µg/mL increases UCFP results by 11% and at 7.5 µg/mL the interference is less than 5%.

Expected Values: Urine ‡ < 11.9 mg/dL < 149.1 mg/day
Spinal Fluid² 15 – 45 mg/dL [150 – 450 mg/L]

This urine reference population consisted of 115 adults (50 females, 65 males). The reference interval was calculated non-parametrically and represents 95% of the population.

Each laboratory should establish its own expected values for urine and cerebrospinal fluid protein as performed on the Dimension Vista® System.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n=5 replicates at the following urine and cerebrospinal fluid protein concentrations are:

UCFP concentration	Acceptable SD Maximum
30 mg/dL [300 mg/L]	2.4 mg/dL [24 mg/L]
100 mg/dL [1000 mg/L]	10.9 mg/dL [109 mg/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Precision ^{3,d}			
Material	Mean mg/dL [mg/L]	Standard Deviation (% CV)	
		Repeatability	Within-Lab
Liquichek™ Urine Chemistry Control			
Level 1	27.8 [278]	0.6 [6] (2)	1.1 [11] (4)
Level 2	74.6 [746]	1.0 [10] (1)	1.4 [14] (2)
Level 3	106.1 [1061]	2.6 [26] (2)	2.7 [27] (3)
Liquichek™ Spinal Fluid Control			
Level 1	38.4 [384]	0.8 [8] (2)	1.1 [11] (3)
Level 2	88.3 [883]	1.2 [12] (1)	1.8 [18] (2)

- d. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.
Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Method Comparison⁴				
Regression Statistics⁵				
Comparative Method	Slope	Intercept mg/dL [mg/L]	Correlation Coefficient	n⁶
Dimension® Rxl System				
Urine	0.92	2.2 [22]	0.999	55
CSF	1.08	-0.8 [-8]	0.999	55

- e. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.

- f. The range of 55 values in the correlation study for urine was 6.3–218.8 mg/dL [63–2188 mg/L] and the range of 55 values for CSF in the correlation study was 7.7–240.6 mg/dL [77–2406 mg/L].

Specificity

Non-Interfering Substances

The UCFP method was evaluated for interference according to CLSI/NCCLS EP7-A2.⁵ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent.

The following substances do not interfere with the UCFP method when present in urine at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at Urinary/Cerebrospinal Fluid Protein concentration of 48 mg/dL [480 mg/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	5 ng/mL	6.15 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Recovery ‡

Recovery of various proteins by the UCFP method is shown below.

Protein	% Recovery
Albumin	89
Gamma Globulin	57
Transferrin	81
Orosomucoid	44
β-Lipoprotein	63
Lysosome	5

Analytical Sensitivity: 5 mg/dL [50 mg/L]

The analytical sensitivity represents the lowest concentration of Urinary/Cerebrospinal Fluid Protein that can be distinguished from zero. This sensitivity is defined as the mean value (n=20) plus two standard deviations of level 1 (0 mg/dL [0 mg/L]), which is system water.

‡ The Dimension Vista® UCFP method (REF K3026) and the Dimension® UCFP (REF DF26) method utilize the same reagents under equivalent reaction conditions. Interfering substances, Expected Values, and Recovery were tested using Dimension® UCFP(REF DF26) and the results are representative of both methods.

Symbols Key: See Adjacent Panel.

Bibliography: See Adjacent Panel.

Dimension Vista®, Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

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1. Fujita Y, Mori I, Kitano S. Color reaction between Pyrogallol Red—Molybdenum(VI) complex and protein. *Bunseki Kagaku* 1983;32:379-386.
2. Burtis CA, Ashwood ER. *Tietz Textbook of Clinical Chemistry*, Philadelphia: WB Saunders Co., 1994:54-55 (techniques and procedures to minimize infections in laboratory workers), 58-69 (specimen collection), 75-76 (separation and storage of specimens), 717-730 (proteins in urine and cerebrospinal fluid).
3. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087, 2004.
4. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087, 2002.
5. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087, 2005.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos	
	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_FIGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare