SIEMENS

Syva[®]

Emit[®] 🔟 LSD Assay

See shaded sections: Updated information from 2017-05 version.



10871618_G



LSD Assay

1 Intended Use

The Emit® II LSD Assay is a homogeneous enzyme immunoassay with a 0.5 ng/mL cutoff. The assay is intended for use in the qualitative and semiquantitative analyses of lysergic acid diethylamide (LSD) in human urine.

The Emit® II LSD Assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS)* is the preferred confirmatory method.¹ Other chemical confirmation methods are available. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.

* Throughout this document, GC/MS refers to GC/MS or GC/MS/MS or both.

2 Summary and Explanation of the Test

Lysergic acid diethylamide (LSD, lysergamide, lysergide, Delysid) is one of the most potent hallucinogenic agents known. LSD produces vivid hallucinations in all of the sensory modalities. Patients experiencing a "bad trip" may exhibit hysterical behavior, hyperactivity, life-threatening hyperthermia, restlessness, anxiety, numbness, vomiting, collapse, coma, catatonia, fever, diarrhea, and blood coagulation disorders. Delayed psychotic reactions can occur weeks after ingestion. Acute symptoms generally last for 6 to 12 hours, gradually declining in intensity.^{2,3}

A common side effect of LSD use is the flashback, a recurrence of hallucinations weeks, months, or even years after the last dose. Fatalities associated with LSD are almost always due to injuries received while under the influence of the drug. Death as a direct result of LSD toxicity is extremely rare, but has been reported.⁴ LSD is taken orally, usually by swallowing drug-impregnated blotter paper. In the 1960s, a typical dose of LSD ranged from 105–300 micrograms, but in recent years, the dosage amount has decreased to 20–80 micrograms.⁵ LSD and its metabolites can be detected in urine for several days following ingestion.³

The Emit® II LSD Assay detects lysergic acid diethylamide in human urine. Positive results for samples containing other compounds structurally unrelated to LSD have not been observed. The cutoff level for distinguishing positive from negative samples is 0.5 ng/mL.

Methods historically used for detecting LSD in biological fluids include radioimmunoassay, GC/MS, thin-layer chromatography, and fluorimetry.

While confirmation techniques other than GC/MS may be adequate for some drugs of abuse, GC/MS is generally accepted as a vigorous confirmation technique for all drugs, since it provides the best level of confidence in the result.

3 Principle

The Emit® II LSD Assay is a homogeneous enzyme immunoassay technique used for the analysis of specific compounds in human urine.⁶ The assay is based on competition between drug in the specimen and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for antibody binding sites. Enzyme activity decreases upon binding to the antibody, so the drug concentration in the specimen can be measured in terms of enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that is measured spectrophotometrically. Endogenous serum G6PDH does not interfere because the coenzyme NAD functions only with the bacterial (*Leuconostoc mesenteroides*) enzyme employed in the assay.

4 Reagents

REF	Product Description	Volume
9L039UL	Emit® II LSD Assay Antibody/Substrate Reagent 1 Mouse monoclonal antibodies reactive to lysergic acid diethylamide (0.15 µg/mL),* G6P (22 mM), NAD (18 mM), preservatives, and stabilizers	26 mL
	Enzyme Reagent 2 lysergic acid diethylamide labeled with bacterial G6PDH (0.18 U/mL),* preservatives, and stabilizers	13 mL

* The antibody titer and enzyme conjugate activity may vary from lot to lot.

Note: Reagents 1 and 2 are sold as a matched set. They should not be interchanged with components of kits with different lot numbers.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions

- · Reagents 1 and 2 contain nonsterile mouse monoclonal antibodies.
- · This kit contains streptomycin sulfate. Please dispose of appropriately.
- Contains sodium azide (<0.1%) as a preservative. Sodium azide can react with copper or lead
 pipes in drain lines to form explosive compounds. Dispose of properly in accordance with
 local regulations.

For in vitro diagnostic use.

Preparation and Storage of Assay Components

Reagents Preparation

The Emit \circledast II LSD Assay reagents are provided in liquid form, ready to use and may be used directly from the refrigerator. Close the reagent bottles when not in use.

Note: Caps must always be replaced on the original containers.

Storage

When not in use, reagents must be stored at 2–8°C (36–46°F), upright and with screw caps tightly closed. Do not freeze reagents. Avoid prolonged exposure to temperatures above 32°C. Improper storage of reagents can affect assay performance.

Stability

If stored as directed, reagents are stable until the expiration date printed on the label. Refer to the application sheet for on instrument stability information.

5 Specimen Collection and Preparation

- Urine specimens may be collected in plastic (i.e., polypropylene, polycarbonate, polyethylene) or glass containers. Some plastics can adsorb certain drugs.
- If not analyzed immediately, specimens may be stored unrefrigerated for up to 7 days following collection. After 7 days, specimens should be stored frozen.
- LSD degrades with prolonged exposure to light. Store specimens in opaque containers or wrap them in opaque paper or aluminum foil and store frozen.⁷
- Frozen specimens must be thawed and mixed thoroughly prior to analysis.
- Specimens should be at 18–25°C (64–77°F) for testing.
- · Specimens with high turbidity should be centrifuged before analysis.
- Urine specimens within the pH range of 4.5-8.0 do not require prior adjustment of pH.
- · When samples are outside the normal pH range, suspect adulteration.
- Adulteration of the urine specimen may cause erroneous results. If adulteration is suspected, obtain another specimen.
- · Human urine specimens should be handled and treated as if they are potentially infectious.

6 Procedure

Materials Provided

Emit® II LSD Assav

Antibody/Substrate Reagent 1 Enzyme Reagent 2

Materials Required But Not Provided

9L009UL	Emit® LSD 0 ng/mL Calibrator
9L109UL	Emit [®] LSD 0.5 ng/mL Calibrator
9L209UL	Emit [®] LSD 1.5 ng/mL Calibrator
9L309UL	Emit [®] LSD 2.5 ng/mL Calibrator
9L419UL	Emit [®] LSD Control Level 1 (0.25 ng/mL)
9L519UL	Emit [®] LSD Control Level 2 (1 ng/mL)

Instruments

Siemens Healthcare Diagnostics provides instructions for using this assay on a number of chemistry analyzers. Contact the Technical Assistance Center in the USA or your local Siemens representative for application sheets.

Analyzers must be capable of maintaining a constant reaction temperature, pipetting specimens/ reagents and measuring enzyme rates precisely, timing the reaction accurately, and mixing reagents thoroughly.

Assay Sequence

To run the assay, see the instrument operator's manual and the application sheets available from Siemens.

Calibration

Note: These reagents are qualified for use with these calibrators only. However, other control material may be used for quality control purposes.

Table 1 —	Emit® LSD	Calibrators fo	r use in	Qualitative or	Semiquantitative	Analysis
-----------	-----------	----------------	----------	----------------	------------------	----------

Desired Cutoff Level (ng/mL)	Additional Recommended Calibrators for Qualitative Analysis to be used as Quality Control (ng/mL)	Required Calibrators for Semiquantitative Analysis (ng/mL)
0.5	0 1.5 2.5	0 0.5 1.5 2.5

Note: For any individual cutoff level, a calibrator is used as either a calibrator or as a control when the assay is used for qualitative analysis. When a calibrator is used as a calibrator for an individual cutoff level, the other level calibrators (above or below it, as listed above) are used as controls.

Qualitative Analysis

Calibrate by running the appropriate Emit® LSD Calibrator Level for the desired cutoff listed in Table 1 in duplicate. Validate the calibration by running controls (see Quality Control). Refer to the Emit® LSD Calibrators instructions for use and the application sheet for additional information and instrument settings. Recalibrate as indicated by control results.

Semiquantitative Analysis

Prepare a calibration curve by running the appropriate Emit® LSD Calibrators listed in Table 1. Validate the calibration by running controls (see Quality Control). Refer to the Emit® LSD Calibrators instructions for use and the application sheet for additional information and instrument settings. Recalibrate as indicated by control results.

Quality Control

Qualitative Analysis

Validate the calibration by assaying controls. Ensure that the result from Emit® LSD Calibrator 0 ng/mL or Emit® LSD Calibrator 1.5 ng/mL or 2.5 ng/mL relates appropriately to the result from the cutoff calibrator chosen from column 1 in Table 1. That is,

- If Emit® LSD Calibrator 0 ng/mL was run, ensure that the result is negative relative to the selected cutoff calibrator level.
- If Emit® LSD Calibrator 1.5 ng/mL or 2.5 ng/mL was run, ensure that the result is positive relative to the selected cutoff calibrator level.

If appropriate, run low and high controls (e.g., the Emit® LSD Control Level 1 and Emit® LSD Control Level 2) with every run. Ensure that the low control tests negative and the high control tests positive relative to the 0.5 ng/mL calibrator (cutoff). Once the calibration is validated, run patient samples.

Semiguantitative Analysis

Validate the calibration curve by assaying commercial controls. Ensure that control results fall within acceptable limits as defined by your laboratory. (See Table 4, Qualitative Analysis Precision, for typical rate separations for the SYVA \circledast -30R Biochemical System).

Once the calibration curve is validated, run urine specimens.

Qualitative and Semiquantitative Analysis

- Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of Quality Control (QC) material with known LSD concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.
- 2. Refer to the instrument operator's manual for appropriate instrument checks.

7 Results

Qualitative Analysis

Refer to Table 1 for the appropriate cutoff Emit® Calibrator. Table 1 contains the concentration of lysergic acid diethylamide present in the selected Emit® Calibrator selected as a cutoff for distinguishing "positive" from "negative" specimens.

Positive Results. A specimen that gives a change in rate value greater than or equal to the Emit® Calibrator cutoff rate value is interpreted as positive. The sample contains LSD.

Negative Results. A specimen that gives a change in rate value less than the Emit® Calibrator cutoff rate value is interpreted as negative: Either the specimen does not contain LSD or LSD is present in concentrations below the cutoff level for this assay.

Semiquantitative Analysis

The semiquantitation of positive results enables the laboratory to determine an appropriate dilution of the specimen for confirmation by GC/MS. Semiquantitation also permits the laboratory to establish quality control procedures and assess control performance. Refer to the Analytical Recovery section for the semiquantitative range.

Using the Emit® II LSD Assay, it is possible to make semiquantitative determinations of LSD. An estimate of relative total drug concentrations may be obtained by running the appropriate Emit® LSD Calibrators: 0 ng/mL, 0.5 ng/mL cutoff, 1.5 ng/mL, 2.5 ng/mL. Refer to the application sheet for instructions.

Immunoassays that produce a single result in the presence of multiple detectable components cannot fully quantitate the concentration of individual components. Interpretation of results must take into account that urine concentrations can vary extensively with fluid intake and other biological variables. A more specific alternative chemical method must be used to obtain a confirmed analytical result (see Section 1, Intended Use).

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Siemens has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. User defined modifications are not supported by Siemens as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Siemens Application Sheets or these instructions for use.

8 Limitations

- · The assay is designed for use with human urine only.
- A positive result from the assay indicates the presence of LSD but does not indicate or measure intoxication.
- · Boric acid is not recommended as a preservative for urine.
- There is a possibility that substances and/or factors not listed (e.g., technical or procedural
 errors) may interfere with the test and cause false results.
- Samples from patients taking chlorpromazine (Thorazine[®]) may produce positive results with this assay.
- Interpretation of results must take into account that urine concentrations can vary extensively
 with fluid intake and other biological variables.
- Immunoassays that produce a single result in the presence of a drug and its metabolites cannot fully quantitate the concentration of individual components.

9 Expected Values

When the Emit® II LSD Assay is used as a qualitative assay, the amount of drugs and metabolites detected by the assay in any given specimen cannot be estimated. The assay results distinguish between positive and negative specimens; positive indicating specimens that contain LSD.

When used semiquantitatively, the assay yields approximate, cumulative concentrations of the metabolites detected by the assay (see Section 7, Results).

10 Specific Performance Characteristics

The data appearing in this section were collected on the SYVA \otimes -30R Biochemical System using the Emit \otimes II LSD Assay and radioimmunoassay (RIA) (Reference Method).

Method Comparison

Qualitative Results

118 samples were analyzed by the ${\sf Emit} \circledast {\sf II}$ LSD Assay and by RIA (Reference Method).

Fifty-Nine (59) samples showed positive results by both methods, while 57 samples showed negative results by both methods. Two samples showed positive by the EmitB II LSD Assay and negative by RIA. Data are summarized in Table 2.

Table 2 — Qualitative Results for the 0.5 ng/mL Cutoff



Semiquantitative Results

Ten (10) of the 57 samples that showed negative results by the Emit® II LSD Assay and RIA qualitatively were also tested semiquantitatively by these methods. All 10 samples showed negative results by both methods.

The 2 samples that showed positive results by the Emit® method and negative results by the RIA method qualitatively were tested semiquantitatively by both methods. Both samples showed positive results by the Emit® method; however, 1 sample showed a negative result and 1 sample showed a positive result by the RIA method. (Both samples produced results very near the cutoff level of each assay.)

Fifty-seven (57) of the 59 samples that showed positive results by the Emit® method and RIA qualitatively were tested semiquantitatively by both methods. Two (2) samples were not included in this comparison because Emit® or RIA semiquantitative results or GC/MS quantitative results could not be obtained. All 57 samples showed positive results by both methods; their concentrations were determined using GC/MS.

Table 3 shows the distribution of the samples according to where their results lie within one of three incremental ranges. Also shown is the sample concentration range, determined by GC/MS, for each incremental range.

Table 3 — Emit \circledast II LSD Assay Comparative Analysis Semiquantitative Results

	Emit®		Emit® RIA	
Range (ng/mL)	Distribution of Samples	GC/MS Concentration Range (ng/mL)	Distribution of Samples	GC/MS Concentration Range (ng/mL)
<0.5	10	N/A	11	0.04*
0.5–1	7	0.02-0.28	11	0.02-1.09
>1	52	0.03-16.86	47	0.18–16.86

* Concentration for the sample >0.5 ng/mL by the Emit ® method and <0.5 ng/mL by RIA.

The Emit $^{\odot}$ II LSD Assay (range 0.25 to 2.5 ng/mL) and RIA (range 0.25 to 1 ng/mL) showed good overall semiquantitative agreement.

Ten (10) samples that showed negative results by the Emit® method also showed negative results by RIA. Seven (7) samples tested by the Emit® method had concentrations in the 0.5 ng/mL to 1 ng/mL range. Of the 7 samples, 6 also had concentrations in the 0.5 ng/mL to 1 ng/mL range when tested by RIA. The remaining 1 of the 7 samples was borderline negative by RIA and borderline positive by the Emit® method. Five (5) of the 11 samples that had concentrations in the 0.5 ng/mL to 1 ng/mL range when tested by RIA vere >1 ng/mL when tested by the Emit® method. The differences in the sample concentrations between the two methods ranged from 0.03 ng/mL to 0.61 ng/mL.

Fifty-two (52) of the positive samples had concentrations greater than 1 ng/mL by the Emit® method; 47 of the 52 samples had concentrations greater than 1 ng/mL by RIA.

Concentrations generated by the Emit® method and RIA showed a similar overall relationship to the corresponding concentrations generated by GC/MS.

Precision

Qualitative precision was determined by assaying 6 samples on 20 days, 2 runs per day in replicates of 6 (N = 240). Semiquantitative precision was determined by assaying 4 samples on 20 days, 2 runs per day in replicates of 4 (N = 160). Precision data were calculated according to the National Committee of Clinical Laboratory Standards (NCCLS) Guideline EP5-T2 (March 1992). Results are summarized in Tables 4 and 5.

Table 4 — Precision: Qualitative Analysis

Calibrator or Control	Mean (mAU/min)	SD	CV (%)
Within-Run			
0 ng/mL Calibrator	151	0.83	0.55
Control Level 1 (0.25 ng/mL)	158	0.89	0.56
0.5 ng/mL Calibrator	164	0.85	0.52
Control Level 2 (1 ng/mL)	176	0.81	0.46
1.5 ng/mL Calibrator	184	0.85	0.46
2.5 ng/mL Calibrator	189	0.92	0.49
Total			
0 ng/mL Calibrator	151	0.94	0.62
Control Level 1 (0.25 ng/mL)	158	0.92	0.59
0.5 ng/mL Calibrator	164	0.94	0.57
Control Level 2 (1 ng/mL)	176	0.88	0.50
1.5 ng/mL Calibrator	184	0.98	0.53
2.5 ng/mL Calibrator	189	1.03	0.54

Table 5 — Precision: Semiquantitative Analysis

Calibrator or Control	Mean (ng/mL) SD		CV (%)	
Within-Run				
Control Level 1 (0.25 ng/mL)	0.25	0.03	13.5	
0.5 ng/mL Calibrator	0.5	0.04	7	
Control Level 2 (1 ng/mL)	1.1	0.03	3.2	
1.5 ng/mL Calibrator	1.5	0.10	6.9	
Total				
Control Level 1 (0.25 ng/mL)	0.25	0.04	15	
0.5 ng/mL Calibrator	0.5	0.04	8.8	
Control Level 2 (1 ng/mL)	1.1	0.04	3.6	
1.5 ng/mL Calibrator	1.5	0.13	8.7	

Analytical Recovery

Qualitative Results

In qualitative spike analysis, the Emit® II LSD Assay correctly identified the mean rate of spiked specimens containing less than the cutoff listed in Table 1 as negative and the mean rate of spiked specimens containing greater than the cutoff listed in Table 1 as positive 84% of the time.

Semiquantitative Results

Negative human urine was spiked with concentrations of LSD at levels throughout the semiquantitative range of 0.25 to 2.5 ng/mL. For each known concentration, drug recovery was calculated using the average concentration obtained by the Emit® II LSD Assay. All spikes were analyzed in triplicate. Semiquantitative results are shown in Table 6.

Table 6 — Emit® II LSD Assay Spike Recovery Semiquantitative Results

Expected LSD Concentration (ng/mL)	Average LSD Concentration by Emit® II LSD Assay (ng/mL)	Recovery (%)
0.28	0.27	96
0.37	0.37	100
0.45	0.42	93
0.58	0.49	84
0.59	0.57	97
0.75	0.76	101
1.03	1.08	105
1.90	2.03	107

Specificity

The Emit® II LSD Assay detects LSD in human urine.

Table 7 lists the concentrations of compounds that produce a result that is approximately equivalent to the 0.5 ng/mL calibrator cutoff.

These concentrations are within the range of levels found in urine following use of the compound or, in the case of metabolites, the parent compound. If a specimen contains more than one compound detected by the assay, lower concentrations than those listed in Table 7 may combine to produce a rate equal to or greater than that of the cutoff calibrator. Data presented are representative of typical performance of this assay.

Table 7 — Concentrations that Produce a Result Approximately Equivalent to the 0.5 ng/mL LSD Cutoff

Compound	Concentration (µg/mL) at the 0.5 ng/mL Cutoff	
D-Amphetamine	500	
M-Chlorophenylpiperazine	33	
Ergonovine	1	
Methadone	400	
D-Methamphetamine	100	
Methysergide	3	
Phencyclidine	30	

Table 8 lists the compounds that produce a negative result by the Emit® II LSD Assay. Specificity testing was performed at the 0.5 ng/mL cutoff, which represents the greatest potential for crossreactivity. Positive results for compounds structurally unrelated to LSD have not been observed.

Table 8 — Concentrations of Compounds Showing a Negative Response

Compound	Concentration Tested (µg/mL) at the 0.5 ng/mL Cutoff	
Benzoylecgonine	1000	
α -Ergocryptine	20	
Ergotamine	100	
Lysergic Acid	100	
Methaqualone	1000	
Morphine	1000	
Oxazepam	250	
Propoxyphene	1000	
Psilocin	100	
Psilocybin	100	
Secobarbital	1000	
Serotonin	1000	
11-nor-∆9-THC-9-COOH	150	
L-Tryptophan	100	

Analytical Sensitivity

The sensitivity level of the Emit® II LSD Assay is 0.13 ng/mL. This level represents the lowest concentration of LSD that can be distinguished from 0 ng/mL with a confidence level of 95%.

11 Bibliography

- 1. Hawks RL, Chiang CN, eds. Urine Testing for Drugs of Abuse. Rockville, MD: National Institute on Drug Abuse (NIDA), Dept of Health and Human Services; 1986. NIDA research monograph-73.
- 2. Abraham HD, Aldridge AM. Adverse consequences of lysergic acid diethylamide. Addiction. 1993:88:1327-1334.
- 3. Baselt RC, Cravey RH. Disposition of Toxic Drugs and Chemicals in Man. 3rd ed. Chicago, IL: Year Book Medical Publishers, Inc; 1989:470–473.
- 4. Fysh RR, Oon MC, Robinson KN, Smith RN, White PC, Whitehouse MJ. A fatal poisoning with LSD. Forensic Sci Int. 1985;28:109–113.
- 5. It never went away: LSD, a sixties drug, attracts young users in the nineties. Washington, DC: Department of Justice, Drug Enforcement Administration, Demand Reduction Section; 1991.
- 6. Henry JB. Clinical Diagnosis and Management by Laboratory Methods. 18th ed. Philadelphia, PA: WB Saunders Co; 1991:866.
- 7. Schwartz RH. LSD. Its rise, fall, and renewed popularity among high school students. Pediatr Clin North Am. 1995;42:403-413.

12 Symbols Key

CON

2	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
\triangle	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
666	Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
∇	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos</n></n></n>
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Límite de temperatura
∏ i]	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
NON	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
CE	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
ONTENTS	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel
	2015-03 EFIGS

For technical assistance, call Siemens Healthcare Diagnostics: 1-800-227-8994 in the USA 1-800-264-0083 in Canada

Outside the USA and Canada, call your local Siemens representative.

Syva[®], Syva[®], and Emit[®] are trademarks of Siemens Healthcare Diagnostics.

Thorazine® is a registered trademark of SmithKline Beecham.

© 2011 Siemens Healthcare Diagnostics All rights reserved.

Siemens Healthcare Diagnostics Inc. 500 GBC Drive Newark, DE 19714 USA

Global Siemens Headquarters Siemens AG Wittelsbacherplatz 2 80333 Muenchen Germany

Global Siemens Healthcare Headquarters Siemens AG Siemens AG Healthcare Sector Henkestrasse 127 91052 Erlangen Germany Phone: +49 9131 84-0 siemens.com/healthcare

Global Division Siemens Healthcare Diagnostics Inc. 511 Benedict Avenue Tarrytown, NY 10591 USA siemens.com/healthcare

Printed in USA 2019-05 10871618 US G