Dimension[®] clinical chemistry system

Heterogeneous Immunoassay Module

Flex[®] reagent cartridge

TPSA

See shaded sections: Updated information from 2015-04 version.

Issue Date 2019-07-29

Total Prostate Specific Antigen

CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

Warning: The concentration of TPSA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the PSA assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining PSA levels serially is changed, additional sequential testing should be carried out. Prior to changing assays, the laboratory MUST confirm baseline values for patients being serially monitored.

Intended Use: The TPSA method for the Dimension[®] clinical chemistry system with the heterogeneous immunoassay module is an *in vitro* diagnostic test intended to quantitatively measure total prostate specific antigen (PSA) in human serum and plasma:

- As an aid in the detection of prostate cancer when used in conjunction with digital rectal exam (DRE) in men 50 years or older. Prostate biopsy is required for diagnosis of cancer.
- As an aid in the management (monitoring) of prostate cancer patients.

Summary: Prostate specific antigen (PSA) is a serine protease of approximately 30000 Daltons produced by the epithelial cells of the prostate gland.^{1,2} The level of PSA in serum and other tissues is normally very low. In malignant prostate disease (prostatic adenocarcinoma) and in non-malignant disorders such as benign prostate hypertrophy (BPH) and prostatitis, the serum level of PSA may become elevated.³


In serum, PSA exists primarily as three forms: complexed with either a 1-antichymotrypsin (ACT) or a 2-macroglobulin and free.^{4,5} The PSA protein associated with a 2-macroglobulin is encapsulated and unavailable for measurement by current immunoassay systems. The Dimension[®] TPSA assay measures both the free and the ACT bound components of serum PSA.

The specificity of PSA to prostate tissue makes it a significant marker in the early detection and management of prostate diseases.

Prostate cancer is the most common type of cancer found in men in the United States and the second leading cause of male cancer mortality, accounting for more than 30,000 deaths in 1999.⁶ Prior to the use of PSA for early detection of prostate cancer, the traditional method of digital rectal examination (DRE) detected considerably fewer tumors.^{3,6} The most sensitive method for early detection of prostate cancer uses both DRE and PSA. The American Cancer Society and the American Urological Association (AUA) recommend that early detection of prostate cancer should be offered to asymptomatic men 50 years of age or older with an estimated life expectancy of more than 10 years. An abnormal DRE and/or an elevated PSA may suggest the presence of prostate cancer; however a prostate biopsy is required for final diagnosis.⁶

PSA testing is also accepted as an adjunctive test in the management of prostate cancer.^{7,8} Serum levels of PSA are most useful when sequential values are obtained and monitored over time. After complete removal of the prostate gland (radical prostatectomy), PSA levels should decline to a very low or nondetectable level. A rise of the serum PSA level in prostatectomy patients indicates residual prostate tissue; recurrence or metastasis of the disease.⁹ Serum PSA levels during radiation treatment should decline and remain at baseline while the patient is in remission.¹⁰

Principles of Procedure: The TPSA method is a one step enzyme immunoassay based on the "sandwich" principle. Sample is incubated with chromium dioxide particles and β-galactosidase conjugate reagent; each coated with monoclonal antibodies recognizing different binding sites on PSA. A chrome particle/PSA^a/conjugate sandwich forms during the incubation period. The incubation conditions are optimized to produce equivalent immunoreactions with both free PSA and PSA-ACT. Unbound conjugate is removed by magnetic separation and washing. The sandwich-bound β-galactosidase is combined with the chromogenic substrate chlorophenol red-β-d-galactopyranoside (CPRG). Hydrolysis of CPRG releases a chromophore, chlorophenol red (CPR). The color change measured at 577 nm due to formation of CPR is directly proportional to the concentration of PSA present in the patient sample.

a. The PSA term used here includes the free PSA and PSA-ACT complexes.

Reagents

Wells ^b	Form	Ingredient	Concentration ^c	Source
1	Liquid	PSA Ab-β-galactosidase ^d	0.09 mg/mL	Mouse, monoclonal
2	Liquid	Chrome diluent	4.8 mg/mL	
3	Tablet ^e	Antibody-CrO ₂ ^d	1.8 mg/mL	Mouse, monoclonal
4, 5, 6	Tablets ^e	CPRG	10.3 mg/mL	
7	Liquid	CPRG diluent	42 mg/mL	

b. Wells are numbered consecutively from the wide end of the cartridge.

c. Nominal value in hydrated cartridge.

d. Antibody titer and conjugate activity may vary from lot to lot.

e. Tablets contain excipients, buffers, and stabilizers.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Contains sodium azide (<0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the Dimension[®] System.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 15 days for wells 1 – 3 and 7
5 days for wells 4 – 6

Specimen Collection and Handling: Serum and lithium heparin plasma can be collected by normal procedures.¹¹ Follow the instructions provided with your specimen collection device for use and processing.¹²

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation.¹³ Clotting time may be increased due to thrombolytic or anticoagulant therapy.

Serum and plasma specimens should be separated from cells within 3 – 4 hours after venipuncture.¹³ Samples should be kept at 4 °C and analyzed within 8 hours. For longer storage, samples may be frozen at -20 °C or colder. For extended storage samples are stable for at least four months if stored at -80 °C.¹⁴

An erroneously elevated PSA level can be observed if the serum specimen from a patient is collected following digital rectal examination (DRE), needle biopsy or transurethral resection.¹⁵

Procedure

Materials Provided

TPSA Flex[®] reagent cartridge, Cat. No. RF451

Materials Required But Not Provided

Reaction Vessels, Cat. No. RXV1A
Chemistry Wash, Cat. No. RD701
Sample Probe Cleaner, Cat. No. RD703
T/F PSA Calibrator Cat. No. RC452
Quality Control Materials

Test Steps

Sampling¹, reagent delivery, mixing, separation, processing and printing of results are automatically performed by the Dimension[®] system with the heterogeneous immunoassay module. For details of this processing, refer to your Dimension[®] Operator's Guide.

f. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Reaction Vessel	
Sample Size	40 µL
Antibody-CrO ₂	30 µL
Antibody-β-galactosidase	50 µL
Incubation Temperature	42 °C Dimension [®] RxL, Xpand [®] , RxL Max [®] , Xpand [®] Plus, EXL [™]
	37 °C Dimension [®] EXL [™] with LM Dimension [®] EXL [™] 200
Incubation Period	9 minutes
Reaction Cuvette	
Transfer Volume	60 µL
CPRG Reagent Volume	150 µL
Temperature	37.0 °C
Reaction Time	5 minutes
Wavelength	577 and 700 nm
Type of Measurement	Bichromatic rate

Calibration	
Assay Range	0.13 – 100.00 ng/mL [µg/L] ^a
Calibration Material	T/F PSA Calibrator, Cat. No. RC452 ^h
Calibration Scheme	Level 1, n = 5 Level 2, n = 2 Level 4, n = 3 Level 5, n = 2 Level 6, n = 3
Units	ng/mL [µg/L]
Typical Calibration Levels	0.0, 4.0, 20.0, 50.0, 108.0 ng/mL [µg/L]
Calibration Frequency	Every 90 days for any lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -1000.0 C ₁ 3000.0 C ₂ -2.0 C ₃ 200.0 C ₄ 0.5

g. Système International d'Unités [SI Units] are in brackets.

h. The T/F PSA calibrator has 6 levels. To calibrate TPSA, use levels: 1, 2, 4, 5, 6. Do not use vial labeled 3.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known PSA concentrations. Unless addressed by your internal laboratory procedures, do not report patient results if quality control is outside acceptance limits.

Results: The instrument automatically calculates and prints the concentration of total PSA in ng/mL [µg/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.13 – 100.00 ng/mL [µg/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 100.00 ng/mL [µg/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water to obtain results within assay range. Enter dilution factor. Reassay. Resulting readout is corrected by dilution.

Autodilution: Refer to your Dimension® Operator's Guide. Recommended Auto Dilute volume is 2 µL.

Results less than 0.13 ng/mL [µg/L] should be reported as "less than 0.13 ng/mL [µg/L]".

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies.¹⁶ Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.¹⁷

PSA levels may be lower in patients who receive hormonal therapy and may not adequately reflect the presence of residual or recurrent disease.⁹

A system malfunction may exist if the following 5-test precision is observed:

PSA Concentration	SD
4.0 ng/mL [µg/L]	>0.22 ng/mL [µg/L]
50.0 ng/mL [µg/L]	>2.63 ng/mL [µg/L]

Expected Values in Detection of Prostate Cancer

The following table summarizes the results of a retrospective study using samples from 1204 men aged 50 years and older who were referred to a urologist for determination of the presence of prostate cancer. These samples were collected from five clinical sites across the United States. All of these men had a biopsy performed. Samples were tested for PSA using the Dimension® TPSA method on the Dimension® RxL system. The PSA results and corresponding DRE and biopsy results are shown below. The study demonstrated that TPSA testing when used in conjunction with DRE, was more effective in detecting prostate cancer than DRE alone.

Distribution of results

	Biopsy Result		
	Benign	Malignant	Total
DRE +	126	116	242
DRE -	663	299	962
Mean PSA (ng/mL)	6.68	10.07	7.73
Total	789	415	1204
PSA >4.0, DRE -	509	269	778
PSA >4.0, DRE +	63	97	160
Total	572	366	938
PSA ≤4.0, DRE -	154	30	184
PSA ≤4.0, DRE +	63	19	82
Total	217	49	266

The Predictive Power (PPV) was estimated as the probability of having a positive biopsy given a positive result for DRE, PSA, and PSA with DRE. The results are as follows:

Method	PPV%	95% Confidence Interval
DRE+ only	48.0	41.7 – 54.0
PSA >4.0 only	39.0	35.9 – 42.1
PSA >4.0 or DRE+	37.8	34.8 – 40.7
PSA >4.0 and DRE+	60.6	50.9 – 80.0
PSA >4.0 and DRE -	34.6	29.2 – 40.0

Serum PSA concentrations, regardless of the value, should not be interpreted as definitive evidence for the presence or absence of prostate cancer. Prostate biopsy is required for the diagnosis of cancer.

Expected values in management of prostate cancer patients – The following table summarizes the distribution of TPSA values determined on a Dimension® RxL system in 393 healthy and 1213 diseased subjects. These patients were referred to a urologist for determination of the presence of prostate cancer.

	Number of Subjects	Distribution of PSA Values, Percent (%)				
		0 – 4	4.1 – 10	10.1 – 30	30.1 – 60	>60
		(ng/mL) [µg/L]	(ng/mL) [µg/L]	(ng/mL) [µg/L]	(ng/mL) [µg/L]	(ng/mL) [µg/L]
Healthy Subjects						
Males 50 – 59 yrs	160	92.0	8.0	0.0	0.0	0.0
Males 60 – 69 yrs	130	89.0	9.0	2.0	0.0	0.0
Males 70+ yrs	103	81	19	0.0	0.0	0.0
Prostate Diseases						
BPH	373	30.9	56.8	11.8	0.5	0.0
Prostatitis	89	28.6	58.7	12.7	0.0	0.0
Prostatic Intraepithelial Neoplasm (PIN)	264	25.0	59.0	16.0	0.0	0.0
Suspicious of Cancer	63	30.1	52.9	17.0	0.0	0.0
Prostate Cancer	424	11.8	54.0	27.6	3.3	3.3

Each laboratory should establish its own expected values for PSA as performed on the Dimension® system.

When changing PSA assays in the course of monitoring a patient, additional sequential testing should be carried out to confirm baseline values.

Specific Performance Characteristics

All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® System Operator's Guide). Data was collected on a Dimension® RxL clinical chemistry system for assessment of non-interfering substances, hook effect, recovery and analytical sensitivity. Data was collected on a Dimension® EXL™ with LM system for assessment of method comparison, precision, specificity, functional sensitivity, limit of detection and limit of blank. Equivalent performance has been demonstrated on all Dimension® systems. The following data represent typical performance for the Dimension® system.

Method Comparison

Method comparison studies were conducted with the Dimension® TPSA Flex® reagent cartridge at three sites. Two-hundred and fifty-five (255) patient samples ranging from 0.15 to 99.7 ng/mL [µg/L] were randomly assayed on the Dimension® EXL™ with LM system vs. the Dimension® EXL system. The data was analyzed using Deming analysis and the pooled data is summarized below.

Regression Statistics

Comparative Method	Slope	95% CI	Intercept ng/mL [µg/L]	95% CI
TPSA on Dimension® EXL	1.05	1.03 to 1.06	0.02	-0.08 to 0.11

Precision

Precision testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A2)¹⁸. Specimens at each level were analyzed in singlet, twice per day, two runs per day, for 20 days. The study was conducted at three sites. Data from a single site and the pooled data from all three sites are presented in the following tables.

Material	Precision from one site		Total SD (%CV)
	Mean ng/mL [µg/L]	Within-Run SD (%CV)	
Plasma pool	0.09	0.01 (13.8)	0.02 (18.0)
Serum pool	0.59	0.02 (2.5)	0.02 (3.8)
Liquichek™ 2	2.08	0.03 (1.6)	0.06 (2.8)
Serum Pool	3.99	0.06 (1.4)	0.12 (2.9)
Liquichek™ 3	19.52	0.37 (1.9)	0.51 (2.6)
Internal QC	74.20	1.27 (1.7)	1.84 (2.5)

Material	Pooled precision from three sites		Total SD (%CV)
	Mean ng/mL [µg/L]	Within-Run SD (%CV)	
Plasma pool	0.09	0.01 (15.6)	0.02 (25.6)
Serum pool	0.61	0.03 (4.8)	0.04 (6.7)
Liquichek™ 2	2.13	0.04 (2.0)	0.09 (4.4)
Serum Pool	4.08	0.11 (2.8)	0.18 (4.0)
Liquichek™ 3	19.65	0.42 (2.1)	0.62 (3.2)
Internal QC	72.51	1.37 (1.9)	2.34 (3.2)

Liquichek™ Immunoassay Plus Control is a trademark of Bio-Rad Laboratories, Irvine, CA 92618.

Specificity

The two monoclonal antibodies used in the Dimension® TPSA test recognize PSA and PSA-ACT on an equimolar basis in the range of 0-100% free PSA. Equimolarity was assessed by testing various mixtures of complexed PSA (cPSA) and free PSA (fPSA) at set total PSA concentrations of 0.15, 1.0, 2.0, 4.0, 10.0, 20.0 and 50.0 ng/mL [µg/L]. An equimolar assay should yield a consistent total PSA result with mixtures of fPSA and cPSA. Test results with the Dimension® TPSA assay indicate that the total concentration of PSA remained constant while varying the amounts of cPSA and fPSA at seven different ratios, such as 0% free/ 100% complexed and 10% free/ 90% complexed PSA. Linear regression slopes for the different total PSA concentrations (observed total PSA assay result versus varying ratios of fPSA: cPSA) were less than 0.1, which indicates equimolarity.

Non-Interfering Substances

The following substances do not interfere with the TPSA method when present in serum at the concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 10% at a PSA level of 4.5 ng/mL [µg/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1322 µmol/L
Albumin	6.8 g/dL	68 g/L
Amikacin	15 mg/dL	256 µmol/L
Aminoglutethimide	6 ng/mL	26 nmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic acid	3 mg/dL	170.3 µmol/L
Bilirubin	60 mg/dL	1026 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cimetidine	10 mg/dL	397 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Creatinine	30 mg/dL	2652 µmol/L
Cyclophosphamide	25 mg/dL	896 µmol/L
Dextran 75	2500 mg/dL	25 g/L
Diazepam	2 mg/dL	70 µmol/L
Diethylstilbestrol	0.02 mg/dL	700 nmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Doxorubicin-HCl	7 mg/dL	121 µmol/L
Erythromycin	20 mg/dL	272 µmol/L
Estramustine phosphate	20 mg/dL	343 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Finasteride	0.2 mg/dL	5.4 µmol/L
Flutamide	1 mg/dL	36 µmol/L
Furosemide	2 mg/dL	60 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Goserelin acetate	0.01 mg/dL	79 nmol/L
Hemoglobin	1000 mg/dL	0.62 mmol/L (monomer)
Heparin	8000 U/L	8000 U/L
Human plasma kallikrein	100 µg/mL	1 µmol/L
Ibuprofen	40 mg/dL	1939 µmol/L
Immunoglobulin G	6 g/dL	60 g/L
Ketoconazole	7 mg/dL	132 µmol/L
Leuprolide acetate	10 mg/dL	86 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lipemia	3000 mg/dL	33.9 mmol/L (triglyceride)
Lithium	3.5 mg/dL	5.04 mmol/L
Megestrol acetate	2 mg/dL	52 µmol/L
Methotrexate	300 mg/dL	7 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Prostatic Acid Phosphatase	1000 ng/mL	1000 µg/L
Protein (low)	4 g/dL	40 g/L
Protein (high)	12 g/dL	120 g/L
Rheumatoid factor	571 IU/mL	571 IU/mL
Salicylic acid	50 mg/dL	3.62 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric acid	20 mg/dL	1.2 mmol/L
Valproic acid	50 mg/dL	3.467 mmol/L

Hook Effect

One step sandwich immunometric assays are susceptible to a high-dose "hook effect", where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.¹⁹ The TPSA method shows no hook effect up to 15000 ng/mL [µg/L].

Recovery

Multiple dilutions of serum and plasma samples with PSA values of 1.54 – 22.40 ng/mL [µg/L] were made with water.

Sample TPSA concentrations were measured and the percent recovery calculated. Recovery ranged from 96.0 to 103% with an average recovery of 101%.

Analytical Sensitivity: 0.05 ng/mL [µg/L]

The analytical sensitivity represents the lowest concentration of PSA that can be distinguished from zero. This sensitivity is defined as the concentration at two standard deviations above the level 1(0 ng/mL [µg/L]) calibrator (n = 20).

Functional Sensitivity: 0.13 ng/mL [µg/L]

The Functional Sensitivity represents the lowest concentration with an observed 20% coefficient of variation.

Limit of Detection and Limit of Blank

The limit of detection (LoD) for TPSA is 0.13 ng/mL [µg/L], determined consistent with CLSI guideline EP17-A and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 150 determinations, with 5 blank and 5 low level samples. The limit of blank (LoB) is 0.09 ng/mL [µg/L]. LoD is the lowest concentration of TPSA that can be detected reliably. LoB is the highest concentration that is likely to be observed for a blank sample.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex®, EXL™, Xpand®, Xpand® Plus, and RxL Max® are trademarks of Siemens Healthcare Diagnostics.

©2009 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system Heterogeneous Immunoassay Module

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Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Lilja H. A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. *J Clin Invest* 1985;76:1899-903.
2. Watt KWK, Lee P-J, Tinkulu T, Chan W-P, Loo R. Human prostate-specific antigen: structural and functional similarities with serine protease. *Proc Natl Acad Sci*. 1986; 83:3166-3170.
3. Coley CM, Barry MJ, Fleming C. et al: Early detection of prostate cancer: I. Prior probability and effectiveness of tests. *ANN Intern Med* 1997; 126:394-406.
4. Stenman V-H, Leinonen J, Alfthan H. et al. A complex between prostate-specific antigen and α 1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: Assay of the complex improves clinical sensitivity for cancer. *Cancer Res* 1991;51:222-6.
5. Lilja H, Christensson A, Dahlen V. et al. Prostate-specific antigen in human serum occurs predominantly in complex with α 1-antichymotrypsin. *Clin Chem* 1991; 37:1618-25.
6. Carroll P. et al: Prostate-Specific Antigen best practice policy – Part I: Early detection and diagnosis of prostate cancer. *Urology* 2001; 57: 217-224.
7. Rogers E, Ohori M, Kassabian VS. et al: Salvage radical prostatectomy: Outcome measured by serum prostate specific antigen levels. *J Urol* 1995; 153: 104-110.
8. Pound CR, Partin AW, Eisenberger MA. et al: Natural history of progression after PSA elevation following radical prostatectomy. *JAMA* 1999; 281: 1591-1597.
9. Morgan WR, Zincke H. et al. Prostate specific antigen values after radical retropubic prostatectomy for adenocarcinoma of the prostate: Impact of adjuvant treatment (hormonal and radiation). *J Urol* 1991; 145:319.
10. Lee W, Hanlon AL, Hanks GE.: Prostate specific antigen nadir level and disease-free survival. *J.Urol* 1996; 156: 450-453.
11. Tietz NW. Textbook of Clinical Chemistry, W.B. Saunders Co. Philadelphia, PA 1999: pp. 42-72.
12. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
13. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
14. Jung K. et al. Preanalytical Determinants of Total and Free Prostate-Specific Antigen and Their Ratio: Blood Collection and Storage Conditions. *Clin Chem* 1998; 44:685-688.
15. Crawford ED, Schutz MJ, Clejan S. et al. The effect of digital rectal examination of prostate specific antigen levels. *JAMA* 1992; 267:2227.
16. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using $F(ab')_2$ conjugate and polyclonal mouse IgG. *Clin Chem* 1992; 38:1737-42.
17. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45:7942-956.
18. Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
19. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook” effect in two-site immunometric assays. *Anal Biochem* 1982; 127:308-15.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contenuto sufficiente per <n> saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

0197