Dimension® clinical chemistry system

Flex® reagent cartridge

ALDL

See shaded sections: Updated information from 2017-08 version.

Issue Date 2019-04-22

Automated LDL Cholesterol

Intended Use: The ALDL method for the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of low-density lipoprotein cholesterol (LDL-C) in human serum and plasma. LDL-C measurements are used in the diagnosis and treatment of lipid disorders such as diabetes mellitus, atherosclerosis, and various liver and renal diseases.

Summary: Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The phospholipid, free cholesterol and protein constitute the outer surface of the lipoprotein particle, while the inner core contains mostly esterified cholesterol and triglycerides. These particles serve to solubilize and transport cholesterol and triglycerides in the bloodstream.

The relative proportions of protein and lipid determine the density of these lipoproteins and provide a basis on which to begin their classification.¹ These classes are: chylomicrons, very-low-density lipoprotein (VLDL), low-density lipoprotein (LDL), intermediate density lipoprotein (IDL), high-density lipoprotein (HDL) and lipoprotein (a) (Lp(a)). LDL is the main cholesterol-containing particle in plasma. When present in excessive amounts, LDL-C can be deposited in the arterial wall resulting in atherosclerosis.²

Clinical studies have shown that the different lipoprotein classes have very distinct and varied effects on coronary artery disease (CAD) risk. Additionally, numerous studies all point to LDL cholesterol as a key factor in the development of atherosclerosis and CAD. For this reason, the Third Report of the National Cholesterol Education Program (NCEP) Expert panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III- ATP III) identified elevated LDL-C as the primary target of cholesterol-lowering therapy. As a result, the cutpoints for initiating treatment are stated in terms of LDL-C concentration.³

Methods for LDL-C measurement assume that total cholesterol is composed primarily of cholesterol in VLDL, IDL, LDL, HDL and Lp(a). LDL-C can be measured using both indirect and direct methods. The Friedewald equation developed in 1972 is the most frequently used indirect method for estimating LDL-C concentration. Using this equation, LDL-C concentration is calculated as follows:

$$[LDL-C] = [Total\ Chol] - [HDL\ Chol] - [Triglyceride]/5$$


All concentrations are in mg/dL. The factor [Triglyceride]/5 is an estimate of VLDL cholesterol concentration and is based on the average ratio of triglyceride to cholesterol in VLDL. In practice, the Friedewald calculation works reasonably well. However, it should not be used with samples that have triglyceride concentrations above 400 mg/dL, when chylomicrons are present (i.e., non-fasting specimens) or in patients with dysbetalipoproteinemia (Type III hyperlipoproteinemia).⁴ At high triglyceride concentrations, LDL-C concentrations are underestimated.

Until recently the only means of consistently measuring LDL-C concentrations accurately was to perform beta-quantification; an expensive, time consuming and labor intensive approach that most clinical laboratories are unable to perform. The Automated Low Density Lipoprotein (ALDL) method is a direct assay not dependent on the Friedewald calculation and is referenced to the beta-quantification determination of LDL-C concentration.

Principles of Procedure:

The ALDL Cholesterol assay is a homogeneous method for directly measuring LDL-C levels in human serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

The method is in a two reagent format and depends on the properties of detergent 1 which solubilizes only non- LDL particles. cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non-color forming reaction. Detergent 2 solubilizes the remaining LDL particles. The soluble LDL-C is then oxidized by the action of cholesterol esterase and cholesterol oxidase forming cholestenone and hydrogen peroxide (H₂O₂). The enzymatic action of peroxidase on H₂O₂ produces color in the presence of N,N-bis(4-sulfobutyl)-m-toluidine, disodium salt (DSBmT) and 4-aminoantipyrine (4-AA) that is measured using a bichromatic (540, 700 nm) endpoint technique. The color produced is directly proportional to the amount of LDL-C present in the sample.

Reagents

Wells ^a	Form	Ingredient	Concentration	Source
1, 2, 3 (Reagent 1)	Liquid	MES Buffer, Detergent 1 Cholesterol Esterase Cholesterol Oxidase Peroxidase 4-aminoantipyrine (4-AA) Ascorbic acid oxidase Preservative	pH 6.3	Cellulomonas sp. Pseudomonas sp. Horseradish
4, 5, 6 (Reagent 2)	Liquid	MES Buffer Detergent 2 DSBmT ^b Preservative	pH 6.3	Curcubita sp.

a. Wells are numbered consecutively from the wide end of the cartridge.

b. N,N-bis(4-sulfobutyl)-m-toluidine, disodium salt

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; Handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8° C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 5 days for wells 1 – 6

Specimen Collection and Handling: Blood should be collected after a 12-hour fast by normal procedures.

Follow the instructions provided with your specimen collection device for use and processing.⁵

Complete clot formation should take place before centrifugation.^{6,7}

Serum, EDTA-treated or heparinized (lithium or sodium heparin) plasma are the recommended specimens. Serum or plasma samples should be removed from cells within 3 hours of venipuncture. Serum or plasma may be refrigerated at 2 – 8 °C for up to 3 days if not tested within 24 hours. For longer storage samples may be frozen at -20 °C for several weeks or at -70 °C or lower for longer periods.⁴

EDTA plasma results should be multiplied by 1.03 to provide serum-equivalent results.⁸

Corvac® and SST® collection tubes do not affect the ALDL method.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.

SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure:

Materials Provided

ALDL Flex® reagent cartridge, Cat. No. DF131

Materials Required But Not Provided

ALDL Calibrator, Cat. No. DC131

Quality Control Materials

Test Steps

Sampling⁹, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume.

Test Conditions

Sample Size	3 µL
Reagent 1 Volume	300 µL
Reagent 2 Volume	100 µL
Temperature	37 °C
Wavelength	540 and 700 nm
Type of Measurement	Bichromatic Endpoint

Calibration

Assay range	5 – 300 mg/dL [0.13 – 7.8 mmol/L] ^d
Calibration Material	ALDL Calibrator, Cat. No. DC131
Calibration Scheme	3 Levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.0259) = [mmol/L]
Typical Calibration Levels	0, 130, 315 mg/dL [0, 3.4, 8.1 mmol/L]
Calibration Frequency	Every 60 days for any one lot
A new calibration is required	- For each lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 1.44 C ₁ 0.98

d. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known low-density lipoprotein concentrations. The National Cholesterol Education Program recommends controls that span the medical decision points and are traceable to the National Reference System (NRS/CHOL) reference materials and methods.⁸

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of LDL-C in mg/dL [mmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5 – 300 mg/dL [0.13 – 7.8 mmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 300 mg/dL [7.8 mmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water to obtain a result within assay range. Enter the dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): If using the auto-dilution feature, results above 300 mg/dL [7.8 mmol/L] will automatically be repeated. The autodilution volume is 2 µL.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in ALDL results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

A system malfunction may exist if the following 5 test precision is observed:

Concentration	SD
130 mg/dL [3.4 mmol/L]	>2.0 mg/dL [0.05 mmol/L]
315 mg/dL [8.1 mmol/L]	>5.0 mg/dL [0.13 mmol/L]

Interfering Substances

Bilirubin (unconjugated) of 80 mg/dL [1368 µmol/L]* will decrease an ALDL result of 124 mg/dL [3.2 mmol/L] by 10%.

Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] will decrease an ALDL result of 122 mg/dL [3.2 mmol/L] by 19%.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values: The National Cholesterol Education Program Adult Treatment Panel III (NCEP- ATP III) 3 provides the following classifications of LDL-C concentrations:

Category	LDL -C	
	mg/dL	[mmol/L]
Optimal	<100	<2.6
Near optimal/above optimal	100-129	2.6-3.3
Borderline High	130-159	3.4-4.1
High	160-189	4.1-4.9
Very High	≥190	≥4.9

Each laboratory should establish its own reference interval for LDL-C as performed on the Dimension® system.

Specific Performance Characteristics*

Material	Precision ^g		
	Mean	Standard Deviation (% CV)	
	mg/dL [mmol/L]	Within-run	Total
HDL Plus® QC			
Level 1	170 [4.4]	1.60 [0.04] (0.94)	4.34 [0.11] (2.55)
Level 2	120 [3.1]	0.87 [0.02] (0.72)	3.38 [0.09] (2.83)
Level 3	53 [1.4]	0.55 [0.01] (1.03)	1.04 [0.03] (1.95)
Serum Pool 1	106 [2.7]	1.42 [0.04] (1.34)	2.72 [0.07] (2.57)
Serum Pool 2	163 [4.2]	2.66 [0.07] (1.63)	3.62 [0.09] (2.22)

e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed on Dimension® (Refer to your Dimension® Operator's Guide).

f. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999).

g. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

HDL Plus® is a registered trademark of Quantimetrix Corporation, Redondo Beach, CA.

Method Comparison

Regression Statistics^h

Comparative Method	Slope	Intercept mg/dL [mmol/L]	Correlation Coefficient	n ⁱ
N-geneous™ LDL Cholesterol Reagent ^l	0.95	4.7 [0.12]	0.997	122
Beta-Quantification ^k	1.01	3.3 [0.08]	0.982	49

h. Model equation for regression statistics is: Result of Dimension® analyzer = [Slope x Comparative method result] + Intercept.

i. The range of LDL-C values in the correlation study was 35 to 273 mg/dL [0.90 to 7.06 mmol/L] for the N-geneous™ Cholesterol Reagent^l and 70 to 246 mg/dL [1.81 to 6.36 mmol/L] for Beta-Quantification. Correlation study was performed using serum samples.

j. N-geneous™ LDL Cholesterol Reagent^l is manufactured by Genzyme Corporation, Cambridge, MA 02139-1562. Testing was performed on the Beckman CX9 chemistry analyzer manufactured by Beckman Coulter Inc., Brea, CA 92821.

k. Beta-Quantification Reference Method for LDL-C is a three-step procedure⁹ involving ultracentrifugation, precipitation with heparin-manganese reagent and quantification with the Abell-Kendall reference method for cholesterol. This study was performed by the Core Laboratory for Clinical Studies at Washington University, St. Louis, Missouri, 63110.

Specificity

HIL Interference

The ALDL method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration [SI Units]	LDL Cholesterol Conc. [SI Units]	Bias (%) ^l
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	123 mg/dL [3.19 mmol/L]	<10
Bilirubin (unconjugated)	60 mg/dL [1026 µmol/L]	124 mg/dL [3.21 mmol/L]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	122 mg/dL [3.16 mmol/L]	<10

l. Analyte results should not be corrected based on this bias.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Non-Interfering Substances

Interference from the following substances at the concentrations indicated when added to a 110 mg/dL [2.8 mmol/L] LDL-C serum pool respectively, is less than 10%.

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Albumin	6.8 g/dL	68 g/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic Acid	5 mg/dL	284 µmol/L
Atorvastatin	3.6 µg/mL	2.98 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol	500 mg/dL	13.0 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Clofibrate	40 mg/dL	1648.1 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Fenofibrate	4.6 mg/dL	127.4 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gemfibrozil	12.4 mg/dL	495.3 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin (lithium)	8 U/mL	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
IgG	5 g/dL	50 g/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium	3.5 mg/dL	5.04 mmol/L
Lovastatin	2 µg/mL	4.94 µmol/L
Niacin	3.1 mg/dL	251.8 µmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Pravastatin	10.3 µg/mL	24.3 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein: Total	3.6 g/dL	36 g/L
Protein: Total	11.8 g/dL	118 g/L
Rheumatoid factors	500 IU/mL	500 IU/mL
Salicylic Acid	50 mg/dL	3.62 mmol/L
Simvastatin	6.9 µg/mL	16.5 µmol/L
Theophylline	25 mg/dL	1388 µmol/L
Triglyceride (endogenous)	1000 mg/dL	11.3 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Cross-Reactivity

Cross-reactivity of HDL-C was evaluated by adding known amounts of HDL-C to a human serum pool containing 103 mg/dL [2.7 mmol/L] of LDL-C. Percent cross-reactivity was calculated as follows:

$$\% \text{ cross reactivity} = \frac{[\text{measured LDL-C}] - [\text{control LDL-C}]}{\text{HDL-C}} \times 100$$

HDL-C Concentration	% Cross-reactivity
100 mg/dL [2.6 mmol/L]	2

Analytical Sensitivity

The sensitivity of the ALDL method is 5 mg/dL [0.13 mmol/L] and represents the lowest concentration of LDL-C that can be distinguished from zero. This sensitivity is defined as the concentration at two standard deviations above the mean (n=20) of Level 1 ALDL Calibrator (0 mg/dL) [0 mmol/L].

Certification

The method has been evaluated by and has met the certification acceptance criteria of the Cholesterol Reference Method Laboratory Network (CRMLN).

Symbols Key: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics

All rights reserved.

Dimension® clinical chemistry system

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Gotto AM. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988;23: Suppl. 1,4.
2. Rifai N, Warnick GR, Dominiczak MH. Handbook of lipoprotein testing. Washington: AACC Press, 1997, pp. 14-15, 145, 273.
3. National Cholesterol Education Program: Third Report of the Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). Journal of American Medical Association, May 16, 2001, p. 3.
4. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry-3rd edition, W.B. Saunders Co., Philadelphia, PA, 1999, pp. 42-72 (specimen collection and processing), pp. 843-845 (Friedewald Equation) and p. 849 (specimen storage).
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 (ISBN 1-56238-555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 (ISBN 1-56238-519-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
8. NIH Publication No. 95 - 3044: NCEP Program: Recommendations on Lipoprotein Measurement - From the Working Group on Lipoprotein Measurement, Sept. 1995, p. 13 (Executive Summary), p. 47 (Recommendations for Manufacturers).

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare