Dimension® clinical chemistry system

Flex® reagent cartridge

ALB

See shaded sections: Updated information from 2015-01 version.

Issue Date 2019-04-22

Albumin


Intended Use: The ALB method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of albumin in human serum and plasma.

Summary: Measurements of albumin are used in the diagnosis and treatment of numerous diseases involving primarily the liver or kidneys.

Albumin is the protein of the highest concentration in plasma. Albumin is formed exclusively in the liver and serves as a transport and binding protein for calcium, fatty acids, bilirubin, hormones, vitamins, trace elements and drugs. It is also of prime importance in maintaining the colloidal osmotic pressure in both the vascular and extravascular spaces.¹ Decreased serum albumin concentration can result from liver disease. It can also result from kidney disease, which allows albumin to escape into the urine. Decreased serum albumin may also be explained by malnutrition or a low protein diet.²

The albumin method is an adaptation of the bromocresol purple (BCP) dye-binding method reported by Carter³ and Louderback, et al.⁴ Because of an enhanced specificity of BCP for albumin, this method is not subject to globulin interference.⁵ Multiple wavelength blanking increases sensitivity and minimizes spectral interference from lipemia.

Principles of Procedure: In the presence of a solubilizing agent, BCP binds to albumin at pH 4.9. The amount of albumin-BCP complex is directly proportional to the albumin concentration. The complex absorbs at 600 nm and is measured using a polychromatic (600, 540, 700 nm) endpoint technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 6	Liquid	BCP dye Acetate Buffer Surfactant Microbial inhibitor	2.7 x 10 ⁻⁴ M

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 2-chloroacetamide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.^{6,7}

Follow the instructions provided with your specimen collection device for use and processing.⁸

Complete clot formation should take place before centrifugation.⁹

Specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.⁹

Heparin and EDTA, in the concentrations routinely found in commercially available blood collection tubes (14.3 U/mL), do not interfere with the ALB method. Lithium heparin (280 U/mL [280,000 U/L]^c) lowers the ALB result by 2.4 g/dL [24 g/L] at an albumin concentration of 4.1 g/dL [41 g/L].

Corvac® and SST® collection tubes do not affect the ALB method.

c. Système International d'Unités [SI Units] are in brackets.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.

SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

ALB Flex® reagent cartridge, Cat. No. DF13

Materials Required But Not Provided

Total Protein/Albumin Calibrator, Cat. No. DC31

Quality Control Materials

Test Steps

Sampling^d, reagent delivery, mixing, and processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume; precise container filling is not required.

Test Conditions

Sample Size	5 µL
Reagent 1 Volume	125 µL
Diluent Volume	370 µL
Temperature	37 °C
Wavelength	540, 600 and 700 nm
Type of Measurement	Polychromatic endpoint

Calibration

Assay Range	0.6 – 8.0 g/dL [6 – 80 g/L]
Calibration Material	Purified human albumin or secondary calibrators such as Total Protein/Albumin Calibrator Cat. No. DC31
Calibration Scheme	3 levels, n = 3
Units	g/dL [g/L] (g/dL x 10) = [g/L]
Typical Calibration Levels	0.5, 4.5, 8.3 g/dL [5, 45, 83 g/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	- For each lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -1.060 C ₁ 0.023

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known albumin concentration.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of albumin in g/dL [g/L] using the calculation scheme illustrated in the Principles of Operation Section of your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.6 – 8.0 g/dL [6 – 80 g/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 8.0 g/dL [80 g/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water to obtain result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

Results less than 0.6 g/dL [6 g/L] should be reported as "less than 0.6 g/dL [6 g/L]".

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in albumin results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
0.69 g/dL [7 g/L]	> 0.10 g/dL [1.0 g/L]
3.84 g/dL [38 g/L]	> 0.12 g/dL [1.2 g/L]

Interfering Substances

CMPF (3-carboxy-4-methyl-5-propyl-2-furanpropanoic acid) present in sera of patients with renal failure has been reported to give falsely low albumin values.^{10, 11}

Lipemia (Intralipid®) at 1000 mg/dL [11.3 mmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values:¹² 3.4 – 5.0 g/dL [34 – 50 g/L]

This reference population consisted of the following: 240 males
240 females

The reference interval was calculated non-parametrically and represents the central 95% of the population.
Posture affects serum protein concentration; it is lower when a subject is in the supine position.¹³
Each laboratory should establish its own reference interval for albumin as performed on the Dimension® system.

Specific Performance Characteristics^a

Material	Precision ^d		
	Mean g/dL [g/L]	Standard Deviation (% CV)	
		Within-run	Total
Multiquel® Control			
Normal	4.23 [42]	0.056 [0.56] (1.3)	0.072 [0.72] (1.7)
Abnormal	3.21 [32]	0.050 [0.50] (1.6)	0.073 [0.73] (2.3)

e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (Refer to your Dimension Operator's Guide).
f. Specimens at each level were analyzed in triplicate for 20 runs. The within-run and total standard deviations were calculated by the analysis of variance method.
Multiquel® is a registered trademark of Bio-Rad Laboratories, Inc., Irvine, CA 92618, USA.

Comparative Method	Method Comparison		
	Regression Statistics ^a		
	Slope	Intercept g/dL [g/L]	Correlation Coefficient
ALB method on the aca® discrete clinical analyzer	1.060	0.13 [1.3]	0.991

g. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

Specificity

HIL Interference
The ALB method was evaluated for interference from hemoglobin, icterus, and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration [SI Units]	Albumin Conc. g/dL [g/L]	Bias (%) ⁱ
Hemoglobin (hemolysate)	1000 mg/dL	2.3	<10
	[0.62 mmol/L] (monomer)	[23.0]	
Bilirubin (unconjugated)	80 mg/dL	2.3	<10
	[1368 µmol/L]	[23.0]	
Lipemia (Intralipid®)	600 mg/dL	2.3	<10
	[6.78 mmol/L]	[23.0]	
	1000 mg/dL	2.3	
	[11.3 mmol/L]	[23.0]	i

h. Analyte results should not be corrected based on this bias.
i. The interference testing at this level tripped a test report message; therefore the magnitude of the interference could not be determined.

Non-Interfering Substances

The following substances have no measurable effect on the ALB method at the concentrations indicated:

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	20 ng/mL	25.6 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	800 mg/dL	174 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	16 µg/mL	29.4 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	80 µg/mL	344 µmol/L
Phenytoin	30 µg/mL	119 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Salicylate	100 mg/dL	7.24 mmol/L
Theophylline	100 µg/mL	555 µmol/L
Total Protein	12 g/L	12 g/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 0.6 g/dL [6 g/L]
The analytical sensitivity represents the lowest concentration of albumin that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of level 1 (0.0 g/dL [0 g/L]), which is system water.
Symbols Key: See adjacent panel.
Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.
©2008 Siemens Healthcare Diagnostics
All rights reserved

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 547-548, (clinical significance).
2. Whicher J, Spence C. When is serum albumin worth measuring? Ann Clin Biochem 1987; 24:572-80.
3. Carter P. ultramicroestimation of human serum albumin: binding cationic dye 5,5'-dibromo-o-cresolsulfonphthalein, Microchem J 1970; 15:531-539.
4. Louderback A, Measley EH, Taylor NA. A new dye-binder technic using bromocresol purple for Determination of albumin in serum, Clin Chem 1968; 14:793-794.
5. Laskey FD, Li ZMC, Shaver DD, Savory MG, Willey DG, Mikolaj BJ, Lantry CL. Evaluation of a Bromocresol purple method for the determination of albumin adapted to the Du Pont aca® discreteClinical analyzer, Clin Biochem 1985; 18:290-296.
6. Tietz NW. Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1986, pp 478-497 (specimen collection and storage recommendations).
7. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 1898 USA, 2003.
8. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection: Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
9. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 (ISBN 1-56238 555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
10. Mabuchi H, Nakahashi H. Underestimation of serum albumin by the bromocresol purple method and a major endogenous ligand in uremia. Clin Chim Acta 1987; 167: 89-96.
11. Maguire GA, Price CP. Bromocresol purple method for serum albumin gives falsely low values in patients with renal insufficiency. Clin Chim Acta 1986; 155: 83-88.
12. Willey DA, Savory J, Lasky F. An Evaluation of a Revised Albumin Method for the aca® discrete clinical analyzer, Du Pont Company, Wilmington, DE, August 1982.
13. Humphrey KR, Gruemer HD, Lott JA. Impact of posture on the reference range for serum proteins and calcium, Clin Chem 1977; 23:1343-1345.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
EXP CCYY-MM-DD	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Límite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
CONTENTS	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel

2006-06_EFIGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare