Dimension® clinical chemistry system

Flex® reagent cartridge

PCHE

See shaded sections: Updated information from 2016-02 version.


Issue Date 2019-04-22

Pseudocholinesterase

Intended Use: The PCHE method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of pseudocholinesterase activity in human serum and heparinized plasma. PCHE measurements may be used in the diagnosis and treatment of cholinesterase inhibition disorders.

Summary: The pseudocholinesterase method is based on the coupled oxidation-reduction indicator reaction described by Gal and Roth.¹

Principles of Procedure: Butyrylthiocholine (BTC) is hydrolyzed by PCHE, releasing thiocholine which directly reduces the blue dye, 2,6-dichlorophenol-indophenol (DIP) to the colorless form. The change in absorbance at 600 nm due to the disappearance of DIP is directly proportional to the PCHE activity and is measured using a bichromatic (600, 700) rate technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 3	Tablet ^c	DIP	
		Buffer	0.06 µmol/L
4 – 6	Tablet ^c	BTC	1.0 µmol/L

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

c. Tablet contains excipients.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the instrument.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.²

Blood collection tubes containing potassium oxalate or sodium fluoride should not be used.

Corvac® and SST® collection tubes do not affect the PCHE method.

Follow the instructions provided with your specimen collection device for use and processing.³

Complete clot formation should take place before centrifugation.

Separated serum/plasma specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.⁴

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.

SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

PCHE Flex® reagent cartridge, Cat. No. DF51

Materials Required But Not Provided

PCHE Verifier, Cat. No. DC28

Quality Control Materials

Test Steps

Sampling,⁴ reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to the Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Size	2 µL
Reagent 1 Volume	300 µL
Reagent 2 Volume	50 µL
Diluent Volume	148 µL
Temperature	37 °C
Wavelength	600 and 700 nm
Type of Measurement	Bichromatic rate

Verification

Assay Range (37 °C)	0 – 14 U/mL [0 – 14000 U/L] ^e
Verification Material	PCHE Verifier, Cat. No. DC28
Verification Scheme	3 levels, n = 3
Units	U/mL [U/L] (U/mL x 1000) = [U/L]
Typical Verification Levels	3, 6, 10 U/mL [3000, 6000, 10000 U/L]
Verification Frequency	Every 3 months for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C_0 0.000 C_1 0.0464
e. Système International d Unités [SI Units] are in brackets.	

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known pseudocholinesterase activity. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the activity of pseudocholinesterase in U/mL [U/L] using the calculation scheme illustrated in your Dimension® Operator's Guide. A change of 0.2 milliabsorbance units (mA) per minute corresponds to pseudocholinesterase activity of 1 U/mL at 37 °C.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0 – 14 U/mL [0 – 14000 U/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 14 U/mL [14000 U/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilutions with Enzyme Diluent (Cat. No. 790035901) to obtain result within the assay range. Enter dilution factor.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions.

Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Activity	SD
4.0 U/mL [4000 U/L]	>0.1 U/mL [100 U/L]
10.1 U/mL [10100 U/L]	>0.3 U/mL [300 U/L]

Interfering Substances

Since BTC has been reported nonreactive toward erythrocytic (true) cholinesterase, partially hemolyzed specimens should not cause interference.⁵

Expected Values:⁶ 7 – 19 U/mL [7000 – 19000 U/L]

(Generally values below the reference interval are clinically significant.)

This reference population consisted of the following:

- 241 subjects (approximately evenly distributed between males and females)
- Age range of 20 to 65 years

This reference interval was calculated non-parametrically and represents the central 95% of the population.

Each laboratory should establish its own reference interval for PCHE as performed on the Dimension® system.

Specific Performance Characteristics¹

Precision^a

Material	Mean U/mL [U/L]	Standard Deviation (% CV)	
		Within-run	Between-day
Fisher Control			
Serum Unassayed			
Level 1	10.6 [10600]	0.2 [200] (2.4)	0.3 [300] (3.1)
Level 2	6.8 [6800]	0.1 [100] (1.7)	0.2 [200] (3.6)

f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (see Operator's Guide).

g. Specimens at each level were analyzed in triplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

Method Comparison

Regression Statistics^h

Comparative Method	Slope	Intercept U/mL [U/L]	Correlation Coefficient	n
PCHE on the aca® discrete clinical analyzer	1.12	0.2 [200]	0.977	102 ⁱ

h. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

i. Range of samples: 0.2 – 22.5 U/mL

Specificity

HIL Interference

The PCHE method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	PCHE Activity U/mL [U/L]	Bias ^j %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	5 [5000]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	5 [5000]	<10
Lipemia (Intralipid®)	3000 mg/dL [34.29 mmol/L]	5 [5000]	<10

j. Analyte results should not be corrected based on this bias.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Non-Interfering Substances

The following substances do not interfere with the PCHE method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at the pseudocholinesterase activity of 14 U/mL [14000 U/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	284 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	252 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	421 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	3000 mg/dL	33.9 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 0.8 U/mL [800 U/L]

The analytical sensitivity represents the lowest activity of pseudocholinesterase that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water.

Symbols Key: See adjacent panel

Bibliography: See adjacent panel

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

- Gal J, and Roth E. Spectrophotometric Methods for Determination of Cholinesterase Activity, *Clinica Chimica Acta* Vol. 2 (1957) 316-326.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Szasz G. Comparison Between Thiocolinesters and o-Nitrophenylbutyrate as Substrates in the Assay of Serum Cholinesterase, *Clin Chem* 1968;14:646.
- Westgard JO, Hund MR, Lahmeyer BL. Estimates of Normal Ranges for 15 Determinations on the Du Pont Automatic Clinical Analyzer, Dept. Of Medicine and Clinical Laboratories, University of Wisconsin Center for Health Science, Madison, WI.

Symbols Key Symbol Schlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbrikante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare