Dimension® clinical chemistry system

Flex® reagent cartridge

THEO

See shaded sections: Updated information from 2015-04 version

Issue Date 2019-04-22

Theophylline

Intended Use: The THEO method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of theophylline in human serum and plasma.

Summary: Theophylline measurements may be used in the diagnosis and treatment of theophylline overdose and in therapeutic drug monitoring. Theophylline is a methylated xanthine, 1,3-dimethylxanthine. It is structurally related to purines and uric acid, as well as to xanthine itself. The most commonly used compounds are aminophylline, the double salt of theophylline and ethylenediamine. About 10% is excreted unchanged in the urine and the remaining 90% of the drug is converted to other compounds before it is eliminated from the body. The biologic half-life of theophylline varies from about 3.5 hours in young children to 8 – 9 hours in most adults. It is substantially prolonged in the presence of liver disease and/or cardiac decompensation.^{1,2,3}

Principles of Procedure: The THEO method is a homogeneous particle enhanced turbidimetric inhibition immunoassay (PETINIA) technique which uses a synthetic particle-theophylline conjugate (PR) and theophylline-specific monoclonal antibody (Ab). Theophylline present in the sample competes with the particles for the antibody, thereby decreasing the rate of aggregation. Hence, the rate of aggregation is inversely proportional to the concentration of theophylline in the sample. The rate of aggregation is measured using a turbidimetric rate at 340 nm.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2	Liquid	Antibody	0.070 g/L ^{c,d}	Mouse, monoclonal
4, 5	Liquid	Particle Reagent	8.0 g/L ^{c,d}	
7	Liquid	Buffer	150 mM	

- Wells are numbered consecutively from the wide end of the cartridge.
- Nominal value per test at manufacture.
- Antibody titer and conjugate activity may vary from lot to lot.
- Contains buffers, stabilizers and preservatives.

Risk and Safety:

H412
P273, P501
Harmful to aquatic life with long lasting effects.

Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 2-methyl-4-isothiazolin-3-one

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Contains human source material. Each donor unit used in the preparation of this product was tested by FDA-approved methods for the presence of antibodies to Human Immunodeficiency Virus type 1 (HIV-1) and type 2 (HIV-2), as well as for Hepatitis B surface Antigen (HBsAg) and antibody to Hepatitis C Virus (HCV), and found to be negative (not repeatedly reactive). Because no testing can offer complete assurance that these or other infectious agents are absent, this material should be handled using good laboratory practice to avoid skin contact and ingestion.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 5
30 days for well 7

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁴

Follow the instructions provided with your specimen collection device for use and processing.⁵

Complete clot formation should take place before centrifugation.⁶

Specimens should be free of particulate matter.

Separated specimens are stable for 8 hours at room temperature and 7 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder for up to 3 months.^{6,7}

Collection tubes containing EDTA, lithium heparin, potassium oxalate or sodium fluoride do not affect the THEO method.

Corvac® and SST® collection tubes do not affect the THEO method.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.
SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

THEO Flex® reagent cartridge, Cat. No. DF71

Materials Provided

Drug Calibrator, Cat. No. DC22B
Quality Control Materials

Materials Required But Not Provided

Test Steps

Sampling,⁸ reagent delivery, mixing, and processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

e. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume. Precise container filling is not required.

Test Conditions

Sample Volume	4 µL
Reagent 1 Volume	305 µL
Reagent 2 Volume	40 µL
Reagent 3 Volume	40 µL
Diluent Volume	111 µL
Temperature	37 °C
Wavelength	340 nm
Type of Measurement	Turbidimetric rate

Calibration

Assay Range 2.0 – 40.0 µg/mL [11 – 222 µmol/L]^f
Calibration Material Drug Calibrator, Cat. No. DC22B
Calibration Scheme 5 levels, n = 2
Units µg/mL [µmol/L]

(µg/mL x 5.55) = [µmol/L]
Typical Calibration Levels 0.0, 5.0, 10.0, 20.0, 40.0 µg/mL
[0, 28, 56, 111, 222 µmol/L]

Calibration Frequency Every 2 months for any one lot
A new calibration is required

Assigned Coefficients
C₀ 330.0
C₁ -350.0
C₂ -2.00
C₃ 14.00
C₄ 0.500

f. Système International d' Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known theophylline concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of theophylline in µg/mL [µmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 2.0 – 40.0 µg/mL [11 – 222 µmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 40.0 µg/mL [222 µmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with theophylline-free serum or plasma or Reagent grade water to obtain result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

THEO results less than 2.0 µg/mL [11 µmol/L] should be reported as "less than 2.0 µg/mL [11 µmol/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
10 µg/mL [56 µmol/L]	>0.6 µg/mL [3.3 µmol/L]
20 µg/mL [111 µmol/L]	>1.0 µg/mL [5.5 µmol/L]

Interfering Substances

A theophylline metabolite, 1,3-dimethyl uric acid, is usually undetectable in samples from patients receiving theophylline. However, it can reach detectable levels in uremic patients. Theophylline values will be increased by 1.0 µg/mL [5.5 µmol/L] in the presence of 10 µg/mL [51.0 µmol/L] of 1,3-dimethyl uric acid.

Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] and above tripped an error flag on this method, so the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Therapeutic Range: Therapeutic theophylline concentrations vary significantly, depending on the individual. A range of 10 – 20 µg/mL [56 – 111 µmol/L] includes effective serum concentrations for many patients; however, some individuals are best treated outside this range. Concentrations greater than 20.0 µg/mL [111 µmol/L] are often associated with toxic symptoms.^{1,2,3}

The physician must determine the appropriate therapeutic range for each patient.

Specific Performance Characteristics^a

Material	Mean µg/mL [µmol/L]	Precision ^b	
		Standard Deviation (% CV) Within-run	Between-day
Control Gilford TDM			
Level I	29.1 [161]	0.4 [2.2] (1.4)	0.9 [5.0] (3.1)
Level II	14.1 [78]	0.2 [1.1] (1.4)	0.5 [2.8] (3.5)
Level III	7.6 [42]	0.1 [0.5] (1.3)	0.2 [1.1] (2.6)

g. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

h. Specimens at each level were analyzed in triplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

Method Comparison

Comparative Method	Regression Statistics ^d			
	Slope	Intercept µg/mL [µmol/L]	Correlation Coefficient	n
THEO method on the aca® discrete clinical analyzer	0.95	0.1 [0.5]	0.986	124

j. Model equation for regression statistics is: [results of-Dimension® system] = slope x [comparative method results] + intercept.

Specificity

HIL Interference

The THEO method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	THEO Concentration µg/mL [µmol/L]	Bias ^f %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	17.8 [89]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	19.3 [106]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	18.7 [104]	<10

j. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the THEO method when present in the amounts indicated, which approximate maximum physiological concentrations. Systematic inaccuracies (bias) due to these substances are less than 1.0 µg/mL [5 µmol/L] at theophylline concentrations of 10.0 µg/mL [55 µmol/L].

Substance	Test Concentration	SI Units
Caffeine	30 µg/mL	154.5 µmol/L
Theobromine	10 µg/mL	55.5 µmol/L
p-Xanthine	10 µg/mL	56.8 µmol/L
3-Methyl-Xanthine	10 µg/mL	60.2 µmol/L
1-Methyluric Acid	10 µg/mL	54.5 µmol/L
8-Chlorotheophylline	50 µg/mL	232.5 µmol/L

The following substances do not interfere with the THEO method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at the THEO concentration of 22.3 µg/mL [124 µmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.024 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	284 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	5 ng/mL	6.15 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 µg/mL	25.9 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	500 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylate	60 mg/dL	4.34 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 1.0 µg/mL [6 µmol/L]

This analytical sensitivity represents the lowest concentration of theophylline that can be distinguished from zero. The sensitivity is defined as the mean value (n = 20) plus two standard deviations of the 0.00 µg/mL [0 µmol/L] level 1 Drug Calibrator.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex® and aca® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Jacobs M, Senior R, Kessler G. Clinical experience with theophylline. Relationships between dosage, serum concentration and toxicity. *J AM Med Assoc* 1976; 235:1983-1986.
2. Finn AL, Taylor WJ. *Individualizing Drug Therapy, Practical Applications of Drug Monitoring*. New York: Gross, Townsend, Frank, Inc., 1981: 31-65.
3. Burtis CA, Ashwood ER, Bruns DE. *Tietz Textbook of Clinical Chemistry and Molecular Biology*, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 1261 (clinical significance), pp. 2314 (reference values).
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
7. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 3rd ed. Washington, DC: AACC Press, 2007: p 856.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFRS

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare