

Emit® II Plus Buprenorphine Assay



10872256_B



Buprenorphine Assay

1 Intended Use

The Emit® II Plus Buprenorphine Assay is a homogeneous enzyme immunoassay with a 5 ng/mL cutoff. The assay is intended for use in laboratories for the qualitative and/or semiquantitative analyses of buprenorphine in human urine. Emit® II Plus assays are designed for use with a number of chemistry analyzers.

The semiquantitative mode is for the purpose of enabling laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as Liquid Chromatography/Mass Spectrometry (LC/MS) or permitting laboratories to establish quality control procedures.

The Emit® II Plus Buprenorphine Assay provides only a preliminary analytical test result. A more specific alternative chemical method(s) must be used to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) or LC/MS are the preferred confirmatory methods.^{1,2} Other chemical confirmation methods are available. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.

2 Summary and Explanation of the Test

Buprenorphine is a semi-synthetic opioid analgesic derived from thebaine, a component of the opium poppy *Papaver somniferum*. Buprenorphine structurally resembles morphine, but has both antagonist and agonist properties.³ Buprenorphine is a schedule III drug. The Food and Drug Administration (FDA) has approved the use of Suboxone®, which contains buprenorphine for treatment of opiate dependency in the US. Under the US Drug Abuse Treatment Act of 2000 (DATA) buprenorphine can be prescribed in a physician's office for treatment of opiate dependency.

It has been shown that buprenorphine has abuse potential and may itself cause dependency. It produces typical opioid effects and side effects such as euphoria and respiratory depression.

Buprenorphine is metabolized in the human liver primarily by N-dealkylation to pharmacologically active norbuprenorphine, which along with the parent compound is conjugated to form buprenorphine glucuronide and norbuprenorphine glucuronide.⁴

The Emit® II Plus Buprenorphine Assay detects buprenorphine and norbuprenorphine in human urine and gives a positive result if these drugs are present at concentrations equal to or greater than the cutoff. Positive results for specimens containing structurally related opioid compounds have not been observed.

While confirmation techniques other than GC/MS or LC/MS may be adequate for some drugs of abuse, GC/MS or LC/MS is generally accepted as a vigorous confirmation technique for all drugs, since it provides the best level of confidence in the result.

3 Principle

The Emit® II Plus Buprenorphine Assay is a homogeneous enzyme immunoassay technique used for the analysis of specific compounds in human urine. The assay is based on competition between drug in the specimen and drug labeled with the recombinant glucose-6-phosphate dehydrogenase (G6PDH) for antibody binding sites. Enzyme activity decreases upon binding to the antibody, so the drug concentration in the specimen can be measured in terms of enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH in the presence of glucose-6-phosphate (G6P), resulting in an absorbance change that is measured spectrophotometrically. Endogenous serum G6PDH does not interfere because the coenzyme NAD functions only with the bacterial (*Leuconostoc mesenteroides*) enzyme employed in the assay.

4 Reagents

| REF | Product Description | Volume |
|---------|---|---------|
| 9S039UL | Emit® II Plus Buprenorphine Assay | |
| 9S029UL | | |
| 9S129UL | | |
| | Antibody/Substrate Reagent 1 | |
| | Mouse monoclonal antibodies to buprenorphine (0.53 µg/mL)*, | 28 mL |
| | NAD (6.9 mM), G6P (10.9 mM), bovine serum albumin, preservatives, | 115 mL |
| | and stabilizers | 1000 mL |
| | Enzyme Reagent 2 | |
| | Norbuprenorphine labeled with bacterial rG6PDH (0.50 µg/mL)*, | 14 mL |
| | HEPES buffer, bovine serum albumin, preservatives, | 50 mL |
| | and stabilizers | 500 mL |

**The antibody titer and enzyme conjugate activity may vary from lot to lot.*

Note: Reagents 1 and 2 are provided as a matched set. They should not be interchanged with components of kits with different lot numbers.

Risk and Safety

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

For *in vitro* diagnostic use.

Preparation of Assay Components

The Emit® II Plus Buprenorphine Assay reagents are provided liquid, ready to use and may be used directly from the refrigerator. Cap the reagent bottles when not in use.

Note: Caps must always be replaced on the original containers.

Storage

When not in use, reagents must be stored at 2–8°C (36–46°F), upright and with screw caps tightly closed. If stored as directed, reagents are stable until the expiration date printed on the label. Refer to the application sheet for instrument stability information. Do not freeze reagents. Avoid prolonged exposure to temperatures above 32°C. **Improper storage of reagents can affect assay performance.**

The purpose of specimen handling and storage information is to provide guidance to users; however, users may validate their own procedures for handling and storing patient samples.

5 Specimen Collection and Preparation

- Urine specimens may be collected in plastic (i.e., polypropylene, polycarbonate, polyethylene) or glass containers. Some plastics can adsorb certain drugs.
- If not analyzed immediately, specimens may be stored refrigerated or unrefrigerated for up to 5 days. After 5 days specimens should be stored frozen at -20°C.^{5,6}
- Frozen specimens must be thawed and mixed thoroughly prior to analysis.
- Specimens with high turbidity should be centrifuged before analysis.
- Urine specimens should be within the pH range of 3.0–11.0.
- Adulteration of the urine specimen may cause erroneous results. If adulteration is suspected, obtain another specimen.
- Human urine specimens should be handled and treated as if they are potentially infectious.
- Repeat freeze/thaw cycles should be avoided.

6 Procedure

Materials Provided

Emit® II Plus Buprenorphine Assay
Reagent 1
Reagent 2

Materials Required But Not Provided

Emit® Calibrator/Control
9A509UL Level 0 (0 ng/mL)

Emit® II Plus Specialty Drug Calibrators/Controls

9S529UL Level 1 (2.5 ng/mL)
9S549UL Level 2 (5 ng/mL)
9S569UL Level 3 (15 ng/mL)
9S589UL Level 4 (25 ng/mL)

Emit® II Plus Specialty Drug Controls

9S609UL Negative (3 ng/mL)
9S619UL Positive (7 ng/mL)

Instruments

Siemens Healthcare Diagnostics provides instructions for using this assay on a number of chemistry analyzers. Contact the Technical Solutions Center in the USA or your local Siemens representative for application sheets.

Analyzers must be capable of maintaining a constant reaction temperature, pipette specimens/reagents, mix thoroughly, measure enzyme rates precisely and time the reaction accurately.

Daily Maintenance

Refer to the instrument Operator’s Guide for maintenance instructions.

Assay Sequence

To run the assay, see the instrument Operator’s Guide and the application sheets available from Siemens.

Calibration

Calibrate the assay using the calibrators listed in Table 1.

Note: These reagents are qualified for use with these calibrators only.

Table 1 — Emit® II Plus Specialty Drug Calibrators/Controls for use in Qualitative or Semiquantitative Analysis

| Cutoff Level (ng/mL) | Required Calibrator/Control for Qualitative Analysis (ng/mL) | Required Calibrators/Controls for Semiquantitative Analysis (ng/mL) |
|----------------------|--|--|
| 5 | Emit® II Plus Specialty Drug Calibrator/Control Level 2 (5) | Emit® Calibrator/Control Level 0 (0) Emit® II Plus Specialty Drug Calibrators/Controls Level 1 (2.5) Level 2 (5) Level 3 (15) Level 4 (25) |

Qualitative Analysis

Calibrate by running the Emit® II Plus Specialty Drug Calibrator/Control Level 2 in duplicate. Validate the calibration by running controls (see Quality Control). Refer to the Emit® II Plus Specialty Drug Calibrators/Controls Instructions for Use and the specific application sheet for additional information. Recalibrate as indicated by control results.

Semiquantitative Analysis

Prepare a calibration curve by running the appropriate Emit® II Plus Specialty Drug Calibrators/Controls listed in Table 1 in duplicate. Validate the calibration by running controls (see Quality Control). Refer to the Emit® II Plus Specialty Drug Calibrators/Controls Instructions for Use and specific application sheet for additional information. Recalibrate as indicated by control results.

Quality Control

Validate the calibration by assaying the controls listed in table 2. Other material may be used for quality control purposes.

Table 2 — Controls for use in Qualitative or Semiquantitative Analysis

| Cutoff Level | Recommended Controls for Qualitative Analysis | Controls to be used for Semiquantitative Analysis |
|--------------|--|---|
| 5 ng/mL | Emit® Calibrator/Control Level 0 (Negative) Emit® II Plus Specialty Drug Calibrator/Control Level 4 (Positive) or Emit® II Plus Specialty Drug Controls may be used | Emit® II Plus Specialty Drug Control Negative Control Positive |

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of Quality Control (QC) material with known buprenorphine concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Refer to the instrument Operator’s Guide for appropriate instrument checks.

Qualitative Analysis

Validate the calibration by assaying controls. Ensure that the control results relates appropriately to the cutoff calibrator result. That is,

- If the Emit® Calibrator/Control Level 0 (0 ng/mL) or the Emit® II Plus Specialty Drug Control Negative was run; ensure that the result is negative relative to the cutoff calibrator level.
- If the Emit® II Plus Specialty Drug Calibrator/Control Level 4 (25 ng/mL) or Emit® II Plus Specialty Drug Control Positive was run; ensure that the result is positive relative to the cutoff calibrator level.

Once the calibration is validated, run urine specimens.

Semiquantitative Analysis

Validate the calibration curve by assaying commercial controls. Ensure that control results fall within acceptable limits as defined by your laboratory.

Once the calibration curve is validated, run urine specimens.

7 Results

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Siemens has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. User defined modifications are not supported by Siemens as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Siemens Application Sheets or these Instructions for Use.

Qualitative Analysis

The Emit® II Plus Specialty Drug Calibrator/Control level 2 (cutoff), which contains a concentration of 5 ng/mL buprenorphine, is used as a reference for distinguishing “positive” from “negative” specimens.

Positive Results A specimen that gives a change in rate value greater than or equal to the Emit® II Plus Specialty Drug Calibrator/Control cutoff rate value is interpreted as positive. The specimen is presumptive positive buprenorphine.

Negative Results A specimen that gives a change in rate value less than the Emit® II Plus Specialty Drug Calibrator/Control cutoff rate value is interpreted as negative; either the specimen does not contain buprenorphine or buprenorphine is present in concentrations below the cutoff level for this assay.

Semiquantitative Analysis

The semiquantitation of positive results enables the laboratory to determine an appropriate dilution of the specimen for confirmation by GC/MS or LC/MS. Semiquantitation also permits the laboratory to establish quality control procedures and assess control performance. Refer to the Recovery section for the semiquantitative range.

Using the Emit® II Plus Buprenorphine Assay, it is possible to make semiquantitative determinations of buprenorphine. An estimate of relative total drug concentrations may be obtained by running the Emit® Calibrator/Control Level 0 (0 ng/mL) and the appropriate Emit® II Plus Specialty Drug Calibrators/Controls Levels 1 (2.5 ng/mL), 2 (5 ng/mL), 3 (15 ng/mL), 4 (25 ng/mL).

8 Limitations

- The assay is designed for use with human urine only.
- A positive result from the assay indicates the presence of buprenorphine but does not indicate or measure intoxication.
- Boric acid is not recommended as a preservative for urine.
- There is a possibility that substances and/or factors not listed (e.g., technical or procedural errors) may interfere with the test and cause false results.
- Interpretation of results must take into account that urine concentrations can vary extensively with fluid intake and other biological variables.
- Immunoassays that produce a single result in the presence of a drug and its metabolites cannot fully quantitate the concentration of individual components.

9 Expected Values

When the Emit® II Plus Buprenorphine Assay is used as a qualitative assay, the amount of drugs and metabolites detected by the assay in any given specimen cannot be estimated. The assay results distinguish between positive and negative specimens—positive indicating specimens that contain buprenorphine.

When used semiquantitatively, the assay yields approximate, cumulative concentrations of buprenorphine and the metabolites detected by the assay (see Specificity in Section 10).

10 Specific Performance Characteristics

The data that appear in this section were collected from the Viva-E® using the Emit® II Plus Buprenorphine Assay.

Method Comparison

Qualitative and Semiquantitative Results

A total of one-hundred twenty seven (127) unaltered buprenorphine samples were analyzed using the Emit® II Plus Buprenorphine Assay and the reference method Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS). Both methods used a 5 ng/mL cutoff. Thirty-three (33) samples were within +/- 50% of the cutoff by LC/MS/MS.

Sixty-five (65) samples showed positive results by both methods, while fifty-four (54) samples showed negative results by both methods. Seven specimens showed negative results by LC/MS/MS but positive by Emit® II Plus Buprenorphine Assay and one positive results by LC/MS/MS was negative by Emit® II Plus Buprenorphine Assay method. All the discordant samples were within +/- 25% of the cutoff by both the Emit® II Plus Buprenorphine Assay and LC/MS/MS. Data are summarized in Tables 3 and 4.

Table 3 — Qualitative and Semiquantitative Accuracy Summary

| | | LC/MS/MS | | | | |
|--------------------------|-----|---|--|--|--|-----|
| | | LOW NEG <50% below the cutoff (<2.5 ng/mL) | NEG Within 50% below the cutoff (2.5 – 4.9 ng/mL) | POS Within 50% above the cutoff (5.0 – 7.5 ng/mL) | HIGH POS >50% above the cutoff (>7.5 ng/mL) | |
| Qualitative Summary | | | | | | |
| Emit® | POS | 0 | 7 | 16 | 49 | 90% |
| | NEG | 45 | 9 | 1 | 0 | 98% |
| Semiquantitative Summary | | | | | | |
| Emit® | POS | 0 | 7 | 16 | 49 | 90% |
| | NEG | 45 | 9 | 1 | 0 | 98% |

Table 4 — Discordant Result Summary

| Sample # | Emit® II Plus Buprenorphine (ng/mL) | LC/MS/MS | | | Emit® vs LC/MS/MS (Pos/Neg) |
|----------|---|--------------------------|-----------------------------|--|-----------------------------------|
| | | Buprenorphine (ng/mL) | Norbuprenorphine (ng/mL) | Buprenorphine + Norbuprenorphine (ng/mL) | |
| 77 | 5.5 | 0 | 4.60 | 4.60 | +/- |
| 190 | 5.6 | 0 | 3.92 | 3.92 | +/- |
| 193 | 6.1 | 0 | 4.97 | 4.97 | +/- |
| 195 | 5.4 | 0 | 4.06 | 4.06 | +/- |
| 226 | 5.8 | 0 | 4.21 | 4.21 | +/- |
| 250 | 5.1 | 0 | 4.13 | 4.13 | +/- |
| 316 | 5.7 | 0 | 3.86 | 3.86 | +/- |
| 338 | 4.3 | 5.12 | 0 | 5.12 | -/+ |

Precision

Repeatability and within lab precision was determined by assaying urine pools spiked with buprenorphine at nine different levels. The testing sequence for each level consisted of two replicates, twice a day, for twenty days (n = 80). Precision data were calculated according to the Clinical and Laboratory Standards Institute (CLSI) Guideline EP5-A2.7 Results are summarized in Tables 5 and 6.

Table 5 — Precision: Qualitative Analysis

| Urine Pool (ng/mL) | % of Cutoff | # of Results | Repeatability Results | Within-Laboratory Results |
|-----------------------|----------------|--------------|-----------------------------|------------------------------|
| 0 | -100% | 80 | 80 Negative | 80 Negative |
| 2.50 | -50% | 80 | 80 Negative | 80 Negative |
| 3 | -40% | 80 | 80 Negative | 80 Negative |
| 3.75 | -25% | 80 | 80 Negative | 80 Negative |
| 5 | cutoff | 80 | 25 Negative/ 55 Positive | 25 Negative/ 55 Positive |
| 6.25 | +25% | 80 | 80 Positive | 80 Positive |
| 7 | +40% | 80 | 80 Positive | 80 Positive |
| 7.50 | +50% | 80 | 80 Positive | 80 Positive |
| 10 | 100% | 80 | 80 Positive | 80 Positive |

Table 6 — Precision: Semiquantitative Analysis

| Urine Pool (ng/mL) | % of Cutoff | # of Results | Repeatability Results | Within-Laboratory Results |
|-----------------------|----------------|--------------|-----------------------------|------------------------------|
| 0 | -100% | 80 | 80 Negative | 80 Negative |
| 2.50 | -50% | 80 | 80 Negative | 80 Negative |
| 3 | -40% | 80 | 80 Negative | 80 Negative |
| 3.75 | -25% | 80 | 80 Negative | 80 Negative |
| 5 | cutoff | 80 | 25 Negative/ 55 Positive | 25 Negative/ 55 Positive |
| 6.25 | +25% | 80 | 80 Positive | 80 Positive |
| 7 | +40% | 80 | 80 Positive | 80 Positive |
| 7.50 | +50% | 80 | 80 Positive | 80 Positive |
| 10 | 100% | 80 | 80 Positive | 80 Positive |

An additional precision study was performed using urine pools prepared by spiking Buprenorphine into drug-free human urine at one concentration level relative to the 5 ng/mL cutoff: -75% below the cutoff and +75% above the cutoff. The studies were performed on the Viva-E® analyzer. The samples were analyzed in duplicate, 40 times for a total of 80 replicates over one day.

Table 7 — Precision: Qualitative Analysis

| Urine Pool (ng/mL) | % of Cutoff | # of Results | Repeatability Results | Within-Laboratory Results |
|-----------------------|----------------|--------------|-----------------------|------------------------------|
| 1.25 | -75% | 80 | 80 Negative | 80 Negative |
| 8.75 | +75% | 80 | 80 Positive | 80 Positive |

Table 8 — Precision: Semi-Quantitative Analysis

| Urine Pool (ng/mL) | % of Cutoff | # of Results | Repeatability Results | Within-Laboratory Results |
|-----------------------|----------------|--------------|-----------------------|------------------------------|
| 1.25 | -75% | 80 | 80 Negative | 80 Negative |
| 8.75 | +75% | 80 | 80 Positive | 80 Positive |

Recovery

Recovery of buprenorphine samples were prepared by spiking known levels of buprenorphine into negative urine pools. Each spiked sample was analyzed in replicates of five using the Emit® II Plus Buprenorphine Assay. Results are shown in Table 9.

Table 9 — Results of Recovery

| Expected Concentration (ng/mL) | Mean Buprenorphine Concentration by Emit® II Plus Buprenorphine Assay (ng/mL) | % Recovery |
|-----------------------------------|---|------------|
| 2 | 2.1 | 105.0 |
| 3 | 3.1 | 103.3 |
| 4 | 3.9 | 97.5 |
| 5 | 5.0 | 100.0 |
| 8 | 7.7 | 96.3 |
| 12 | 11.1 | 92.5 |
| 18 | 17.7 | 98.3 |
| 22 | 21.0 | 95.5 |
| 25 | 23.9 | 95.6 |

Specificity

Buprenorphine Metabolites

The buprenorphine metabolites norbuprenorphine, buprenorphine glucuronide and norbuprenorphine glucuronide were spiked into a negative urine pool at the levels shown in Table 10 and run at n=5 replicates. The samples were assayed and the mean recovery results were determined.

Table 10 — Buprenorphine Metabolite Recovery

| Compound | Concentration Tested (ng/mL) | Mean Recovery (ng/mL) | % Cross-reactivity |
|------------------------------|---------------------------------|--------------------------|--------------------|
| Norbuprenorphine | 5 | 4.6 | 92.6 |
| Buprenorphine Glucuronide | 1000 | 0.9 | 0.1 |
| Norbuprenorphine Glucuronide | 1000 | 1.2 | 0.1 |

Structurally Related Compounds

The cross-reactivity of the Emit® II Plus Buprenorphine assay was evaluated against structurally related drugs. The structurally related drugs listed below were spiked into negative urine pools and % cross-reactivity was calculated. The compounds at stated concentrations produced a negative response relative to a 5 ng/mL cutoff.

Table 11 — Structurally Related Compounds

| Compound | Concentration Tested (ng/mL) | % Cross-Reactivity |
|------------------------|---------------------------------|--------------------|
| 6-Acetylcodeine | 100,000 | <0.01 |
| 6-Acetylmorphine | 100,000 | <0.01 |
| Codeine | 100,000 | <0.01 |
| Dextromethorphan | 100,000 | <0.01 |
| Dihydrocodeine | 100,000 | <0.01 |
| Ethyl Morphine | 100,000 | <0.01 |
| Heroin | 100,000 | <0.01 |
| Hydrocodone | 100,000 | <0.01 |
| Hydromorphone | 100,000 | <0.01 |
| Levorphanol | 100,000 | <0.01 |
| Morphine | 100,000 | <0.01 |
| Morphine 3-glucuronide | 100,000 | <0.01 |
| Morphine 6-glucuronide | 100,000 | <0.01 |
| Nalorphine | 100,000 | <0.01 |
| Naloxone | 100,000 | <0.01 |
| Naltrexone | 100,000 | <0.01 |
| Norcodeine | 100,000 | <0.01 |
| Normorphine | 100,000 | <0.01 |
| Noroxycodone | 100,000 | <0.01 |
| Noroxymorphone | 100,000 | <0.01 |
| Oxycodone | 100,000 | <0.01 |
| Oxymorphone | 100,000 | <0.01 |

Structurally Unrelated Compounds

The interference of structurally unrelated compounds and common over the counter drugs was evaluated qualitatively and semiquantitatively at the concentration listed below. The compounds were spiked into two levels of controls at +/- 40% (3 ng/mL and 7 ng/mL) of the cutoff concentration. The protocol used follows the CLSI Guideline, EP7-A2.⁸

At the stated concentration, the sample did not give a false response relative to the 5 ng/mL cutoff.

Table 12 — Structurally Unrelated Compounds

| Compound | Concentration | | |
|-----------------------------|-------------------|-----------------------------|------------------------------|
| | Tested (µg/mL) | -40% of Cutoff (3 ng/mL) | + 40% of Cutoff (7 ng/mL) |
| 10, 11-Dihydrocarbamazepine | 85 | Negative | Positive |
| Acetaminophen | 1000 | Negative | Positive |
| Acetylsalicylic Acid | 1500 | Negative | Positive |
| Amitriptyline | 100 | Negative | Positive |
| Amoxicillin | 500 | Negative | Positive |
| AZT (Zidovudine) | 2000 | Negative | Positive |
| Benzoylcegonine | 1000 | Negative | Positive |
| Brompheniramine | 75 | Negative | Positive |
| Caffeine | 1000 | Negative | Positive |
| Captopril | 500 | Negative | Positive |
| Chlordiazepoxide | 100 | Negative | Positive |
| Chlorpromazine | 10 | Negative | Positive |
| Cimetidine | 1000 | Negative | Positive |
| Clomipramine | 2.5 | Negative | Positive |
| Clonidine | 1000 | Negative | Positive |
| Cyclobenzaprine | 125 | Negative | Positive |
| d-Amphetamine | 700 | Negative | Positive |
| Desipramine | 800 | Negative | Positive |
| Diazepam | 100 | Negative | Positive |
| Digoxin | 0.01 | Negative | Positive |
| Diphenhydramine | 1000 | Negative | Positive |
| d-Methamphetamine | 500 | Negative | Positive |
| Doxepine | 100 | Negative | Positive |
| EDDP | 1000 | Negative | Positive |
| EMDP | 100 | Negative | Positive |
| Enalapril | 500 | Negative | Positive |
| Fluoxetine | 500 | Negative | Positive |
| Glutethimide | 500 | Negative | Positive |
| Haloperidol | 100 | Negative | Positive |
| Hydroxyzine | 500 | Negative | Positive |
| Ibuprophen | 1000 | Negative | Positive |
| Imipramine | 200 | Negative | Positive |
| Ketamine | 100 | Negative | Positive |
| Ketorolac Tromethamine | 400 | Negative | Positive |
| LAAM (L-α-Acetylmethadol) | 25 | Negative | Positive |
| L-Cotinine | 100 | Negative | Positive |
| Levofloxacin | 100 | Negative | Positive |
| Levothyroxine (L-Thyroxine) | 50 | Negative | Positive |
| Lidocaine | 1000 | Negative | Positive |
| Lormetazepam | 1.0 | Negative | Positive |
| LSD | 10 | Negative | Positive |
| MDMA (Ecstasy) | 1000 | Negative | Positive |
| Meperidine | 800 | Negative | Positive |
| Methadone | 500 | Negative | Positive |
| Methaqualone | 600 | Negative | Positive |
| NAPA | 400 | Negative | Positive |
| Naproxen | 1000 | Negative | Positive |
| Nicotinic Acid | 500 | Negative | Positive |
| Nifedipine | 500 | Negative | Positive |
| Nordiazepam | 100 | Negative | Positive |
| Nortryptiline | 250 | Negative | Positive |
| Oxazepam | 300 | Negative | Positive |
| Perphenazine | 150 | Negative | Positive |
| Phencyclidine | 900 | Negative | Positive |
| Phenobarbital | 500 | Negative | Positive |
| Phenelzine | 100 | Negative | Positive |
| Phenytoin | 1000 | Negative | Positive |
| Procainamide | 1000 | Negative | Positive |
| Procyclidine | 800 | Negative | Positive |
| Promethazine | 100 | Negative | Positive |
| Propoxyphene | 1000 | Negative | Positive |
| Protriptyline | 200 | Negative | Positive |
| Pseudoephedrine | 1000 | Negative | Positive |
| Quinacrine | 900 | Negative | Positive |
| Ranitidine | 1000 | Negative | Positive |
| Ritalin | 1000 | Negative | Positive |
| Salicylic Acid | 500 | Negative | Positive |

Table 12 — Structurally Unrelated Compounds (continued)

| Compound | Concentration | | |
|-----------------------|-------------------|-----------------------------|------------------------------|
| | Tested (µg/mL) | -40% of Cutoff (3 ng/mL) | + 40% of Cutoff (7 ng/mL) |
| Scopolamine | 500 | Negative | Positive |
| Secobarbital | 1000 | Negative | Positive |
| Tapentadol | 100 | Negative | Positive |
| THC | 100 | Negative | Positive |
| Thioridazine | 100 | Negative | Positive |
| Tramadol | 1000 | Negative | Positive |
| Trazodone | 5 | Negative | Positive |
| Trimethoprim | 1000 | Negative | Positive |
| Triprolidine (zymine) | 50 | Negative | Positive |
| Tyramine | 100 | Negative | Positive |
| Verapamil | 500 | Negative | Positive |
| Zolpidem | 100 | Negative | Positive |

Endogenous Substances

The endogenous substances were evaluated qualitatively and semiquantitatively at the concentrations listed below. The compounds were spiked into two levels of controls at +/- 40% (3 ng/mL and 7 ng/mL) of the cutoff concentration. The protocol used follows the CLSI Guideline, EP7-A2.⁸

At the stated concentration, the sample did not give a false response relative to the 5 ng/mL cutoff.

Table 13 — Endogenous Substances: Semiquantitative Results

| Substance | Concentration Tested | -40% of Cutoff (3 ng/mL) | + 40% of Cutoff (7 ng/mL) |
|------------------------|-------------------------|-----------------------------|------------------------------|
| Acetone | 1.0 g/dL | Negative | Positive |
| Ascorbic Acid | 1.5 g/dL | Negative | Positive |
| Conjugated Bilirubin | 2.0 mg/dL | Negative | Positive |
| Unconjugated Bilirubin | 2.0 mg/dL | Negative | Positive |
| Creatinine | 0.5 g/dL | Negative | Positive |
| Ethanol | 1.0 g/dL | Negative | Positive |
| IgG | 0.5 g/dL | Negative | Positive |
| Glucose | 2.0 g/dL | Negative | Positive |
| Hemoglobin | 115 mg/dL | Negative | Positive |
| Human Serum Albumin | 0.5 g/dL | Negative | Positive |
| Oxalic Acid | 0.1 g/dL | Negative | Positive |
| Riboflavin | 7.5 mg/dL | Negative | Positive |
| Sodium Chloride | 6.0 g/dL | Negative | Positive |
| Urea | 6.0 g/dL | Negative | Positive |
| Sodium Azide | 1% | Negative | Positive |
| Sodium fluoride | 1% | Negative | Positive |
| Galactose | 1.0 g/dL | Negative | Positive |


















Specific Gravity and pH

Negative urine pools with specific gravity values ranging from 1.002–1.035 and pH values ranging from 3.0–11.0 were tested in the presence of two levels of controls at +/- 40% (3 ng/mL and 7 ng/mL) of the cutoff concentration and no interference was observed.

12 Bibliography

1. Kronstrand R, Selden TG, and Josefsson M. Analysis of Buprenorphine, Norbuprenorphine and Their Glucuronides in Urine by Liquid Chromatography-Mass Spectrometry, *JAT*. 2003;27:464-470.
2. Fuller DC. A Simple Gas Chromatography-Mass Spectrometry Procedure for the Simultaneous Determination of Buprenorphine and Norbuprenorphine in Human Urine, *JAT*. 2008;32:626-630.
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, Biomedical Publications, 8th Ed. 2008:730-35.
4. Cone EJ. et al, The Metabolism and Excretion of Buprenorphine in Humans, *Drug Metabolism and Disposition*. 1984;12:577:581.
5. ARUP Consult Physicians Guide <http://ltd.aruplab.com/Tests/Pub/2010092>
6. Huang W, Moody DE, McCance-Katz EF. The In Vivo Glucuronidation of Buprenorphine and Norbuprenorphine Determined by Liquid Chromatography-Electrospray Ionization-Tandem Mass Spectrometry, *Therapeutic Drug Monitoring*. 2006;28(2):245-251.
7. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
8. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP7-A2 [ISBN 1-56238-584-4]. CLSI, Wayne, PA USA, 2005.

13 Symbols Key

| | |
|--|--|
|  | Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar |
|  | Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad |
|  | Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote |
|  | Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo |
|  | Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso |
|  | Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante |
|  | Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea |
|  | Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos |
|  | In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro |
|  | Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Límite de temperatura |
|  | Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso |
|  | Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril |
|  | CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE |
|  | Contents / Inhalt / Contenu / Contenido / Contenido |
|  | Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución |
|  | Level / Konzentration / Niveau / Livello / Nivel |
|  | Prescription Device (US Only) / Verschreibungspflichtiges Medizinprodukt (nur USA) / Dispositif soumis à prescription (États-Unis uniquement) / Dispositivo su prescrizione (solo USA) / Dispositivo de prescripción médica (solo EE. UU.) |

2015-10_EFIGS

For technical assistance, call Siemens Healthcare Diagnostics:

1-800-227-8994 in the USA

1-800-264-0083 in Canada

Outside the USA and Canada, call your local Siemens representative.

Notice: Adulteration of reagents, use of instruments without appropriate capabilities, or other failure to follow instructions as set forth in this labeling can affect performance characteristics and stated or implied label claims.

U.S. Patent Numbers 7,220,842 B2 and 7,863,427 B2

 Syva® and Emit® are trademarks of Siemens Healthcare Diagnostics.

© 2015 Siemens Healthcare Diagnostics

All rights reserved.

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare



Printed in USA
2019-05
10872256_B