# IMMULITE®

## AFP

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**SIEMENS**

# IMMULITE®/IMMULITE 1000 AFP

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of alpha-fetoprotein (AFP) in either of two contexts: (a) serial measurements in human serum to aid in the management of patients with nonseminomatous testicular cancer; or (b) measurements in maternal serum and amniotic fluid during gestational weeks 15 through 20 — used in conjunction with ultrasonography or amniography — to aid in detection of fetal open neural tube defects.

Catalog Number: **LKAP1** (100 tests)

Test Code: **AF** Color: **Light Gray**

**Caution:** In the United States, federal law restricts this device to sale by or on the order of a physician.

The concentration of AFP in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. **The results reported by the laboratory to the physician must include the identity of the assay used. Values obtained with different AFP assays cannot be used interchangeably.** Before changing assays, the laboratory must: (a) for cancer management — confirm baseline values for patients being serially monitored; (b) for prenatal testing — establish a range of normal values for the new assay based on normal sera and amniotic fluids from pregnant women with confirmed gestational age.

## Summary and Explanation

Alpha-fetoprotein (AFP) is a single-chain glycoprotein with a molecular mass of approximately 70,000 daltons. AFP shares considerable sequence homology with albumin, and is produced by the fetus primarily in cells of the yolk sac, gastrointestinal tract and liver. AFP appears

as a major serum protein in the fetus, but its concentration decreases rapidly toward birth.<sup>1,2,3</sup> The reappearance of elevated AFP concentrations in adult serum has been observed not only during pregnancy, but also in conjunction with several benign and malignant diseases.

### Testicular Cancer

Elevated levels of AFP have been observed not only in patients with nonseminomatous testicular cancer, but also in patients with other malignancies such as hepatocellular carcinoma, ovarian cancer, gastrointestinal cancer and pulmonary cancer.<sup>8-15</sup> Serum AFP is frequently elevated in benign hepatic conditions such as acute viral hepatitis, chronic active hepatitis and cirrhosis. Conditions of pregnancy, ataxia telangiectasia and hereditary tyrosinemia have also presented with elevated concentrations of AFP.<sup>8-15</sup>

Seminomas, in pure form, do not present with elevated concentrations of AFP. However, elevated concentrations of serum AFP have been observed in patients diagnosed with seminomatous testicular cancer accompanied by nonseminomatous metastases.<sup>9,16,18,19</sup> During chemotherapy, patients with advanced seminoma and hepatic dysfunction have also presented with elevated serum AFP concentrations.<sup>20</sup> The interpretation of elevated AFP concentrations in patients with seminoma requires special consideration and should assist the clinician in the selection of appropriate therapy.<sup>8,15,21</sup>

The clinical utility of AFP measurement as an aid in the management of patients with nonseminomatous testicular cancer is well documented.<sup>9,16,17,18,22</sup> AFP measurement has found clinical application as an aid in assessing the extent of disease.<sup>18,22-26</sup>

Serial measurements of serum AFP have been shown to reflect the effectiveness of therapeutic regimens in patients with nonseminomatous testicular tumors.<sup>9,15,17,26,27</sup> Post-surgical determinations of AFP are particularly valuable. The presence of residual tumor is strongly suggested if post-operative AFP concentrations fail to return to normal.<sup>9,15,28,29</sup> The accurate interpretation of post-surgical changes in AFP

concentration requires consideration of its metabolic decay rate.<sup>21,22,24,25</sup> When utilizing AFP for monitoring therapy or disease recurrence during chemotherapy, it should be noted that levels often fall rapidly during chemotherapy, reaching normal levels while tumor masses are still evident.<sup>17,21</sup> In such instances, completion of the planned therapy has been recommended.<sup>21</sup>

Following therapy or surgery, serial measurements of AFP have also proved clinically useful when monitoring for progression or recurrence of disease in patients with nonseminomatous testicular cancer. It has been reported that AFP levels frequently rise during disease progression and fall during disease remission.<sup>9,17,18</sup> Elevated AFP levels have frequently been observed to accompany tumor recurrence before progressive disease is clinically evident.<sup>9,18</sup>

### **Fetal Open Neural Tube Defects**

AFP is detectable not only in fetal serum, but also in amniotic fluid and maternal serum. A concentration gradient exists such that when the fetal serum AFP level is 2000 kIU/mL, the amniotic fluid AFP (AFAFP) level is 20 kIU/mL, and the maternal serum AFP (MSAFP) is 0.02 kIU/mL. In normal pregnancy, the fetal serum AFP concentration peaks at 14 weeks gestation.<sup>34</sup> The AFAFP concentration peaks at about 12 weeks and the MSAFP peaks at approximately 28–32 weeks gestation.<sup>36</sup> The fall in AFAFP concentration reflects the fall in fetal serum AFP concentration which results from increased fetal size and fluid volume.<sup>34</sup> Elevated levels of MSAFP and AFAFP may occur most often due to multiple pregnancy and due to incorrect gestational age.

Measurement of AFP concentrations is clinically valuable in screening for open NTDs and other fetal abnormalities;<sup>35</sup> pregnancies associated with open NTDs present with elevated AFP levels. Excess AFP gains access to amniotic fluid, and to a lesser extent to the maternal serum, by transudation across the exposed surface of the fetus or across damaged glomeruli.<sup>35,37</sup> These conditions are found in open NTDs including open spina bifida and anencephaly, omphalocele, and congenital nephrosis.<sup>32,38</sup> Additional causes of elevated AFP concentrations

including both maternal and fetal sources are impending spontaneous abortion, fetal distress or death, oligohydramnios, toxemia, gastroschisis, Meckel's syndrome, sacrococcygeal teratoma, Turner's syndrome and maternal hepatic and oncologic disorders.<sup>35</sup>

Recommended protocols for open NTD screening have been published.<sup>33,35</sup> The cutoff levels for maternal serum and amniotic fluid can be chosen to optimize the needs of the populations being tested based upon varying prevalence of open NTDs. Cutoffs commonly utilize multiples of the median of 2.0 or 2.5 for MSAFP and AFAFP testing. The optimal time for screening MSAFP is between the 16th and 18th weeks of pregnancy, although screening is still effective before or after this period. Elevated AFP concentrations may be subjected to a repeat sampling and analysis to exclude transient rises.

More commonly, ultrasonography is employed to rule out multiple pregnancies and to confirm gestational age. Ultrasonography may also identify signs of open NTDs, particularly anencephaly which is a large, easily visualized lesion. If correction for gestational age or multiple pregnancy does not result in an AFP concentration within the normal range, then diagnostic ultrasonography and/or amniotic fluid sampling is indicated. The greatest diagnostic power can be achieved by combining biochemical analysis of amniotic fluid and diagnostic ultrasonography in cases of a positive MSAFP screen.<sup>35</sup>

Elevated MSAFP results are not diagnostic for NTDs and should not be considered a cause for termination of pregnancy. An overlap exists in the distributions of AFP concentrations from pregnancies with and without open NTDs. Closed NTDs, for example, are not usually associated with increased MSAFP or AFAFP concentrations. Thus, further testing is required to define fetal status. In light of these considerations and the multiple causes for elevated AFP concentrations, all clinical information should be evaluated and confirmatory tests performed wherever possible before reaching a diagnosis.

AFP can be measured by several immunologic techniques, depending on the degree of sensitivity desired. Radial immunodiffusion, countercurrent immunoelectrophoresis, and rocket immunoelectrophoresis are three techniques well suited for research applications. Enzyme-linked immunosorbent assays and radioimmunoassays of both competitive and non-competitive designs have been successfully employed clinically both for maternal serum and amniotic fluid measurements.

**Note:** IMMULITE AFP Physician Brochure (Cat. #ZS1105) and Patient Brochure (Cat. #ZS1106), explaining the use of AFP prenatal testing to aid in the detection of fetal open NTD are available by calling Siemens Healthcare Diagnostics Customer Services 1-800-372-1782 or your National Distributor.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 AFP is a solid-phase, two-site sequential chemiluminescent immunometric assay.

**Incubation Cycles:** 2 × 30 minutes

### Specimen Collection

**Serum:** Collect blood by venipuncture<sup>31</sup> into plain tubes and separate the serum from the cells as soon as possible. Specimens must be obtained before amniocentesis to obtain a valid specimen.

Lipemic samples should be clarified with the help of an ultracentrifuge before assay of serum.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 AFP has not been tested with all possible variations of tube types.

**Amniotic Fluid:** Collect amniotic fluid by amniocentesis into plain tubes. Samples should be obtained by aseptic transabdominal amniocentesis performed by an experienced obstetrician during the second trimester of pregnancy in women with confirmed gestational age. Centrifuge the specimen, retaining a portion of the clear supernatant. Inspect both supernatant and sediment for signs of blood or hemoglobin, as contamination by even trace amounts of fetal material will raise the apparent AFP concentration of the sample, rendering it unsuitable for analysis. The origin of the fetal material should be determined by a test for fetal hemoglobin. If fetal contamination has occurred and the AFP concentration is elevated, an additional specimen should be obtained after 7 to 10 days for evaluation. Amniotic fluid contamination by maternal serum may reflect accurate AFP levels provided the degree of contamination is not sufficient to dilute the sample. Henceforth in this package insert, *amniotic fluid* refers to the clear supernatant obtained from amniotic fluid by centrifugation.

**Timing: It is essential to know the gestational age to evaluate AFP results. The recommended time for collection is 16 to 18 weeks for serum, 16 to 20 weeks for amniotic fluid. Serum samples must be collected before amniocentesis since this procedure may lead to spuriously elevated maternal serum levels persisting for 2 to 3 weeks.**

### Volume Required

**Serum:** 10 µL (sample cup must contain at least 100 µL more than the total volume required).

**Amniotic Fluid:** 10 µL of prediluted amniotic fluid specimen (sample cup must contain at least 100 µL more than the total volume required).

All amniotic fluid samples must first be diluted 1-in-101 using AFP Sample Diluent, e.g. by adding 10 µL of amniotic fluid to 1.0 mL of AFP Sample Diluent. Results of the diluted sample must be multiplied by the dilution factor to obtain the final AFP concentration.

### Storage

**Serum:** 3 days at 2–8°C. Freeze at –20°C if not assayed within 3 days.

**Amniotic Fluid:** Amniotic fluid samples should be stored at –20°C. Aliquot if necessary to avoid repeated freezing and thawing. Allow the sample to come to room temperature (15–28°C) before assay, and mix by *gentle* swirling or inversion. Do not attempt to thaw specimens by heating them in a waterbath. If specimens are to be mailed, samples should be packed in dry ice if the time in transit exceeds 72 hours, or if elevated temperatures are a concern, as in warm climates or during the summer. If a repeat analysis is required, the original type of specimen should be taken to maintain consistency of results.

### Warnings and Precautions

For *in vitro* diagnostic use.

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight (see insert).

**Water:** Use distilled or deionized water.

### Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

#### AFP Test Units (LAP1)

Each barcode-labeled unit contains one bead coated with murine monoclonal anti-AFP. Stable at 2–8°C until expiration date.

**LKAP1:** 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

#### AFP Reagent Wedges (LAPA, LAPB)

Two barcode-labeled wedges.

**LAPA:** one wedge (7.5 mL) of a protein buffer/nonhuman serum matrix.

**LAPB:** one wedge (7.5 mL) of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal rabbit anti-AFP in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.

**LKAP1:** 1 set

#### AFP Adjustors (LAPL, LAPH)

Two vials (Low and High), 2 mL each, of AFP in a bovine serum matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LKAP1:** 1 set

### Kit Components Supplied Separately

#### AFP Sample Diluent (LAPZ)

For the manual dilution of high serum samples and amniotic fluid samples. 25 mL of an AFP-free bovine serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**Analysis of amniotic fluid requires a 1-in-101 dilution of the sample (manual dilution with AFP Sample Diluent).**

**Results of the diluted sample must be multiplied by the dilution factor to obtain the final AFP concentration.**

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)  
**LSCC:** Sample Cup Caps (optional)

Also Required  
 Sample transfer pipets; distilled or deionized water

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE and IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:**  
 2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of AFP.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

### AFP Values in Testicular Cancer Patients

In a study involving two clinical sites, 119 serum samples from men in apparent good health (median age: 61, central 95%: 27 to 79 years) were processed by the IMMULITE AFP assay. The results ranged from 0.5 to 5.5 IU/mL, with a median of 1.6 IU/mL and a 99th percentile of 5 IU/mL.

The study also included men with testicular cancer; patients with other malignancies (of liver, bladder, kidney, pancreas, lung, prostate and colon); patients with nonmalignant conditions

(such as cirrhosis, hepatitis B and C, ulcerative colitis, emphysema, colon and rectal polyps); and a few women in apparent good health. The distribution of IMMULITE AFP results is tabulated below (with the total number for each group in parentheses).

<b>IU/mL:</b>	<b>&lt; 5</b>	<b>5–15</b>	<b>15–100</b>	<b>&gt; 100</b>
<b>Males</b>				
Healthy Males (119)				
	118	1	—	—
Seminomatous Testicular Cancer (6)				
	6	—	—	—
Nonseminomatous Testicular Cancer (60)				
	14	8	15	23
Liver Cancer (10)				
	3	—	2	5
Other Malignant Diseases (40)				
	36	1	—	3
Cirrhosis (4)				
	3	1	—	—
Hepatitis (24)				
	19	4	1	—
Other Nonmalignant Diseases (6)				
	5	—	—	1
<b>Females</b>				
Healthy Females (29)				
	29	—	—	—
Malignant Diseases (20)				
	18	—	1	1
Nonmalignant Diseases (16)				
	15	—	1	—

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Patients with nonseminomatous testicular cancer can be expected to have a distribution of AFP values both within and above the reference range for apparently healthy adult male subjects. In pure form, seminomas do not present with elevated serum AFP levels. However, elevated AFP levels have been observed in patients diagnosed with seminomas accompanied by metastases of nonseminomatous testicular cancer.<sup>8</sup>

A significant increase of AFP levels in patients considered free of metastatic tumor may indicate the development of metastasis. Elevated levels after surgery may indicate incomplete removal of the tumor or the presence of metastases.

Elevated levels of serum AFP are associated with benign liver conditions such as hepatitis and cirrhosis. Most (95%) of the patients with these benign diseases have AFP levels lower than 200 ng/mL (165 IU/mL).<sup>8-15</sup>

### AFP Values in Maternal Serum and Amniotic Fluid

Due to potential variation in testing at different laboratories, it is recommended that a particular testing center determine its own set of median AFP values for weeks 15 to 20 of gestation, measured in the population to be screened. Cutoff values commonly utilize multiples of the medians (MoM) of 2.0 or 2.5 for maternal serum and amniotic fluid testing. Each AFP test result can then be expressed as a multiple of the unaffected population median value. This is obtained by dividing the AFP value by the median value for its corresponding gestational week. Gestational weeks are defined as completed gestational weeks; e.g., 16 weeks, 6 days would be considered the 16th week. It has been recommended that median and MoM values determined for each gestational week be based upon at least 100 maternal sera and 50 amniotic fluids from unaffected singleton pregnancies with confirmed gestational age.

Provided below are medians for *maternal serum* samples, generated using the IMMULITE and IMMULITE 2000 AFP assays, calculated by a weighted log-linear regression from data collected from unaffected, singleton pregnancies at three clinical sites in the United States:

Gestational Week	No. of Specimens	Medians IU/mL*	Multiples of Regressed Medians (IU/mL)		
			2.0	2.5	3.0
15	370	24.9	49.8	62.3	74.7
16	605	28.5	57.0	71.3	85.5
17	569	32.6	65.2	81.5	97.8
18	431	37.2	74.4	93.0	111.6
19	221	42.5	85.0	106.3	127.5
20	91	48.6	97.2	121.5	145.8

\*Regressed

Provided below are medians for *amniotic fluid* samples, calculated by a weighted log-linear regression from data collected from unaffected, singleton pregnancies at two clinical sites in the United States:

Gestational Week	No. of Specimens	Medians kIU/mL*	Multiples of Regressed Medians (kIU/mL)		
			2.0	2.5	3.0
15	76	13.0	26.0	32.5	39.0
16	89	10.7	21.4	26.8	32.1
17	53	8.73	17.5	21.8	26.2
18	54	7.14	14.3	17.9	21.4
19	46	5.84	11.7	14.6	17.5
20	23	4.78	9.56	12.0	14.3

\*Regressed

### Limitations

**Diagnosis:** The occurrence of elevated serum AFP levels in conditions other than nonseminomatous testicular cancer precludes the use of AFP measurements in the diagnosis of nonseminomatous testicular cancer.

**Screening:** AFP measurements cannot be recommended as a screening procedure to detect cancer in the general population. Elevated concentrations of serum AFP have been observed not only in patients with nonseminomatous testicular cancer but also in malignant conditions such as hepatocellular carcinoma, ovarian cancer, and gastrointestinal and pulmonary cancer. Benign hepatic conditions such as acute viral hepatitis, chronic active hepatitis and cirrhosis may present with elevated concentrations of serum AFP. Elevated AFP concentrations have also been observed in pregnancy, ataxia telangiectasia and hereditary tyrosinemia.

**Prenatal Testing: A reliable AFP evaluation for prenatal testing requires precise determination of the gestational age. Underestimation of the gestational age may lead to false positive determination, while overestimation of gestational age may result in a false negative interpretation. When gestational age is uncertain, confirmation with ultrasonography is indicated.**

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing

interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See tables and graphs for data representative of the assay's performance. All data, except comparisons between the IMMULITE and IMMULITE 1000 analyzers, were generated on the IMMULITE analyzer. Results are expressed in IU/mL. Unless otherwise noted, all were generated on serum samples collected from testicular cancer patients.

### Conversion Factor:

IU/mL × 1.21 → ng/mL

**Calibration Range:** up to 300 IU/mL (363 ng/mL) (WHO 1st IS 72/225)

**Analytical Sensitivity:** 0.2 IU/mL (0.24 ng/mL)

**High-dose Hook:** None up to 450,000 IU/mL

**Precision:** Seven samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Serum and amniotic fluid samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Serum samples spiked 1-in-20 with three AFP solutions (9.8, 31 and 67 IU/mL) were assayed. Amniotic samples spiked 1-in-20 with three high amniotic fluid samples (10,000, 20,000 and 36,000 IU/mL) were also assayed. (See "Recovery" tables for representative data.)

**Specificity:** The assay is highly specific for AFP. (See "Specificity" table.)

**Bilirubin (unconjugated):** A small but (by *t*-test) statistically significant effect. (See "Bilirubin" table.)

**Biotin:** Specimens that contain biotin at a concentration of 3500 ng/mL demonstrate a less than or equal to 10% change in results.

**Hemolysis:** Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

## Method Comparison – Testicular

**Cancer Studies:** The assay was compared to a legally marketed assay (Kit C), at two clinical sites. At the first (in the northwestern United States), a total of 264 specimens were evaluated, including samples from male patients with **nonseminomatous testicular cancer** and other malignant and nonmalignant conditions, as well as samples from a few female patients. The results were compared qualitatively, relative to cutoffs based on the stated 99th percentiles for healthy males — namely, 5 IU/mL for IMMULITE, and 8.9 ng/mL (7.36 IU/mL) for Kit C.

Kit C	IMMULITE AFP		Relative Sensitivity (95% CI)	Relative Specificity (95% CI)
	Positive	Negative		
Positive	74	4	94.9%	97.3%
Negative	5	181	(87.4–98.6%)	(93.8–99.1%)

Agreement: 96.6%

At the second site (in the southern United States), a total of 213 specimens were studied, including samples from male patients with seminomatous and nonseminomatous **testicular cancer**, and other malignant conditions, as well as samples from a few female patients.

Kit C	IMMULITE AFP		Relative Sensitivity (95% CI)	Relative Specificity (95%CI)
	Positive	Negative		
Positive	62	3	95.4%	97.3%
Negative	4	144	(87.1–99.0%)	(93.2–99.3%)

Agreement: 96.7%

The 424 results from the two sites which were within range by both assays were also compared by linear regression analysis:

(IML) = 0.83 (Kit C) – 0.17 IU/mL  
 r = 0.99  
 n = 424

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.81	–0.65
Upper	0.84	0.31

### Method Comparison – Neural Tube

**Defect Studies:** In two separate clinical studies conducted in the United States, IMMULITE AFP results were compared to two legally marketed assays (Kit A and Kit B) in a linear regression for **maternal serum** samples, in the range from nondetectable to 300 IU/mL.

(IML) = 0.83 (Kit B) + 6.51 IU/mL  
 r = 0.97  
 n = 1006

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.82	5.93
Upper	0.85	7.09

(IML) = 0.87 (Kit A) + 2.72 IU/mL  
 r = 0.96  
 n = 346

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.84	1.64
Upper	0.89	3.79

In one of the studies above, IMMULITE AFP results were also compared to Kit B in a linear regression for **amniotic fluid** samples, in the range from nondetectable to 173 kIU/mL\*.

(IML) = 0.78 (Kit B) + 0.17 kIU/mL  
 r = 0.99  
 n = 185

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.76	–0.19
Upper	0.79	0.53

\* Amniotic fluid samples were diluted 1-in-101 off-line before they were tested by the IMMULITE instrument.

The assay was also compared to IMMULITE 2000 AFP on amniotic fluid samples, in the range from approximately 3 to 20 kIU/mL\*. (See “Method Comparison 1” graph.)

By linear regression:

(IML 2000) = 1.03 (IML) + 0.52 kIU/mL  
 r = 0.96  
 n = 46

Means:  
 10.0 kIU/mL (IML)  
 10.8 kIU/mL (IML 2000)

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.93	–0.50
Upper	1.12	1.54

\* Amniotic fluid samples were diluted 1-in-101 off-line before they were tested by the IMMULITE instrument.

The assay was also compared to IMMULITE 2000 AFP on maternal serum samples, in the range from approximately 10 to 120 IU/mL. By linear regression:

(IML 2000) = 1.01 (IML) + 0.154 IU/mL  
 r = 0.982  
 n = 346

Means:  
 33.8 IU/mL (IML)  
 34.3 IU/mL (IML 2000)

95% Confidence Interval (CI)	Slope	Intercept
Lower	0.99	–0.60
Upper	1.03	0.91

### Clinical Sensitivity for Maternal Serum, n = 13:

Gestational Week	% > 2.0 MoM	% > 2.5 MoM	% > 3.0 MoM
15–20	92.3%	69.2%	69.2%
95% CI for All Samples	64.0%–99.8%	38.6%–90.9%	38.6%–90.9%

### Clinical Specificity for Maternal Serum:

Gestational Week	n	% ≤ 2.0 MoM	% ≤ 2.5 MoM	% ≤ 3.0 MoM
15	173	96.0%	98.8%	98.8%
16	411	98.1%	99.3%	99.5%
17	372	96.5%	99.7%	100%
18	204	95.1%	99.0%	100%
19	108	94.4%	99.1%	100%
20	50	100%	100%	100%
15–20	1318	96.7%	99.3%	99.7%
95% CI for All Samples		95.5%–97.6%	98.7%–99.7%	99.2%–99.9%

### Clinical Sensitivity for Amniotic Fluid, n = 10:

Gestational Week	% > 2.0 MoM	% > 2.5 MoM	% > 3.0 MoM
15–20	90.0%	90.0%	90.0%
95% CI for All Samples	55.5%–99.8%	55.5%–99.8%	55.5%–99.8%

### Clinical Specificity for Amniotic Fluid:

Gestational Week	n	% ≤ 2.0 MoM	% ≤ 2.5 MoM	% ≤ 3.0 MoM
15	23	100%	100%	100%
16	39	97.4%	100%	100%
17	25	100%	100%	100%
18	34	94.1%	97.1%	100%
19	33	100%	100%	100%
20	13	100%	100%	100%
15–20	167	98.2%	99.4%	100%
95% CI for All Samples		94.8%–99.6%	96.7%–100%	97.8%–100%

### Method Comparison with Maternal Serum – IMMULITE 1000 vs. IMMULITE:

A total of 110 maternal serum specimens were tested by the IMMULITE/ IMMULITE 1000 AFP assays on the two systems. (See “Method Comparison 2” graph.) By linear regression:

$$(IML\ 1000) = 0.99 (IML) + 0.13\ IU/mL$$

$$r = 0.996$$

Means:

2.88 IU/mL (IML 1000)

2.77 IU/mL (IML)

### Method Comparison with Amniotic Fluid – IMMULITE 1000 vs. IMMULITE:

A total of 91 amniotic fluid specimens were tested by the IMMULITE/ IMMULITE 1000 AFP assays on the two systems. (See “Method Comparison 3” graph.) By linear regression:

$$(IML\ 1000) = 1.00 (IML) - 0.05\ kIU/mL^*$$

$$r = 0.991$$

Means:

11.5 kIU/mL (IML 1000)

11.5 kIU/mL (IML)

\*Amniotic fluid samples were diluted 1-in-101 off-line before they were tested by the IMMULITE and IMMULITE 1000 instruments.

## References

### Testicular Cancer

- Herberman, ed. Immunodiagnosis of cancer. 1979:101.
- Int J Cancer 1971;7:218.
- Cancer Res 1972;32:979.
- Scand J Clin Lab Investig 1956;8:174.
- Acta Unio Internationalis Contra Cancrum 1963;19:80.
- Vopr Med Khim 1964;10:90.
- Int J Cancer 1968;3:364.
- Hum Pathol 1979;10:557.
- Lancet 1976;2:433.
- Adv Cancer Res 1971;14:295.
- Proc Natl Acad Sci USA 1973;70:526.
- Kirkpatrick, ed. Alpha-fetoprotein. 1981:115.
- Cancer 1979;44:984.
- Cancer Res 1975;35:991.
- Cancer 1974;34:1510.
- Med Intelligence 1976;295:1237.
- Cancer 1976;37:215.
- J Urol 1977;118:994.
- Cancer 1978;42:2768.
- J Clin Oncol 1986;4:41.
- Rose, ed. Manual of clinical laboratory immunology. 1986:810.
- Cancer 1981;47:328.
- Urol Clin North Am 1977;4:393.
- Ravitch, ed. Current problems in surgery. 1978:1.
- Cancer 1980;45:1755.
- Pavone-Macaluso, ed. Testicular cancer and other tumors. 1983:63.
- J Urol 1978;119:759.
- Kirkpatrick, ed. Alpha-fetoprotein. 1981:135.
- N Eng J Med 1977;296:693.
- Eur J Clin Chem Biochem 1993;31:517.
- NCCLS H3-A3, 1991.

### Fetal Open Neural Tube Defects

- Brock DJH. Prenatal diagnosis – chemical methods. Br Med Bull 1976;32:16.
- Crandall BF. Alpha-fetoprotein: a review. CRC Crit Rev Clin Lab Sci 1981;Sept:127-85.
- Gitlin D. Normal biology of AFP. Ann NY Acad Sci 1975:259:17-28.
- Haddow JE, et al. Fetal disorders associated with elevated MSAFP values. Foundation Blood Res 1990:1.
- Leek AE, Chard T. Proceedings of colloquium on alpha-fetoprotein. In: Masseyeff R, editor. Paris: L'Institut National de la Sante et de la Recherche Medicale 563, (Nice) 1974.
- Second Report of the Collaborative Study on Alpha-Fetoprotein in Relation to Neural Tube Defects. Amniotic fluid alpha-fetoprotein measurement in antenatal diagnosis of

anencephaly and open spina bifida in early pregnancy. Lancet 1979;ii:651. 38) Seppälä M. Fetal pathophysiology of human  $\alpha$ -fetoprotein. Ann N Y Acad Sci 1975;259:59.

## Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

## Tables and Graphs

### Precision (IU/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	0.70	0.047	6.7%	0.059	8.4%
2	2.63	0.195	7.4%	0.217	8.3%
3	12.9	0.63	4.9%	1.04	8.1%
4	32.3	1.46	4.5%	2.65	8.2%
5	45.1	2.48	5.5%	3.69	8.2%
6	60.9	2.81	4.6%	4.22	6.9%
7	187	8.64	4.6%	10.9	5.8%

### Linearity (IU/mL) - Serum

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	32 in 32 <sup>5</sup>	160	—	—
	16 in 32	80	80	100%
	8 in 32	40	40	100%
	4 in 32	21	20	105%
	2 in 32	11	10	110%
	1 in 32	5.3	5.0	106%
2	64 in 64	184	—	—
	32 in 64	91	92	99%
	16 in 64	48	46	104%
	8 in 64	24	23	104%
	4 in 64	12	11	109%
	2 in 64	6.2	5.7	109%
	1 in 64	3.3	2.9	114%
3	64 in 64	231	—	—
	32 in 64	118	115	103%
	16 in 64	57	58	98%
	8 in 64	29	29	100%
	4 in 64	14	14	100%
	2 in 64	7.3	7.2	101%
	1 in 64	3.5	3.6	97%
4	64 in 64	235	—	—
	32 in 64	114	117	97%
	16 in 64	61	59	103%
	8 in 64	30	29	103%
	4 in 64	16	15	107%
	2 in 64	7.7	7.3	105%
	1 in 64	3.9	3.7	105%

### Linearity (IU/mL) - Amniotic Fluid

	Total Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	1 in 100 <sup>5</sup>	138	—	—
	1 in 200	66	69	96%
	1 in 400	35	35	100%
	1 in 800	14	17	82%
	1 in 1600	9.0	8.6	105%
	1 in 3200	4.3	4.3	99%
2	1 in 100	175	—	—
	1 in 200	81	88	92%
	1 in 400	46	44	105%
	1 in 800	23	22	105%
	1 in 1600	12	11	109%
	1 in 3200	6.5	5.5	119%
3	1 in 100	281	—	—
	1 in 200	119	141	85%
	1 in 400	60	70	86%
	1 in 800	31	35	89%
	1 in 1600	16	18	89%
	1 in 3200	8.4	8.8	95%

### Recovery (IU/mL) - Serum

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	10.6	—	—
	A	11.4	10.6	108%
	B	12.3	11.6	106%
	C	13.7	13.4	102%
2	—	59	—	—
	A	55	57	96%
	B	58	58	100%
	C	58	59	98%
3	—	117	—	—
	A	112	112	100%
	B	109	113	96%
	C	113	115	98%

### Recovery (IU/mL) - Amniotic Fluid

	Low Amniotic Sample <sup>1</sup>	Spiking High Amniotic Sample <sup>2</sup>	Observed <sup>3</sup>	Expected <sup>4</sup>	%O/E <sup>5</sup>
1	—	—	244	—	—
		A	813	732	111%
		B	1371	1232	111%
		C	1943	2032	96%
2	—	—	4483	—	—
		A	4519	4759	95%
		B	6678	5259	127%
		C	4896	6059	81%
3	—	—	9050	—	—
		A	9054	9098	100%
		B	10,042	9598	105%
		C	10,640	10,398	102%

### Specificity


Compound <sup>1</sup>	Amount Added <sup>2</sup>	% Cross reactivity <sup>3</sup>
Human IgG	10 mg/mL	ND
Human Hemoglobin	192 mg/dL	ND
Human Transferrin	400 mg/dL	ND
Human Serum Albumin	60 mg/mL	ND
Cisplatin	100 µg/mL	ND
Cyclophosphamide	1000 µg/mL	ND
Doxorubicin	100 µg/mL	ND
5-Fluorouracil	1000 µg/mL	ND
Mitomycin C	100 µg/mL	ND
Vincristine	1000 ng/mL	ND

ND: not detectable<sup>4</sup>

### Bilirubin



	Bilirubin (unconjugated) <sup>1</sup>				
	100 mg/L		200 mg/L		
	Expected <sup>2</sup>	Observed <sup>3</sup>	%O/E <sup>4</sup>	Observed	%O/E
1	5.3	4.9	92%	5.1	96%
2	5.7	5.3	93%	5.6	98%
3	26	24	94%	25	97%
4	48	46	95%	46	95%
5	49	46	93%	46	93%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2019 Siemens Healthcare Diagnostics. All rights reserved.