Dimension Vista® System

Flex® reagent cartridge

BENZ

See shaded sections: Updated information from 2019-05 version.

Issue Date 2019-10-29

Urine Benzodiazepines Screen


Intended Use: The BENZ method is an *in vitro* diagnostic test for the qualitative and semi-quantitative determination of benzodiazepines in human urine on the Dimension Vista® System. Measurements obtained with the BENZ method are used in the diagnosis and treatment of benzodiazepines use or overdose.

The BENZ method provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2} Other chemical confirmation methods are available. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary: Benzodiazepines are sedative-hypnotic drugs that are structurally similar and include widely used drugs such as chlordiazepoxide, diazepam, and oxazepam. The different benzodiazepines are absorbed at different rates, and the timing of their psychoactive effects varies with the absorption rate. Benzodiazepines are usually taken orally and are metabolized in the liver. Some benzodiazepine metabolites are pharmacologically active.¹ Benzodiazepines potentiate the effect of other central nervous system depressants, such as ethyl alcohol.¹

Principles of Procedure: Syva® Emit® II Plus¹ reagents are contained in the BENZ Flex® reagent cartridge. The BENZ cartridge is used on the Dimension Vista® System to provide an automated application of this assay. It is based on the competition for antibody binding sites between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH).

Matched lots of polyclonal antibody reactive to lorazepam and lorazepam labeled with glucose-6-phosphate dehydrogenase are used in this Syva® Emit® II Plus methodology. The reaction sequence is as shown:

Where: Ab = antibody reactive to diazepam
BENZ = benzodiazepines
BENZ-G6PDH = diazepam-glucose-6-phosphate dehydrogenase conjugate

The concentration of drug in the sample determines the amount of BENZ-glucose-6-phosphate dehydrogenase (BENZ-G6PDH) conjugate that is bound to the antibody. The unbound conjugate catalyzes the oxidation of glucose-6-phosphate, with the simultaneous reduction of NAD⁺ to NADH, more rapidly than does the bound conjugate. The rate of increasing absorbance at 340 nm due to the increase in NADH is related to the concentration of drug in the sample.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1–6	Liquid	Diazepam enzyme conjugates, ^b bovine serum albumin	c, d	
7–12	Liquid	Antibody ^b reactive to diazepam, bovine serum albumin, NAD ⁺ , G6P	c, d	Sheep

- Wells are numbered consecutively from the wide end of the cartridge.
- Antibody and conjugate are a matched pair.
- Antibody titer and conjugate activity may vary from lot to lot.
- Contain buffers, stabilizers and preservatives.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2–8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 2 days for wells 1–12

Specimen Collection and Handling:^{1,2,3,4} Normal procedures for collecting urine may be used for specimens to be analyzed with the BENZ Flex® reagent cartridge.

- Urine specimens should be collected in clean, unbreakable, and leak proof containers. Freshly voided urine specimens should be used. If not analyzed immediately, specimens should be stored refrigerated for less than 24 hours. Specimens should be frozen if storage longer than 24 hours is required.
- Fresh urine specimens do not require any special pretreatment.
- Specimens may be encountered that display turbidity. It is recommended that such specimens be centrifuged before analysis.
- No additives or preservatives are required.
- Boric acid should not be used as a preservative.
- Specimens should be within the pH range of 5–8. Specimens with a pH outside this range should be adjusted to this range by the addition of 1N HCl or 1N NaOH before analysis.
- Specimens should be at a temperature of 20–25 °C (68–77 °F) before analysis.
- Adulteration of the urine specimen may cause erroneous results. If adulteration is suspected, obtain a fresh specimen.
- Human urine specimens should be handled and treated as if they are potentially infectious.

Procedure

Materials Provided

BENZ Flex® reagent cartridge, Cat. No. K5097

Materials Required But Not Provided

UDAT CAL, Cat. No. KC510
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

Note

Sample cups must be filled consistently with 1.0 mL of sample, control, or calibrator using an automatic pipette. Plastic transfer pipettes should not be used. This helps maintain a constant cup surface to sample volume ratio and minimizes drug loss from the sample due to adherence of drugs to the cup surface.

Test Conditions

Sample Volume (delivered to the cuvette)	4 µL
Reagent 1 (Antibody, NAD, G6P) Volume	98 µL
Reagent 2 (Enzyme Conjugate) Volume	42 µL
Temperature	37.0 °C
Reaction Time	5.3 minutes
Wavelength	340 and 600 nm
Type of Measurement	Rate

Calibration

Cutoff Concentration	200 ng/mL
Calibration Material	0 ng/mL lorazepam Drugs of Abuse Calibrator (UDAT CAL) Level A 100 ng/mL lorazepam level prepared automatically using Drugs of Abuse Calibrator (UDAT CAL) Levels A and B 200 ng/mL lorazepam level prepared automatically using Drugs of Abuse Calibrator (UDAT CAL) Levels A and B 300 ng/mL lorazepam level prepared automatically using Drugs of Abuse Calibrator (UDAT CAL) Levels A and B 1000 ng/mL lorazepam Drugs of Abuse Calibrator (UDAT CAL) Level B
Typical Calibration Levels	0, 100, 200, 300, and 1000 ng/mL Validate the calibration by assaying a positive and negative control.
Calibration range	0–1000 ng/mL
Calibration Scheme	5 levels n=5
Calibration Frequency	Every 30 days for any one lot Calibration interval may be extended based on acceptable verification of calibration

A new calibration is required:

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

Quality Control

At least once each day of use, analyze a positive and a negative control relative to the cutoff concentration using a suitable urine based control material. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Users should follow the appropriate federal, state, or local guidelines concerning the assaying of quality controls.

Results:

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Qualitative Mode: The results from a BENZ analysis can be reported qualitatively as **negative or positive** relative to the cutoff for BENZ. **Positive** results are only preliminary and suggest that the sample is likely to contain benzodiazepines. **Negative** results indicate that either the sample does not contain benzodiazepines or that benzodiazepines are present in concentrations below the cutoff level.

Semi-quantitative Mode: Results of a BENZ analysis can be reported semi-quantitatively with a numerical unit concentration (ng/mL) result followed by a **negative or positive** designation. The concentration result in ng/mL represents approximate cumulative concentration of drugs detected by the BENZ analysis. The semi-quantitation numerical result may be used for quality control, calibration or for instrument troubleshooting. Semi-quantitative results in ng/mL are equivalent to drug concentration only when the sample is a standard containing pure lormetazepam. BENZ results greater than or equal to the cutoff concentration should be reported as positive. BENZ results below the cutoff concentration should be reported as negative. Results below the cutoff concentration indicate that either no benzodiazepines are present in the sample or benzodiazepines are present at concentrations below the cutoff level.

Analytical Measurement Range (AMR): 30–900 ng/mL

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Autodilution (AD): This feature is only available for samples tested in the semi-quantitative mode. The autodilute sample volume is 12 µL (dilution factor = 5). Refer to your Dimension Vista® Operator's Guide.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in BENZ results. Refer to your Dimension Vista® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

- A positive result suggests the likely presence of benzodiazepines but does not indicate or measure intoxication.
- The presence of benzodiazepines in urine is only an indication of recent exposure to or use of benzodiazepines.
- The psychological and physiological effects of benzodiazepines do not necessarily correlate with urinary concentration.
- A positive BENZ result suggests the likely presence of drug and its metabolites. The BENZ method cannot fully quantitate the concentration of individual components.
- Interpretation of results must take into account that urine concentrations vary extensively with fluid intake, and other biological variables.
- There is a possibility that other substances and/or factors not listed above may interfere with the test and cause false results, e.g., technical or procedural errors.
- NEGATIVE results for specimens with concentrations below the assay range may be accompanied by an "assay range" or by a "below assay range" message. These results should be reported as NEGATIVE.
- POSITIVE results for specimens with concentrations above the assay range may be accompanied by an "assay range" or by an "above assay range" message. These results should be reported as POSITIVE.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n=5 replicates at the following lormetazepam concentrations are:

Concentration	SD
200 ng/mL lormetazepam	46 ng/mL

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Precision^{5,a}

Material ^f	Mean ng/mL	Repeatability		Within-Lab	
		SD ng/mL	CV %	SD ng/mL	CV %
Level 1 (150 ng/mL)	150	12.6	8	16.3	11
Cutoff Level (200 ng/mL)	225	10.9	5	16.8	7
Level 2 (250 ng/mL)	267	14.8	6	20.0	7

e. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

f. Syva® Emit® Calibrators and in house controls.

Summary of Clinical Results

One hundred native urine specimens were tested with the BENZ Flex® cartridge on the Dimension Vista® System (cutoff=200 ng/mL) and with the BENZ Flex® cartridge on the Dimension® clinical chemistry system (cutoff=200 ng/mL). All positive specimens were also analyzed by GC/MS.

BENZ Flex® Reagent Cartridge on Dimension® system (cutoff 200 ng/mL)

BENZ Flex® Reagent Cartridge on Dimension Vista® System (cutoff 200 ng/mL)

	+	-
+	55	0
-	0	45

GC/MS (cutoff 200 ng/mL)

BENZ Flex® Reagent Cartridge on Dimension Vista® System (cutoff 200 ng/mL)

	+	-
+	56	12*
-	0	16

*Discrepant (ng/mL):

GC/MS	Compound	Dimension Vista® System
117	Nordiazepam	243
130	Nordiazepam	296
137	Oxazepam	316
135	Nordiazepam	317
148	Oxazepam	424
145	Nordiazepam	368
156	Nordiazepam	463
165	Nordiazepam	539
174	Nordiazepam	581
97	Oxazepam	>900
188	Nordiazepam	742
*	*190 ng/mL oxazepam	>900
	*162 ng/mL desmethyldiazepam	
	*88 ng/mL temazepam	

Specificity⁶

Non Interfering Substances

Each of the following compounds when added to urine containing lormetazepam at +/- 25% concentration of the cutoff do not yield a false response relative to the 200 ng/mL cutoff:

Compound	Concentration
Acetone	1.0 g/dL
Ascorbic Acid	1.5 g/dL
Bilirubin	2.0 mg/dL
Creatinine	0.5 g/dL
Ethanol	1.0 g/dL
Gamma Globulin	0.5 g/dL
Glucose	2.0 g/dL
Hemoglobin	115 mg/dL
Human Serum Albumin	0.5 g/dL
Oxalic Acid	0.1 g/dL
Riboflavin	7.5 mg/dL
Sodium Chloride	6.0 g/dL
Urea	6.0 g/dL

‡ Each of the following compounds was added to drug free urine and gave negative BENZ results at the concentration listed for the 200 ng/mL cutoff.

Compound	Concentration
Acetaminophen	1000 µg/mL
α-Acetyl-N,N-dinormethadol (dinor LAAM)	25 µg/mL
l-α-Acetylmethadol (LAAM)	25 µg/mL
N-Acetylprocainamide (NAPA)	400 µg/mL
Acetylsalicylic Acid	1000 µg/mL
Amitriptyline	1000 µg/mL
d-Amphetamine	1000 µg/mL
Benzoylcegonine	1000 µg/mL
Buprenorphine	100 µg/mL
Caffeine	1000 µg/mL
Cimetidine	1000 µg/mL
Clomipramine	2.5 µg/mL
Clonidine	1000 µg/mL
Codeine	500 µg/mL
Cotinine	100 µg/mL
Cyclobenzaprine	28 µg/mL
Desipramine	800 µg/mL
Diphenhydramine	1000 µg/mL
Doxepin	10 µg/mL
2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	1000 µg/mL
Fluoxetine	500 µg/mL
Glutethimide	500 µg/mL
Ibuprofen	1000 µg/mL
Ketamine	100 µg/mL
Ketorolac Tromethamine	1000 µg/mL
LSD	10 ng/mL
Meperidine	1000 µg/mL
d-Methamphetamine	2 µg/mL
Methaqualone	1500 µg/mL
Morphine	1000 µg/mL
Naproxen	1000 µg/mL
Nortriptyline	750 µg/mL
Phencyclidine	1000 µg/mL
Phenytoin	1000 µg/mL
Promethazine	1000 µg/mL
Propoxyphene	1000 µg/mL
Ranitidine	900 µg/mL
Scopolamine	500 µg/mL

Secobarbital	1000 µg/mL
11-nor- Δ^8 -THC-9-COOH	0.2 µg/mL
Thioridazine	100 µg/mL
Tramadol	100 µg/mL
Tyramine	100 µg/mL
Zidovudine (AZT)	2 mg/mL
Zolpidem	100 µg/mL

Cross-reactivity ‡

The table below gives the compounds this assay is designed to detect and the levels at which the compounds have been found to give a response approximately equivalent to that of the 200 ng/mL lorazepam cutoff. Each concentration represents the reactivity level for the stated compound when it is added to a negative urine specimen. If a sample contains more than one compound detected by the assay, lower concentrations than those listed below may combine to produce a rate approximately equivalent to or greater than that of the cutoff calibrator.

Compound	Concentration (ng/mL)
Alprazolam	65
7-aminoclonazepam	5300
7-aminoflunitrazepam	930
Bromazepam	630
Chlordiazepoxide	3300
Clobazam	260
Clonazepam	210
Clorazepate	g
Clotiazepam	380
Demoxepam	1600
<i>N</i> -Desalkylflurazepam	130
<i>N</i> -Desmethyldiazepam	110
Diazepam	70
Estazolam	90
Flunitrazepam	140
Flurazepam	190
Halazepam	110
α -Hydroxyalprazolam	100
α -Hydroxyalprazolam Glucuronide	110 ^h
1- <i>N</i> -Hydroxyethylflurazepam	150
α -Hydroxytriazolam	130
Ketazolam	100
Lorazepam	600
Lorazepam glucuronide	>20000 ^h
Medazepam	150
Midazolam	130
Nitrazepam	78
Norchlordiazepoxide	4500
Oxazepam	250
Oxazepam glucuronide	>30000 ^h
Prazepam	90
Temazepam	140
Temazepam Glucuronide	>20000 ^h
Tetrazepam	70
Triazolam	130

- g. Clorazepate degrades rapidly in stomach acid to nordiazepam. Nordiazepam hydroxylates to oxazepam.
h. Glucuronide metabolite of α -Hydroxyalprazolam cross-reacts with this assay. Other glucuronide metabolites such as Lorazepam, Oxazepam, and Temazepam cross react to a limited extent. The cross-reactivity of other glucuronide metabolites with this assay is not known.

Analytical Sensitivity

The sensitivity of the BENZ method is 30 ng/mL for the 200 ng/mL cutoff and it represents the lowest concentration of BENZ that can be distinguished from zero. Sensitivity is defined as the concentration at two standard deviations above 0 ng/mL using Drugs of Abuse Calibrator Level A (n=20).

‡ The Dimension Vista® BENZ method (REF K5097) and the Dimension® BENZ (REF DF97A) method utilize the same reagents under equivalent reaction conditions. Interfering substances and Cross-reactivity were tested using Dimension® BENZ (REF DF97A) and the results are representative of both methods.

Symbols Key: See Adjacent Panel.

Bibliography: See Adjacent Panel.

Dimension Vista®, Dimension®, Flex®, Emit® and Syva® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

1. Emit® II Plus Benzodiazepines Assay Instructions for Use, 2000.
2. Hawks RL, Chiang CN. Urine testing for drugs of abuse, NIDA Research Monograph 73, U.S. Government Printing Office, Washington, DC 20402, 1986:24-29.
3. Burtis CA, Ashwood ER. *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Philadelphia, PA 1999, Third edition pp 39@41 (Biological Hazards), pp 1191-1193 (Drugs of Abuse).
4. Baselt RC, Cravey RH. Disposition of toxic drugs and chemicals in man. Chicago: Year Book Publishers, Inc., 2000: 245 -248.
5. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
6. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare