

## EMIT® *tox*™ Acetaminophen

## Application Sheet

Shading indicates technical content that differs from the previous version.

### Emit® *tox*™ Acetaminophen Application Sheet

**For the AU400®, AU600®, AU640®,  
AU680®, AU2700®, AU5400®  
Clinical Chemistry Systems**

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the instrument manual for additional instructions.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthineers to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

### Reagents

These reagents are qualified for use with the Calibrator(s) listed below only.

| <u>Assay</u>                                 | <u>Catalog Number</u> |
|--|-----------------------|
| Emit® <i>tox</i> ™ Acetaminophen Assay       | 7A319UL               |
| Emit® <i>tox</i> ™ Acetaminophen Assay       | OSR7A229              |
| <u>Calibrators</u>                           | <u>Catalog Number</u> |
| Emit® <i>tox</i> ™ Acetaminophen Calibrators | 7A409UL               |

### Storage

Reagents which are in use may be stored on board the analyzer for up to 4 weeks or as long as the control results fall within acceptable limits.

### Instrument

### Calibration

Prepare a calibration curve whenever a new lot of reagent is used or as indicated by control results. Calibrate by placing the appropriate calibrators in the assigned positions in the calibration rack (yellow rack). Run a reagent blank (blue rack) with each calibration curve.

### Instrument Settings

See page 2.

### Results

Results are reported in µg/mL [µmol/L].

# EMIT<sup>®</sup> tox™ Acetaminophen

## Application Sheet

### Instrument

#### Instrument Settings

##### General Screen

Reagent ID: 517

| General                  |              | LIH      | ISE            | Range |
|--------------------------|--------------|----------|----------------|-------|
| Test Name                | User Defined |          | <div>◀ ▶</div> |       |
| Type                     | Serum        |          | Operation Yes  |       |
| Sample Volume            | 3.0          | uL       | Dilution       | 0 uL  |
| Reagents R1 Volume       | 150          | uL       | Dilution       | 0 uL  |
| Reagents R2 Volume       | 75           | uL       | Dilution       | 0 uL  |
| Wavelength Pri           | 340          |          | Sec            | 410   |
| Method                   | RATE         |          |                |       |
| Reaction Slope           | +            |          |                |       |
| Measuring Point 1 First  | 15           |          | Last           | 23    |
| Measuring Point 2 First  | NA           |          | Last           | NA    |
| Linearity                | 100          | %        |                |       |
| No-Lag-Time              |              |          |                |       |
| Pre-Dilution Rate        |              |          |                |       |
| Min OD                   | Max OD       |          |                |       |
| Reagent OD Limit         | L -2.0000    | H 2.5000 |                |       |
| Dynamic Range            | L 0.3        | H 200    |                |       |
| Correlation Factor       | A 1.000000   | B 0.0    |                |       |
| Onboard Stability Period | User Defined |          |                |       |

#### Range Screen

| General         |              | LIH          | ISE            | Range |
|-----------------|--------------|--------------|----------------|-------|
| Test Name       | User Defined |              | <div>◀ ▶</div> |       |
| Type            | Serum        |              |                |       |
| Value Flag      | Value        |              | Level L        | NA    |
| Level H         | NA           |              |                |       |
| Normal Ranges   | Age L        | Age H        |                |       |
| Sex             | Year         | Month        | Year           | Month |
| 1               |              |              |                |       |
| 2               |              |              |                |       |
| 3               |              |              |                |       |
| 4               |              |              |                |       |
| 5               |              |              |                |       |
| 6               |              |              |                |       |
| 7 None Selected |              |              |                |       |
| 8 Out of Range  |              |              |                |       |
| L               | H            |              |                |       |
| Panic Value     | User Defined | User Defined | Unit           | µg/mL |
| Decimal Places  | 1            |              |                |       |

### Calibration Parameters Screen

| General           |                                      | ISE                          |              |
|-------------------|--------------------------------------|------------------------------|--------------|
| Test Name         | User Defined                         | ◀ ▶                          | Type Serum   |
| Calibration Type  | 5AB                                  | Formula                      | EIA Type 1   |
|                   | Counts                               | 1                            | Process      |
|                   | CONC                                 |                              |              |
| Point 1           | Cal No                               | OD                           | CONC         |
| Point 2           | User Def                             | User Def                     | 10.0         |
| Point 3           | User Def                             | User Def                     | 25.0         |
| Point 4           | User Def                             | User Def                     | 50.0         |
| Point 5           | User Def                             | User Def                     | 100.0        |
| Point 6           | User Def                             | User Def                     | 200.0        |
| Point 7           |                                      |                              |              |
| 1-Point Cal Point |                                      | Factor/OD-L                  | Factor/OD-H  |
|                   |                                      | -2.000000                    | 2.500000     |
|                   |                                      | -2.000000                    | 2.500000     |
|                   |                                      | -2.000000                    | 2.500000     |
|                   |                                      | -2.000000                    | 2.500000     |
|                   |                                      | -2.000000                    | 2.500000     |
|                   |                                      |                              |              |
|                   |                                      |                              |              |
| 1-Point Cal Point | <input type="checkbox"/> With CONC-0 | Slope Check                  | +            |
| MB Type Factor    | NA                                   | Advanced Calibration         | User Def     |
|                   |                                      | Calibration Stability Period | User Defined |

### Performance

#### Method Comparison

Clinical specimens were tested using each Emit® tox™ Acetaminophen Assay on the AU600 analyzer and on the Syva®-30R analyzer. The results from the AU600 are as follows:

|                         |       |
|-------------------------|-------|
| Slope                   | 0.98  |
| Intercept (µg/mL)       | -1.34 |
| Correlation Coefficient | 0.992 |
| Number of Samples       | 50    |

#### Precision

Within run precision was calculated according to NCCLS Guideline EP5-A by running 2 replicates of each control level twice a day for 20 days (N=80). Total precision was also calculated from these data and presented in µg/mL.

|      | Within Run Precision |         |         | Total Precision |         |         |
|------|----------------------|---------|---------|-----------------|---------|---------|
|      | Level 1              | Level 2 | Level 3 | Level 1         | Level 2 | Level 3 |
| Mean | 15.8                 | 45.2    | 135.3   | 15.8            | 45.2    | 135.3   |
| SD   | 0.6                  | 1.1     | 5.1     | 0.7             | 1.9     | 6.8     |
| % CV | 3.6                  | 2.5     | 3.8     | 4.4             | 4.1     | 5.0     |

**Endogenous Substances**

No clinically significant interference has been found in samples to which 800 mg/dL hemoglobin, 30 mg/dL bilirubin or 750 mg/dL triglycerides were added to simulate hemolytic, icteric, or lipemic samples.

**Analytical Sensitivity**

The sensitivity level of the Emit® *tox*™ Acetaminophen Assay on the AU600 is 0.3 µg/mL acetaminophen. This level represents the lowest concentration of acetaminophen that can be distinguished from 0 ng/mL with a confidence level of 95%.










**NOTE:** Performance on the AU400, AU600, AU640, AU680, AU2700, and AU5400 series analyzers has been shown to be equivalent.

# EMIT® tox™ Acetaminophen

## Application Sheet



















### Definition of Symbols

The following symbols may appear on the product labeling:

| Symbol  | Symbol Title   | Source               | Symbol  | Symbol Title  | Source               |
|---|--|----------------------|---|---|----------------------|
|    | Manufacturer   | 5.1.1 <sup>a</sup>   |    | Authorized representative in the European Community | 5.1.2 <sup>a</sup>   |
|    | Use-by date  | 5.1.4 <sup>a</sup>   |    | Authorized representative in Switzerland            | Proprietary          |
|    | Catalog number   | 5.1.6 <sup>a</sup>   |    | Batch code  | 5.1.5 <sup>a</sup>   |
|    | Consult Instructions for Use                                       | 5.4.3 <sup>a</sup>   |    | Contains sufficient for <n> tests                   | 5.5.5 <sup>a</sup>   |
|  | Internet URL address to access the electronic instructions for use | Proprietary          |  | Version of Instructions for Use                     | Proprietary          |
|  | In vitro diagnostic medical device                                 | 5.5.1 <sup>a</sup>   |  | Revision  | Proprietary          |
|  | Prescription device (US only)                                      | FDA <sup>c</sup>     |  | Unique Device Identifier                            | 5.7.10 <sup>b</sup>  |
|  | CE Marking with Notified Body                                      | EU IVDR <sup>d</sup> |  | CE Marking  | EU IVDR <sup>d</sup> |
|  | Temperature limit  | 5.3.7 <sup>a</sup>   |  | Keep away from sunlight                             | 5.3.2 <sup>a</sup>   |
|  | Upper limit of temperature   | 5.3.6 <sup>a</sup>   |  | Lower limit of temperature                          | 5.3.5 <sup>a</sup>   |
|  | Do not re-use  | 5.4.2 <sup>a</sup>   |  | Do not freeze                                       | Proprietary          |

# EMIT® *tox*™ Acetaminophen

## Application Sheet

| Symbol  | Symbol Title                     | Source               | Symbol  | Symbol Title                  | Source               |
|---|----------------------------------|----------------------|---|-------------------------------|----------------------|
|    | Recycle                          | 1135 <sup>a</sup>    |    | This way up                   | 0623 <sup>a</sup>    |
|    | Biological risks                 | 5.4.1 <sup>a</sup>   |    | Caution                       | 5.4.4 <sup>a</sup>   |
|    | Common Units                     | Proprietary          |    | International System of Units | Proprietary          |
|    | Document face up <sup>f</sup>    | 1952 <sup>a</sup>    |    | Date format (year-month-day)  | N/A                  |
|   | Non-sterile                      | Proprietary          |    | Date format (year-month)      | N/A                  |
|  | Reconstitution volume            | Proprietary          |  | Contents                      | Proprietary          |
|  | For forensic/toxicology use only | Proprietary          |  | Level                         | Proprietary          |
|  | Dropper                          | Proprietary          |  | Cassette                      | Proprietary          |
|  | Not for self-testing             | EU IVDR <sup>d</sup> |  | Not for near-patient testing  | EU IVDR <sup>d</sup> |

<sup>a</sup> International Standard Organization (ISO). ISO 15223-1 Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied.

<sup>b</sup> ISO 15223-1:2020-04.

<sup>c</sup> Federal Register. Vol. 81, No 115. Wednesday, June 15, 2016. Rules and Regulations: 38911.

<sup>d</sup> IVDR REGULATION (EU) 2017/746

<sup>e</sup> International Standard Organization (ISO). ISO 7000 Graphical symbols for use on equipment.

<sup>f</sup> Indicates Assay-eNote.

**Syva®**

## **EMIT® *tox*™ Acetaminophen**

## **Application Sheet**

Emit®, *tox*™ and Syva® are trademarks of Siemens Healthineers.

AU®, AU400®, AU600®, AU640®, AU680®, AU2700®, and AU5400® are registered trademarks of Beckman Coulter, Inc.

**For technical assistance:**  
**Beckman Coulter customers, contact the**  
**Customer Technical Support Center at**  
**1-800-854-3633 (USA & Canada)**

**In other countries, please contact your local**  
**Beckman Coulter representative.**

**Siemens Healthineers customers,**  
**contact the Technical Solutions Center**  
**at 1-800-227-8994 In the USA.**

### **Technical Assistance**

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.  
[siemens-healthineers.com](http://siemens-healthineers.com)

© 2010–2024 Siemens Healthineers  
All rights reserved.