

EMIT® II Plus Ecstasy Assay

Application Sheet

Shading indicates technical content that differs from the previous version.

Emit® II Plus Ecstasy Assay Application Sheet

**For the AU400®, AU600®, AU640®,
AU680®, AU2700®, AU5400®
Clinical Chemistry Systems**

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the instrument manual for additional instructions.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthineers to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

Reagents

These reagents are qualified for use with the Calibrators/Controls listed below only. Other material may be used however, for quality control purposes.

Assay	Catalog Number	
	28 mL Kit	1000 mL Kit
Emit® II Plus Ecstasy Assay	9X029UL	9X129UL
Emit® II Plus Ecstasy Assay	OSR9X229	

Calibrator/Control

Catalog Number

Emit® Calibrator/Control Level 0	9A509UL
Emit® II Plus 6-AM/Ecstasy Calibrator/Control Level 1	9R529UL
Emit® II Plus 6-AM/Ecstasy Calibrator/Control Level 2	9R549UL
Emit® II Plus 6-AM/Ecstasy Calibrator/Control Level 3	9R569UL
Emit® II Plus 6-AM/Ecstasy Calibrator/Control Level 4	9R589UL

Storage

Reagents which are in use may be stored on board the analyzer for up to 4 weeks or as long as the control results fall within acceptable limits.

Qualitative

Calibration

Calibrate by running the Level for the appropriate cutoff listed in Table 1. Select the appropriate Option to use based on how results should be flagged (see Results section below).

Run a reagent blank (blue rack) and place the appropriate calibrator in a white rack for Option 1 or in its assigned position in a yellow rack for Option 2 or 3. Recalibrate as indicated by control results.

Results

Results are reported based on the Option chosen.

Option

- Results are printed as OD values. Positive samples **are NOT** flagged. Operator must compare the sample response to the cutoff calibrator response and determine if the sample is positive or negative.
- The cutoff is normalized to 100. Positive samples are ≥ 100 and are flagged with a (P).
- The cutoff is normalized to 0. Positive samples are ≥ 0 and are flagged with a (P).

EMIT® II Plus Ecstasy Assay

Application Sheet

Table 1

Assay	Level 2 ng/mL	Level 3 ng/mL
Ecstasy	300	500

Instrument Settings

See page 2.

Qualitative

Instrument Settings

General Screen

Reagent ID: 525

General		LIH	ISE	Range				
Test Name	User Defined				Type	Urine	Operation	Yes
Sample	Volume	\$	uL	Dilution	0	uL	Pre-Dilution Rate	1
Reagents	R1 Volume	130	uL	Dilution	0	uL	Min OD	Max OD
	R2 Volume	55	uL	Dilution	0	uL	L -2.000	H 2.500
Wavelength	Pri	340		Sec	410		Reagent OD Limit	
Method	FIXED				First L	-2.000	First H	2.500
Reaction Slope	+				Last L	-2.000	Last H	2.500
Measuring Point 1	First	13		Last	15		Dynamic Range	
Measuring Point 2	First			Last			L User Defined	H User Defined
Linearity			%			Correlation Factor		
No-Lag-Time						A	1	B †
						Onboard Stability Period	User Defined	

\$ – Sample Volume: Enter 7.0 uL for the 300 cutoff; 5.0 uL for the 500 cutoff

† – Correlation Factor B: For Results Option 1 enter 0; for Option 2 enter 0; for Option 3 enter -100

Range Screen

General LIH ISE Range									
Test Name		User Defined		◀ ▶		Type		Urine	
Value Flag		Flag		Level L		-9999		Level H	
Normal Ranges		Age L		Age H					
	Sex	Year	Month	Year	Month	L		H	
<input type="checkbox"/> 1									
<input type="checkbox"/> 2									
<input type="checkbox"/> 3									
<input type="checkbox"/> 4									
<input type="checkbox"/> 5									
<input type="checkbox"/> 6									
7 None Selected									
8 Out of Range									
Panic Value		L		H		Unit		Decimal Places	
								0	

†† – Level High: For Results Option 1 enter 9999; for Option 2 enter 100; for Option 3 enter 0

Calibration Specific Screen — Results Option 1

General ISE									
Test Name		User Defined		◀ ▶		Type		Urine	
Calibration Type		MB		Formula		Y=AX+B		Counts	
								2	
								Process	
								OD	
	Cal No	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 1									
Point 2									
Point 3									
Point 4									
Point 5									
Point 6									
Point 7									
1-Point Cal Point		<input type="checkbox"/> With CONC-0		Slope Check				Advanced Calibration	
								User Def	
MB Type Factor		1000.000		Calibration Stability Period		User Defined			

EMIT® II Plus Ecstasy Assay

Application Sheet

Calibration Specific Screen — Results Option 2 or 3

General ISE					
Test Name	User Defined		Type		Urine
Calibration Type	AB		Formula	Y=AX+B	
	Counts	2		Process	OD
	Cal No	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1	User Def	User Def	100	-9999	9999
Point 2					
Point 3					
Point 4					
Point 5					
Point 6					
Point 7					
1-Point Cal Point	<input type="checkbox"/>	<input type="checkbox"/> With CONC-0		Slope Check	<input type="checkbox"/>
Advanced Calibration			User Def		
MB Type Factor			Calibration Stability Period	User Defined	

Semi-Quantitative

Calibration

Prepare a calibration curve by running the Calibrator/ Controls in the order listed in Table 2. Recalibrate as indicated by control results. Run a reagent blank (blue rack) daily to ensure consistent day to day control results.

Results

Results are reported as ng/mL.

Table 2

Assay	Level 0	Level 1 ng/mL	Level 2 ng/mL	Level 3 ng/mL	Level 4 ng/mL
Ecstasy 300	0	150	300	500	NA
Ecstasy 500	NA	150	300	500	1000

Instrument Settings

General Screen

General		LIH	ISE	Range
Test Name	User Defined		◀ ▶	Type Urine Operation Yes
Sample Volume	§ uL		Dilution 0 uL	Pre-Dilution Rate 1
Reagents R1 Volume	130 uL		Dilution 0 uL	Min OD Max OD
R2 Volume	55 uL		Dilution 10 uL	L -2.000 H 2.500
Wavelength Pri	340		Sec 410	Reagent OD Limit
Method	FIXED			First L -2.000 First H 2.500
Reaction Slope	+			Last L -2.000 Last H 2.500
Measuring Point 1 First	13		Last 15	Dynamic Range
Measuring Point 2 First			Last	L 75 H †
Linearity			%	Correlation Factor
No-Lag-Time				A 1 B 0
				Onboard Stability Period User Defined

§ – Sample Volume: Enter 7.0 uL for the 300 cutoff; 5.0 uL for the 500 cutoff

† – Dynamic Range High: Enter 500 for the 300 cutoff; 1000 for the 500 cutoff

Range Screen

General		LIH		ISE		Range	
Test Name		User Defined		◀ ▶		Type Urine	
Value Flag		Value		Level L		User Defined	
Normal Ranges		Age L		Age H		Level H †	
	Sex	Year	Month	Year	Month	L	H
<input type="checkbox"/> 1							
<input type="checkbox"/> 2							
<input type="checkbox"/> 3							
<input type="checkbox"/> 4							
<input type="checkbox"/> 5							
<input type="checkbox"/> 6							
7 None Selected							
8 Out of Range							
Panic Value		L		H		Unit	ng/mL
						Decimal Places	0

† – Level High: Enter the concentration of the cutoff calibrator from Table 1

Calibration Specific Screen

General		ISE	
Test Name		User Defined	
◀ ▶		Type Urine	
Calibration Type		4AB	
Formula		Polygonal	
Counts		2	
Process		CONC	
	Cal No	OD	CONC
Point 1	User Def	User Def	†
Point 2	User Def	User Def	†
Point 3	User Def	User Def	†
Point 4	User Def	User Def	†
Point 5			
Point 6			
Point 7			
1-Point Cal Point	<input type="checkbox"/> With CONC-0	Slope Check	+
MB Type Factor		Advanced Calibration	
		User Def	
Calibration Stability Period		User Defined	

† – CONC: Enter the Calibrator/Control concentration from Table 2 for the desired cutoff

Performance

Method Comparison

Clinical urine specimens were tested using the Emit® II Plus Ecstasy Assay on the AU600 analyzer and on the SYVA®-30R analyzer. Discordant samples were not observed at the 500 ng/mL cutoff. The discordant specimen at the 300 ng/mL cutoff contained 265 ng/mL methylenedioxymethamphetamine (MDMA) by GC/MS analysis. The results from the AU600 are listed below along with the percent agreement with the SYVA®-30R analyzer.

Assay	Positive	Negative	% Agreement
Ecstasy 300	56	52	98
Ecstasy 500	51	48	100

Precision

Within run precision was calculated according to NCCLS Guideline EP5-A by running 2 replicates of the cutoff Calibrator/Control and positive and negative controls twice a day for 20 days (N=80). Total precision was also calculated from these data. The following data are presented in mAU/min.

Ecstasy (300 ng/mL cutoff)

	Within Run Precision			Total Precision		
	Cutoff Cal.	Control 75%	Control 125%	Cutoff Cal.	Control 75%	Control 125%
Mean	465	435	500	465	435	500
SD	4.2	5.6	5.6	4.9	6.7	7.5
CV%	0.9	1.3	1.1	1.1	1.5	1.5

Ecstasy (500 ng/mL cutoff)

	Within Run Precision			Total Precision		
	Cutoff Cal.	Control 75%	Control 125%	Cutoff Cal.	Control 75%	Control 125%
Mean	585	539	612	585	539	612
SD	2.7	3.1	2.9	3.4	3.3	3.2
CV%	0.8	1.0	0.8	0.8	1.0	0.8

Analytical Recovery

Negative human urine specimens were spiked with concentrations of MDMA.

Qualitative analysis of the specimens spiked with drug concentrations lower than the cutoff concentration were correctly identified as negative 100% of the time. Specimens spiked with drug concentrations greater than the cutoff were correctly identified as positive 100% of the time.

Results from semi-quantitative analysis of the specimens are listed below.

Ecstasy (300 ng/mL cutoff)

Concentration (ng/mL)	EIA Type 1 Mean (ng/mL)	Polygonal Mean (ng/mL)
100	107	78
250	257	260
375	393	402
450	448	455

Ecstasy (500 ng/mL cutoff)

Concentration (ng/mL)	EIA Type 1 Mean (ng/mL)	Polygonal Mean (ng/mL)
250	259	267
375	404	413
450	485	468
550	613	644
750	791	871







NOTE: Performance on the AU400®, AU600®, AU640®, AU680®, AU2700®, and AU5400® series analyzers has been shown to be equivalent.

EMIT® II Plus Ecstasy Assay

Application Sheet



















Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Manufacturer	5.1.1 ^a		Authorized representative in the European Community	5.1.2 ^a
	Use-by date	5.1.4 ^a		Authorized representative in Switzerland	Proprietary
	Catalog number	5.1.6 ^a		Batch code	5.1.5 ^a
	Consult Instructions for Use	5.4.3 ^a		Contains sufficient for <n> tests	5.5.5 ^a
	Internet URL address to access the electronic instructions for use	Proprietary		Version of Instructions for Use	Proprietary
	<i>In vitro</i> diagnostic medical device	5.5.1 ^a		Revision	Proprietary
RxOnly	Prescription device (US only)	FDA ^c		Unique Device Identifier	5.7.10 ^b
	CE Marking with Notified Body	EU IVDR ^d		CE Marking	EU IVDR ^d
	Temperature limit	5.3.7 ^a		Keep away from sunlight	5.3.2 ^a
	Upper limit of temperature	5.3.6 ^a		Lower limit of temperature	5.3.5 ^a
	Do not re-use	5.4.2 ^a		Do not freeze	Proprietary

EMIT® II Plus Ecstasy Assay

Application Sheet

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Recycle	1135 ^a		This way up	0623 ^a
	Biological risks	5.4.1 ^a		Caution	5.4.4 ^a
	Common Units	Proprietary		International System of Units	Proprietary
	Document face up ^f	1952 ^a		Date format (year-month-day)	N/A
	Non-sterile	Proprietary		Date format (year-month)	N/A
	Reconstitution volume	Proprietary		Contents	Proprietary
	For forensic/toxicology use only	Proprietary		Level	Proprietary
	Dropper	Proprietary		Cassette	Proprietary
	Not for self-testing	EU IVDR ^d		Not for near-patient testing	EU IVDR ^d

^a International Standard Organization (ISO). ISO 15223-1 Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied.

^b ISO 15223-1:2020-04.

^c Federal Register. Vol. 81, No 115. Wednesday, June 15, 2016. Rules and Regulations: 38911.

^d IVDR REGULATION (EU) 2017/746

^e International Standard Organization (ISO). ISO 7000 Graphical symbols for use on equipment.

^f Indicates Assay-eNote.

Syva®**EMIT® II Plus Ecstasy Assay****Application Sheet**

Emit® and Syva® are trademarks of Siemens Healthineers.

AU®, AU400®, AU600®, AU640®, AU680®, AU2700®, and AU5400® are registered trademarks of Beckman Coulter, Inc.

For technical assistance:

**Beckman Coulter customers, contact the
Customer Technical Support Center at
1-800-854-3633 (USA & Canada)**

**In other countries, please contact your local
Beckman Coulter representative.**

**Siemens Healthineers customers,
contact the Technical Solutions Center
at 1-800-227-8994 In the USA.**

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.
[siemens-healthineers.com](https://www.siemens-healthineers.com)

© 2010–2024 Siemens Healthineers

All rights reserved.