

EMIT® II Plus Ethyl Alcohol Assay

Application Sheet

Shading indicates technical content that differs from the previous version.

Emit® II Plus Ethyl Alcohol Application Sheet

For the AU480® and AU5800® Clinical Chemistry Systems

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the instrument manual for additional instructions.

Results of this test should always be interpreted in conjunction with the patient's medical history clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthineers to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

Reagents

These reagents are qualified for use with the Calibrators/ Controls listed below only. Other material may be used however, for quality control purposes.

Assay	Catalog Number		
	28 mL Kit	115 mL Kit	1000 mL Kit
Emit® II Plus Ethyl Alcohol Assay	9K039UL	9K309UL	9K409UL
Emit® II Plus Ethyl Alcohol Assay	OSR9K229		

Emit® Ethyl Alcohol Calibrators

Negative	9K029UL
100 mg/dL	9K059UL

Emit® Ethyl Alcohol Controls

Low	9K049UL
High	9K079UL

Storage

Reagents which are in use may be stored on board the analyzer for up to 4 weeks or as long as the control results fall within acceptable limits.

Instrument

Calibration

Calibrate by running the negative and 100 mg/dL calibrators whenever a new lot of reagents is used or as indicated by control results.

Instrument Settings

See page 2.

Results

Results are reported in mg/dL or g/L.

NOTE: To convert mg/dL to g/L ethyl alcohol, multiply by 0.01.

Instrument Settings

General Screen

Reagent ID: 524

General		LIH		ISE		HbA1c		Calculated Tests		Range	
Test Name	Test # . User Defined		<	>	Type	Ser or Urine		Operation		Yes	
Sample Volume	2.8 uL		Dilution		0 uL		OD Limit				
Pre-Dilution Rate	1		Min OD		-2.0000		Max OD		2.5000		
Reagent Volume	R1(R1-1) 147 uL		Dilution		0 uL		Reagent OD limit				
			First Low		-2.0000		High		2.5000		
			Last Low		-2.0000		High		2.5000		
	R2(R2-2) 63 uL		Dilution		0 uL		Dynamic Range				
			Low		10.0		High		600.0		
Common Reagent	Type			Name			Correlation Factor	A	1.00000		B
Wave Length	Pri.	340 nm		Sec.	410 nm		Factor for Maker	A			B
Method											
Reaction Slope	+		Onboard Stability Period		User Defined		Day		User Defined		Hour
Measuring Point-1	First	12		Last	16						
Measuring Point-2	First			Last							
Linearity Limit											
Lag Time Check											

Range Screen

General		LIH		ISE		HbA1c		Calculated Tests		Range	
Test Name	Test # . User Defined		<	>	Type	Ser or Urine					
Value/Flag Level	Value										
	Low		High		Panic Value		Low		High		
	User Defined		User Defined		User Defined		User Defined		User Defined		
Specific Ranges											
	Sex	From Year	Month	To Year	Month	Low	High				
<input type="checkbox"/> 1											
<input type="checkbox"/> 2											
<input type="checkbox"/> 3											
<input type="checkbox"/> 4											
<input type="checkbox"/> 5											
<input type="checkbox"/> 6											
<input type="checkbox"/> 7	No Demographics										
<input type="checkbox"/> 8	Not within expected values										
Unit	*		Decimal Places		1						

*mg/dL or g/L

Calibration Specific Screen

Calibrators		Calibration Specific		STAT Table Calibration	
General		ISE			
Test Name		Test # . User Defined		Type Ser/Urine <input type="checkbox"/> Use Serum Cal	
Calibration Type		2AB		Formula Polygonal Counts 2	
<Calibration Parameters>		Slope Check +			
OD Range					
Calibrator	OD	Conc	Low	High	
Point-1	Cal # . User Defined	0.0	-2.0000	2.5000	
Point-2	Cal # . User Defined	*	-2.0000	2.5000	
Point-3					
Point-4					
Point-5					
Point-6					
Point-7					
Point-8					
Point-9					
Point-10					
Allowable Range Check					
<input type="checkbox"/> Reagent Blank					
<input type="checkbox"/> Calibration					
Advanced Calibration Operation Yes					
Interval (RB/ACAL) Lot/Lot					
<input type="checkbox"/> Lot Calibration					
<Point Cal. For Master Curve> No. of Correction Points <input type="checkbox"/> Use Master Curve <input type="checkbox"/>					
OD Range					
Calibrator	OD	Conc	Low	High	
Point-1					
Point-2					
Stability					
Reagent Blank <input type="checkbox"/> Day <input type="checkbox"/> Hour					
Calibration <input type="checkbox"/> Day <input type="checkbox"/> Hour					
MB Type Factor <input type="checkbox"/> 1-Point Calibration Point None with Conc-0					

*100.0 for mg/dL or 1.0 for g/L

Performance

Method Comparison

Clinical urine and serum specimens were tested using the Emit® II Plus Ethyl Alcohol Assay on the AU600® analyzer and on the SYVA®-30R analyzer.

Ethyl Alcohol Urine

Slope	0.97
Intercept	2.9
Correlation Coefficient	1.00
Number of Samples	50

Ethyl Alcohol Serum

Slope	1.05
Intercept	-9.6
Correlation Coefficient	0.993
Number of Samples	50

Precision

Within run precision was calculated according to NCCLS Guideline EP5-A by running 2 replicates of the 100 mg/dL calibrator with positive and negative controls twice a day for 20 days (N=80). Total precision was also calculated from these data.

Ethyl Alcohol

	Within Run Precision			Total Precision		
	100 mg/dL	40 mg/dL	300 mg/dL	100 mg/dL	40 mg/dL	300 mg/dL
Mean	100	40	300	100	40	296
SD	0.9	0.5	2.5	2.4	1.7	5.6
CV%	0.8	1.1	0.8	2.4	4.1	1.9

Analytical Recovery

Negative human urine and serum were spiked with ethyl alcohol at concentrations throughout the assay range. Recovery results on the AU600® are listed below.

Ethyl Alcohol: Urine

Concentration (mg/dL)	Mean (mg/dL)
25	26.7
80	82.4
200	208
400	370

Ethyl Alcohol: Serum

Concentration (mg/dL)	Mean (mg/dL)
25	27.5
80	80.1
200	190
400	375

Analytical Sensitivity

The sensitivity level of the Emit® II Plus Ethyl Alcohol Assay on the AU600® is <10 mg/dL ethyl alcohol. This level represents the lowest concentration of ethyl alcohol that can be distinguished from zero mg/dL with a confidence level of 95%.









NOTE: Performance on the AU480®, AU5800® and AU600® series analyzers has been shown to be equivalent.

EMIT® II Plus Ethyl Alcohol Assay

Application Sheet



















Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Manufacturer	5.1.1 ^a		Authorized representative in the European Community	5.1.2 ^a
	Use-by date	5.1.4 ^a		Authorized representative in Switzerland	Proprietary
	Catalog number	5.1.6 ^a		Batch code	5.1.5 ^a
	Consult Instructions for Use	5.4.3 ^a		Contains sufficient for <n> tests	5.5.5 ^a
	Internet URL address to access the electronic instructions for use	Proprietary		Version of Instructions for Use	Proprietary
	<i>In vitro</i> diagnostic medical device	5.5.1 ^a		Revision	Proprietary
RxOnly	Prescription device (US only)	FDA ^c		Unique Device Identifier	5.7.10 ^b
	CE Marking with Notified Body	EU IVDR ^d		CE Marking	EU IVDR ^d
	Temperature limit	5.3.7 ^a		Keep away from sunlight	5.3.2 ^a
	Upper limit of temperature	5.3.6 ^a		Lower limit of temperature	5.3.5 ^a
	Do not re-use	5.4.2 ^a		Do not freeze	Proprietary

EMIT® II Plus Ethyl Alcohol Assay

Application Sheet

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Recycle	1135 ^a		This way up	0623 ^a
	Biological risks	5.4.1 ^a		Caution	5.4.4 ^a
	Common Units	Proprietary		International System of Units	Proprietary
	Document face up ^f	1952 ^a		YYYY-MM-DD Date format (year-month-day)	N/A
	Non-sterile	Proprietary		YYYY-MM Date format (year-month)	N/A
	Reconstitution volume	Proprietary		CONTENTS Contents	Proprietary
	For forensic/toxicology use only	Proprietary		LEVEL Level	Proprietary
	Dropper	Proprietary		CASSETTE Cassette	Proprietary
	Not for self-testing	EU IVDR ^d		Not for near-patient testing	EU IVDR ^d

^a International Standard Organization (ISO). ISO 15223-1 Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied.

^b ISO 15223-1:2020-04.

^c Federal Register. Vol. 81, No 115. Wednesday, June 15, 2016. Rules and Regulations: 38911.

^d IVDR REGULATION (EU) 2017/746

^e International Standard Organization (ISO). ISO 7000 Graphical symbols for use on equipment.

^f Indicates Assay-eNote.

Syva®

EMIT® II Plus Ethyl Alcohol Assay

Application Sheet

Emit® and Syva® are trademarks of Siemens Healthineers.

AU®, AU480®, AU600®, and AU5800® are registered trademarks of Beckman Coulter, Inc.

For technical assistance:

**Beckman Coulter customers, contact the
Customer Technical Support Center at
1-800-854-3633 (USA & Canada)**

**In other countries, please contact your local
Beckman Coulter representative.**

**Siemens Healthineers customers, contact the Technical Solutions
Center at 1-800-227-8994 in the USA.**

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.
siemens-healthineers.com

© 2009–2024 Siemens Healthineers

All rights reserved.