## **Syva**<sup>®</sup>

## **Specific Gravity Validity Test**

## Specific Gravity Validity Test Application Sheet

## For the AU480<sup>®</sup> and AU5800<sup>®</sup> Clinical Chemistry Systems

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the instrument manual for additional instructions.

The parameters defined in this application sheet have been developed by Siemens Healthcare Diagnostics to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

### Reagents

These reagents are qualified for use with the Calibrators listed below only. Other material may be used for quality control purposes.

Test	Catalog Number		
	<u>Small</u>	Large	
Syva® Specific Gravity Validity Test	3T899UL	3T699UL	
Calibrator	Cata	alog Number	
Syva® Specific Gravity Validity Calibrator 1.0030	;	3T619UL	
Syva® Specific Gravity Validity Calibrator 1.0200	;	3T629UL	

#### Storage

Reagents which are in use may be stored on board the analyzer for up to 4 weeks or as long as the control results fall within acceptable limits.

Instrument

#### Calibration

Calibrate by running the Syva® Specific Gravity Calibrators listed in Table 1 in the calibration rack (yellow rack). Run a reagent blank (blue rack). Recalibrate whenever a new lot of reagent is used or as indicated by control results.

#### Table 1

Test	Calibrators		
Specific Gravity	1.0030	1.0200	

#### Instrument Settings

See page 2.

#### Results

Results are reported in specific gravity units.

Note: Specimen results <1.0010 should be reported as <1.0010.

## **Syva**®

## **Specific Gravity Validity Test**

## **Application Sheet**

## Instrument Settings

#### **General Screen**

General	LIH	ISE	HbA1c Calculated Tests Range
Test Name Use	er Defined <	> Type Urir	e Operation Yes
Sample Volume Pre-Dilution Rat	3.2 uL e 1	Dilution 0 uL	OD Limit Min OD <u>-2.0000</u> Max OD <u>2.5000</u>
Reagent Volume	R1(R1-1) 157 uL	Dilution 0 uL	Reagent OD limit First Low <u>-2.0000</u> High <u>2.5000</u> Last Low <u>2.0000</u> High <u>2.5000</u>
	R2(R2-2) 0 uL	Dilution 0 uL	2000 -2.0000 High 2.5000
Common Reagent Wave Length Method Reaction Slope Measuring Point-1	Type Pri. 600 nm END1 + First 0	Name Secnm Last	Dynamic RangeLow1.0010High1.0250Correlation FactorA1.0000B0.0000Factor for MakerABBOnboard Stability PeriodUser DefinedDayUser DefinedHour
Measuring Point-2 Linearity Limit Lag Time Check	First %	Last	

## **Range Screen**

General LIH ISE	Hb	A1c	Calculated Tests	Range
Test Name Test # . User Defined < > Typ	Urine Urine			
Level			Panic Value	
	Low	High	Low	High
Specific Ranges	User Defined	User Defined	User Defined	User Defined
From To	Low	Lliab		
	LOW			
7 No Demographics	-99999	99999		
8 Not within expected values				
Unit Decimal Places 4				

FOR DISTRIBUTION WITHIN THE US ONLY

## **Syva**®

## **Specific Gravity Validity Test**

## **Application Sheet**

### Calibration Parameters

Calibrators Calibration Specific STAT Table Calibration						
General	General ISE					
Test Nam	Test Name Test # . User Defined < > Type Urine Use Serum Cal					
Calibratio	n Type 🛛 🗛		Formula	Y = AX +	B Counts	S 2
<calibrati< td=""><td>on Parameters&gt;</td><td></td><td></td><td><b>–</b> .</td><td>-</td><td>Slope Check None</td></calibrati<>	on Parameters>			<b>–</b> .	-	Slope Check None
			_	Factor	Range	
	Calibrator	OD	Conc	Low	High	1
Point-1	Cal # . User Defined		1.0030	-9999999	9999999	
Point-2	Cal # . User Defined		1.0200	-9999999	9999999	Allowable Range Check
Point-3						
Point-4						
Point-5						
Point-6						Advanced Calibration
Point-7						
Point-8						Interval (RB/ACAL) Lot/Lot
Point-9						
Point-10	For Mostor Cumuos		Correction	Deinte	Lies Master C	
<point cal.="" curve="" for="" master=""> No. of Correction Points Use Master Curve</point>						
Factor Range						
	Calibrator	OD	Conc	Low	High	Stability User Defined
Point-1						Reagent Blank
Point-2						Calibration Day Hour
МВ Туре	Factor 1	-Point	Calibratio	on Point No	one with Con	ic-0

Performance

### **Method Comparison**

Sixty-five (65) human urine specimens were analyzed using the Specific Gravity Validity Test on the AU600® and the Hitachi 717 analyzers. Results are shown in the tables below.

Test	≥1.0030	<1.0030	% Agreement
SG 1.0030 Cutoff	47	16	96.9

Test	≥1.0200	<1.0200	% Agreement
SG 1.0200 Cutoff	16	49	100

#### Precision

Within run precision was calculated according to NCCLS Guideline EP5–A by running 2 replicates of each urine sample twice a day for 10 days (N=40). Total precision

was also calculated from these data. The following data are presented in specific gravity units.

	Within Run Precision			
	Level 1	Level 2	Level 3	Level 4
Mean	1.0005	1.0028	1.0204	1.0245
SD	0.0003	0.0003	0.0002	0.0004
% CV	0.03	0.03	0.02	0.04

	Total Precision			
	Level 1	Level 2	Level 3	Level 4
Mean	1.0005	1.0028	1.0204	1.0245
SD	0.0005	0.0003	0.0006	0.0007
% CV	0.05	0.03	0.05	0.07

FOR DISTRIBUTION WITHIN THE US ONLY

## **Syva**<sup>®</sup>

## **Specific Gravity Validity Test**

## **Application Sheet**

#### Linearity

Linearity was tested for the Syva® Specific Gravity Validity Test by preparing 8 samples throughout the assay range. Five replicates of each sample were analyzed on the AU600® analyzer. Linearity was acceptable across the assay range. Results are shown below.

Refractometer	AU600®
1.0000	0.9988
1.0010	1.0007
1.0020	1.0019
1.0030	1.0026
1.0040	1.0034
1.0150	1.0157
1.0200	1.0206
1.0250	1.0250

**NOTE:** Performance on the AU480®, AU5800® and AU600® series analyzers has been shown to be equivalent.

Syva® is a trademark of Siemens Healthcare Diagnostics.

AU®, AU480®, AU600®, and AU5800® are registered trademarks of Beckman Coulter, Inc.

#### For technical assistance, call Siemens Healthcare Diagnostics: 1-800-227-8994 in the USA

© 2009 Siemens Healthcare Diagnostics All rights reserved.

#### Rev. AU4SPGR.4

Siemens Healthcare Diagnostics Inc. 500 GBC Drive Newark, DE 19714 USA

Global Siemens Headquarters Siemens AG Wittelsbacherplatz 2 80333 Muenchen Germany

Global Siemens Healthcare Headquarters Siemens AG Healthcare Sector Henkestrasse 127 91052 Erlangen Germany Phone: +49 9131 84-0 siemens.com/healthcare Global Division Siemens Healthcare Diagnostics Inc. 511 Benedict Avenue Tarrytown, NY 10591 USA siemens.com/healthcare

FOR DISTRIBUTION WITHIN THE US ONLY