

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

Shading indicates technical content that differs from the previous version.

### Emit® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid Application Sheet

**For the AU400®, AU600®, AU640®,  
AU680®, AU2700®, AU5400®  
Clinical Chemistry Systems**

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the instrument manual for additional instructions.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthineers to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

### Reagents

These reagents are qualified for use with the Calibrators/Controls listed below only. Other material may be used however, for quality control purposes.

Assay	Catalog Number		
	28 mL Kit	115 mL Kit	1000 mL Kit
Emit® II Plus Cannabinoid Assay	9N039UL	9N029UL	9N129UL
Emit® II Plus Cannabinoid Assay	OSR9N229		

### Calibrator/Control

### Catalog Number

Emit® Calibrator/Control Level 0	9A509UL
Emit® Calibrator/Control Level 2	9A549UL
Emit® Calibrator/Control Level 3	9A569UL
Emit® Calibrator/Control Level 4	9A589UL
Emit® Calibrator/Control Level 5	9A609UL

### Storage

Reagents which are in use may be stored on board the analyzer for up to 4 weeks or as long as the control results fall within acceptable limits.

### Qualitative

### Calibration

Calibrate by running the Emit® Calibrator/Control Level 4. Select the appropriate Option to use based on how results should be flagged (see Results section below).

Run a reagent blank (blue rack) and place the appropriate calibrator in a white rack for Option 1 or in its assigned position in a yellow rack for Option 2 or 3. Recalibrate as indicated by control results.

### Results

Results are reported based on the Option chosen.

### Option

1. Results are printed as OD values. Positive samples **are NOT** flagged. Operator must compare the sample response to the cutoff calibrator response and determine if the sample is positive or negative.
2. The cutoff is normalized to 100. Positive samples are  $\geq 100$  and are flagged with a (P).
3. The cutoff is normalized to 0. Positive samples are  $\geq 0$  and are flagged with a (P).

### Instrument Settings

See page 2.

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

### Qualitative

#### Instrument Settings

#### General Screen

Reagent ID: 531

General		LIH	ISE	Range
Test Name	User Defined		◀ ▶	Type <b>Urine</b> Operation <b>Yes</b>
Sample Volume	3.0	uL	Dilution 0	uL Pre-Dilution Rate 1
Reagents R1 Volume	150	uL	Dilution 0	uL Min OD Max OD
R2 Volume	64	uL	Dilution 0	uL L -2.000 H 2.500
Wavelength Pri	340		Sec 410	
Method	FIXED			Reagent OD Limit First L -2.000 First H 2.500
Reaction Slope	+			Last L -2.000 Last H 2.500
Measuring Point 1 First	15		Last 17	Dynamic Range L User Defined H User Defined
Measuring Point 2 First			Last	Correlation Factor
Linearity		%		A 1 B †
No-Lag-Time				Onboard Stability Period User Defined

† – Correlation Factor B: For Results Option 1 enter 0; for Option 2 enter 0; for Option 3 enter -100

#### Range Screen

General		LIH	ISE	Range
Test Name	User Defined		◀ ▶	Type <b>Urine</b>
Value Flag	Flag		Level L -9999	Level H ††
Normal Ranges	Age L	Age H		
	Sex	Year	Month	Year
1				
2				
3				
4				
5				
6				
7 None Selected				
8 Out of Range				
	L	H		
Panic Value			Unit	Decimal Places 0

†† – Level High: For Results Option 1 enter 9999; for Option 2 enter 100; for Option 3 enter 0

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

### Calibration Specific Screen — Results Option 1

General   ISE					
Test Name	User Defined		Type		Urine
Calibration Type	MB		Formula	Y=AX+B	
	Counts	2		Process	OD
	Cal No	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1					
Point 2					
Point 3					
Point 4					
Point 5					
Point 6					
Point 7					
1-Point Cal Point	<input type="checkbox"/>	With CONC-0		Slope Check	
Advanced Calibration			User Def		
MB Type Factor	1000.000		Calibration Stability Period	User Defined	

### Calibration Specific Screen — Results Option 2 or 3

General   ISE					
Test Name	User Defined		Type		Urine
Calibration Type	AB		Formula	Y=AX+B	
	Counts	2		Process	OD
	Cal No	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1	User Def	User Def	100	-9999	9999
Point 2					
Point 3					
Point 4					
Point 5					
Point 6					
Point 7					
1-Point Cal Point	<input type="checkbox"/>	With CONC-0		Slope Check	None
Advanced Calibration			User Def		
MB Type Factor			Calibration Stability Period	User Defined	

**Syva®**

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

### Semi-Quantitative

### Calibration

Prepare a calibration curve by running the Calibrator/  
Controls in the order listed in Table 1. Recalibrate as  
indicated by control results. Run a reagent blank (blue  
rack) daily to ensure consistent day to day control results.

**Table 1**

Assay	Level 2 ng/mL	Level 3 ng/mL	Level 4 ng/mL	Level 5 ng/mL
Cannabinoid	20	50	100	200

### Results

Results are reported as ng/mL.

### Instrument Settings

#### General Screen

General		LIH	ISE	Range
Test Name	User Defined		◀ ▶	Type <b>Urine</b> Operation <b>Yes</b>
Sample	Volume	3.0	uL	Dilution 0 uL
Reagents	R1 Volume	150	uL	Dilution 0 uL
	R2 Volume	64	uL	Dilution 0 uL
Wavelength	Pri	340		Sec 410
Method	FIXED			
Reaction Slope	+			
Measuring Point 1	First	15		Last 17
Measuring Point 2	First			Last
Linearity				
No-Lag-Time				
Pre-Dilution Rate 1				
Min OD Max OD				
Reagent OD Limit				
L -2.000 H 2.500				
First L -2.000 First H 2.500				
Last L -2.000 Last H 2.500				
Dynamic Range				
L 28 H 200				
Correlation Factor				
A 1 B 0				
Onboard Stability Period User Defined				

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

### Range Screen

General   LIH   ISE   <b>Range</b>									
Test Name		<b>User Defined</b>		◀ ▶		Type		<b>Urine</b>	
Value Flag		<b>Value</b>		Level L		<b>-9999</b>		Level H	
Normal Ranges		Age L		Age H					
	Sex	Year	Month	Year	Month	L		H	
<input type="checkbox"/> 1									
<input type="checkbox"/> 2									
<input type="checkbox"/> 3									
<input type="checkbox"/> 4									
<input type="checkbox"/> 5									
<input type="checkbox"/> 6									
7 None Selected									
8 Out of Range									
Panic Value		L		H		Unit		<b>ng/mL</b>	
						Decimal Places		<b>0</b>	

### Calibration Specific Screen

General   ISE									
Test Name		<b>User Defined</b>		◀ ▶		Type		<b>Urine</b>	
Calibration Type		<b>4AB</b>		Formula		<b>Polygonal</b>		Counts	
								<b>2</b>	
								Process	
								<b>CONC</b>	
	Cal No	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 1	<b>User Def</b>	<b>User Def</b>	<b>20</b>	<b>-2.000</b>	<b>2.500</b>				
Point 2	<b>User Def</b>	<b>User Def</b>	<b>50</b>	<b>-2.000</b>	<b>2.500</b>				
Point 3	<b>User Def</b>	<b>User Def</b>	<b>100</b>	<b>-2.000</b>	<b>2.500</b>				
Point 4	<b>User Def</b>	<b>User Def</b>	<b>200</b>	<b>-2.000</b>	<b>2.500</b>				
Point 5									
Point 6									
Point 7									
1-Point Cal Point		<input type="checkbox"/> With CONC-0		Slope Check		<b>+</b>		Advanced Calibration	
								<b>User Def</b>	
MB Type Factor				Calibration Stability Period		<b>User Defined</b>			

## EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

### Performance

#### Method Comparison

Clinical urine specimens were tested using the Emit® II Plus Cannabinoid assay on the AU400 analyzer and using the corresponding Emit® II Plus assay on the SYVA®-30R analyzer. Specimens positive by either method contained 11-nor- $\Delta^9$ -THC-9-COOH. The results from the AU400 are listed below along with the percent agreement with the SYVA®-30R.

		SYVA®-30R	
		Positive	Negative
AU400	Negative	0	51
	Positive	46	3
		Agreement = 97%	

#### Precision

Within run precision was calculated according to NCCLS Guideline EP5-A by running 2 replicates of the cutoff Calibrator/Control and positive and negative controls twice a day for 20 days (N=80). Total precision was also calculated from these data.

#### Qualitative (mA/min)

	Repeatability			Within-Lab		
	Cutoff	Control	Control	Cutoff	Control	Control
	Cal.	75%	125%	Cal.	75%	125%
Mean	396	366	430	396	366	430
SD	3.6	3.3	3.8	9.4	10.7	8.5
CV%	0.92	0.91	0.89	2.4	2.9	2.0

#### Semi-Quantitative (ng/mL)

	Repeatability			Within-Lab		
	Cutoff	Control	Control	Cutoff	Control	Control
	Cal.	75%	125%	Cal.	75%	125%
Mean	95	72	138	95	72	138
SD	3.5	2.6	5.7	7.4	7.6	10.7
CV%	3.7	3.6	4.1	7.8	10.6	7.8

#### Analytical Recovery

Negative human urine specimens were spiked with concentrations of 11-nor- $\Delta^9$ -THC-9-COOH.

Qualitative analyses of the specimens spiked with drug concentrations lower than the cutoff concentration were correctly identified as negative 100% of the time. Specimens spiked with drug concentrations greater than the cutoff were correctly identified as positive 100% of the time.

Results from semi-quantitative analysis of the specimens are listed below.

#### Semi-Quantitative

Concentration (ng/mL)	Mean (ng/mL)
40	38
50	49
75	76
125	163
150	184
180	196
205	201

#### Analytical Sensitivity

The sensitivity level of the Emit® II Plus Semi-Quantitative Cannabinoid 100 Assay on the AU400 is 28 ng/mL 11-nor- $\Delta^9$ -THC-9-COOH. This level represents the lowest concentration of 11-nor- $\Delta^9$ -THC-9-COOH that can be distinguished from zero ng/mL with a confidence level of 95%.

**NOTE:** Performance on the AU400, AU600, AU640, AU680, AU2700, and AU5400 series analyzers has been shown to be equivalent.

#### Specificity Exception

Compounds listed in the following table produce a positive result at the 50 ng/mL cutoff at the concentration listed

Compound	Concentration (ng/mL)
11-Hydroxy- $\Delta^8$ -THC	233













**NOTE:** Performance on the AU480®, AU5800® and AU400® series analyzers has been shown to be equivalent.

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet



















### Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Manufacturer	5.1.1 <sup>a</sup>		Authorized representative in the European Community	5.1.2 <sup>a</sup>
	Use-by date	5.1.4 <sup>a</sup>		Authorized representative in Switzerland	Proprietary
	Catalog number	5.1.6 <sup>a</sup>		Batch code	5.1.5 <sup>a</sup>
	Consult Instructions for Use	5.4.3 <sup>a</sup>		Contains sufficient for <n> tests	5.5.5 <sup>a</sup>
	Internet URL address to access the electronic instructions for use	Proprietary		Version of Instructions for Use	Proprietary
	<i>In vitro</i> diagnostic medical device	5.5.1 <sup>a</sup>		Revision	Proprietary
<b>RxOnly</b>	Prescription device (US only)	FDA <sup>c</sup>		Unique Device Identifier	5.7.10 <sup>b</sup>
	CE Marking with Notified Body	EU IVDR <sup>d</sup>		CE Marking	EU IVDR <sup>d</sup>
	Temperature limit	5.3.7 <sup>a</sup>		Keep away from sunlight	5.3.2 <sup>a</sup>
	Upper limit of temperature	5.3.6 <sup>a</sup>		Lower limit of temperature	5.3.5 <sup>a</sup>
	Do not re-use	5.4.2 <sup>a</sup>		Do not freeze	Proprietary

# EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid

## Application Sheet

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Recycle	1135 <sup>e</sup>		This way up	0623 <sup>e</sup>
	Biological risks	5.4.1 <sup>a</sup>		Caution	5.4.4 <sup>a</sup>
	Common Units	Proprietary		International System of Units	Proprietary
	Document face up <sup>f</sup>	1952 <sup>e</sup>		Date format (year-month-day)	N/A
	Non-sterile	Proprietary		Date format (year-month)	N/A
	Reconstitution volume	Proprietary		Contents	Proprietary
	For forensic/toxicology use only	Proprietary		Level	Proprietary
	Dropper	Proprietary		Cassette	Proprietary
	Not for self-testing	EU IVDR <sup>d</sup>		Not for near-patient testing	EU IVDR <sup>d</sup>

<sup>a</sup> International Standard Organization (ISO). ISO 15223-1 Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied.

<sup>b</sup> ISO 15223-1:2020-04.

<sup>c</sup> Federal Register. Vol. 81, No 115. Wednesday, June 15, 2016. Rules and Regulations: 38911.

<sup>d</sup> IVDR REGULATION (EU) 2017/746

<sup>e</sup> International Standard Organization (ISO). ISO 7000 Graphical symbols for use on equipment.

<sup>f</sup> Indicates Assay-eNote.



**Syva®**

## **EMIT® II Plus Cannabinoid 100/ Semi-Quantitative Cannabinoid**

## **Application Sheet**

Emit® and Syva® are trademarks of Siemens Healthineers.

AU®, AU400®, AU600®, AU640®, AU680®, AU2700®, and AU5400® are registered trademarks of Beckman Coulter, Inc.

### **For technical assistance:**

**Beckman Coulter customers, contact the  
Customer Technical Support Center at  
1-800-854-3633 (USA & Canada)**

**In other countries, please contact your local  
Beckman Coulter representative.**

**Siemens Healthineers customers,  
contact the Technical Solutions Center  
at 1-800-227-8994 In the USA.**

### **Technical Assistance**

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.  
[siemens-healthineers.com](https://www.siemens-healthineers.com)

© 2010–2024 Siemens Healthineers

All rights reserved.