

EMIT[®] II Plus 6-Acetylmorphine Assay

Shading indicates technical content that differs from the previous version.

Emit[®] II Plus 6-Acetylmorphine Assay Application Sheet

For the AU480[®] and AU5800[®] Clinical Chemistry Systems

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the instrument manual for additional instructions.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthineers to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

Reagents

These reagents are qualified for use with the Calibrators/ Controls listed below only. Other material may be used, however, for quality control purposes.

Emit II Plus 6-Acetylmorphine Assay	Catalog Number		
28 mL Kit	9R039UL		
115 mL Kit	9R029UL		
1000 mL Kit	9R129UL		

Application Sheet

Calibrator/Control	Catalog Number
Emit Calibrator/Control Level 0	9A509UL
Emit II Plus 6-AM/Ecstasy Calibrator/Control Level 1	9R529UL
Emit II Plus 6-AM/Ecstasy Calibrator/Control Level 2	9R549UL
Emit II Plus 6-AM/Ecstasy Calibrator/Control Level 3	9R569UL
Emit II Plus 6-AM/Ecstasy Calibrator/Control Level 4	9R589UL

Preparation, Storage and Stability of Reagent 2

R2-9R039UL: Add 13 mL of distilled or deionized water

R2-9R029UL: Add 45 mL of distilled or deionized water

R2-9R129UL: Add 450 mL of distilled or deionized water

After preparation, allow reagent to equilibrate on the analyzer for 1 hour.

Table 1 – Preparation, Storage and Stability of Reagent 2

Component	Storage Temperature	Reconstitution Volume	Unopened Stability	Reagent Stability After Reconstitution
Reagent 2 9R039UL	2–8°C	13 mL	Exp. date	30 days
Reagent 2 9R029UL	2–8°C	45 mL	Exp. date	30 days
Reagent 2 9R129UL	2–8°C	450 mL	Exp. date	30 days

Storage

Reagents which are in use may be stored on board the analyzer for up to 4 weeks.

Calibration

Results

EMIT[®] II Plus 6-Acetylmorphine Assay

Qualitative

Calibrate by running the Level for the cutoff listed in

flagged (see Results section below).

indicated by control results.

Table 2. Recalibrate with each new hydration. Select the appropriate Option to use based on how results should be

Run a reagent blank (blue rack) and place the appropriate

calibrator in a white rack for Option 1 or in its assigned position in a yellow rack for Option 2 or 3. Recalibrate as

Results are reported based on the Option chosen.

Option

- 1. Results are printed as OD values. Positive samples are NOT flagged. Operator must compare the sample response to the cutoff calibrator response and determine if the sample is positive or negative.
- 2. The cutoff is normalized to 100. Positive samples are \geq 100 and are flagged with a (P).
- 3. The cutoff is normalized to 0. Positive samples are ≥ 0 and are flagged with a (P).

Table 2

	Level 2
Assay	ng/mL
6-Acetylmorphine	10

Instrument Settings

See page 2.

Qualitative

Instrument Settings

General LIH ISE	HbA1c Calculated Tests Range
Test Name Test # . User Defined < > Type	Urine Operation Yes
Sample Volume7.5uLDilution0Pre-Dilution Rate1	uL OD Limit Min OD <u>-2.0000</u> Max OD <u>2.5000</u>
Reagent Volume R1(R1-1) 80 uL Dilution 0	uL Reagent OD limit First Low <u>-2.0000</u> High <u>2.5000</u> Last Low <u>-2.0000</u> High <u>2.5000</u>
R2(R2-2) 40 uL Dilution 0]uL
	Dynamic Range Low User Def High User Def
Common Reagent Type Name	Correlation Factor A 1 B †
Wave Length Pri. 340 nm Sec. 410	nm Factor for Maker A B
Method FIXED	
Reaction Slope +	Onboard Stability Period User Day User Hour
Measuring Point-1 First 13 Last 15	
Measuring Point-2 First Last	
Linearity Limit %	
Lag Time Check	



Application Sheet





EMIT® II Plus 6-Acetylmorphine Assay

Application Sheet

Range Screen

General LIH ISE		A1c	Calculated Tests	Range
Test Name Test # . User Defined > Type	De Urine			
Value/Flag Flag			Panic Value	
	Low	High	Low	High
Specific Ranges	-999999	†		
From To				
Sex Year Month Year Month	Low	High	1	
]		
] []	
7 No Demographics] []	
8 Not within expected values]]	
Unit Decimal Places 0				

† – Level High: For Results Option 1 enter 9999999; for Option 2 enter 100; for Option 3 enter 0

Calibration Specific Screen – Results Option 1

Calibrators Calibration Specific STAT Table Calibration
General ISE
Test Name Test # . User Defined < > Type Urine
Calibration Type MB Formula Y = AX + B Counts 2
<calibration parameters=""> Slope Check None</calibration>
Factor Range
Calibrator OD Conc Low High
Point-1
Point-2 Allowable Range Check
Point-3 Reagent Blank
Point-4 Calibration
Point-5
Point-6 Advanced Calibration
Point-7 Operation Yes
Point-8 Interval (RB/ACAL)
Point-9
Point-10
<point cal.="" curve="" for="" master=""> No. of Correction Points Use Master Curve</point>
OD Range
Calibrator OD Conc Low High Stability User Defined
Point-1 Reagent Blank Day Hou
Point-2 Calibration Day Hou
MB Type Factor 1000.000 1-Point Calibration Point with Conc-0



EMIT® II Plus 6-Acetylmorphine Assay

Application Sheet

Calibrators Calibration Specific STAT Table Calibration General ISE						
Test Nam	e Test # . User De	fined	< >	Тур	De Urine	Use Serum Cal
Calibratio	n Type AB	•	Formula	Y = AX +	B Counts	s 2 Slope Check None
				Factor	Range	
	Calibrator	OD	Conc	Low	High	
Point-1	Cal # . User Defined		100	-999999	999999	
Point-2 Point-3 Point-4 Point-5 Point-6 Point-7 Point-8 Point-9 Point-10 <point cal.<="" td=""><td>For Master Curve></td><td>• No. of</td><td>Correction</td><td>Points</td><td>Use Master C</td><td>Allowable Range Check Reagent Blank Calibration Advanced Calibration Operation Interval (RB/ACAL) Lot/Lot Urve</td></point>	For Master Curve>	• No. of	Correction	Points	Use Master C	Allowable Range Check Reagent Blank Calibration Advanced Calibration Operation Interval (RB/ACAL) Lot/Lot Urve
				OD Ra	ange	
Point-1 Point-2 MB Type	Calibrator Factor	OD 1-Point	Conc Calibratio	Low	High with Cor	Stability Reagent Blank Calibration nc-0

Calibration Specific Screen – Results Option 2 or 3



EMIT® II Plus 6-Acetylmorphine Assay

Application Sheet

Semi-Quantitative

Calibration

Prepare a calibration curve by running the Calibrator/ Controls in the order listed in Table 3. After each hydration run controls. Recalibrate as indicated by control results. Run a reagent blank (blue rack) daily to ensure consistent day to day control results.

Table 3

Assay	Level 1	Level 2	Level 3	Level 4
	ng/mL	ng/mL	ng/mL	ng/mL
6-Acetylmorphine	5	10	15	20

Results

Results are reported as ng/mL.

Instrument Settings

General Screen

General	LIH	ISE	HbA1c Calculated Tests Range
Test Name Test # .	User Defined <	> Type Urin	e Operation Yes
Sample Volume Pre-Dilution Rat	7.5 uL	Dilution ⁰ uL	OD Limit Min OD <u>-2.0000</u> Max OD <u>2.5000</u>
Reagent Volume	R1(R1-1) 80 uL	Dilution 0 uL	Reagent OD limit
			First Low <u>-2.0000</u> High <u>2.5000</u>
	R2(R2-2) 40 uL	Dilution 0 uL	Last Low <u>-2.0000</u> High <u>2.5000</u>
			Dynamic Range Low 2.9 High 20
Common Reagent	Туре	Name	Correlation Factor A 1 B 0
Wave Length	Pri. 340 nm	Sec. 410 nm	Factor for Maker A B
Method	FIXED		
Reaction Slope	+		Onboard Stability Period User Day Defined Hour
Measuring Point-1	First 13	Last 15	Defined
Measuring Point-2	First	Last	
Linearity Limit	%		
Lag Time Check			



EMIT® II Plus 6-Acetylmorphine Assay

Application Sheet

Range Screen

General LIH ISE	Hb	A1c	Calculated Tests	Range			
Test Name Test # . User Defined < > Type Urine							
Value/Flag Value			Panic Value				
	Low	High	Low	High			
Specific Ranges	User Defined	10					
From To							
Sex Year Month Year Month	Low	High	٦				
7 No Demographics							
8 Not within expected values							
Unit ng/mL Decimal Places 1							

Calibration Specific Screen

Calibrators Calibration Specific STAT Table Calibration General ISE							
Test Name Test # . User Defined < > Type Urine Use Serum Cal							
Calibratio	Calibration Type 4AB Formula † Counts 2						
<calibrati< td=""><td>on Parameters></td><td></td><td></td><td></td><td></td><td>Slope Check +</td></calibrati<>	on Parameters>					Slope Check +	
				OD Ra	inge		
	Calibrator	OD	Conc	Low	High		
Point-1	Cal # . User Defined		5	-2.0000	2.5000	-	
Point-2	Cal # . User Defined		10	-2.0000	2.5000	Allowable Range Check	
Point-3	Cal # . User Defined		15	-2.0000	2.5000	Reagent Blank	
Point-4	Cal # . User Defined		20	-2.0000	2.5000	☐ Calibration	
Point-5							
Point-6						Advanced Calibration	
Point-7						Operation Yes	
Point-8							
Point-9						Interval (RB/ACAL) Lot/Lot	
Point-10						Lot Calibration	
<point cal.="" curve="" for="" master=""> No. of Correction Points Use Master Curve</point>							
OD Range							
	Calibrator	OD	Conc	Low	High	Stability User Defined	
Point-1						Reagent Blank	
Point-2						Calibration	
MB Type Factor 1-Point Calibration Point None with Conc-0							

† – Formula: enter EIA Type 1 or Polygonal



EMIT[®] II Plus 6-Acetylmorphine Assay

Application Sheet

Performance

Data was collected on the AU400® analyzer. Performance on the AU400, AU480, and AU5800 series analyzers has been shown to be equivalent.

Method Comparison

Clinical urine specimens were tested using the Emit® II Plus 6-Acetylmorphine Assay on the AU400 analyzer and on the Viva-E® analyzer. The results are listed below along with the percent agreement with the Viva-E® analyzer.

Assay	Positive	Negative	% Agreement	
6-Acetylmorphine	50	55	100	

Precision

Repeatability was calculated according to CLSI Guideline EP5-A2 by running 2 replicates of a cutoff Calibrator/ Control and positive and negative controls twice a day for 20 days (N=80). Within-lab was also calculated from these data. The following data are presented in mAU/min.

6-Acetylmorphine

	R	epeatabilit	y	Within-Lab			
	Control (Control		Control	Control	
	Cutoff	75%	125%	Cutoff	75%	125%	
Mean	575	526	590	575	526	590	
SD	4.8	4.1	5.4	15.2	17.5	15.3	
CV%	0.8	0.8	0.9	2.7	3.3	2.6	

Analytical Recovery

Negative human urine specimens were spiked with concentrations of 6-Acetylmorphine.

Qualitative analysis of the specimens spiked with drug concentrations lower than the cutoff concentration were correctly identified as negative 100% of the time. Specimens spiked with drug concentrations greater than the cutoff were correctly identified as positive 100% of the time.

Results from semi-quantitative analysis of the specimens are listed below.

6-Acetylmorphine

Concentration (ng/mL)	EIA Type 1 Mean (ng/mL)	Polygonal Mean (ng/mL)
5.0	5.6	5.5
7.5	8.2	8.5
12.5	13.5	13.5
15.0	16.3	16.2
17.5	18.9	18.9
20.0	21.2	21.3

Syva®

EMIT® II Plus 6-Acetylmorphine Assay

Application Sheet

Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Manufacturer	5.1.1ª	EC REP	Authorized representative in the European Community	5.1.2ª
	Use-by date	5.1.4ª	CH REP	Authorized representative in Switzerland	Proprietary
REF	Catalog number	5.1.6ª	LOT	Batch code	5.1.5ª
ĺ	Consult Instructions for Use	5.4.3ª	Σ	Contains sufficient for <n> tests</n>	5.5.5ª
[] i]	Internet URL address to access the electronic instructions for use	Proprietary	Rev. XX	Version of Instructions for Use	Proprietary
IVD	In vitro diagnostic medical device	5.5.1ª	Rev. REVISION	Revision	Proprietary
RxOnly	Prescription device (US only)	FDA°	UDI	Unique Device Identifier	5.7.10 ^b
C xxxx	CE Marking with Notified Body	EU IVDR ^d	CE	CE Marking	EU IVDR ^d
X	Temperature limit	5.3.7ª	×	Keep away from sunlight	5.3.2ª
X	Upper limit of temperature	5.3.6ª	X	Lower limit of temperature	5.3.5ª
(Do not re-use	5.4.2ª		Do not freeze	Proprietary



EMIT® II Plus 6-Acetylmorphine Assay

Application Sheet

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Recycle	1135°	<u> </u>	This way up	0623°
Ś	Biological risks	5.4.1ª	\triangle	Caution	5.4.4ª
UNITS C	Common Units	Proprietary	UNITS SI	International System of Units	Proprietary
	Document face up ^r	1952°	YYYY-MM-DD	Date format (year-month-day)	N/A
NON	Non-sterile	Proprietary	ҮҮҮҮ-ММ	Date format (year-month)	N/A
\longrightarrow	Reconstitution volume	Proprietary	CONTENTS	Contents	Proprietary
FORENSIC TOXICOLOGY	For forensic/toxicology use only	Proprietary	LEVEL	Level	Proprietary
DROPPER	Dropper	Proprietary	CASSETTE	Cassette	Proprietary
	Not for self-testing	EU IVDR ^d		Not for near-patient testing	EU IVDR ^d

^a International Standard Organization (ISO). ISO 15223-1 Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied.

^b ISO 15223-1:2020-04.

^c Federal Register. Vol. 81, No 115. Wednesday, June 15, 2016. Rules and Regulations: 38911.

^d IVDR REGULATION (EU) 2017/746

^e International Standard Organization (ISO). ISO 7000 Graphical symbols for use on equipment.

f Indicates Assay-eNote.



EMIT[®] II Plus 6-Acetylmorphine Assay

Application Sheet

For technical assistance:

Beckman Coulter customers, contact the Customer Technical Support Center at 1-800-854-3633 (USA & Canada)

In other countries, please contact your local Beckman Coulter representative.

Siemens Healthineers customers, contact the Technical Solutions Center at 1-800-227-8994 in the USA

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established. siemens-healthineers.com

Emit®, Viva-E® and Syva® are trademarks of Siemens Healthineers.

AU®, AU400®, AU480®, and AU5800® are registered trademarks of Beckman Coulter, Inc.

© 2011–2024 Siemens Healthineers All rights reserved.

Siemens Healthcare Diagnostics Inc. 500 GBC Drive Newark, DE 19714 USA

CE 0197 963q_EN Rev. 06, 2024-04