

ADVIA® Chemistry XPT

Systems

Phenytoin_2 (PHNY_2)

Current Revision and Date ^a	Rev. D, 2020-02
Product Name	ADVIA® Chemistry Phenytoin_2 (PHNY_2) Reagents REF 10377506
Systems	ADVIA Chemistry XPT System
Materials Required but Not Provided	ADVIA Chemistry Drug Calibrator I REF 10376770 Reagent container adapters Commercially available controls
Specimen Types	Human serum, plasma
Assay Principle	Enzyme Multiplied Immunoassay Technique (EMIT)
Assay Range ^b	Serum: 1.0–(34.6–40.0) µg/mL (4.0–(137.0–158.4 µmol/L)) Plasma: 1.0–(34.6–40.0) µg/mL (4.0–(137.0–158.4 µmol/L))
Reagent Storage	2–8°C
Reagent On-System Stability	28 days
Reagent Code	74727

^a In Rev. C or later, a vertical bar in the margin indicates a technical update to the previous version.

^b The phenytoin concentration in the ADVIA Chemistry Drug Calibrator I - Level 5 varies from 34.6–40.0 µg/mL (137.0–158.4 µmol/L).

Intended Use

For *in vitro* diagnostic use in the quantitative analysis of phenytoin in human serum or plasma on ADVIA® Chemistry XPT systems.

Summary and Explanation

Monitoring phenytoin concentrations in serum, along with careful clinical assessment, is the most effective means of improving seizure control, reducing the risk of toxicity, and minimizing the need for additional anticonvulsant medication, based on the following reasons:^{1,2}

- Serum phenytoin concentrations correlate better with pharmacologic activity than does dosage because of individual differences in absorption, metabolism, disease states, concomitant medication, and compliance. Serum concentration monitoring helps physicians individualize dosage regimens.
- The hepatic enzyme system for metabolizing phenytoin can become saturated within the drug's therapeutic range. When this occurs, small dosage alterations can lead to unexpected drug accumulation and clinical toxicity.
- Phenytoin is safe and effective only in a narrow range of serum concentrations.

Techniques historically used to monitor serum phenytoin concentrations include chromatographic assays, and immunoassay.¹⁻⁴

Principles of the Procedure

The ADVIA Chemistry Phenytoin_2 (PHNY_2) assay is a homogeneous immunoassay that is used for the quantitative analysis of phenytoin (free and protein-bound) in serum or plasma.^{5,6} The ADVIA Chemistry PHNY_2 assay uses the Syva® Emit® 2000 Phenytoin reagent in ADVIA Chemistry containers.

The ADVIA Chemistry PHNY_2 assay is based on competition for antibody binding sites between phenytoin in the sample and phenytoin labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH). Upon binding to the antibody, enzyme activity decreases. As a result, the phenytoin concentration in the samples can be measured in terms of enzyme activity.

Active enzyme converts oxidized nicotinamide adenine dinucleotide (NAD) to reduced nicotinamide adenine dinucleotide (NADH), resulting in an absorbance change that is measured spectrophotometrically. Endogenous serum G6PDH does not interfere because the coenzyme functions only with the bacterial (*Leuconostoc mesenteroides*) enzyme used in the assay.

Reagents

Reagent	Description	Storage	Reagent Stability
REF 10377506	ADVIA Chemistry Phenytoin_2 (PHNY_2) Reagents		
Phenytoin_2 Reagent 1 PHNY_2 R1	11.4 mL in 20-mL containers Mouse monoclonal antibodies reactive to phenytoin (53.4 µg/mL) ^a Glucose-6-Phosphate (G6P) (22 mmol/L) Nicotinamide adenine dinucleotide (NAD) (18 mmol/L) Methylisothiazolinone (MIT) (0.1%) Bovine serum albumin, preservatives and stabilizers	2–8°C	Unopened: Stable until the expiration date on product. On-system: 28 days
Phenytoin_2 Reagent 2 PHNY_2 R2	7.0 mL in 20-mL containers Phenytoin labeled with bacterial G6PDH (0.24 U/mL) ^a Tris buffer Methylisothiazolinone (MIT) (0.1%) Bovine serum albumin, preservatives and stabilizers	2–8°C	Unopened: Stable until the expiration date on product. On-system: 28 days

^a The antibody titer and enzyme conjugate activity may vary from lot to lot.

Warnings and Precautions

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare.



Caution

This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner, and in compliance with prevailing regulatory requirements.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

For *in vitro* diagnostic use.

Preparing Reagents

All reagents are liquid and ready to use.

Before use, gently invert the capped reagent to disrupt bubbles and ensure homogeneity. If bubbles still exist or foam is present, use a clean transfer pipette to aspirate them from the reagent container prior to use.

Storing and Stability

Unopened reagents are stable until the expiration date printed on the product label when stored at 2–8°C.

Do not freeze reagents. Avoid prolonged exposure to temperatures above 32°C.

Specimen Collection and Handling

Siemens Healthcare Diagnostics validated serum and plasma (lithium heparin, sodium heparin, potassium EDTA, citrate, and sodium fluoride/potassium oxalate anticoagulants) for the ADVIA Chemistry PHNY_2 assay.

The purpose of handling and storage information is to provide guidance to users. It is the responsibility of the individual laboratory to use all available references and/or its own studies when establishing alternate stability criteria to meet specific needs.

Collecting the Specimen

Follow these guidelines for specimens used for this assay:

- Do not use whole blood for this assay.
- Some sample dilution may occur when samples are collected in tubes that contain citrate anticoagulant. Consider the amount of dilution and the need to correct for the dilution when interpreting assay results for these samples.
- Certain pharmacokinetic factors influence the correct time of sample collection after the last drug dose. These factors include dosage form, mode of administration, concomitant drug therapy, and biological variations affecting drug disposition.^{1,2}
- Draw a sample within 2 to 4 hours after an intravenous loading dose and, at steady state, collect a specimen representing the trough level just before the next scheduled dose.¹
- Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁷ Follow the instructions provided with your specimen collection device for use and processing.⁸
- Complete clot formation should take place before centrifugation.
- Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of 2 hours from the time of collection.⁹
- Specimens should be free of particulate matter.
- Specimens should be as fresh as possible.
- Handle all human serum and plasma samples as potentially infectious. Handle and dispose of human samples in accordance with established good laboratory practices and in compliance with all local and regulatory requirements.¹⁰⁻¹²

Storing the Specimen

Store the serum and plasma refrigerated at 2–8°C.

Specimens may be stored refrigerated at 2–8°C for 1 month, or stored frozen for up to 3 months.¹³

Transporting the Specimen

When transporting specimens, maintain the specimen temperature at 2–8°C.¹³

Procedure

Materials Provided

The following materials are provided:

Item	Contents	Number of Tests
REF 10377506	Reagent 1: 4 × 20-mL containers Reagent 2: 4 × 20-mL containers	4 × 100

Materials Required but Not Provided

The following materials are required to perform this assay, but are not provided:

Item	Description
REF 10376770	ADVIA Chemistry Drug Calibrator I
REF 10316975	20-mL reagent container adapter for 40-mL slot
REF 10723030	20-mL reagent container adapter for 70-mL slot
	Commercially available control materials

Assay Procedure

Sampling, reagent delivery, mixing, and processing are automatically performed by the ADVIA Chemistry system.

For detailed information on performing the procedure, refer to the system operating instructions.

Preparing the System

For detailed information on preparing the system, refer to the system operating instructions.

Preparing the Samples

Before placing samples on the system, ensure that samples have the following characteristics:

- Samples are free of fibrin or other particulate matter.
- Samples are free of bubbles.

On-System Stability

The ADVIA Chemistry PHNY_2 reagents are stable on the system for 28 days.

Do not use reagents beyond the expiration date.

Performing Calibration

To calibrate the ADVIA Chemistry PHNY_2 assay, use ADVIA Chemistry Drug Calibrator I, REF 10376770.

Enter the lot-specific calibrator values that are provided with each lot of calibrator. Perform the calibration as described in the calibrator instructions for use.

Calibration Frequency

Calibrate the assay every 28 days.

Calibrate the assay after the following events:

- When the reagent lot number changes
- When a reagent pack is replaced by a new reagent pack with the same lot number, and the previous reagent pack was recalibrated during use
- After replacing critical optical or hydraulic components
- When indicated by quality control procedures

Individual laboratory quality control programs and procedures may require more frequent calibration.

Reagent Blank (RBL) Frequency

The ADVIA Chemistry system measures the RBL during assay calibration.

Run an additional RBL on the same reagent pack every day.

Run an additional RBL when a reagent pack is replaced by a new reagent pack with the same lot number, and an additional reagent blank was run during use.

Note Use purified, distilled, or deionized water, or ADVIA Chemistry Drug Calibrator I - Level 1 as the sample for the RBL in the ADVIA Chemistry PHNY_2 assay.

For more information regarding running daily reagent blanks for multiple standard methods on the ADVIA Chemistry systems, refer to the Customer Bulletin entitled: *Performing Reagent Blank (RBL) on Multi-Standard (MSTD) Assays* (PN 10818532, latest revision).

Performing Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze 2 levels (low and high) of a commercially available quality control (QC) material with known phenytoin concentrations.

A satisfactory level of performance is achieved when the analyte values obtained are within the expected control range for the system or within your range, as determined by an appropriate internal laboratory quality control scheme.

The actual frequency of control in a laboratory is based on many factors, such as workflow, system experience, and government regulation. Each laboratory should evaluate the controls based on the frequency established by their laboratory guidelines.

Also, assay controls under the following conditions:

- Whenever you use a new reagent lot
- Following any system maintenance, cleaning, or troubleshooting procedure
- After performing a new calibration or an additional reagent blank

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Taking Corrective Action

If the quality control results do not fall within the expected control range or within the laboratory's established values, do not report results. Take the following actions:

1. Determine and correct the cause of the unacceptable control results:
 - a. Verify that the assay was performed according to the instructions for use.
 - b. Verify that the materials are not expired.
 - c. Verify that required maintenance was performed.
 - d. Rerun the assay with fresh quality control samples, and confirm that quality control results are within acceptable limits before running patient samples.
 - e. If the quality control results are not within acceptable limits, recalibrate the assay, and repeat the prior step.
 - f. If necessary, contact your local technical support provider or distributor for assistance.
2. After corrective action is complete, repeat required testing of patient samples before reporting results.

Perform corrective actions in accordance with your established laboratory protocol.

Results

Calculation of Results

The system calculates and reports results based on the absorbance measurements of the test sample during the test, and of the calibrator(s) from calibration.

The instrument calculates the concentration of phenytoin in $\mu\text{g/mL}$ (common units) or $\mu\text{mol/L}$ (SI units).

Conversion factor: $\mu\text{g/mL} \times 3.96 = \mu\text{mol/L}$

Interpretation of Results

Results of this assay should always be interpreted in conjunction with the patient's medical history, clinical presentation, and other findings.

Limitations

A number of substances cause physiological changes in serum or plasma analyte concentrations. A comprehensive discussion of possible interfering substances, their serum or plasma concentrations, and their possible physiological involvements is beyond the scope of this document. Consult the listed reference for specific details on known potential interfering substances.¹⁴

As with any chemical reaction, you must be alert to the possible effect on results of unknown interferences from medications or endogenous substances. The laboratory and physician must evaluate all patient results in light of the total clinical status of the patient.

Therapeutic Range

The therapeutic range for phenytoin is listed in the following table.

Therapeutic Range	10–20 $\mu\text{g/mL}$ (39.6–79.2 $\mu\text{mol/L}$) ^{1,2,5,15}
Toxic Range	> 20 $\mu\text{g/mL}$ (> 79.2 $\mu\text{mol/L}$) ^{1,5,15}

The factors that can influence the relationship between phenytoin serum or plasma concentrations and clinical response include the type and severity of seizures, age, general state of health, and use of other drugs.

The concentration of phenytoin in serum or plasma depends on the following parameters when considering how to interpret results:^{1,2}

- The time of the last drug dose
- The mode of administration
- Concomitant drug therapy
- Sample condition
- Time of sample collection
- Individual variations in absorption, distribution, biotransformation, and excretion

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation, and other findings.

Siemens provides this information for reference. As with all *in vitro* diagnostic assays, each laboratory should determine its own reference ranges for the diagnostic evaluation of patient results. Consider these values as a guideline only.

Performance Characteristics

Analytical Measuring Range

The analytical range is from 1.0 µg/mL (4.0 µmol/L) to the ADVIA Chemistry Drug Calibrator I - Level 5, which varies from 34.6–40.0 µg/mL (137.0–158.4 µmol/L).

Results that are below the low end of the assay range are flagged < **Conc Range**. You should report the test result as < 1.0 µg/mL (< 4.0 µmol/L).

Results that are above the high end of the assay range are flagged > **Conc Range**.

Extended Measuring Range

Siemens has validated an automatic rerun condition for this assay that extends the reportable range to 103.8 µg/mL (411.0 µmol/L). You may configure the system to trigger automatic reruns. Rerun results will be flagged **Autorepeat**.

Sensitivity

The ADVIA Chemistry PHNY_2 assay performance at low levels was analyzed as described in CLSI protocol EP17-A2, and the limit of blank (LoB) and limit of detection (LoD) were determined.¹⁶

The LoB is the highest measurement result that is likely to be observed on a blank sample. The LoB for the ADVIA Chemistry PHNY_2 assay is 0.5 µg/mL (2.0 µmol/L).

The LoD is the smallest amount that this assay can reliably detect to determine the presence or absence of an analyte. The LoD for the ADVIA Chemistry PHNY_2 assay is 1.0 µg/mL (4.0 µmol/L).

The LoB and LoD values are determined with proportions of false positives (α) less than 5% and false negatives (β) less than 5%, based on 180 determinations with 120 blank and 60 low-level sample replicates.

Precision

The precision of the ADVIA Chemistry PHNY_2 assay was analyzed as described in the CLSI document EP05-A2.¹⁷ Each sample was assayed 3 times per run, 2 runs per day, for at least 20 days.

Precision: Common Units

Specimen Type	N	Mean (µg/mL)	Repeatability (Within-Run)		Between-Run		Between-Day		Within-Lab (Total)	
			SD ^a (µg/mL)	CV ^b (%)	SD (µg/mL)	CV (%)	SD (µg/mL)	CV (%)	SD (µg/mL)	CV (%)
Serum Control 1	120	3.7	0.22	6.0	0.00	0.0	0.17	4.5	0.28	7.5
Serum Control 2	120	8.9	0.49	5.5	0.34	3.8	0.00	0.0	0.59	6.7
Serum Control 3	120	17.8	0.66	3.7	0.29	1.6	0.49	2.8	0.88	4.9
Serum Pool 1	120	10.1	0.59	5.9	0.17	1.7	0.34	3.4	0.71	7.0
Serum Pool 2	120	19.5	1.11	5.7	0.66	3.4	0.12	0.6	1.30	6.6

^a SD = standard deviation

^b CV = coefficient of variation

Precision: SI Units

Specimen Type	N	Mean (µmol/L)	Repeatability (Within-Run)		Between-Run		Between-Day		Within-Lab (Total)	
			SD ^a (µmol/L)	CV ^b (%)	SD (µmol/L)	CV (%)	SD (µmol/L)	CV (%)	SD (µmol/L)	CV (%)
Serum Control 1	120	14.7	0.87	6.0	0.00	0.0	0.67	4.5	1.11	7.5
Serum Control 2	120	35.2	1.94	5.5	1.35	3.8	0.00	0.0	2.34	6.7
Serum Control 3	120	70.5	2.61	3.7	1.15	1.6	1.94	2.8	3.48	4.9
Serum Pool 1	120	40.0	2.34	5.9	0.67	1.7	1.35	3.4	2.81	7.0
Serum Pool 2	120	77.2	4.40	5.7	2.61	3.4	0.48	0.6	5.15	6.6

^a SD = standard deviation

^b CV = coefficient of variation

Actual results will vary depending on the study design, and on the sample and sample population used. Results obtained at individual laboratories may vary from the data provided.

Accuracy / Method Comparison

The performance of the ADVIA Chemistry PHNY_2 assay (y) was compared with the performance of the comparison assay on the indicated system (x).

Specimen Type	Comparison Assay (x)	N	r	Regression Equation	Sy.x	Sample Range
Serum	ADVIA 2400 PHNY_2	101	0.994	y = 0.97x - 0.33 µg/mL y = 0.97x - 1.31 µmol/L	0.86 µg/mL 3.41 µmol/L	1.3–31.1 µg/mL 5.1–123.2 µmol/L
Serum	ADVIA 1800 PHNY_2	100	0.991	y = 0.98x - 0.02 µg/mL y = 0.98x - 0.08 µmol/L	1.00 µg/mL 3.96 µmol/L	1.2–32.2 µg/mL 4.8–127.5 µmol/L
Plasma (Lithium Heparin)	ADVIA 1800 PHNY_2	57	0.987	y = 1.00x + 0.08 µg/mL y = 1.00x + 0.32 µmol/L	1.59 µg/mL 6.30 µmol/L	2.3–37.6 µg/mL 9.1–148.9 µmol/L

Specimen Type	Comparison Assay (x)	N	r	Regression Equation	Sy.x	Sample Range
Plasma ^a (Lithium Heparin) ADVIA 1650/1800	ADVIA 1650/1800 PHNY_2 - Serum	58	0.99	y = 0.98x + 0.31 µg/mL y = 0.98x + 1.23 µmol/L	0.91 µg/mL 3.60 µmol/L	1.6–32.0 µg/mL 6.3–126.7 µmol/L
Plasma (Sodium Heparin)	ADVIA 1800 PHNY_2	53	0.990	y = 1.01x - 0.71 µg/mL y = 1.01x - 2.81 µmol/L	1.45 µg/mL 5.74 µmol/L	2.9–38.6 µg/mL 11.5–152.9 µmol/L
Plasma ^a (Sodium Heparin) ADVIA 1650/1800	ADVIA 1650/1800 PHNY_2 - Serum	58	0.99	y = 0.98x + 0.21 µg/mL y = 0.98x + 0.83 µmol/L	0.93 µg/mL 3.68 µmol/L	1.6–32.0 µg/mL 6.3–126.7 µmol/L
Plasma (Potassium EDTA)	ADVIA 1800 PHNY_2	58	0.987	y = 0.96x - 0.20 µg/mL y = 0.96x - 0.79 µmol/L	1.61 µg/mL 6.38 µmol/L	3.6–39.6 µg/mL 14.3–156.8 µmol/L
Plasma ^a (Potassium EDTA) ADVIA 1650/1800	ADVIA 1650/1800 PHNY_2 - Serum	58	0.99	y = 0.99x + 0.32 µg/mL y = 0.99x + 1.27 µmol/L	0.78 µg/mL 3.09 µmol/L	1.6–32.0 µg/mL 6.3–126.7 µmol/L
Plasma (Citrate)	ADVIA 1800 PHNY_2	49	0.989	y = 0.97x - 0.23 µg/mL y = 0.97x - 0.91 µmol/L	1.41 µg/mL 5.58 µmol/L	2.8–38.4 µg/mL 11.1–152.1 µmol/L
Plasma ^a (Citrate) ADVIA 1650/1800	ADVIA 1650/1800 PHNY_2 - Serum	58	0.99	y = 1.00x + 0.11 µg/mL y = 1.00x + 0.44 µmol/L	0.85 µg/mL 3.37 µmol/L	1.6–32.0 µg/mL 6.3–126.7 µmol/L
Plasma (NaFl-KOx) ^b	ADVIA 1800 PHNY_2	50	0.985	y = 1.00x - 0.32 µg/mL y = 1.00x - 1.27 µmol/L	1.47 µg/mL 5.82 µmol/L	3.7–38.7 µg/mL 14.7–153.3 µmol/L
Plasma ^a (NaFl-KOx) ^b ADVIA 1650/1800	ADVIA 1650/1800 PHNY_2 - Serum	58	0.99	y = 0.99x + 0.36 µg/mL y = 0.99x + 1.43 µmol/L	0.87 µg/mL 3.45 µmol/L	1.6–32.0 µg/mL 6.3–126.7 µmol/L

^a Matrix comparison. Correlations between serum and plasma samples on ADVIA 1650/1800 Chemistry systems are provided for reference.

^b sodium fluoride/potassium oxalate

The correlation of the assay may vary depending on the study design, comparable method, and sample population. Results obtained at individual laboratories may vary from the data provided.

Interferences

Interferent	Interferent Level	Phenytoin Sample Concentration	Interference
Bilirubin (conjugated)	60 mg/dL (1026 µmol/L)	9.1 µg/mL (36.0 µmol/L)	NSI ^a
	60 mg/dL (1026 µmol/L)	18.4 µg/mL (72.9 µmol/L)	NSI
Bilirubin (unconjugated)	60 mg/dL (1026 µmol/L)	8.9 µg/mL (35.2 µmol/L)	NSI
	60 mg/dL (1026 µmol/L)	18.6 µg/mL (73.7 µmol/L)	NSI

Interferent	Interferent Level	Phenytoin Sample Concentration	Interference
Hemolysis (hemoglobin)	1000 mg/dL (10.0 g/L)	8.9 µg/mL (35.2 µmol/L)	NSI
	1000 mg/dL (10.0 g/L)	17.7 µg/mL (70.1 µmol/L)	NSI
Lipemia ^b (triglycerides from Intralipid)	1000 mg/dL (11.3 mmol/L)	9.2 µg/mL (36.4 µmol/L)	NSI
	1000 mg/dL (11.3 mmol/L)	17.6 µg/mL (69.7 µmol/L)	NSI

^a NSI = No significant interference. A percentage effect $\geq 10\%$ is considered a significant interference.

^b SI units calculated as triolein

Note There is poor correlation between turbidity and triglyceride concentration in a lipemic sample.¹⁸

Actual results will vary depending on the study design, the levels of the potential interferences tested, and the samples used. Results obtained at individual laboratories may vary from the data provided.

Specificity

The ADVIA Chemistry PHNY_2 assay measures the total (protein-bound plus unbound) phenytoin concentration in serum or plasma. The specificity of the method was evaluated by assaying compounds whose chemical structure or concurrent usage could potentially interfere with the ADVIA Chemistry PHNY_2 assay. This data was collected on an automated chemistry system using method parameters equivalent to those used on the ADVIA Chemistry XPT system.⁶

The following compounds did not interfere with the ADVIA Chemistry PHNY_2 assay when tested in the presence of 10 µg/mL phenytoin. Levels tested were at or above maximum pharmacological concentrations.

Compound Tested	Concentration Tested (µg/mL)
Amitriptyline	25
Amobarbital	75
Carbamazepine	500
Carbamazepine-10,11-epoxide	500
Chlordiazepoxide	60
Chlorpromazine	8
Clorazepate	500
Diazepam	60
Ethosuximide	500
Ethotoin	200
5-Ethyl-5-phenylhydantoin	200
Glutethimide	200
5-(p-Hydroxyphenyl)-5-phenylhydantoin	50

Compound Tested	Concentration Tested (µg/mL)
5-(p-Hydroxyphenyl)-5-phenylhydantoin glucuronide	1000
Imipramine	5
Mephenytoin	35
Mephobarbital	500
Methsuximide	150
Pentobarbital	100
Phenobarbital	500
Phensuximide	500
2-Phenyl-2-ethyl-malondiamide (PEMA)	500
Primidone	200
Promethazine	10
Secobarbital	25
Valproic Acid	1000

The average crossreactivity of fosphenytoin (concentration tested 15–100 µg/mL), as measured with PHNY_2 method performed on the ADVIA 1650 Chemistry system, was found to be 3.5%.¹⁹

Standardization

The ADVIA Chemistry PHNY_2 assay is traceable to an internal standard manufactured using United States Pharmacopoeia (USP) material. Assigned values of ADVIA Chemistry Drug Calibrator I are traceable to this standardization.

Technical Assistance

For customer support, please contact your local technical support provider or distributor.
siemens.com/healthcare




















References

1. Finn AL, Olanow CW: Phenytoin: Therapeutic use and serum concentration monitoring. In: Taylor WJ, Finn AL eds: *Individualizing Drug Therapy: Practical Applications of Drug Monitoring. Vol 2.* New York: Gross, Townsend, Frank, Inc; 1981:63–85.
2. Winter ME, Tozer TN. Phenytoin. In: Evans WE, Schentag JJ, Jusko WJ eds. *Principles of Therapeutic Drug Monitoring. 2nd ed.* Spokane, Washington: Applied Therapeutics, Inc; 1987:493–539.
3. Leppik IE, Oles KS, Sheehan ML, et al. Phenytoin and phenobarbital concentrations in serum: a comparison of Ames Seralyzer with GLC, TDX, and EMIT. *Ther Drug Monit.* 1989;11:73–78.
4. Rainbow SJ, Dawson CM, Tickner TR. Direct serum injection high-performance liquid chromatographic method for the simultaneous determination of phenobarbital, carbamazepine, and phenytoin. *J Chromatogr.* 1990;527:389–396.
5. Pincus MR, Abraham NZ Jr. Toxicology and therapeutic drug monitoring. In: McPherson R, Pincus MR eds: *Henry's Clinical Diagnosis and Management by Laboratory Methods. 21st ed.* Philadelphia, PA: Saunders Elsevier; 2007:297–325.

6. Siemens Syva Emit 2000 Phenytoin Assay [package insert]. Glasgow, DE: Siemens Healthcare Diagnostics; 2006.
7. Clinical and Laboratory Standards Institute (formerly NCCLS). *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Guideline—Sixth Edition*. CLSI document GP41-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
8. Clinical and Laboratory Standards Institute (formerly NCCLS). *Tubes and Additives for Venous Blood Specimen Collection: Approved Standard; Approved Guideline—Sixth Edition*. CLSI document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
9. Clinical and Laboratory Standards Institute (formerly NCCLS). *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Fourth Edition*. CLSI document GP44-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
10. Centers for Disease Control and Prevention. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*. 1988; 37:377–82, 387–8.
11. Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
12. Federal Occupational Safety and Health Administration. *Bloodborne Pathogens Standard*. 29 CFR 1910.1030.
13. Kumps AH. Therapeutic drug monitoring: A comprehensive and critical review of analytical methods for anticonvulsive drugs. *J Neurol*. 1982;228:1–16.
14. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Washington, DC: AACC Press; 2000.
15. Levine M, McCollom R, Chang T, et al. *The Drug Monit*. 1988;10:50–57.
16. Clinical and Laboratory Standards Institute (formerly NCCLS). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
17. Clinical and Laboratory Standards Institute (formerly NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI document EP05-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
18. Bornhorst JA, Roberts RF, Roberts WL. Assay-specific differences in lipemic interference in native and Intralipid-supplemented samples. *Clin Chem*. 2004; 50:2197-201.
19. Data on file at Siemens Healthcare Diagnostics.

Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Definition	Symbol	Definition
	<i>In vitro</i> diagnostic medical device	 REF	Catalog number
	Legal manufacturer		Authorized Representative in the European Community
	CE Mark		CE Mark with identification number of notified body
	Consult instructions for use		Biological risk
	Keep away from sunlight and heat		Temperature limitation
	Lower limit of temperature		Upper limit of temperature
	Do not freeze (> 0°C)		Up
	Use by		Contains sufficient for (n) tests
	Recycle		Printed with soy ink
Rev.	Revision	YYYY-MM-DD	Date format (year-month-day)
	Batch code	RxOnly	Prescription Device (US only)

Trademarks

ADVIA, Syva, and Emit are trademarks of Siemens Healthcare Diagnostics.

Intralipid is a trademark of Fresenius Kabi AB.

© 2014–2020 Siemens Healthcare Diagnostics. All rights reserved.

 Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591 USA

Global Siemens Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

