3gAllergy™ Specific IgE Universal Kit

For use on IMMULITE® 2000 systems

(IMMULITE 2000 and IMMULITE 2000 XPi automated immunoassay analyzers)

SIEMENS

IMMULITE® 2000 3gAllergy™ Specific IgE Universal Kit

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of allergen-specific IgE in human serum, as an aid in the clinical diagnosis of IgE-mediated allergic disorders.

Catalog Number: **L2KUN6** (600 tests)

Test Code: **SPE** Color: **Light Gray**

Summary and Explanation

Many allergies are mediated by immunoglobulins of the IgE class. In sensitized individuals suffering from this immediate (atopic or anaphylactic) type of allergy, IgE molecules act as points of contact between the allergen and specialized cells that release histamine and other agents upon exposure to the allergen; this initiates the events which we recognize as allergic reactions.^{5,9} When evaluated in the light of other clinical and laboratory findings, *in vitro* allergen-specific IgE tests can help the physician identify the allergen (or allergens) to which an individual is sensitive.

Principle of the Procedure

IMMULITE 2000 3gAllergy™ Specific IgE is a solid-phase, two-step, chemiluminescent immunoassay that exploits liquid phase kinetics in a bead format.^{10,11} It represents a significant advance over conventional methods relying on allergens attached to a solid-phase support, such as a paper disk.

The solid phase (bead) is coated with anti-ligand. The liquid phase consists of alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgE antibody in a human/nonhuman serum buffer matrix. Individual ligand-labeled allergens and panels in liquid phase are not supplied in the kit but are required for the run.

In the first cycle, the patient sample and the ligand-labeled specific allergen are incubated together with the coated bead for 30 minutes. During this time, specific IgE in the sample binds to the ligand-labeled allergen, which in turn binds to the anti-ligand on the bead. Unbound sample is then removed by centrifugal washes.

In the second cycle, the enzyme conjugated monoclonal murine anti-human IgE antibody is added to the original reaction tube for additional 30 minutes incubation. The enzyme conjugated monoclonal murine anti-human IgE antibody binds to immobilized IgE. The unbound enzyme conjugate is removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 2 × 30 minutes
Time to first result: 65 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators. IMMULITE 2000 3gAllergy™ Specific IgE has not been tested with all possible variations of tube types.

Volume Required: 50 µL serum

Storage: 7 days at 2–8°C*, or 6 months at –20°C.¹³

* Data on file at Siemens Healthcare Diagnostics.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.³⁷⁻³⁹

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: 2-methyl-2H-isothiazol-3-one, nonylphenol, ethoxylated; 3gAllergy Specific IgE Reagent Wedge

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Caution: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Safety data sheets (MSDS) available on www.siemens.com/diagnostics.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source

materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

3gAllergy™ Specific IgE Bead Pack (L2UN12)

With barcode. 200 beads, coated with anti-ligand. Stable at 2–8°C until expiration date.

L2KUN6: 3 packs

3gAllergy™ Specific IgE Reagent Wedge (L2UNA6)

With barcode. 30 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgE antibody in a human/nonhuman serum buffer matrix, dispensed equally into chambers B and C. Stable at 2–8°C until expiration date.

L2KUN6: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

3gAllergy™ Specific IgE Adjustors (L2UNJ3, L2UNJ4)

Two vials (Low and High), 2.0 mL each, of human IgE in a nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KUN6: 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

3gAllergy™ Specific IgE Adjustor Antibody (L2UN51)

Two tubes, 2.75 mL each of liquid, ready-to-use ligand-labeled polyclonal goat anti-human IgE antibody, with preservative. Stable at 2–8°C until expiration date. This reagent is placed in the Allergen Holder Wedge when running the Specific IgE Adjustors.

L2KUN6: 2 sets

3gAllergy™ Specific IgE (SPE) Universal Kit Controls (L2UNC1, L2UNC2)

Two vials, 2 mL each of human IgE in a nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KUN6: 2 sets

Refer to the control insert for concentration levels.

Before use, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

3gAllergy™ Specific IgE Control Antibody (L2UNS2)

Two tubes, 2.75 mL each of liquid, ready-to-use ligand-labeled polyclonal goat anti-human IgE antibody, with preservative. Stable at 2–8°C until expiration date. This reagent is placed in the Allergen Holder Wedge when running the IMMULITE 2000 IgE Controls.

L2KUN6: 2 sets

Kit Components Supplied Separately

3gAllergy™ Specific IgE Sample Diluent (L2UNZ)

For on-board dilution of samples. One vial, concentrated (ready-to-use), human serum albumin matrix, with preservative.

Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C. Dispose of in accordance with applicable laws.

L2UNZ: 25 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2UNZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2AW1-3: Allergen Holder Wedges (barcoded)

L2AW1, 400930-02: serially coded 1–33

L2AW2, 400930-03: serially coded 34–66

L2AW3, 400930-04: serially coded 67–99

L2ATC: Allergen Tube Caps

L2ATS2: Allergen Tube Septa

Also Available

MC6LCM, DC1LCM, DC2LCM,

L2SNCCM: Human serum-based allergen-specific IgE controls

Also Required

Distilled or deionized water; test tubes; controls.

3gAllergy™ Specific Allergens and Mixed Allergen Panels

Individual allergens are packaged and sold in 40 and 20 test modules containing 2.75 mL and 1.75 mL, respectively.

Mixed allergen panels are packaged and sold in 40 test modules only containing 2.75 mL each.

Each allergen tube contains specific allergens or a mix of allergens in a protein-based buffer matrix, with preservative.

Do *not* freeze. Store refrigerated: stable at 2–8°C until the expiration date marked on the label or 90 days onboard instrument.

Do *not* use if signs of microbial contamination such as a cloudy appearance have been observed.

Individual allergens and mixed allergen panels are intended for use with the IMMULITE 2000 systems analyzers. For a full list and catalog numbers, please refer to the 3gAllergy™ Menu.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Allergen Loading

- 1 Select an open position on the Reagent Carousel through the software.
- 2 Replace caps on Allergen Tubes with septa. Do not invert the Allergen Tube once the septum is installed.
- 3 Place Allergen Tubes containing IMMULITE 2000 specific allergens, specific allergen panels, Specific IgE Adjustor Antibody, and/or Specific IgE Control Antibody in the Allergen Holder Wedge, with the barcodes facing the open side of the wedge.
- 4 Close the wedge and scan the allergen barcodes with the hand-held imaging scanner.
- 5 Once scanning is complete, load the Allergen Holder Wedge into the reagent carousel.
- 6 Repeat this procedure to load subsequent Allergen Holder Wedges.

The Allergen Holder Wedge **must** be scanned prior to installation into the instrument to ensure correct instrument operation. Removing or replacing any vials within the allergen wedge will require rescanning of the wedge with the barcode scanner to update the allergen information.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: To monitor adjustment, use the controls supplied with the kit. Specific allergen controls are also available, see "Kit Components Supplied Separately" section. The specific allergen controls (DC1LCM, DC2LCM, MC6LCM and L2SNCCM) should be tested to monitor allergen performance.

Interpretation of Results

Individual Allergen Results:

The class number is an indication of the amount of endogenous IgE specific for the selected allergen. Quantitative values (kU/L) and interpretation of class results for two scoring systems (standard and extended) are provided in the tables below.¹⁶

The **standard** classification system utilizes the following class cutoffs:

Class	kU/L	Reactivity for Individual/ Panel Allergen(s)
0*	< 0.10	Absent or ND†
	0.10–0.34	Very Low
I	0.35–0.69	Low
II	0.70–3.49	Moderate
III	3.50–17.49	High
IV	17.5–52.49	Very High
V	52.5–99.99	
VI	≥ 100	

* Class 0 in the standard system signifies: not detectable by second-generation assays.

† ND: not detectable by IMMULITE 2000 3gAllergy.

The **extended** classification system utilizes the following class cutoffs.

Class	kU/L	Reactivity for Individual/ Panel Allergen(s)
0	< 0.10	Absent or ND†
0/1	0.10–0.24	Very Low
I	0.25–0.39	Low
II	0.40–1.29	Moderate
III	1.30–3.89	High
IV	3.90–14.99	Very High
V	15.00–24.99	
VI	≥ 25	

† ND: not detectable by IMMULITE 2000 3gAllergy.

The choice of classification systems can be made by the user within the IMMULITE 2000 operational software.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Mixed Allergen Panel Results

A *positive* result (see Interpretation of Results for standard and extended classification systems) with a Mixed Allergen Panel indicates that antibodies to one or more of the component allergens in the panel are present in elevated amounts in the patient serum sample. To identify the allergen-specific IgE, the sample should be retested with individual allergens corresponding to components of the panel.

A *negative* result (< 0.10 kU/L) indicates the absence or nondetectable levels of IgE specific for the panel's component allergens.

Panel value results cannot be compared with quantitative results based on testing for individual allergens, nor can they be considered the cumulative total of individual allergen results.

Limitations

A definitive clinical diagnosis should not be made solely on the basis of an *in vitro* allergen-specific IgE result. Diagnosis should be made by the physician only after all clinical and laboratory findings have been evaluated.^{2,4}

In vitro allergen-specific IgE results should not be used as a definitive guide to selecting an initial dose for immunotherapy. A skin test with the proposed initial dilution of the allergenic extract should be performed first to demonstrate the patient's ability to tolerate the dose.

In food allergy, circulating IgE antibodies may remain undetectable if directed towards allergens which are revealed or altered during processing or digestion and which therefore do not exist in the original food for which the patient is tested.^{1,3}

Identical results for different allergens may not be associated with clinically equivalent manifestations, due to differences in IgE-binding capacity.⁶

The user should be aware of the possibility of clinical crossreactivity within an allergen family.^{7,8}

In drug or occupational allergy, a negative result may be observed in patients who are hypersensitive to drug or occupational allergens in the following circumstances:

- The symptoms are mediated without IgE involvement.
- The sample is collected less than 2 weeks after the allergic reaction. The test should be repeated after 2 weeks to confirm results.
- The sample is collected a long period of time after the last allergic reaction occurred. It has been shown that the IgE antibody concentration decreases with time.

The following special considerations apply to latex allergy testing:

- The possibility of clinical crossreactivity exists between latex and certain foods including avocado, banana, chestnut, and kiwi.¹²
- Since the latex assay measures allergen-specific IgE, type IV delayed reaction or irritation from latex will not be detected.

Class 0 results for insect venoms indicate absent or very low levels of circulating venom-specific IgE antibodies. Such results do not preclude existence of current or future clinical hypersensitivity to insect sting.

Do not use this assay in patients taking biotin supplements.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in kU/L. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

IMMULITE 2000 3gAllergy™ Specific IgE has been the subject of a number of published studies.^{14,15}

Reportable Range: 0.10–100 kU/L (2nd WHO IRP 75/502 human serum IgE)

Analytical Sensitivity: Limit of Blank (highest value expected for a sample with

no analyte; determined in accordance with CLSI EP17-A³⁵): 0.03 kU/L

Limit of Detection (lowest detectable concentration; determined in accordance with CLSI EP17-A³⁵): 0.10 kU/L

Functional Sensitivity: (concentration with 20% coefficient of variation (CV) determined in accordance with CLSI EP5-A2³⁶): 0.20 kU/L

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Specificity: The antibodies are highly specific for human IgE and exhibit no crossreactivity to other human Immunoglobulin classes.

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin:

Biotin Test Level (ng/mL)	Wheat* Concentration (kU/L)	
	0.74	3.9
Bias (%)		
5.0	-4	-6
9.0	-9	-13
19.0	-18	-22
32.0	-37	-39
100	-60	-66
1500	-89	-92

*Wheat allergen data is considered representative of the assay's performance.

Specimens that contain biotin at a concentration of 5 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely depressed results for patient samples.

The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy

adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.⁴⁰ Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.⁴¹

Hemolysis: Presence of hemoglobin in concentrations up to 500 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to AlaSTAT Microplate Allergen-Specific IgE, a standard scoring assay, on 7520 samples.

IMMULITE 2000	6			1	3	31	42	303
	5				11	76	33	35
	4			3	178	168	26	10
	3		1	292	604	49	2	2
	2	16	229	886	98			
	1	108	172	94	4			
	0	3912	118	11	2			
		0	1	2	3	4	5	6
AlaSTAT Microplate								

Total Agreement: 97%
Relative Sensitivity: 96%
Relative Specificity: 97%
Total Identical: 81%
Class 1 Identical: 99%
Class 2 Identical: 100%

Means of Class Scores:
1.29 (AlaSTAT)
1.37 (IML 2000)

The same data was analyzed in a regression analysis with the following results:

(IML 2000) = 0.989 (AlaSTAT) + 1.21 kU/L
 $r = 0.87$

Means:
8.53 kU/L (AlaSTAT)
9.64 kU/L (IMMULITE 2000)

Clinical Performance: Samples were assayed in singlicate across 24 specific allergens. Clinical diagnosis was based on patient donor history and/or skin test data.¹⁷⁻³⁴

Results of Clinical Study*, n = 4214

IMMULITE 2000	Clinical Diagnosis	
	Clinical	Normal
Positive	595	101
Negative	537	2981

Agreement: 84.9%
95% Confidence Interval: 84–86%

Specificity: 96.7%
95% Confidence Interval: 96–97%

Sensitivity: 52.6%
95% Confidence Interval: 50–55%

Allergens Included: Common Sagebrush (W43), Wingscale (W75), Red Maple (T27), White Hickory (T41), Sweet Gum (T211), Red Cedar (T219), American Cockroach (I206), Bayberry/Sweet Gale (T218), Baccharis (W67), Dog Fennel (W46), Locust Tree (T208), Live Oak (T103), *Hormodendrum hordei* (M45), *Stemphylium solani* (M88), Red Mulberry (T71), Privet (T210), White Bald Cypress (T37), Egg (F245), Cashew (F202), Clam (F207), Oyster (F290), Pistachio (F203), Scallop (F338), and Walnut (F256).

*Data on file at Siemens Healthcare Diagnostics.

References

- 1) Aas K. The diagnosis of hypersensitivity to ingested foods. *Clin Allergy* 1978;8:39–50.
- 2) Barbee RA, et al. Longitudinal changes in allergen skin test reactivity in a community population sample. *J Allergy Clin Immunol* 1987;79:16–24.
- 3) Bleumink E. Food allergy: the chemical nature of the substance eliciting symptoms. *World Rev Nutr Diet* 1970;12: 505–70.
- 4) Bloch K, Salvaggio J. Use and interpretation of diagnostic immunologic laboratory tests. *JAMA* 1982; 246:2734–58.
- 5) Halpern GM. Markers of human allergic disease. *J Clin Immunoassay* 1983;6:131–8.
- 6) Lichtenstein LM, et al. IgE antibody measurements in ragweed hay fever; relationship to clinical severity and the results of immunotherapy. *J Clin Invest* 1973; 52:472–82.
- 7) Lowenstein H. Cross reactions among pollen antigens. *Allergy* 1980;35:198–200.
- 8) Weber RW, Nelson HS. Pollen allergens and their interrelationships. *Clin Rev Allergy* 1985;3: 291–318.
- 9) Wide L, Bennich H, Johansson SGO. Diagnosis of allergy by an *in vitro* test for allergen antibodies. *Lancet* 1967;2:1105–7.
- 10) El Shami AS, Alaba O. Liquid-phase *in vitro* allergen-specific IgE assay with *in situ* immobilization. *Adv Biosci* 1989;74:191–201.

- 11) Alaba O, El Shami AS. Evaluation of non-specific IgE binding: comparison of two *in vitro* allergen assays. *Adv Biosci* 1989;74:203–14.
- 12) Pecquet C. IgE-mediated allergy to latex in 80 patients. Presented at the XVth European Congress of Allergology and Clinical Immunology, Paris, 12 May 1992.
- 13) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:358.
- 14) Li TM, Chuang T, Tse S, Hovanec-Burns D, El Shami AS. Development and validation of a third generation allergen-specific IgE assay on the continuous random access IMMULITE 2000 analyzer. *Ann Clin Lab Sci* 2004;34(1):67–74. Available via www.AnnClinLabSci.org.
- 15) Guilloux L, Hamberger C. Évaluation du dosage des IgE spécifiques sur l'IMMULITE® 2000 DPC. *Immuno-analyse & Biologie Spécialisée*. 2004;19(1):71–80. Available via www.sciencedirect.com.
- 16) Hoffman DR. Comparison of methods of performing the radioallergosorbent test: Phadebas, Fada-Nalebuff, and Hoffman protocols. *Ann Allergy*. 1980 Dec;45(6):343–6.
- 17) Bernstein IL, et al. Allergy Diagnostic Testing: An Updated Practice Parameter. *Ann Allergy Asthma Immunol* 2008;100(3) Sup 3:S1–S148.
- 18) Weber RW. Cross-reactivity of pollen allergens: impact on allergen immunotherapy. *Ann Allergy Asthma Immunol* 2007;99:203–212.
- 19) Weber RW. Cross-reactivity of Pollen Allergens. *Curr Allergy Asthma Rep* 2004;4:401–408.
- 20) White JF, Bernstein DL. Key pollen allergens in North America. *Ann Allergy Asthma Immunol* 2003;91:425–35.
- 21) Compes E, et al. Hypersensitivity to black locust (*Robinia pseudoacacia*) pollen: “allergy mirages”. *Ann Allergy Asthma Immunol* 2006; 96(4):586–92.
- 22) Willison L, et al. Pistachio vicilin, Pis v 3, is immunoglobulin E-reactive and cross-reactive with the homologous cashew allergen, Ana o 1. *Clin Exp Allergy* 2008;38:1229–38.
- 23) Clark A, Ewan P. Good prognosis, clinical features, and circumstances of peanut and tree nut reactions in children treated by a specialist allergy center. *J Allergy Clin Immunol* 2008;122(2):286–9.
- 24) Barre et al. Vicilin allergens of peanut and tree nuts (walnut, hazelnut, and cashew nut) share structurally related IgE-binding epitopes. *Mol Immunol* 2008;45:1231–40.
- 25) Maloney J, Rudengren M, Ahlstedt S, Bock S, Sampson H. The use of serum-specific IgE measurements for diagnosis of peanut, tree nut, and see allergy. *J Allergy Clin Immunol* 2008;122(1):145–151.
- 26) Goetz D, Whisman B, Goetz. Cross-reactivity among edible nuts: double immunodiffusion, crossed immunoelectrophoresis, and human specific IgE serologic surveys. *Ann Allergy Asthma Immunology* 2005;95:45–52.
- 27) <http://www.allergen.org/Allergen.aspx> International Union of Immunological Societies – Allergen Nomenclature Sub-Committee.
- 28) Leung P, et al. IgE reactivity against a cross-reactive allergen in crustacean and mollusca: evidence for tropomyosin as the common allergen. *J Allergy Clin Immunol* 1996;98(5):

954-61. 29) Aalberse R, Akkerdaas J, van Ree R. Cross-reactivity of IgE antibodies to allergens. *Allergy* 2001;56:478-90. 30) Sidenius K, Hallas T, Poulsen L, Mosbech H. Allergen cross-reactivity between house-dust mites and other invertebrates. *Allergy* 2001;56:723-33. 31) Ferreira F, Hawranek P, Gruber N, Wopfner N. Allergic cross-reactivity: from gene to the clinic. *Allergy* 2004;59:243-67. 32) Zhang Y, Hiroaki M, Morita E. Case report: cross-reactivity among shrimp, crab and scallops in a patient with a seafood allergy. *J Dermatology* 2006;3:174-77. 33) Purohit A, et al. role of tropomyosin as a cross-reacting allergen in sensitization to cockroach in patients from Martinique (French Caribbean Island) with a respiratory allergy to mite and a food allergy to crab and shrimp. *Europ Ann Allergy and Clin Immunol* 2007;39:85-88. 34) Nakamura A, et al. Effect of maillard reaction on allergenicity of scallop tropomyosin. *J Agric Food Chem* 2005;53:7559-64. 35) CLSI. Protocols for the Determination of Limits of Detection and Quantitation; Approved Guideline. CLSI document EP17-A Vol. 24 (No. 34). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004. 36) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI Document EP5-A2. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004. 37) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8. 38) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 39) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030. 40) Grimsey P, et al. 2017. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *Int. J. Pharmacokinet.* 2(4), 247-256. 41) Piketty ML, et al. 2017. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 55(6):817-825.

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (kU/L)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	%CV ⁵	SD	%CV	
1	0.5	0.05	10	0.06		12
2	0.7	0.07	10	0.08		11
3	1.1	0.04	3.6	0.05		4.5
4	1.1	0.08	7.3	0.08		7.3
5	1.5	0.08	5.3	0.09		6.0
6	1.6	0.15	9.4	0.16		10
7	1.8	0.10	5.6	0.10		5.6
8	1.9	0.14	7.4	0.16		8.4
9	2.1	0.12	5.7	0.14		6.7
10	2.4	0.10	4.2	0.14		5.8
11	3.0	0.15	5.0	0.17		5.7
12	3.4	0.12	3.5	0.16		4.7
13	3.6	0.20	5.6	0.24		6.7
14	4.1	0.28	6.8	0.32		7.8
15	6.1	0.28	4.6	0.36		5.9
16	8.9	0.44	4.9	0.56		6.3
17	9.5	0.53	5.6	0.65		6.8
18	10	0.36	3.6	0.50		5.0
19	12	0.56	4.7	0.76		6.3
20	15	0.67	4.5	0.67		4.5
21	19	1.7	8.9	1.9		10
22	43	1.9	4.4	2.6		6.0

Technical Assistance

Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

Linearity (kU/L)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	32 in 32 ⁵	96.6	—	—
	16 in 32	53.8	48.3	111
	8 in 32	26.6	24.2	110
	4 in 32	12.1	12.1	100
	2 in 32	5.74	6.04	95
	1 in 32	3.04	3.02	101
2	32 in 32	53.7	—	—
	16 in 32	30.1	26.8	112
	8 in 32	13.8	13.4	103
	4 in 32	6.81	6.71	102
	2 in 32	3.40	3.35	101
	1 in 32	1.70	1.68	101
3	32 in 32	36.4	—	—
	16 in 32	17.3	18.2	94
	8 in 32	9.48	9.11	104
	4 in 32	4.78	4.55	106
	2 in 32	2.48	2.28	109
	1 in 32	1.21	1.14	106

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IMMULITE and 3gAllergy™ are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.
CONTROL AB