Dimension Vista® System

Dimension Vista® 500, Dimension Vista® 1000T, Dimension Vista® 1500, Dimension Vista® 3000T Systems

Flex® reagent cartridge

CEA

See shaded sections: Updated information from 2019-05 version.

Issue Date 2020-08

Carcinoembryonic Antigen

Warning: The concentration of CEA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CEA assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining CEA levels serially is changed, additional sequential testing should be carried out. Prior to changing assays, the laboratory MUST confirm baseline values for patients being serially monitored.

Intended Use: The CEA method is an *in vitro* diagnostic test for the quantitative measurement of carcinoembryonic antigen in human serum and sodium or lithium heparinized plasma on the Dimension Vista® System.

Measurements of carcinoembryonic antigen are used as an aid in the management of cancer patients in whom changing CEA concentrations have been observed.

Summary: Carcinoembryonic antigen (CEA) was first isolated from a colon carcinoma specimen in 1965.¹ The antigen is identical to the CEA produced by digestive tract tissue of the normal fetus. CEA is a heterogeneous cell surface glycoprotein (molecular mass 150 – 300 kDa, 45 – 55% carbohydrate) belonging to the immunoglobulin “gene superfamily”.² The heterogeneous nature of CEA accounts for the differences in diagnostic cutoff values reported for different immunoassays. It is therefore essential that serial measurements employ the same immunoassay.³

CEA is present at low concentrations in normal adult serum. Elevated levels of serum CEA are seen in smokers, in benign liver disease (hepatitis, cirrhosis), benign gastric and intestinal disease (polyposis, ulcerative colitis), benign breast disease, pulmonary infection, emphysema and renal failure. Abnormally elevated levels of serum CEA are seen in cancers of the colon, rectum, lung, breast, liver, pancreas, stomach and ovary.⁴ Abnormal serum CEA levels have been found to be a better indicator for monitoring colorectal cancer than for other cancers.⁵

Globally, cancer of the colon and rectum is the third highest occurring cancer in males, the fourth highest occurring cancer in females and has a high morbidity and mortality. It is most curable when detected at an early stage. Clinically CEA levels are not useful in screening since fewer than 25% of patients with colon cancer have elevated CEA levels on diagnosis.^{6,7} After surgical resection of the colon, CEA is very useful in monitoring the success of treatment and in surveillance for recurrence and metastasis. Serial monitoring for a period of 2 years is the recommended protocol since return to elevated serum CEA levels is often the first abnormality observed.⁸

Principles of Procedure: The CEA method is a homogeneous, sandwich chemiluminescent immunoassay based on LOCI® technology. The LOCI® reagents include two synthetic bead reagents and a biotinylated anti-CEA monoclonal antibody fragment. The first bead reagent (Chemibeads) is coated with an anti-CEA monoclonal antibody and contains a chemiluminescent dye. The second bead reagent (Sensibeads) is coated with streptavidin and contains a photosensitizer dye. Sample is incubated with biotinylated antibody and Chemibeads to form bead-CEA-biotinylated antibody sandwiches. Sensibeads are added and bind to the biotin to form bead-pair immunocomplexes. Illumination of the complex at 680 nm generates singlet oxygen from Sensibeads which diffuses into the Chemibeads, triggering a chemiluminescent reaction. The resulting signal is measured at 612 nm and is a direct function of the CEA concentration in the sample.^{9,10}

Reagents

Wells ^{a,b}	Form	Ingredient	Concentration ^c	Source
1 – 4	Liquid	CEA biotinylated antibody	5 µg/mL	Mouse monoclonal
5 – 8	Liquid	CEA Chemibeads	100 µg/mL	Mouse monoclonal
9 – 12	Liquid	Streptavidin		
		Sensibeads	100 µg/mL	Recombinant <i>E. coli</i>

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Wells 1 – 12 contain buffers, stabilizers and preservatives.

c. Nominal value per well in a cartridge.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Used LOCI® reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open well stability: 7 days for wells 1 – 12

Specimen Collection and Handling: Recommended specimen types: serum and plasma (lithium and sodium heparin).

Samples and controls stabilized with sodium azide cannot be used.

Collect serum and plasma using recommended procedures.¹¹ Follow the instructions provided with your specimen collection device for use and processing.

Samples should be stored at 4 °C and analyzed within one week.¹²

For longer storage, samples may be frozen at -20 °C or colder for 4 months. Specimens should be free of particulate matter. Lipemic or frozen samples, which become turbid after thawing, must be clarified by centrifugation (10 minutes at approximately 15000 x g) prior to testing.

The purpose of specimen storage information is to provide guidance to customers; however, customers may validate their own procedures for storing patient samples.

Procedure

Materials Provided

CEA Flex® reagent cartridge, Cat. No. K6453

Materials Required But Not Provided

LOCI 5 CAL, Cat. No. KC600 or LOCI 6 CAL, Cat. No. KC604
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

	Test Conditions
Sample Volume (delivered to the vessel)	2 µL
Biotinylated Antibody Volume	20 µL
Chemibead Reagent Volume	20 µL
Sensibead Reagent Volume	100 µL
Temperature	37.0 °C
Reaction time	10 minutes
Wavelength	680 nm Illumination 612 nm Emission
Type of Measurement	Chemiluminescence
Calibration	
Calibration Material	LOCI 5 CAL, Cat. No. KC600 or LOCI 6 CAL, Cat. No. KC604
Calibration Scheme	5 levels, n = 3
Units	ng/mL [µg/L] ^d
Typical Calibration Levels	Level 1 (CAL A): 0 ng/mL [µg/L] Level 2 (CAL B): 5 ng/mL [µg/L] Level 3 (CAL C): 100 ng/mL [µg/L] Level 4 (CAL D): 500 ng/mL [µg/L] Level 5 (CAL E): 1050 ng/mL [µg/L]
Calibration Frequency	Every 30 days for any one lot Calibration interval may be extended based on acceptable verification of calibration.

A new calibration is required:

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

d. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze two levels of a Quality Control (QC) material with known carcinoembryonic antigen concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of carcinoembryonic antigen in ng/mL [µg/L] using the calculation scheme described in your Dimension Vista® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.2 – 1000.0 ng/mL [µg/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 1000 ng/mL [µg/L] should be repeated on dilution.

Manual Dilution: Dilute with Laboratory reagent grade water up to 1:10 to obtain results within AMR range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 20 µL for serum/plasma (dilution factor = 10). Refer to your Dimension Vista® Operator's Guide.

- Samples with results less than 0.2 ng/mL [µg/L] will be reported as “less than 0.2 ng/mL [µg/L]” by the instrument.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in carcinoembryonic antigen results. Refer to your Dimension Vista® Operator’s Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies.^{13,14} Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.

Results of this test should be interpreted in light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information. Elevated levels of CEA may occur in non neoplastic conditions; therefore, this test is not intended for diagnosis or screening for cancer.

Expected Values: Non-smokers: 0.0 – 3.0 ng/mL [µg/L]
Smokers: 0.0 – 5.0 ng/mL [µg/L]

The expected values were calculated non-parametrically and represent results determined from a population of healthy adults (n = 347; 96% of 198 non-smokers and 96.6% of 149 smokers).

Distribution of CEA results:

	n	0.0 – 3.0 (%)	3.1 – 5.0 (%)	5.1 – 10.0 (%)	>10.0 (%)
Non-Smokers	198	190 (96.0)	8 (4.0)	0 (0.0)	0 (0.0)
Smokers	149	113 (75.8)	31 (20.8)	5 (3.4)	0 (0.0)
Total	347	303 (87.2)	39 (11.2)	5 (1.4)	0 (0.0)
Colorectal cancer	74*	48 (64.9)	6 (8.1)	7 (9.4)	13 (17.6)

* The 74 colorectal cancer samples were the first samples available for each serial patient. These were not baseline samples.

Each laboratory should establish its own expected values for CEA as performed on the Dimension Vista® System.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n = 5 replicates at the following analyte concentrations are:

CEA concentration	Acceptable SD Maximum
5.2 ng/mL [µg/L]	1 ng/mL [µg/L]
547 ng/mL [µg/L]	42 ng/mL [µg/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

CEA Serial Monitoring Performance

Stored retrospective serial samples from 75 male and female colorectal cancer patients were obtained from a supplier to the medical device industry. One patient was excluded from statistical analysis when it was determined that all inclusion criteria were not met. A minimum of 3 serial samples were obtained from each of the 74 available patients (total number of samples 291).

The reference change value (RCV) was used to determine if a significant change in CEA occurred.¹⁵ For this calculation, the RCV for each assay (Dimension Vista® CEA and the predicate) was derived by taking into account the published biological variation for CEA and the within lab precision of the specific assay.¹⁶ The RCV for the Vista® CEA method was calculated to be 36.2% and that of the predicate to be 36.7%.

Changes in CEA concentrations and in disease status were analyzed on a per visit basis. Patients were categorized as Active/Progressive, Responding, Stable, or No Evidence of Disease (NED) by the attending physicians based on the clinical information (medical imaging, physical examination, and other clinical investigations). All 74 patient sets were analyzed to determine how the change in disease status per sequential pair (n = 217) varied according to RCV. Table 1 shows the distribution of results for the 291 samples when compared to disease status.

Table 1. Per Visit Dimension Vista® CEA Value vs. Disease State

Change in CEA	Change in Disease State				
	Responding n (%)	Stable n (%)	No Evidence of Disease n (%)	Progression n (%)	Total
>36.2% increase	6 (2.8%)	14 (6.5%)	6 (2.8%)	32 (14.8%)	58 (26.7%)
No significant Change	12 (5.5%)	33 (15.2%)	62 (28.6%)	22 (10.1%)	129 (59.5%)
>36.2% decrease	5 (2.3%)	15 (6.9%)	6 (2.8%)	4 (1.8%)	30 (13.8%)
Total	23 (10.6%)	62 (28.6%)	74 (34.1%)	58 (26.7%)	217 (100.0%)

Per visit clinical performance results for the Dimension Vista® CEA test and predicate device are given in Table 2 & 3. In this evaluation, disease status was classified as “Progression” and “No Progression” with “No Progression” consisting of responding, stable, and no evidence of disease. Using these classifications, sensitivity and specificity were determined.

Table 2: Per Visit Dimension Vista® CEA RCV vs. Disease State

	Progression	No-Progression	Total
>36.2% increase	32	26	58
≤36.2% increase	26	133	159
Total	58	159	217

	Estimate	Exact 95% Confidence Interval
% Overall Agreement	76.0%	(69.8% – 81.6%)
% Sensitivity	55.2%	(41.5% – 68.3%)
% Specificity	83.6%	(77.0% – 89.0%)

Table 3: Per Visit Predicate CEA RCV vs. Disease State

	Progression	No-Progression	Total
>36.7% increase	32	32	64
≤36.7% increase	26	127	153
Total	58	159	217

	Estimate	Exact 95% Confidence Interval
% Overall Agreement	73.3%	(66.9% – 79.0%)
% Sensitivity	55.2%	(41.5% – 68.3%)
% Specificity	79.9%	(72.8% – 85.8%)

The per visit concordances were pooled by taking the correlation structure within each patient series into consideration.¹⁷ Efficacy is demonstrated when the sum of sensitivity and specificity is greater than one. Non-parametric estimates for the 95% confidence intervals were derived using a bootstrap resampling technique with 2000 iterations. For the Dimension Vista® CEA method the bootstrap 95% CI was 1.2416 to 1.5227 for the sum of Sensitivity + Specificity, and for the comparative method the bootstrap 95% CI for the sum of the Sensitivity + Specificity was 1.2180 to 1.4771.¹⁸ This demonstrated that both tests are effective.

All specimens were also analyzed to determine agreement between the two assays using their respective RCV. These results are shown in Table 4.

Table 4: Dimension Vista® CEA Comparison to Predicate CEA Assay (on a per visit basis)

	Predicate CEA		Total
	>36.7% increase	≤36.7% increase	
Vista CEA			
>36.2% increase	55	3	58
≤36.2% increase	9	150	159
Total	64	153	217

	Estimate	Exact 95% Confidence Interval
% Overall Agreement	94.5% (205/217)	(90.5% – 97.1%)
% Positive Agreement	85.9% (55/64)	(75.0% – 93.4%)
% Negative Agreement	98.0% (150/153)	(94.4% – 99.6%)

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Material	Precision ^{19a}		
	Mean ng/mL [µg/L]	Standard Deviation (%CV) Repeatability	Within-Lab
Liquichek™ Immunoassay Plus Control			
Level 1	2.1	0.1 (2.9)	0.1 (3.4)
Level 2	26.2	0.6 (2.2)	0.8 (2.9)
Serum pool 1	0.9	0.02 (2.3)	0.02 (2.6)
Serum pool 2	12.8	0.3 (2.6)	0.4 (3.1)
Serum pool 3	67.5	1.3 (1.9)	1.6 (2.4)
Serum pool 4	478.2	6.2 (1.3)	10.0 (2.1)
Serum pool 5	756.4	13.6 (1.8)	24.9 (3.3)
Plasma pool	239.7	5.2 (2.2)	8.6 (3.6)

e. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618.

Method Comparison²⁰

Comparative Method	Slope	Regression Statistics ^f		n
		Intercept ng/mL [µg/L]	Correlation Coefficient	
ACCESS® CEA	1.01	9.01	0.989	141 ^g
ACCESS® CEA	1.04	0.44	0.970	46 ^h

ACCESS® Immunoassay System is a registered trademark of Beckman Coulter, Inc.
f. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.
g. The range of 141 values in the correlation study was 0.8 – 974 ng/mL [µg/L].
h. The range of 46 values in the correlation study was 0.8 – 17.1 ng/mL [µg/L].

Serum and Plasma Equivalency

No clinically significant difference was observed between serum (x) and plasma (y) samples as shown in the following regression analyses.

Sample comparison	Slope	Intercept ng/mL [µg/L]	Correlation Coefficient	N
Lithium heparin versus serum	1.00	1.45	0.997	54 ⁱ
Sodium heparin versus serum	0.99	2.55	0.998	54 ⁱ

i. The range of CEA values in the correlation study was 1.2 to 992.6 ng/mL [µg/L].

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The CEA method was evaluated for interference according to CLSI/NCCLS EP7-A2.²¹

Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Carcinoembryonic antigen ng/mL [µg/L]	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer) 1000 mg/dL [0.62 mmol/L]	5 500	< 10 < 10
Bilirubin (unconjugated)	60 mg/dL [1026 µmol/L]	5 500	< 10 < 10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	5 500	< 10 < 10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	5 500	< 10 < 10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany
* Analyte results should not be corrected based on this bias.

Interfering Substances

The CEA method was evaluated for interference according to CLSI/NCCLS EP7-A2.²¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Analyte Concentration	Biotin Test Level (ng/mL)			
	100	250	500	1200
	% Bias			
5 ng/mL	0.4	-1.2	-2.3	-17.6

- Specimens that contain biotin at a concentration of 500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely depressed results for patient samples.
- The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.²² Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.²³

Non Interfering Substances

The following substances do not interfere with the CEA method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at carcinoembryonic concentration of 1.3 to 541 ng/mL [µg/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Amikacin	8 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Bevacizumab (Avastin®)	125 mg/dL	8.4 µmol/L
Bleomycin	3.3 mg/dL	23.6 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Carboplatin	90.1 mg/dL	2.4 mmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	13 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2.7 mmol/L
Cyclophosphamide	360.4 mg/dL	12.9 mmol/L
Dextran 40	6 g/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Doxorubicin	16.5 mg/dL	284 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
5-Fluorouracil	130.9 mg/dL	10.1 mmol/L
Folinic acid (Leucovorin)	6.5 mg/dL	127 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Hemoglobin	2 g/L	0.12 mmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Irinotecan hydrochloride	100 mg/dL	1.6 mmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Methotrexate	450.5 mg/dL	9.9 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Octreotide (Sandostatin®)	50 mg/dL	504 µmol/L

Oxaliplatin (Eloxatin®)	10 mg/dL	252 µmol/L
Pericillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	8 g/dL	80 g/L
Rheumatoid Factors	500 IU/mL	500 IU/mL
Salicylic Acid	60 mg/dL	4.34 mmol/L
Tamoxifen	6 mg/dL	162 µmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	3060 mg/dL	34.6 mmol/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vinblastine	4 mg/dL	44 µmol/L
Vincristine	0.44 mg/dL	4.8 µmol/L

Hook effect

One step sandwich immunoassays are susceptible to a high-dose “hook effect”, where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.²⁴ The CEA method generated signal high enough to trigger the Above Assay Range flag and thus did not hook back to generate falsely low value with CEA up to 225,000 ng/mL [µg/L].

Cross-reactivity

The following substances were evaluated for cross-reactivity with the CEA method when present in serum and plasma in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 10% at CEA concentration of 5 ng/mL [µg/L].

Substance	Concentration
NCA	500 ng/mL
NCA-2	100 ng/mL

Spiking Recovery

Known amounts of CEA, approximately 5, 15, 75, and 500 ng/mL [µg/L], were added to a human serum pool with a baseline CEA value of 3.4 ng/mL. The sample concentrations were measured and the percent recovery ranged from 94.0% to 100.4% with a mean recovery of 97.1%.

$$\% \text{ Recovery} = \frac{\text{Value obtained} - \text{Baseline}}{\text{Amount Added}} \times 100$$

Dilution Recovery

Five serum samples with CEA values from 157.5 to 751.5 ng/mL [µg/L] were diluted 1:2, 1:3, 1:4, 1:5, and 1:10 with Reagent grade water and assayed for recovery. The recoveries ranged from 98 % to 109 % with a mean of 102.9%.

Limit of Detection

The Limit of Detection (LoD - the lowest concentration that can be detected reliably) for CEA is 0.2 ng/mL [µg/L]. It was determined consistent with CLSI guideline EP17-A.²⁵ and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 15 determinations, with 4 blank and 4 low level samples. The Limit of Blank (LoB) is the highest concentration that is likely to be observed for a blank sample and is 0.12 ng/mL [µg/L].

Symbols Key: See Adjacent Panel.

Bibliography: See Adjacent Panel.

Dimension®, Dimension Vista®, LOCI® and Flex® are trademarks of Siemens Healthcare Diagnostics.

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©2008 Siemens Healthcare Diagnostics

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Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

1. Gold P, Freedman SO. Demonstration of Tumor Specific Antigens in Human Colonic Carcinomata by Immunological Tolerance and Absorption Techniques. *J Exp Med*.1965; 121: 439-462.
2. Chan DW, Booth RA, Diamandis EP. Tumor Markers in Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Fourth Edition. 2006; 768-769.
3. Paus E, Almasbak H et al. A Single-chain Fv-based Immunofluorometric Assay Specific for the CEA Variant NCA-2. *J Immunol Methods*. 2003; 286: 125-139.
4. Tietz Clinical Guide to Laboratory Tests, Fourth Edition. 2006; 220-222.
5. NACB: Practice Guidelines And Recommendations For Use Of Tumor Markers In The Clinic *Quality Requirements [Section 2]*. National Academy of Clinical Biochemistry Guidelines on Quality Requirements for the Use of Tumor Markers. Catharine Sturgeon, Elizabeth Hammond, Soo-Ling Ch'ng, György Sölétormos, Daniel F Hayes.
6. Perkins GL, Slater ED et al. Serum Tumor Markers. *Am Family Physician*.2003; 68: 1075-1082.
7. Goldstein MJ, Mitchell EP. Carcinoembryonic Antigen in the Staging and Follow-up of Patients with Colorectal Cancer. *Cancer Invest*.2005; 23: 338-351.
8. Fakih MG. CEA Monitoring in Colorectal Cancer. *Oncology*. 2006; 20.
9. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al. Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogenous immunoassay method. *Clin Chem* 42:9 1996, 1518-1526.
10. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et al. Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. *Proc.Natl. Acad. Sci. USA*, Vol 91, pp 5426 – 5430, June 1994 *Biochemistry*.
11. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
12. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition. AACC Press Washington DC 1997.
13. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays *Clin Chem* 1999; 45:7; 942-956.
14. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab)₂ conjugate and polyclonal mouse IgG. *Clin Chem* 1992; 38:1737-1742.
15. Fraser Callum G. PhD, Biological Variation: From Principles to Practice, Washington, DC: AACC Press, 2001 [pp. 71 – 78].
16. Ricos C, Alvarez V, Cava F, Garcia-Lario JV, Hernandez A, Jimenez CV, Minchinela J, Perich C, Simon M. "Current databases on biologic variation: pros, cons and progress." *Scand J Clin Lab Invest* 1999; 59: 491-500. <<http://www.westgard.com/biodatabase1.html>>
17. B. Emir, S. Wieand, John Q.S., and S. Cha. "Analysis of repeated markers used to predict progression of cancer". *Statistics in Medicine*. 1998 [2563-2578].
18. Davison AC, Hinkley DV. Bootstrap methods and their application. Cambridge, NY : Cambridge University Press, 1997.
19. Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
20. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
21. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
22. Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet*. 2017;2(4): 247–256.
23. Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med*. 2017;55(6):817-825.
24. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. *Anal Biochem* 1982;127:308-315.
25. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic *in vitro* / Dispositivo medico-diagnostico *in vitro* / Producto sanitario para diagnóstico *in vitro*

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

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All rights reserved.