Dimension® EXL™ integrated chemistry system

LOCI® Module

Dimension® EXL™ 200, Dimension® EXL™ with LM Systems

Flex® reagent cartridge

VB12

See shaded sections: Updated information from 2019-05 version.

Issue Date 2020-08

LOCI Vitamin B12

Intended Use: The VB12 method is an *in vitro* diagnostic test for the quantitative measurement of vitamin B12 (B12) in human serum and plasma on the Dimension® EXL™ with LM integrated chemistry system. Measurements of vitamin B12 may be used in the diagnosis of vitamin B12 deficiency.

Summary: Vitamin B12, or cobalamin, is an essential vitamin with a molecular mass of 1355 Da that is found in a variety of foods such as fish, shellfish, meats and dairy products. Intrinsic factor (IF), transcobalamin II (TCII) and haptocorrin (HC) are binding proteins necessary for the assimilation, transport and delivery of B12 to the blood and body tissues. Vitamin B12 is primarily stored in the liver and released on demand. The body uses B12 very efficiently, reabsorbing B12 from the small intestine and returning it to the liver so little is excreted and nutritional deficiency is extremely rare.¹ Vitamin B12 is necessary for DNA synthesis, normal red blood cell maturation and myelin sheath formation and maintenance. It is a coenzyme in the conversion of methylmalonic acid to succinic acid and in the synthesis of methionine.¹

Vitamin B12 deficiency is one of the causes of megaloblastic anemia, a disease in which red blood cells are larger than normal and the ratio of nucleus size to cell cytoplasm is increased. Since folic acid deficiency can also cause megaloblastic anemia, measurement of serum B12 levels is an important part of the differential diagnosis.¹ Vitamin B12 deficiency also causes macrocytic anemias which are characterized by abnormal red blood cell maturation and early release from the bone marrow. Pernicious anemia is a macrocytic anemia. In this disease, an absence of IF prevents normal absorption of B12. In both megaloblastic anemia caused by B12 deficiency and pernicious anemia, treatment with B12 is the therapeutic course.¹

Vitamin B12 deficiency can also lead to abnormal neurologic and psychiatric symptoms such as ataxia, muscle weakness, dementia, psychosis and mood disturbances.² Many patients show neurological changes without developing macrocytic anemia. Populations at risk for B12 deficiency include strict vegetarians, the elderly and populations with increased B12 requirements associated with pregnancy, thyrotoxicosis, hemolytic anemia, hemorrhage, malignancy and liver or kidney disease.² Early diagnosis of B12 deficiency is crucial because of the latent nature of this disorder and the risk of irreversible neurological damage. Recent studies suggest that in addition to serum B12 levels, folic acid, methylmalonic acid and homocysteine should be measured to improve the specificity of the diagnosis.^{3,4,5} Elevated B12 levels are seen in hematological disorders (chronic myelogenous leukemia, promyelocytic leukemia, polycythemia vera) and in liver disorders (acute hepatitis, cirrhosis, hepatocellular carcinoma).^{2,6}

Principles of Procedure: The vitamin B12 method is a homogeneous, competitive chemiluminescent immunoassay based on LOCI® technology. LOCI® reagents include two synthetic bead reagents and biotinylated intrinsic factor (IF). The first bead reagent (Chemibead) is coated with a B12 derivative and contains a chemiluminescent dye. The second bead reagent (Sensibead) is coated with streptavidin and contains photosensitive dye. The patient sample is pretreated with sodium hydroxide (NaOH) and dithioerythritol (DTE) to release the serum B12 from its carrier proteins. Potassium cyanide (KCN) is added to convert all the forms of B12 into a single, cyanocobalamin form, and dicyanocobinamide is added to keep the B12 from rebinding with the carrier proteins. After the sample pretreatment, the biotinylated IF and Chemibead reagents are added sequentially to the reaction vessel. Vitamin B12 from the sample competes with the B12 Chemibead for a limited amount of biotinylated IF. Sensibead reagent is then added and binds to the biotin to form bead pair immunocomplexes. Illumination of the complex at 680 nm generates single oxygen from the Sensibeads which diffuses to the Chemibeads triggering a chemiluminescent reaction. The resulting signal is measured at 612 nm and is an inverse function of the concentration of vitamin B12 in the sample.^{7,8}

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1	Liquid	VB12 Sensibead reagent	360 µg/mL	Recombinant <i>E. coli</i>
2	Liquid	VB12 Chemibead reagent	200 µg/mL	
3	Empty			
4	Liquid	Biotinylated IF	3 ng/mL	porcine
		Dicyanocobinimide	60 ng/mL	
5	Tablet	Dithioerythritol (DTE)	12.5 mg/mL	
6	Liquid	Sodium Hydroxide (NaOH)	0.75 N	
		Potassium Cyanide (KCN)	3 mM	
7 – 8	Empty			

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

Risk and Safety

H290, H314, H317
P280, P301 + P310 + P331, P303 + P361 + P353 + P310,
P305 + P351 + P338, P390, P501

Danger!

May be corrosive to metals. Causes severe skin burns and eye damage. May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Absorb spillage to prevent material damage. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Sodium hydroxide; 5-chloro-2-methyl-3(2H)-isothiazolone mixture with 2-methyl-3(2H)-isothiazolone

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions:

Used reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion. Reaction vessels are designed for single use only.

For *in vitro* diagnostic use.

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the Dimension® EXL™ system.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 2 and 4 – 6

Specimen Collection and Handling

Recommended specimen types: serum and plasma (lithium or sodium heparin, and EDTA).

Samples and controls stabilized with azide cannot be used.

Samples are stable at room temperature for up to 24 hours. Samples may be stored refrigerated at 2 – 8 °C for up to 48 hours. If testing is delayed beyond 48 hours, samples should be frozen at -20 °C or colder. Mix thoroughly after thawing. Avoid repeated freezing and thawing.⁹

Protect samples from light. Avoid using hemolyzed samples.⁹

Follow the instructions provided with your specimen collection device for use and processing.¹⁰

Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection. Specimens should be free of particulate matter.¹¹

The purpose of specimen handling and storage information is to provide guidance to users; however, users may validate their own procedures for handling and storing patient samples.

Procedure**Materials Provided**

LOCI VB12 Flex® reagent cartridge, Cat. No. RF642

Materials Required But Not Provided

LOCI ANEMIA Calibrator, Cat. No. RC640

Reaction Vessels, Cat. No. RXV1

Reagent Probe Cleaner, Cat. No. RD702

Sample Probe Cleaner, Cat. No. RD703

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension® EXL™ System. For details of this processing, refer to your Dimension® EXL™ Operator's Guide.

Test Conditions

Sample Volume (delivered to the vessel)	12 µL
Extractant Reagent Volume	15 µL
DTE	12 µL
Biotinylated IF/Neutralizer Reagent Volume	50 µL
Chemibead Reagent Volume	15 µL
Sensibead Reagent Volume	50 µL
Temperature	37 °C
Reaction time	32 minutes
Wavelength	Illumination 680 nm, emission 612 nm
Type of measurement	Chemiluminescence

Calibration

Calibration Material LOCI ANEMIA Calibrator, Cat. No. RC640

Calibration Scheme Level 1, n = 5
Levels 2 – 5, n = 3

Units pg/mL [pmol/L]^c
(pg/mL x 0.7378) = [pmol/L]

Typical Calibration Levels
Level 1: 45 pg/mL [33 pmol/L]
Level 2: 200 pg/mL [148 pmol/L]
Level 3: 500 pg/mL [369 pmol/L]
Level 4: 1000 pg/mL [738 pmol/L]
Level 5: 2100 pg/mL [1549 pmol/L]

Calibration Frequency Every 21 days for any one lot

- A new calibration is required:
- For each new lot of Flex® reagent cartridges
 - After major maintenance or service, if indicated by quality control results
 - As indicated in laboratory quality control procedures
 - When required by government regulations

c. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known vitamin B12 concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates the concentration of vitamin B12 in pg/mL [pmol/L] using the calculation scheme described in the Dimension® EXL™ Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 80 – 2000 pg/mL [59 – 1476 pmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 2000 pg/mL [1476 pmol/L] are reported as "Above Assay Range" and should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water to obtain results within the analytical measurement range. Recommended dilution factor = 3. Enter dilution factor on the instrument. Re-assay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 4 µL (dilution factor = 3) for serum/plasma. Refer to your Dimension® system manual.

- Samples with results less than 80 pg/mL [59 pmol/L] should be reported as "less than 80 pg/mL [59 pmol/L]."

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in vitamin B12 results. Refer to your Dimension® EXL™ Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{12,13}

Intrinsic Factor blocking antibodies are present in approximately half of pernicious anemia patients.¹⁴ There is a low frequency possibility that these antibodies may not be completely inactivated during the reaction pretreatment step. If test results are in conflict with the clinical diagnosis, the sample can be tested for intrinsic factor blocking antibodies.

United States Federal law restricts this device to sale by or on the order of a physician or other practitioner licensed by the laws of the state in which he practices, to use or order the use of this device.

Interfering Substances

The VB12 method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁵ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Dextran 40 at 6 g/dL [1500 µmol/L] decreases VB12 results by -14.9% at 200 pg/mL [148 pmol/L] and by -11.8% at 1000 pg/mL [738 pmol/L].

Hemoglobin at 500 mg/dL [0.31 mmol/L] increases VB12 results by 15.2% at 200 pg/mL [148 pmol/L] and less than 10% at 1000 pg/mL [738 pmol/L].

Analyte Concentration	Biotin Test Level (ng/mL)				
	50	100	250	500	1200
	% Bias				
200 pg/ml	1.3	5.9	46	126.1	Above AMR
1000 pg/ml	N/A	2.3	12.4	40.4	Below AMR

• Specimens that contain biotin at a concentration of 100 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely elevated results for patient samples.

• The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.¹⁶ Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.¹⁷

Expected Values:

Region	Range	
	pg/mL	pmol/L
US	193 – 986	142 – 727
Europe	182 – 625	134 – 461

To evaluate differences that may be observed with respect to population and dietary status, each laboratory should establish its own expected values for vitamin B12 as performed on the Dimension® EXL™ system. The values shown in the table above were confirmed using a transference protocol from earlier studies performed using samples from a population of apparently healthy adults obtained in the United States and Germany.

A system malfunction may exist if the acceptable SD maximum is exceeded

VB12 concentration	Acceptable SD Maximum
200 pg/mL [148 pmol/L]	42 pg/mL [31 pmol/L]
1000 pg/mL [738 pmol/L]	104 pg/mL [77 pmol/L]

Specific Performance Characteristics

The following data was collected on the Dimension® EXL™ system and represents typical performance.

Material	Mean pg/mL [pmol/L]	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Precision ^{18,d}			
Bio-Rad Liquichek™ Immunoassay Control			
Level 1	290 [214]	13.2 [9.7] (4.6)	15.0 [11.1] (5.2)
Level 2	498 [367]	11.7 [8.6] (2.3)	18.2 [13.5] (3.7)
Level 3	645 [476]	15.9 [11.7] (2.5)	21.5 [15.9] (3.3)
Serum Pool 1	180 [133]	10.1 [7.4] (5.6)	11.6 [8.6] (6.5)
Plasma Pool (Lithium heparin)	467 [345]	13.1 [9.7] (2.8)	15.8 [11.7] (3.4)
Serum Pool 2	978 [722]	24.9 [18.4] (2.5)	27.7 [20.4] (2.8)
Serum Pool 3	1733 [1279]	27.4 [20.2] (1.6)	35.2 [26.0] (2.0)

d. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618 USA

Comparative Method	Slope	95% CI	Method Comparison ¹⁹ Regression Statistics ^e		Correlation Coefficient	n
			Intercept pg/mL [pmol/L]	95% CI pg/mL [pmol/L]		
Dimension Vista VB12	1.03	1.02 – 1.04	-6.4 [-4.7]	-12.9 – 0.4 [-9.5 – 0.3]	0.997	166 ^f

e. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was Passing-Bablok

f. The range of vitamin B12 values in the method comparison study was 60 – 1966 pg/mL [44 – 1451 pmol/L]. Correlation coefficient was calculated from least squares regression.

Serum and Plasma Equivalence

Comparison testing results on the Dimension® EXL™ system for serum separator tubes (SST), lithium heparin plasma, sodium heparin and EDTA plasma versus serum from a plain red top tube is provided.

Comparative specimen	Slope	Intercept pg/mL [pmol/L]	n
Lithium heparin plasma	1.00	-4.33 [-3.19]	78
Sodium heparin plasma	1.02	-10.73 [-7.92]	78
EDTA plasma	1.00	-7.26 [-5.36]	78
SST serum	1.01	-2.31 [-1.70]	78

Specificity

Hemolysis, Icterus and Lipemia (HIL) interference

The VB12 method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁵ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Vitamin B12 pg/mL [pmol/L]	Bias* %
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L]	200 [148]	15
		1000 [738]	<10
		200 [148]	10
Bilirubin (unconjugated)	300 mg/dL [0.19 mmol/L]	200 [148]	<10
	60 mg/dL [1026 µmol/L]	1000 [738]	<10
Bilirubin (conjugated)	20 mg/dL [342 µmol/L]	200 [148]	<10
		1000 [738]	<10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	200 [148]	<10
		1000 [738]	<10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the VB12 method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances do not exceed 10% at vitamin B12 concentrations of approximately 200 pg/mL [148 pmol/L] and 1000 pg/mL [738 pmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Amikacin	8.0 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5.0 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2.7 mmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L

Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25,000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Rheumatoid Factor	500 IU/mL	500,000 IU/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	3000 mg/dL	33.9 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vancomycin	10 mg/dL	69 µmol/L

Cross Reactivity

Cobinamide at a concentration of 200 ng/mL [192 nmol/L] was evaluated and found to have insignificant cross reactivity with the VB12 assay when present in serum containing approximately 200 pg/mL [148 pmol/L] of vitamin B12.

The percent cross reactivity was calculated as follows:

$$\% \text{ Cross-reactivity} = \frac{[\text{measured VB12}] - [\text{control VB12}]}{[\text{compound}]} \times 100$$

Recovery

Known amounts of vitamin B12 were added to 2 human serum, 2 SST, 2 lithium heparin, 2 sodium heparin and 2 EDTA plasma samples with baseline B12 value between 136 – 197 pg/mL [100 – 145 pmol/L]. The sample concentrations were measured and the percent recovery ranged from 92.0% to 108.9% with a mean recovery of 99.8%.

$$\% \text{ Recovery} = \frac{[\text{Value obtained}] - [\text{Baseline}]}{[\text{Amount Added}]} \times 100$$

Accuracy to WHO Standard 03/178

A study was performed to evaluate the accuracy of the VB12 assay using the World Health Organization Vitamin B12 International Standard 03/178 using two Flex® reagent lots on two instruments. The assay demonstrated a 2.1% difference from the target value of 480 pg/mL [354 pmol/L].

Limit of Detection and Limit of Blank^g

The Limit of Detection (LoD) for VB12 is 50 pg/mL [37 pmol/L], determined consistent with CLSI guideline EP17-A²⁰ and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 168 determinations, with 5 blank and 9 low level samples. The Limit of Blank (LoB) is 30 pg/mL [22 pmol/L].

g. LoD is the lowest concentration of analyte that can be detected reliably. LoB is the highest concentration that is likely to be observed for a blank sample.

Limit of Quantitation

The Limit of Quantitation (LoQ) for VB12 is 80 pg/mL [59 pmol/L], determined consistent with CLSI Guideline EP17-A and an inter-assay imprecision of <20% Coefficient of Variation (CV).

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Dimension Vista®, LOCI® and Flex® are trademarks of Siemens Healthcare Diagnostics.

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© 2014 Siemens Healthcare Diagnostics

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Shenkin A, Baines M, Fell G S et al. Vitamins and Trace Elements. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Fourth Edition, 2006; 1100-1105.
2. Tietz Clinical Guide to Laboratory Tests. Fourth Edition 2006; 1124-1127.
3. Snow CF. *Laboratory Diagnosis of Vitamin B12 and Folate Deficiency*. Arch Intern Med. 1999; 159:1289-1298.
4. Ward PCJ. *Modern Approaches to the Investigation of Vitamin B12 Deficiency*. Clin Lab Med. 2002; 22:435-445.
5. Klee GG. *Cobalamin and Folate Evaluation: Measurement of Methylmalonic Acid and Homocysteine vs Vitamin B12 and Folate*. Clin Chem. 2000;46: 1277-1283.
6. Ermens AAM, Vlasveld LT, Lindemans J. *Significance of Elevated Cobalamin (Vitamin B12) Levels in Blood*. Clin Biochem. 2003; 36:585-590.
7. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al., *Luminescent oxygen channeling assay (LOCI®): sensitive, broadly applicable homogeneous immunoassay method*. Clin Chem. 1996; 42:9:1518-1526.
8. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et al. *Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence*. Proc.Natl.Acad.Sci. USA, 1994; 91:5426-5430.
9. Tietz Clinical Guide to Laboratory Tests. Fourth Edition. Alan H. B. Wu editor, 2006; 410-413, 1124-1127.
10. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
11. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
12. Kricka LJ. *Human anti-animal antibody interference in immunological assays*. Clin Chem. 1999;45: 7:942-956.
13. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab)2 conjugate and polyclonal mouse IgG. Clin Chem 1992; 38: 1737-1742.
14. Klee GG. Cobalamin and Folate Evaluation: Measurement of Methylmalonic Acid and Homocysteine vs Vitamin B12 and Folate. Clin. Chem. 2000; 46:8(B), 1277-1283.
15. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
16. Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. Int. J. Pharmacokinet. 2017;2(4):247-256.
17. Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 2017;55(6):817-825
18. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A (ISBN 1-56238-551-8). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898,USA 2004.
19. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
20. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A (ISBN 1-56238-551-8). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898,USA 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_ERGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

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All rights reserved.