Dimension® EXL™ integrated chemistry system LOCI® Module

Dimension® EXL™ 200, Dimension® EXL™ with LM Systems

Flex® reagent cartridge

VITD

Issue Date 2020-08

See shaded sections: Updated information from 2019-04 version.

LOCI Vitamin D Total Assay

Intended Use: The LOCI Vitamin D Total assay is an *in vitro* diagnostic test for the quantitative measurement of total 25(OH)vitamin D in human serum and plasma on the Dimension® EXL™ integrated chemistry system with LOCI® Module. Measurements of vitamin D are used in the assessment of vitamin D sufficiency.

Summary: Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. Vitamin D is essential for the formation and maintenance of strong, healthy bones. Vitamin D deficiency can result from inadequate exposure to the sun, inadequate alimentary intake, decreased absorption, abnormal metabolism, or vitamin D resistance.¹

While there are many metabolites of vitamin D, the total 25(OH)vitamin D, that is, the sum of 25(OH)vitamin D₂ and 25(OH)vitamin D₃, is the most reliable indicator of vitamin D status.^{2,3} Vitamin D₂ is derived from plant sources, whereas vitamin D₃ is derived primarily from the conversion of 7-dehydrocholesterol in the skin by UV-B radiation from sunlight and secondarily from animal sources.² Whether consumed or produced, both forms of vitamin D (D₂ and D₃) are metabolized by the liver to 25(OH)vitamin D, and then converted in the liver or kidney into the active metabolites 1,25(OH)₂vitamin D₂ and 1,25(OH)₂vitamin D₃.^{2,4}

Vitamin D metabolites are bound to a carrier protein in the plasma and distributed throughout the body. The most reliable clinical indicator of vitamin D status is 25(OH)vitamin D because serum and plasma 25(OH)vitamin D levels reflect the body's storage levels of vitamin D, and 25(OH)vitamin D correlates with the clinical symptoms of vitamin D deficiency.⁵

Principles of Procedure: The LOCI Vitamin D Total assay is a homogeneous competitive chemiluminescent immunoassay based on LOCI® technology. The assay measures the total 25(OH)vitamin D concentration [comprising both 25(OH)vitamin D₂ and 25(OH)vitamin D₃] in both serum and plasma. LOCI Vitamin D Total reagents include a releasing reagent, biotinylated monoclonal antibody, and two synthetic bead reagents. Patient sample is incubated with the releasing reagent to release 25(OH)vitamin D molecules from the vitamin D-binding proteins. The reaction mixture is then incubated with biotinylated antibody to form a 25(OH)vitamin D/biotinylated antibody complex.

Chemibeads containing 25(OH)vitamin D₃ analog and chemiluminescent dye are added to remove the excess free biotinylated antibody. Streptavidin-coated Sensibeads containing a photosensitive dye are added to bind the biotinylated antibody. Aggregates of the Chemibead analog/biotinylated antibody/streptavidin Sensibeads are formed as a result. Illumination of the reaction mixture by light at 680 nm generates singlet oxygen from the Sensibeads, which diffuses into the Chemibeads and triggers a chemiluminescent reaction. The resulting chemiluminescent signal is measured at 612 nm and is inversely proportional to the concentration of total 25(OH)vitamin D in the sample.

Reagents

Wells ^{a,b}	Form	Ingredient	Concentration ^c	Source
1 – 2 (Reagent 4)	Liquid	Releasing Reagent		
3 – 4 (Reagent 1)	Liquid	Biotinylated antibody	0.75 µg/mL	Sheep monoclonal
5 – 6 (Reagent 2)	Liquid	Chemibeads	200 µg/mL	Sheep monoclonal
7 – 8 (Reagent 3)	Liquid	Sensibeads	500 µg/mL	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Contain buffers, stabilizers and preservatives.

c. Nominal value per well in a cartridge.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P363, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Wash contaminated clothing before reuse. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions:

Caution: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Used cuvettes and/or reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C.

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 7 days for wells 1 – 8

Specimen Collection and Handling:

Recommended specimen types: serum, lithium heparin plasma, and EDTA plasma.

Samples and controls stabilized with sodium azide cannot be used.

Collect serum and plasma specimens using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁶ Processed samples may be refrigerated at 2 – 8 °C for 7 days if not tested within 24 hours. For longer storage, samples may be frozen at -20 °C for three months. Samples may be frozen up to 3 times. Thawed samples should be mixed thoroughly before analysis.

Specimens must be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation. If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of particulate matter.⁷

Follow the instructions provided with your specimen collection device for use and processing.⁸ The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Procedure

Materials Provided

VITD Flex® reagent cartridge, Cat. No. RF634

Materials Required But Not Provided

VITD CAL, Cat. No. RC634

Reaction Vessels, Cat. No. RXV1A

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension® EXL™. For details of this processing, refer to your Dimension® EXL™ Operator's Guide.

The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Volume (delivered to the reaction vessel)	8 µL
Reagent 4 Releasing Reagent Volume	18 µL
Reagent 1 Biotinylated Antibody Reagent Volume	24 µL
Reagent 2 Chemibeads Reagent Volume	24 µL
Reagent 3 Sensibeads Reagent Volume	84 µL
Temperature	37 °C
Reaction time	31.0 minutes
Wavelength	Illumination (680 nm), Emission (612 nm)
Type of Measurement	Chemiluminescence

Calibration

Assay Range 5.0 – 150.0 ng/mL [12.5 – 375.0 nmol/L]^d

Calibration Material VITD Calibrator, Cat. No. RC634

Calibration Scheme 5 levels, n = 3

Units ng/mL [nmol/L]
(ng/mL x 2.5) = [nmol/L]

Typical Calibration or Verification Levels
Level 1: 0.0 ng/mL [0.0 nmol/L]
Level 2: 12.0 ng/mL [30.0 nmol/L]
Level 3: 30.0 ng/mL [75.0 nmol/L]
Level 4: 75.0 ng/mL [187.5 nmol/L]
Level 5: 165.0 ng/mL [412.5 nmol/L]

Calibration Frequency Every 7 days for any one lot

A new calibration is required:

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

Assigned Coefficients
C₀: 809.0
C₁: -811.0
C₂: -1.46
C₃: 2.48
C₄: 0.5

d. Système International d'Unités [SI Units] are in brackets.

Standardization

The LOCI Vitamin D Total assay is standardized using internal standards which are traceable to the ID-LC-MS/MS 25(OH)vitamin D Reference Method Procedure (RMP). The ID-LC-MS/MS is traceable to the National Institute of Standards and Technology (NIST) Standard Reference Material (SRM) 2972.

For 40 samples in the range of 7.9 ng/mL [19.8 nmol/L] to 141.9 ng/mL [354.8 nmol/L] the relationship between the LOCI Vitamin D Total assay(y) and the ID-LC-MS/MS 25(OH)vitamin D RMP(x) is described using the Deming regression:

$$\text{LOCI Vitamin D Total assay} = 0.99 (\text{ID-LC-MS/MS}) + 4.6 \text{ ng/mL [11.5 nmol/L]}, r = 1.000$$

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known vitamin D concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of total vitamin D in ng/mL [nmol/L] using the calculation scheme described in your Dimension® EXL™ Operator’s Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5.0 – 150.0 ng/mL [12.5 – 375.0 nmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 150.0 ng/mL [375.0 nmol/L] are reported as “Above Assay Range”. Manual dilution is not required. There is no autodilute feature for VITD.
- Samples with results less than 5.0 ng/mL [12.5 nmol/L] should be reported as “less than 5.0 ng/mL [12.5 nmol/L]”.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in VITD results. Refer to your Dimension® EXL™ Operator’s Guide for the meaning of report flags and comments. Any result containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported.

Do not use hemolyzed samples. Hemoglobin at 500 mg/dL causes an increase in the vitamin D concentration by 11% at 15 ng/mL [37.5 nmol/L], 19% at 31.7 ng/mL [79.2 nmol/L] and 13% at 75.1 ng/mL [187.7 nmol/L].

Do not use lipemic samples. Lipemia (Intralipid) at 500 mg/dL causes a decrease in vitamin D concentration by 10% at 14.9 ng/mL [37.2 nmol/L], 11% at 32.3 ng/mL [80.8 nmol/L] and 8% at 75.9 ng/mL [189.8 nmol/L].

Paricalcitol (Zemlar) at 24 ng/mL has shown a 93.8% cross-reactivity at 28.8 ng/mL [72.0 nmol/L] of vitamin D and 70.8% cross-reactivity at 54.7 ng/mL [136.8 nmol/L] of vitamin D.

Interfering Substances

The VITD assay was evaluated for interference according to CLSI EP07-A2.⁹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. *Bias exceeding 10% is considered interference.

- Total Protein at 19 g/dL causes a decrease in the vitamin D concentration by 6% at 14.3 ng/mL [35.8 nmol/L], 11% at 30.9 ng/mL [77.3 nmol/L] and 12% at 77.3 ng/mL [193.3 nmol/L].
- Unconjugated bilirubin at 80 mg/dL causes an increase in the vitamin D concentration by 13% at 13.5 ng/mL [33.7 nmol/L], 12% at 27.3 ng/mL [68.2 nmol/L] and 6% at 66.7 ng/mL [166 nmol/L].
- Conjugated bilirubin at 40 mg/dL causes an increase in vitamin D concentration by 17% at 14.9 ng/mL [37.2 nmol/L], 12% at 31.7 ng/mL [79.2 nmol/L] and 11% at 75.3 ng/mL [188.2 nmol/L].
- Cholesterol at 350 mg/dL causes a decrease in vitamin D concentration by 10% at 25.2 ng/mL [63.0 nmol/L] and 12% at 64.0 ng/mL [160.0 nmol/L].
- HAMA Anti-Mouse Antibodies (HAMA): Patient samples may contain heterophilic antibodies, of which human anti-mouse antibodies (HAMA) are the most commonly encountered. Heterophilic antibodies such as HAMA may cause falsely decreased or elevated results in immunoassay tests. A test result that is inconsistent with the patient’s clinical presentation and history should be interpreted with caution.

Siemens characterized the impact of HAMA in accordance with CLSI EP07-A2.⁹ A total of 20 samples containing vitamin D at 5.0 ng/mL [12.5 nmol/L] to 87.5 ng/mL [218.7 nmol/L] with varying concentrations of HAMA, from 8 ng/mL to 161,000 ng/mL were tested for HAMA interference. Mean (n=5) percent bias of 2.8% was observed on a sample with a vitamin D concentration of 28.3 ng/mL [70.7 nmol/L] containing 161,000 ng/mL of HAMA. Maximum mean (n=5) percent bias of 10% was observed on a sample with a vitamin D concentration of 5.0 ng/mL [12.5 nmol/L] containing 106 ng/mL of HAMA.

***Analyte results should not be corrected based on observed bias.**
See additional Hemolysis, Icterus, Lipemia (HIL) results under Specificity section.

Analyte Concentration	Biotin Test Level (ng/mL)			
	100	250	500	1200
	% Bias			
15 ng/mL	4.5	6.1	12.1	56.1
30 ng/mL	0.7	1.3	8.8	36.2
75 ng/mL	0.8	2.2	5.5	26.4

- Specimens that contain biotin at a concentration of 250 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely elevated results for patient samples.
- The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.¹⁰ Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.¹¹

Health Based Reference Values:

Health-based reference values based on the Clinical Guidelines Subcommittee of the Endocrine Society Task Force are recommended to replace population based reference values. A review of the available literature¹²⁻¹⁸ suggests that the recommendations for 25(OH)vitamin D levels are summarized in the table below. Consult the listed references for discussion of vitamin D toxicity levels.^{19, 20}

Vitamin D Status	Range, Adult ¹²⁻¹⁵ 25(OH)vitamin D ng/mL [nmol/L]
Deficient	< 20 [< 50]
Insufficient	20 – < 30 [50 – < 75]
Sufficient	30 – 100 [75 – 250]

Observed Values in an Apparently Healthy Population:

Data using the LOCI Vitamin D Total assay were obtained on serum samples collected from 252 adults: 246 adults not taking supplements containing vitamin D, and 6 adults taking supplements containing vitamin D. To represent a broad spectrum of UV light exposure in intended use population, the blood samples tested were collected during different seasons and from subjects residing in diverse U.S. geographical locations. Samples with abnormal values for PTH, calcium, magnesium, phosphorus, and TSH were excluded from this study.

Based on the 95% confidence interval, the following values were established following CLSI guideline EP28-A3c.²¹

To assure proper representation of specific populations, each laboratory should establish its own reference intervals.

The following 25(OH)vitamin D values were obtained:

Median	28.3 ng/mL [70.8 nmol/L]
Observed 95% interval	13.9 – 61.0 ng/mL
(2.5 th and 97.5 th percentiles)	[34.8 – 152.5 nmol/L]

A system malfunction may exist if the maximum 5-test SD precision is exceeded.

Concentration	Maximum 5-test SD
12.0 ng/mL [30.0 nmol/L]	1.5 ng/mL [3.8 nmol/L]
30.0 ng/mL [75.0 nmol/L]	2.0 ng/mL [5.0 nmol/L]

Specific Performance Characteristics

The following data represent typical performance for the VITD assay on the Dimension® EXL™.

Material	Mean ng/mL [nmol/L]	Precision ^{22,a}		Within-Lab Precision	
		Repeatability SD	%CV	SD	%CV
	ng/mL [nmol/L]	ng/mL [nmol/L]		ng/mL [nmol/L]	%CV
Tri-Level Vitamin D Plus ^f					
QC (Low)	18.9 [47.3]	0.58 [1.45]	3.1	1.01 [2.53]	5.4
QC (Level 1)	38.7 [96.8]	1.02 [2.55]	2.6	2.02 [5.05]	5.2
QC (Level 2)	89.6 [224.0]	1.72 [4.29]	1.9	3.67 [9.17]	4.1
Serum 1	8.2 [20.5]	0.46 [1.15]	5.6	0.71 [1.79]	8.7
Serum 2	29.4 [73.5]	0.76 [1.90]	2.6	1.46 [3.65]	5.0
Serum 3	76.5 [191.3]	1.63 [4.08]	2.1	3.11 [7.78]	4.1
Plasma	25.2 [63.0]	0.44 [1.11]	1.8	0.78 [1.96]	3.1

- e. CLSI EP05-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days for a total of 80 replicates.
- f. UTAK Tri-Level Vitamin D Plus Serum Control is a product of UTAK Laboratories, Valencia, CA 91355.

Method Comparison²³

Method comparison study based on CLSI EP09-A3 between the LOCI Vitamin D Total Assay and an ID-LC-MS/MS reference measurement procedure (RMP) developed by University of Ghent generated the following results using Passing-Bablok regression and Pearson correlation. The number and range of vitamin D values on these samples are indicated below.

Comparative Method	Slope ^a	Regression Statistics		Correlation Coefficient	n ^b
		Intercept ^a ng/mL [nmol/L]			
ID-LC-MS/MS 25(OH)vitamin D	1.06	0.4 [1.1]		0.977	163

- g. Slope and intercept were derived using Passing-Bablok analysis:
y = LOCI Vitamin D Total assay; x = ID-LC-MS/MS 25(OH)vitamin D
- h. Vitamin D values ranged from 5.2 – 126.1 ng/mL [13.0 – 315.3 nmol/L] derived from 163 samples, 154 native and 9 contrived samples.

Serum and Plasma Equivalency

No clinically significant difference was observed between serum (x) and plasma (y) as shown in the following ordinary least squares regression.

Sample Comparison	Slope	Intercept ng/mL [nmol/L]	Correlation Coefficient	n	Min ng/mL [nmol/L]	Max ng/mL [nmol/L]
Lithium heparin plasma vs serum (red top)	0.99	0.1 [0.3]	0.992	70	11.5 [28.7]	146.4 [366.0]
EDTA plasma vs serum (red top)	0.98	1.4 [3.4]	0.997	70	11.6 [29.0]	142.5 [356.2]
Serum SST vs serum (red top)	0.99	0.5 [1.2]	0.997	70	11.4 [28.5]	148.6 [371.5]

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The VITD assay was evaluated for interference according to CLSI EP07-A2.⁹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Vitamin D ng/mL [nmol/L]	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer) 250 mg/dL [0.155 mmol/L]	15.0 [37.5]	5
		31.7 [79.2]	9
		75.1 [187.7]	8
		15.0 [37.5]	11
		31.7 [79.2]	19
Bilirubin (unconjugated)	500 mg/dL [0.31 mmol/L]	75.3 [188.7]	13
		13.5 [33.7]	7
		27.3 [68.2]	7
		66.7 [166.7]	8
		13.5 [33.7]	13
Bilirubin (conjugated)	40 mg/dL [684.5 µmol/L]	31.7 [79.2]	12
		66.7 [166.7]	6
		14.9 [37.2]	7
		31.7 [79.2]	4
		75.3 [188.2]	3
Lipemia	20 mg/dL [342 µmol/L]	14.9 [37.2]	17
		31.7 [79.2]	12
		75.3 [188.2]	11
		14.9 [37.2]	-5
		32.3 [80.8]	-6
	40 mg/dL [684 µmol/L]	75.9 [189.8]	-5
		14.9 [37.2]	-10
		32.3 [80.8]	-11
		75.9 [189.8]	-8

* Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances were evaluated for interference using CLSI EP07-A29 and found not to interfere with the VITD assay when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances do not exceed 10% at vitamin D concentrations of 13.9 – 16.9 ng/mL [34.8 – 42.3 nmol/L], 28.4 – 32.0 ng/mL [71.0 – 80.0 nmol/L] and 70.2 – 77.3 ng/mL [175.5 – 193.3 nmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1322 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Bilirubin (unconjugated)	40 mg/dL	684.5 µmol/L
Bilirubin (conjugated)	20 mg/dL	342 µmol/L
Cholesterol	300 mg/dL	7.7 mmol/L
Dextran 40	6000 mg/dL	1.50 mmol/L
Hemoglobin	250 mg/dL	0.155 mmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	30 mg/dL	1455 µmol/L
Lipemia (Intralipid)	300 mg/dL	3 g/L
Protein (Albumin)	5 g/dL	50 g/L
Protein (Total)	15.9 g/dL	159 g/L
Rheumatoid Factor	500 IU/mL	500 IU/mL
Salicylic Acid	65 mg/dL	4.7 mmol/L
Triglycerides	686 mg/dL	7.75 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L

Cross Reactivity

Cross reactivity testing was performed to determine the effect of potential cross-reactants on the LOCI Vitamin D Total assay using CLSI EP07-A2.⁹

Percent cross-reactivity was calculated as follows:

$$\% \text{ Cross-reactivity} = \frac{[\text{measured vitamin D}] - [\text{control vitamin D}]}{[\text{cross-reactant}]} \times 100$$

The following substances were evaluated for cross-reactivity with the VITD assay when present in serum at the vitamin D concentrations listed in the following table.

Cross-reactant	Cross-reactant Concentration (ng/mL)	% Cross-reactivity @ Vitamin D Concentration	% Cross-reactivity @ Vitamin D Concentration
Vitamin D ₂ (ergocalciferol)	1000	0.1 @ 27.5 ng/mL [68.8 nmol/L]	-0.1 @ 52.0 ng/mL [130.0 nmol/L]
Vitamin D ₃ (cholecalciferol)	1000	-0.1 @ 27.5 ng/mL [68.8 nmol/L]	0.1 @ 52.0 ng/mL [130.0 nmol/L]
* 1,25(OH) ₂ vitamin D ₂	0.5	151.9 @ 27.5 ng/mL [68.8 nmol/L]	1.3 @ 52.0 ng/mL [130.0 nmol/L]
* 1,25(OH) ₂ vitamin D ₃	0.5	1.7 @ 27.5 ng/mL [68.8 nmol/L]	-224.5 @ 52.0 ng/mL [130.0 nmol/L]
3-epi-25(OH)vitamin D ₃	100	3.9 @ 27.5 ng/mL [68.8 nmol/L]	2.5 @ 52.0 ng/mL [130.0 nmol/L]
Paricalcitol/Zemplar	24	93.8 @ 28.8 ng/mL [72.0 nmol/L]	70.8 @ 54.7 ng/mL [136.8 nmol/L]
1αOH Vitamin D ₃ (alfacalcidol)	3000	0.1 @ 26.5 ng/mL [66.3 nmol/L]	0.0 @ 51.6 ng/mL [129.0 nmol/L]
24, 25(OH) ₂ vitamin D ₃	60	2.5 @ 26.5 ng/mL [66.3 nmol/L]	1.2 @ 51.6 ng/mL [129.0 nmol/L]
25(OH)vitamin D ₂	30	95.1 @ 28.3 ng/mL [70.8 nmol/L]	94.1 @ 54.7 ng/mL [136.8 nmol/L]
25(OH)vitamin D ₃	30	88.6 @ 28.3 ng/mL [70.8 nmol/L]	90.1 @ 54.7 ng/mL [136.8 nmol/L]

* 1,25(OH)₂vitamin D₂ and 1,25(OH)₂vitamin D₃ were tested at supra-physiological concentrations.²⁴

Limit of Detection¹ and Limit of Blank^{25,j}

The Limit of Detection (LoD) for the VITD assay is 1.3 ng/mL [3.3 nmol/L], determined consistent with CLSI guideline EP17-A2 and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 120 determinations, with 60 blank and 60 low level replicates. The Limit of Blank (LoB) is 0.8 ng/mL [2.0 nmol/L].

i. LoD is the lowest concentration of analyte that can be detected reliably. LoB is the highest concentration that is likely to be observed for a blank sample.

Limit of Quantitation^{25,j}

The limit of quantitation (LoQ) was determined consistent with the CLSI Guideline EP17-A2. The LoQ for the VITD assay is 5.0 ng/mL [12.5 nmol/L] at a total precision of ≤ 20% CV.

j. LoQ is defined as the lowest concentration of 25(OH)vitamin D that can be detected at total precision CV of ≤ 20%.

Linearity²⁵

Linearity across the assay range, 5.0 – 150.0 ng/mL [12.5 – 375.0 nmol/L] was confirmed by testing a sample with a high concentration of vitamin D and serially diluting with a sample with low concentration of vitamin D. Each dilution was assayed in replicates of n=5. Linear regression of observed (y) vs. expected values (x) yielded a slope of 1.02, an intercept of 1.4 ng/mL and a correlation coefficient (r) of 0.998. Percent bias between observed and predicted values was within the allowable 10% deviation from the linear fit.

Recovery

Known amounts of 25(OH)vitamin D₂ (30 ng/mL [75.0 nmol/L]) and 25(OH)vitamin D₃ (30 ng/mL [75.0 nmol/L]) were added to separate serum samples containing vitamin D at 28.3 ng/mL [70.8 nmol/L] and 54.7 ng/mL [136.8 nmol/L].

The calculated percent recoveries are listed below:

Vitamin D at 28.3 ng/mL [70.8 nmol/L]		
Added	Amount added ng/mL [nmol/L]	% Recovery
25(OH)vitamin D ₂	30.0 [75.0]	95
25(OH)vitamin D ₃	30.0 [75.0]	89
Vitamin D at 54.7 ng/mL [136.8 nmol/L]		
Added	Amount added ng/mL [nmol/L]	% Recovery
25(OH)vitamin D ₂	30.0 [75.0]	94
25(OH)vitamin D ₃	30.0 [75.0]	90

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, EXL™, LOCI® and Flex® are trademarks of Siemens Healthcare Diagnostics.

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©2016 Siemens Healthcare Diagnostics

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

- Holick MF. Vitamin D Deficiency. *N Engl J Med*. 2007; 357: 266–81.
- Holick MF. Vitamin D deficiency. *N Engl J Med* 2007 July 19; 357(3): 266–81.
- Lai JK, Lucas RM, Clements MS, Harrison SL, Banks E. Assessing vitamin D status: pitfalls for the unwary. *Mol Nutr Food Res* 2010 August; 54(8): 1062–71.
- Rollins G. Vitamin D testing—what's the right answer? Labs grapple with confusing analytics, evidence. *Clin Lab News* 2009 July; 35(7): 1–9.
- Hart GR. Medical Conditions Associated with Vitamin D Deficiency and the Clinical Consequences. Review Series Volume 1, Sept 2004; 1–8.
- CLSI. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. CLSI document GP-41-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
- CLSI. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. CLSI document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI document GP44-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet*. 2017;2(4): 247–256.
- Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med*. 2017;55(6):817–825
- Holick MF et al. Evaluation, Treatment, and Prevention of Vitamin D Deficiency: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab* 2011; 96 (7): 1911–1930.
- Bischoff-Ferrari HA, Giovannucci E, Willett WC, et al. Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes. *Am J Clin Nutr* 2006; 84:18–28.
- Holick MF. Vitamin D status: measurement, interpretation, and clinical application. *Ann Epidemiol*. 2009; 19(2): 73–78.
- KDOQI Clinical Practice Guidelines for Bone Metabolism and Disease in Children With Chronic Kidney Disease. http://www.kidney.org/PROFESSIONALS/kdoqi/guidelines_pedbone/guide8.htm
- Soubrierbielle JC, Body JJ, Lappe JM, et al. Vitamin D and musculoskeletal health, cardiovascular disease, autoimmunity and cancer: Recommendations for clinical practice. *Autoimmun Rev* 2010; 9: 709–715.
- Dawson-Hughes B, Heaney RP, Holick MF. Estimates of optimal vitamin D status. *Osteoporos Int*. 2005; 16(7): 713–716.
- Vieth R. Why the minimum desirable serum 25-hydroxyvitamin D level should be 75 nmol/L (30 ng/mL). *Best Pract Res Clin Endocrinol Metab*. 2011; 25(4): 681–91.
- Holick MF. The role of vitamin D for bone health and fracture prevention. *Curr Osteoporosis Rep*. 2006; 4(3): 96–102.
- Özkan B, Hatun, Ş, Bereket A. Vitamin D intoxication. *Turk J Pediatr*. 2012; 54: 93–98.
- CLSI. *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. CLSI Document EP28-A3c Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- CLSI. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI document EP05-A2. CLSI, Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- Smith DC, Johnson CS, et al. A Phase I Trial of Calcitriol (1,25 dihydroxycholecalciferol) in patients with Advanced Malignancy. *Clin Cancer Res* June 1999; 5: 1339.
- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic *in vitro* / Dispositivo medico-diagnostico *in vitro* / Producto sanitario para diagnóstico *in vitro*

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Límite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

Prescription Device (US Only) / Verschreibungspflichtiges Medizinprodukt (nur USA) / Dispositif soumis à prescription (États-Unis uniquement) / Dispositivo su prescrizione (solo USA) / Dispositivo de prescripción médica (solo EE. UU.)

2015-11_EFEGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

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All rights reserved.