Dimension® EXL™ integrated chemistry system LOCI® Module

Dimension® EXL™ 200, Dimension® EXL™ with LM Systems

Flex® reagent cartridge

FOLA

See shaded sections: Updated information from 2019-04 version.

Issue Date 2020-08

LOCI Folate

Intended Use: The FOLA method is an *in vitro* diagnostic test for the quantitative measurement of folate in human serum and plasma on the Dimension® EXL™ integrated chemistry system with LOCI® module.

Measurements of folate may be used in the diagnosis and treatment of folate deficiency.

Summary: The terms folate and folic acid are often used interchangeably for a water-soluble B complex vitamin (B9) found in dark green leafy vegetables, fruits, dairy products and cereal. The main function of folate coenzymes in the body is the transfer of one-carbon units in a variety of reactions critical to the synthesis of DNA, RNA and amino acids. In nucleic acid metabolism, folic acid is involved in synthesis of DNA from its precursors and in the synthesis of methionine, which is required for the synthesis of S-adenosylmethionine (SAM), a methyl donor in many biological methylation reactions in DNA and RNA. Folic acid is also involved in amino acid metabolism. The synthesis of methionine from homocysteine requires both folate and vitamin B12 dependent enzymes. A folate deficiency can result in decreased synthesis of methionine and a buildup of homocysteine.^{1,2} Clinical studies indicate that a mild to moderate increase in homocysteine is associated with atherosclerotic vascular disease such as coronary artery disease and stroke.^{3,4,5}

Macrocytic anemia is the major clinical manifestation of folate deficiency. It is characterized by abnormal maturation of red blood cell precursors in the bone marrow, the presence of megaloblasts and decreased red blood cell survival. Both folate and vitamin B12 deficiency can cause macrocytic anemia. Folate supplementation can mask B12 deficiency because the associated anemia responds to folate alone. Misdiagnosis delays treatment of the deficiency allowing irreversible neurological abnormalities to progress. Appropriate treatment depends on the differential diagnosis of the deficiency.^{1,2}

The main causes of folate deficiency are absence of intestinal microorganisms, poor intestinal absorption (surgical resection, celiac disease), increased demands (pregnancy, liver disease, and malignancies), insufficient dietary uptake (alcoholism), anti-folate drugs (methotrexate) and anticonvulsants (carbamazepine, phenobarbital, phenytoin, valproic acid). Although serum folate measurement provides an early index of folate status, red blood cell folate more closely reflects tissue stores and is considered the most reliable indicator of folate status.^{1,6}

Principles of Procedure: The Folate method is a homogeneous, competitive chemiluminescent immunoassay based on LOCI® technology. LOCI® reagents include two synthetic bead reagents and labeled folate binding protein (FBP). The first bead reagent (Chemibeads) is coated with a folic acid derivative and contains a chemiluminescent dye. The second bead reagent (Sensibeads) is coated with streptavidin and contains photosensitive dye. Before the immunological portion of the reaction is initiated, the patient sample is pretreated with Sodium Hydroxide (NaOH) and Dithioerythritol (DTE) to release serum folate from endogenous folate binding protein (FBP) and to maintain 5-methyltetrahydrofolate in its reduced form. After the sample pretreatment, chemibeads and labeled folate binding reagent are added sequentially to the reaction vessel. Folate from the patient sample competes with the folate-chemibead for a limited amount of labeled FBP. Sensibeads are then added and bind to the biotinylated portion of the labeled FBP to form bead pair immunocomplexes. Illumination of the complex by light at 680 nm generates singlet oxygen from the Sensibeads which diffuses to the Chemibeads triggering a chemiluminescent reaction. The resulting signal is measured at 612 nm and is an inverse function of the concentration of folate in the sample.^{7,8}

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 2	Empty			
3	Liquid	Folate Chemibeads	400 µg/mL	
4	Liquid	Streptavidin Sensibeads	500 µg/mL	Recombinant <i>E. coli</i>
5	Tablet ^c	Dithioerythritol (DTE)	39 mg/mL	
6	Tablet ^c	Dithioerythritol (DTE)	39 mg/mL	
7	Liquid	Sodium Hydroxide (NaOH)	0.5 N	
8	Liquid	Folate Binding Protein (FBP)	0.225 µg/mL	bovine
		Biotinylated Antibody	1.8 µg/mL	Mouse, monoclonal

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Tablets contain excipients, buffers, and stabilizers.

Risk and Safety:

H290, H314, H317
P280, P301 + P310 + P331,
P303 + P361 + P353 + P310, P305 + P310, P501

Warning!

May be corrosive to metals. Causes severe skin burns and eye damage. May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor/physician. IF IN EYES: Immediately call a POISON CENTER or doctor/physician. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone; Sodium hydroxide

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Used HM reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion. Reaction vessels are designed for single use only.

For *in vitro* diagnostic use.

Reagent Preparation: Hydration, diluting and mixing are automatically performed by the Dimension® system.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open well stability: 5 days for wells 3 – 8

Specimen Collection and Handling: Recommended specimen types: serum and plasma (lithium or sodium heparin).

Samples and controls stabilized with azide cannot be used.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁹

Samples may be stored refrigerated at 2 – 8 °C for up to 8 hours. If testing is delayed beyond 8 hours, samples should be frozen at -20 °C or colder. Protect samples from light. **Avoid using hemolyzed samples.**¹⁰ Mix thoroughly after thawing. Avoid repeated freezing and thawing.¹¹

Follow the instructions provided with your specimen collection device for use and processing.¹² Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection. Specimens should be free of particulate matter.¹¹

The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Procedure

Materials Provided

FOLA Flex® reagent cartridge, Cat. No. RF644

Materials Required But Not Provided

LOCI Anemia Calibrator, Cat. No. RC640
HM Reaction Vessels, Cat. No. RXV1A
Reagent Probe Cleaner, Cat. No. RD702
Sample Probe Cleaner, Cat. No. RD703
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension® EXL™ System. For details of this processing, refer to your Dimension® EXL™ Operator's Guide.

Test Conditions

Sample Volume (delivered to the vessel)	10 µL
DTE Reagent Volume	10 µL
Extractant Reagent Volume	10 µL
FBP/Bio-Mab Reagent Volume	30 µL
Chemibead Reagent Volume	20 µL
Sensibead Reagent Volume	30 µL
Temperature	37 °C
Reaction time	21 minutes
Illumination	680 nm, emission 612 nm
Type of measurement	Chemiluminescence

Calibration

Assay Range	0.5 – 20.0 ng/mL [1.1 – 45.3 nmol/L] ^d
Calibration Material	LOCI Anemia Calibrator, Cat. No. RC640
Calibration Scheme	5 levels, n = 3
Units	ng/mL [nmol/L] (ng/mL x 2.266) = [nmol/L]
Typical Calibration Levels	Level 1: 0 ng/mL [0 nmol/L] Level 2: 2.5 ng/mL [5.7 nmol/L] Level 3: 5.0 ng/mL [11.3 nmol/L] Level 4: 10.0 ng/mL [22.7 nmol/L] Level 5: 21.0 ng/mL [47.6 nmol/L]
Calibration Frequency	Every 21 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations

d. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known folate concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates the concentration of folate in ng/mL [nmol/L] using the calculation scheme described in the Dimension® EXL™ Operator’s Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.5 – 20.0 ng/mL [1.1 – 45.3 nmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 20.0 ng/mL [45.3 nmol/L] can be reported as > 20.0 ng/mL [45.3 nmol/L] or repeated on dilution.
Manual Dilution: Dilute with Reagent grade water to obtain results within reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.
- Autodilution (AD): The autodilute sample volume is 2 µL (dilution factor = 5) for serum/plasma. Refer to your Dimension® EXL™ Operator’s Guide.
- Samples with results less than 0.5 ng/mL [1.1 nmol/L] should be reported as “less than 0.5 ng/mL [1.1 nmol/L]”.

Limitations of Procedure

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{13,14}

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in folate results. Refer to your Dimension® EXL™ Operator’s Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported.

- **Serum folate measurements should not be made on hemolyzed specimens.**
- Methotrexate and Leucovorin (folinic acid) interfere with the measurement of folate. These chemotherapeutic drugs cross-react with folate binding proteins in the folate assay. Therefore, patients receiving these drugs should not be tested for folate by this method.

Interfering substances

The FOLA method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁵ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

- Cimetidine at 2.0 mg/dL [79.2 µmol/L] increases folate results by 14% at 5.2 ng/mL [11.8 nmol/L] and decreases folate results by 14% at 16.8 ng/mL [38.1 nmol/L].

Analyte Concentration	Biotin Test Level (ng/mL)			
	100	250	500	1200
	% Bias			
5.2 ng/mL	2.8	5.9	35.1	Above AMR
16.8 ng/mL	-1.4	3.6	24.8	Above AMR

- Specimens that contain biotin at a concentration of 250 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely elevated results for patient samples.
- The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.¹⁶ Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.¹⁷

Expected Values: 8.6 – 58.9 ng/mL [19.5 – 133.5 nmol/L]

The reference interval represents the central 95% of results determined non-parametrically from a population of 120 apparently healthy adults (60 males and 60 females, 24 – 67 years of age) and was performed in accordance with CLSI/NCCLS C28-A2.¹⁸ All specimens tested were from US adults, a population typically supplemented with dietary folate.

Because folic acid levels are strongly influenced by diet and dietary supplementation, population-based reference intervals can show marked demographic differences. In the US, for example, fortification of enriched grain products, as required by the Food and Drug Administration (FDA) since the mid-1990s, has led to an estimated doubling of the mean plasma folate level among subjects not using vitamin supplements and a decrease in the prevalence of low folate levels, i.e. levels below 3 ng/mL (7 nmol/L).¹⁹

Each laboratory should establish its own expected values for FOLA as performed on the Dimension® EXL™ system.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n = 5 replicates at the following analyte concentrations are:

FOLA concentration	Acceptable SD Maximum
2.5 ng/mL [5.7 nmol/L]	0.5 ng/mL [1.0 nmol/L]
10.0 ng/mL [22.7 nmol/L]	1.2 ng/mL [2.8 nmol/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® EXL™ integrated chemistry system.

Precision ^{20,e}			
Material	Mean	Standard Deviation (% CV)	
	ng/mL [nmol/L]	Repeatability	Within-Lab
Bio-Rad Liquichek™ Immunoassay Control			
Level 1	2.1 [4.8]	0.09 [0.20] (4.3)	0.16 [0.36] (7.6)
Level 2	6.6 [15.0]	0.27 [0.61] (4.1)	0.36 [0.82] (5.5)
Level 3	9.6 [21.8]	0.30 [0.67] (3.1)	0.49 [1.11] (5.1)
Serum Pool 1	2.3 [5.2]	0.11 [0.25] (4.8)	0.15 [0.34] (6.5)
Serum Pool 2	16.7 [37.8]	0.36 [0.82] (2.2)	0.66 [1.50] (4.0)

- e. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.
Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Method Comparison ²¹ Regression Statistics ^f			
Comparative Method	Slope	Intercept ng/mL [nmol/L]	Correlation Coefficient
Dimension Vista® FOL	1.01	0.05 [0.11]	0.99

- f. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.
- g. The range of folate values in the method comparison study was 0.6 – 19.2 ng/mL [1.4 – 43.5 nmol/L].

Specificity

Hemolysis, Icterus, Lipemia (HIL) interference

The FOLA method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁵ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Folate ng/mL [nmol/L]	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer)	See Limitations	
Bilirubin (unconjugated)	20 mg/dL [342 µmol/L]	5.2 [11.8] 16.8 [38.1]	< 10 < 10
Bilirubin (conjugated)	20 mg/dL [342 µmol/L]	5.2 [11.8] 16.8 [38.1]	< 10 < 10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	5.2 [11.8] 16.8 [38.1]	< 10 < 10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* **Analyte results should not be corrected based on this bias.**

Non Interfering Substances

The following substances do not interfere with the FOLA method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at folate concentrations of 5.2 ng/mL [11.8 nmol/L] and 16.8 ng/mL [38.1 nmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Amikacin	8.0 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5.0 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Creatinine	30 mg/dL	2.65 mmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	3000 mg/dL	33.9 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vancomycin	10 mg/dL	69 µmol/L

Recovery

Serum and plasma samples containing folate levels spanning the analytical measurement range were tested individually and as a 1:1 mixture. The observed values were then compared to the expected values for each mixture.

Sample Type	FOLA ng/mL [nmol/L]	FOLA ng/mL [nmol/L]	1:1 Mixture		FOLA % Recovery
			Expected ng/mL [nmol/L]	Obtained ng/mL [nmol/L]	
Serum	3.0 [6.8]	16.4 [37.2]	9.7 [22.0]	10.1 [22.9]	104%
Serum	2.2 [5.0]	12.5 [28.3]	7.4 [16.8]	6.8 [15.4]	92%
Serum	4.7 [10.7]	14.6 [33.1]	9.7 [22.0]	9.6 [21.8]	99%
Serum	3.2 [7.3]	10.4 [23.6]	6.8 [15.4]	6.7 [15.2]	99%
Plasma	3.2 [7.3]	15.2 [34.4]	9.2 [20.8]	9.1 [20.6]	99%
Plasma	2.9 [6.6]	13.4 [30.4]	8.2 [18.5]	7.9 [17.9]	97%
Plasma	6.8 [15.4]	16.8 [38.1]	11.8 [26.7]	12.0 [27.2]	102%
Mean Recovery:					99%

Limit of Detection: 0.5 ng/mL [1.1 nmol/L]^h

The Limit of Detection (LoD) for FOLA is 0.5 ng/mL [1.1 nmol/L], determined consistent with CLSI guideline EP17-A²² and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 120 determinations, with 5 blank and 5 low level samples.

The Limit of Blank (LoB) is 0.2 ng/mL [0.5 nmol/L].

h. LoD is the lowest concentration of analyte that can be detected reliably. LoB is the highest concentration that is likely to be observed in a blank sample.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® EXL™, Dimension Vista®, LOCI® and Flex® are trademarks of Siemens Healthcare Diagnostics.

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© 2014 Siemens Healthcare Diagnostics

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía:

- Shenkin A, Baines M, Fell G S et. al. Vitamins and Trace Elements. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Fourth Edition, 2006; 1100-1105.
- Tietz Clinical Guide to Laboratory Tests. Fourth Edition 2006; 1124-1127.
- Snow CF. *Laboratory Diagnosis of Vitamin B12 and Folate Deficiency*. Arch Intern Med. 1999; 159:1289-1298.
- Ward PCJ. *Modern Approaches to the Investigation of Vitamin B12 Deficiency*. Clin Lab Med. 2002; 22:435-445.
- Klee GG. *Cobalamin and Folate Evaluation: Measurement of Methylmalonic Acid and Homocysteine vs Vitamin B12 and Folate*. Clin Chem. 2000; 46:1277-1283.
- Ermens AAM, Vlasveld LT, Lindemans J. *Significance of Elevated Cobalamin (Vitamin B12) Levels in Blood*. Clin Biochem. 2003; 36:585-590.
- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al., *Luminescent oxygen channeling assay (LOC®): sensitive, broadly applicable homogeneous immunoassay method*. Clin Chem. 1996; 42:9:1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al. *Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence*. Proc.Natl.Acad.Sci. 2004; 91:5426-5430.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Tietz Clinical Guide to Laboratory Tests. Fourth Edition. Alan H. B. Wu editor, 2006; 410-413, 1124-1127.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Kricka LJ. *Human anti-animal antibody interference in immunological assays*. Clin Chem. 1999; 45:7:942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. Clin Chem 1992; 38: 1737-1742.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4): 247–256.
- Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 2017;55(6):817-825.
- Clinical and Laboratory Standards Institute/NCCLS. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition*. CLSI/NCCLS document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Pfeiffer C, Johnson C, Jain R, et al. Trends in blood folate and vitamin B-12 concentrations in the United States 1988-2004. *American Journal of Clinical Nutrition*, 2007; 86; 718-727
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542 9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical Laboratory Standards Institute/NCCLS. *Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contenuto sufficiente per <n> saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFBS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

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All rights reserved.