Dimension® EXL™ integrated chemistry system

LOCI® Module

Flex® reagent cartridge

NTP

See shaded sections: Updated information from 2019-05 version.

Issue Date 2020-08

N-terminal Pro-Brain Natriuretic Peptide

Intended Use: The NTP method is an *in vitro* diagnostic assay for the quantitative measurement of N-terminal pro-brain natriuretic peptide (NT-proBNP) in human serum and plasma on the Dimension® EXL™ integrated chemistry system with LOCI® Module. In individuals suspected of having congestive heart failure (CHF), measurements of NT-proBNP are used as an aid in the diagnosis and assessment of severity. The test is further indicated for the risk stratification of patients with acute coronary syndrome and heart failure.

Summary: Left ventricular dysfunction can occur as part of coronary heart disease, arterial hypertension, valvular disease, and primary myocardial disease. If the left ventricular dysfunction remains untreated and is progressive, the potential for mortality is high, e.g., due to sudden cardiac death.

Chronic cardiac insufficiency is a clinical syndrome caused by impairment of the cardiac pumping function. Based on the symptoms, the severity of cardiac insufficiency is classified in stages I-IV by the New York Heart Association (NYHA). Clinical tests and imaging procedures are used to diagnose left ventricular dysfunction.¹

The significance of natriuretic peptides in the control of cardiovascular system function has been established. Initial studies reveal that natriuretic peptides can be used for diagnosing clinical problems associated with left ventricular dysfunction.² The following natriuretic peptides have been reported in the literature: atrial natriuretic peptide (ANP), brain natriuretic peptide (BNP), and C-type natriuretic peptide (CNP).^{3,4}

ANP and BNP have natriuretic and diuretic properties. As antagonists of the renin-angiotensin-aldosterone system, they influence the electrolyte and fluid balance in an organism.^{5,6} In subjects with left ventricular dysfunction, serum and plasma concentrations of BNP increase, as do the concentrations of the biologically inactive prohormone, proBNP. ProBNP consists of 108 amino acids. It is secreted mainly by the left ventricle of the heart and, in this process, is cleaved into physiologically active BNP (77 – 108), and the N-terminal fragment NT-proBNP (1 – 76).⁴

Studies indicate that NT-proBNP can be used in diagnostic and prognostic applications.⁷⁻¹³ The concentration of NT-proBNP in plasma indicates the prognosis for the left ventricular dysfunction. It is also useful in assigning symptoms to cardiac or non-cardiac causes.

Determination of NT-proBNP helps to identify subjects with left ventricular dysfunction. Changes in NT-proBNP concentration can be used to evaluate the success of treatment in patients with left ventricular dysfunction. There are indications that NT-proBNP, due to its functions, is suitable for use in assessing vascular remodeling, and therefore contributes to the establishment of individualized rehabilitation procedures.^{10,11}

Fisher et al. found that congestive heart failure patients with elevated NT-proBNP values had a one year mortality rate of 53% compared to 11% in patients with lower values.¹² In the Gusto IV study which involved more than 6800 patients it was shown that NT-proBNP was the strongest independent predictor of one year mortality in patients with acute coronary syndrome.¹³

The European Society of Cardiology Task Force for the Diagnosis and Treatment of Chronic Heart Failure recommends in their guidelines that natriuretic peptides including NT-proBNP "may be most useful clinically as a rule out test due to consistent and very high negative predictive values."

NT-proBNP levels are increased in patients with unstable angina and following myocardial infarction.^{14,15} Although not diagnostic for these conditions, several studies indicate NT-proBNP measurements provide prognostic information for the short- and long-term risk stratification of patients with unstable angina or myocardial infarction.^{15,16,17}

Principles of Procedure: The NTP method is a one-step sandwich chemiluminescent immunoassay based on LOCI® technology. LOCI® reagents include two synthetic bead reagents and a biotinylated monoclonal antibody fragment which recognize an epitope located in the N-terminal part of proBNP. The first bead reagent (Sensibeads) is coated with streptavidin and contains photosensitive dye. The second bead reagent (Chemibeads) is coated with a second antibody specific for a second independent epitope on NT-proBNP and contains chemiluminescent dye. Sample is incubated with Chemibeads and biotinylated antibody to form a particle/NT-proBNP/biotinylated antibody sandwich. Sensibeads then are added and bind to the biotin to form a bead-aggregated immunocomplex. Illumination of the complex by light at 680 nm generates singlet oxygen from Sensibeads, which diffuses to the Chemibeads and triggers a chemiluminescent reaction. The resulting chemiluminescent signal is measured at 612 nm and is directly proportional to the concentration of NT-proBNP in the sample.^{18,19,20}

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 2	Liquid	Biotinylated antibody	8 µg/mL	Sheep, monoclonal
3 – 4	Liquid	NT-proBNP Chemibeads	150 µg/mL	Sheep, monoclonal
5 – 6	Liquid	Streptavidin Sensibeads	1400 µg/mL	Recombinant <i>E. coli</i>
7 – 8	Liquid	Assay buffer		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!
May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Used HM reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open well stability: 3 days for wells 1 – 8

Specimen Collection and Handling: Recommended specimen types: Serum or plasma samples (lithium heparin, sodium heparin, and EDTA).

Samples and controls stabilized with azide cannot be used.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.^{21,22,23} Each laboratory must determine its appropriate sample timing/collection protocol.

Follow the instructions provided with your specimen collection device for use and processing.²³

Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.²⁴ Specimens should be free of particulate matter.

If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of interference by particulate matter.^{21,24}

Samples are stable for three days when stored at 20 – 25 °C or at 2 – 8 °C or for 12 months when frozen at -20 °C or below.²⁵ Repetitive freezing and thawing of specimens should be avoided.

Ensure that patient samples, calibrators, and controls are equilibrated at ambient temperature (22 – 28 °C) before testing. Samples containing precipitates must be centrifuged before performing the assay.

Procedure

Materials Provided

NTP Flex® reagent cartridge, Cat. No. RF623

Materials Required But Not Provided

LOCI NTP CAL, Cat. No. RC623
Dimension® Sample Diluent, Cat. No. 791092901
HM reaction vessels Cat. No. RXV1A
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension® EXL™ System. For details of this processing, refer to your Dimension® EXL™ Operator's Guide.

Test Conditions

Sample Volume (delivered to the HM reaction vessel)	8 µL
Chemibead Reagent Volume	20 µL
Biotinylated Antibody Reagent Volume	20 µL
Sensibead Volume	15 µL
Assay Buffer Volume	100 µL
Temperature	37.0 °C
Reaction Time	14 minutes
Wavelength	680 and 612 nm
Type of Measurement	Chemiluminescence

Calibration	
Assay Range	5 – 35000 pg/mL
Calibration Material	NTP Calibrator, Cat. No. RC623
Calibration Scheme	5 levels, n = 3
Units	pg/mL
Typical Calibration Levels	Level 1: 0 pg/mL Level 2: 250 pg/mL Level 3: 1500 pg/mL Level 4: 12000 pg/mL Level 5: 36750 pg/mL
Calibration Frequency	Every 30 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known NT-proBNP concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of NT-proBNP in pg/mL using the calculation scheme described in your Dimension® EXL™ Operator’s Guide.

Results of this test should always be interpreted in conjunction with the patient’s medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5 – 35000 pg/mL

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 35000 pg/mL can be reported as > 35000 pg/mL or can be manually diluted.
Manual Dilution: Dilute with Sample Diluent (Cat. No. 791092901) to obtain results within the reportable range. The recommended dilution factor is 2. Enter dilution factor on instrument. Reassay. Resulting Readout is corrected for dilution. Dilutions of 1:5 have demonstrated maximum deviations of 20% from the theoretical value.
- Samples with results less than 5 pg/mL should be reported as less than 5 pg/mL.

Limitations of Procedure

Concentrations of natriuretic peptides may be elevated in patients with acute myocardial infarction, patients that are candidates for renal dialysis, and patients that have undergone renal dialysis.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed.

A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{26, 27}

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in NT-proBNP results. Refer to your Dimension® EXL™ Operator’s Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported.

Testing specimens from renal dysfunction patients taking biotin may lead to false negative results. Therefore, do not use this device in patients with renal impairment (eGFR <60), unless it is confirmed that the patient is not taking biotin. Patients taking more than 20 mg/day of biotin may have falsely negative results. Do not use this test in patients taking more than 20 mg/day of biotin.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n = 5 replicates at the following analyte concentrations are:

NT-proBNP concentration	Acceptable SD Maximum
152 pg/mL	13.6 pg/mL
446 pg/mL	44.7 pg/mL

A system malfunction may exist if the acceptable SD maximum is exceeded.

Expected Values

NT-proBNP concentrations in the Reference Group are shown in the following tables. Each laboratory should establish its own reference interval for the NTP method as determined on the Dimension® EXL™ System. The recommended medical decision thresholds, by age group, are:

Patients < 75 years:	125 pg/mL
Patients ≥ 75 years:	450 pg/mL

Reference Study Group†

NT-proBNP concentrations were determined in 318 individuals without congestive heart failure (163 women and 155 men) using the Dimension Vista® PBNP method. This population included apparently healthy individuals and individuals with diabetes, hypertension, and pulmonary disease. The statistics for NT-proBNP concentrations in the reference study group are shown in the following table.

Reference Study Group

All	< 55 yrs	55–64 yrs	65–74 yrs	≥ 75 yrs
Mean	45.9	54.6	175.5	451.3
SD	53.6	62.3	197.6	987.1
Median	27.2	28.8	106.2	173.5
95th Percentile	131.2	206.6	880.8	1615.7
% < 125 pg/mL	94%	80%	61%	-
% < 450 pg/mL	-	-	-	79%
N	163	15	18	122

Males	< 55 yrs	55–64 yrs	65–74 yrs	≥ 75 yrs
Mean	39.0	56.7	219.4	520.8
SD	57.8	58.2	242.4	1080.9
Median	16.6	25.4	123.6	173.5
95th Percentile	168.8	154.7	880.8	1969.7
% < 125 pg/mL	93%	83%	55%	-
% < 450 pg/mL	-	-	-	73%
N	76	6	11	62

Females	< 55 yrs	55–64 yrs	65–74 yrs	≥ 75 yrs
Mean	51.9	53.2	106.6	379.5
SD	49.2	68.3	60.8	883.1
Median	39.4	28.8	87.7	167.3
95th Percentile	124.0	206.6	215.2	1453.2
% < 125 pg/mL	95%	78%	71%	-
% < 450 pg/mL	-	-	-	85%
N	87	9	7	60

Disease Study Group†

Blood samples were obtained from 269 patients diagnosed with congestive heart failure (CHF). NT-proBNP concentrations were determined using the Dimension Vista® PBNP. The population included 78 women and 191 men. The descriptive statistics and New York Heart Association (NYHA) functional classes are provided below. Each laboratory should establish a reference range that represents the patient population that is to be evaluated. In addition, laboratories should be aware of their institution’s current practice for the evaluation of CHF.

CHF Population – All	< 55 yrs	55–64 yrs	65–74 yrs	≥ 75 yrs
Mean	6131.8	4455.5	6168.1	8142.7
SD	13470.3	9296.0	10950.5	12699.1
Median	1550.9	1639.5	2541.3	3780.5
95th Percentile	28099.6	19481.5	17768.6	26923.1
% > 125 pg/mL	89%	92%	95%	-
% > 450 pg/mL	-	-	-	94%
N	72	72	59	66

CHF Population – Males	< 55 yrs	55–64 yrs	65–74 yrs	≥ 75 yrs
Mean	7771.4	4954	6537.2	8827.3
SD	15797.2	10536.7	12383.9	12856.6
Median	2526.1	1699.8	2541.3	5529.2
95th Percentile	31237.1	19481.5	17768.6	18634.8
% > 125 pg/mL	94%	94%	93%	-
% > 450 pg/mL	-	-	-	96%
N	48	51	43	49

CHF Population – Females	< 55 yrs	55–64 yrs	65–74 yrs	≥ 75 yrs
Mean	2852.6	3244.9	5176.1	6169.7
SD	5774.0	5200.7	5727.8	12397.4
Median	452.3	1516.9	2270.8	1125.3
95th Percentile	15569.3	7903.1	16889.9	45824.4
% > 125 pg/mL	79%	86%	100%	-
% > 450 pg/mL	-	-	-	88%
N	24	21	16	17

CHF Population – All	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	6184.5	3941.9	3648.8	6643.1	11089.2
SD	11737.8	8461.3	4897.0	13567.0	14428.2
Median	2282.2	1161.0	1091.7	2541.3	4537.7
5th Percentile	101.1	100.4	46.6	143.1	74.6
95th Percentile	22331.9	14184.1	14238.5	17497.1	28099.6
% > Cutoff	92%	93%	90%	94%	88%
Minimum	5.1	52.8	24.5	5.1	39.6
Maximum	91446.4	52327.0	18634.8	91446.4	78045.2
N	269	45	60	124	40

CHF Population – Males	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	7012.1	4567.5	4402.2	7298.4	14351.2
SD	12991.9	9204.3	5167.4	14789.5	17258.9
Median	3162.6	1639.5	2230.9	3442.2	6015.9
5th Percentile	126.1	136.3	46.5	126.1	101.1
95th Percentile	24004.7	19481.5	14238.5	17497.1	31237.1
% > Cutoff	94%	97%	93%	95%	91%
Minimum	5.1	72.5	44.7	5.1	74.6
Maximum	91446.4	52327.0	18634.8	91446.4	78045.2
N	191	37	40	91	23


CHF Population – Females

NYHA Functional Class					
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	4157.8	1048.3	2142.1	4836.0	6676.0
SD	7561.7	1636.1	4008.9	9366.5	7864.9
Median	950.2	225.4	705.3	1286.7	3780.5
5th Percentile	84.3	52.8	24.5	143.1	39.6
95th Percentile	22331.9	4960.2	11020.8	26923.1	22997.0
% > Cutoff	87%	75%	85%	94%	82%
Minimum	24.5	52.8	24.5	126.8	39.6
Maximum	45824.4	4960.2	15893.8	45824.4	22997.0
N	78	8	20	33	17


These results show that there is a relationship between the severity of the clinical signs and symptoms of CHF and the median NT-proBNP concentration, demonstrating that the Dimension Vista® PBNP method can be used as an aid in the diagnosis of all degrees of CHF severity including asymptomatic patients.

Interpretation of Results‡

The Receiver Operator Characteristics (ROC) Curve presents the clinical sensitivity and specificity at various cutoffs for the 269 patients diagnosed with CHF and 318 subjects without CHF. The ROC curve for the Dimension Vista® assay is shown below. The area under ROC curve (AUC) for the Dimension Vista® PBNP assay is 0.928 with a 95% confidence interval of 0.907 to 0.949.

A box and whiskers plot of the clinical study population is presented below. Recommended clinical thresholds are 125 pg/mL for patients younger than 75 years and 450 pg/mL for patients 75 years and older. Six disease group samples with values above the assay range are not displayed in the plot.

Sensitivity and Specificity vs. Age and Gender‡

The tables below show the clinical sensitivity and specificity of the Dimension Vista® PBNP assay using a cutoff of 125 pg/mL for patients younger than 75 years and 450 pg/mL for patients 75 years or older.

Males		
	< 75 yrs	≥ 75 yrs
Sensitivity	94% (133/142)	96% (47/49)
95% Confidence Interval	90 – 98	90 – 100
Specificity	87% (81/93)	73% (45/62)
95% Confidence Interval	80 – 94	61 – 84

Females		
	< 75 yrs	≥ 75 yrs
Sensitivity	89% (54/61)	88% (15/17)
95% Confidence Interval	81 – 97	73 – 100
Specificity	92% (95/103)	85% (51/60)
95% Confidence Interval	87 – 97	76 – 94

Age-matched ROC analysis‡

An age-matched ROC analysis of the clinical data was performed via the weighted method described in Kondratovich (2002).²⁸ The resulting AUC is 0.929 with a 95% confidence interval of 0.909 to 0.950.

‡ The Dimension Vista® PBNP method (REF K6423A) and the Dimension® NTP (REF RF623) method utilize the same reagents under equivalent reaction conditions and the results are representative of both methods.

Specific Performance Characteristics

The following data represent typical performance for the Dimension EXL™ System.

Material	Mean pg/mL	Precision ^{29, c}	
		Standard Deviation (%CV)	
		Repeatability	Within-Lab
Serum Pool 1	148	2.7 (1.9)	3.8 (2.5)
Serum Pool 2	434	6.9 (1.6)	10.2 (2.4)
Serum Pool 3	946	17.4 (1.8)	20.9 (2.2)
Serum Pool 4	4315	96.6 (2.2)	144.2 (3.3)
Internal QC	30643	1732.5 (5.7)	2761.3 (9.01)

c. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Method Comparison³⁰

Serum and Plasma Equivalency^d

Comparison testing results on the Dimension® EXL™ System for lithium heparin, sodium heparin and EDTA plasma samples versus serum samples is provided as follows:

Comparative Specimen	Slope	Intercept pg/mL	Correlation Coefficient	n
Li Heparin	0.99	-25.9	0.9932	50
Na Heparin	0.98	-43.3	0.9919	53
EDTA	0.96	-0.4	0.9988	57

d. Model equation by Passing-Bablok regression is: [plasma results] = slope x [serum results] + intercept.

Regression Statistics^{e, f}

Comparative Method	Slope	Intercept pg/mL	Correlation Coefficient	n
Dimension Vista® PBNP (monoclonal)	0.950	-1.303	0.998	109 ^g

e. CLSI/NCCLS EP9-A2 was used. The method used to fit the regression line was Passing-Bablok Statistics.
f. Model equation by Passing-Bablok linear regression is: [results for Dimension® EXL™ system] = slope x [comparative method results] + intercept.
g. The range of 109 values in the correlation study was 36 – 33493 pg/mL.

Specificity

HIL interference

The NTP method was evaluated for interference to hemoglobin, icterus and lipemia (HIL) according to CLSI/NCCLS EP7-A2.³¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	NT-proBNP pg/mL	Bias* %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] ^h	120	< 10
		1400	< 10
Bilirubin (unconjugated)	60 mg/dL [1026 µmol/L]	120	< 10
		1400	< 10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	120	< 10
		1400	< 10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	120	< 10
		1400	< 10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

*Analyte results should not be corrected based on this bias.

h. Système International d'Unités [SI units] are in brackets.

Interfering Substances

The NTP method was evaluated for interference according to CLSI/NCCLS EP7-A2.³¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Analyte Concentration	Biotin Test Level (ng/mL)			
	100	250	500	1200
	% Bias			
120 pg/mL	-1.6	-4	-16.3	-74.8
1400 pg/mL	-2.5	-6.4	-17.1	-71.9

* Specimens that contain biotin at a concentration of 250 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely depressed results for patient samples.

* The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.³² Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.³³

Non Interfering Substances

The following substances do not interfere with the NTP method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at NT-proBNP concentration of 120 and 1400 pg/mL.

Compound	Concentration	SI Units
Abciximab	21 µg/mL	0.4 µmol/L
Acetaminophen	20 mg/dL	1324 µmol/L
Allopurinol	2.5 mg/dL	184 µmol/L
Amikacin	8 mg/dL	137 µmol/L
Amiodarone	20 µg/mL	31 µmol/L
Amiodipine Besylate	4 µg/mL	7.1 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Atenolol	1.0 mg/dL	37.6 µmol/L
Atorvastatin	32 mg/dL	572.9 µmol/L
Caffeine	6 mg/dL	308 µmol/L

Captopril	5.0 mg/dL	230 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Cinnarizine	3.0 mg/dL	81.4 µmol/L
Clopidogrel bisulfate	30 µg/mL	71.4 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Cyclosporine A	4000 ng/mL	3.3 µmol/L
Dextran 40	3000 mg/dL	750 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Diclofenac	60 µg/mL	188.6 µmol/L
Digitoxin	60 ng/mL	78.6 nmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Diltiazem	120 µg/mL	266.1 µmol/L
Dipyridamole	30 µg/mL	59.4 µmol/L
Dopamine	65 mg/dL	3.4 mmol/L
Enalapril maleate	16 µg/mL	32.5 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Hydralazine	20 µg/mL	101.7 µmol/L
Hydrochlorothiazide	20 µg/mL	67.2 µmol/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Indomethacin	16 µg/mL	36.9 µmol/L
Isosorbide dinitrate	6.0 mg/dL	254 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lisinopril	16 µg/mL	36.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Lovastatin	16 µg/mL	39.6 µmol/L
L-thyroxine	60 µg/dL	0.77 µmol/L
Methyldopa	2.5 mg/dL	118 µmol/L
Milrinone lactate	2.4 µg/mL	8 µmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nifedipine	6.0 mg/dL	173 µmol/L
Nitrofurantoin	40 µg/mL	168 µmol/L
Nitroglycerin	0.16 µg/mL	0.577 µmol/L
Oxazepam	12 µg/mL	41.8 µmol/L
Oxytetracycline	100 µg/mL	217.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Probenecid	200 µg/mL	700.9 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Propranolol	0.15 mg/dL	5.1 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12.0 g/dL	120 g/L
Quinidine	20 µg/mL	61.6 µmol/L
Rheumatoid Factor	500 IU/mL	500 IU/mL
Salicylic Acid	60 mg/dL	4.34 mmol/L
Simvastatin	32 µg/mL	76.5 µmol/L
Sulfamethoxazole	320 µg/mL	1263.3 µmol/L
Theophylline	4 mg/dL	222 µmol/L
Trimethoprim	64 µg/mL	220.7 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Verapamil	16 mg/dL	351.6 µmol/L
Warfarin	4 ug/dL	13 µmol/L

Hook Effect

One step sandwich immunoassays are susceptible to a high dose "hook effect", where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.³⁴ Such samples must be diluted and reassayed prior to reporting the results. The NTP method shows no hook effect up to 400,000 pg/mL.

Cross-reactivity

The following substances have no significant cross-reactivity (less than 1%) at the concentrations indicated when added to samples containing 0 and approximately 125 pg/mL NT-proBNP:

Substance	Concentration
ANP ²⁸	3.1 µg/mL
NT-proANP ₁₋₃₀ (preproANP ₂₆₋₅₅)	3.5 µg/mL
NT-proANP ₃₁₋₆₇ (preproANP ₅₆₋₉₂)	1.0 ng/mL
NT-proANP ₇₉₋₉₈ (preproANP ₁₀₄₋₁₂₃)	1.0 ng/mL
BNP ₃₂ (Natreacor®)	3.5 µg/mL
CNP ₂₂	2.2 µg/mL
DNP	1.0 ng/mL
VNP	1.0 ng/mL
Adrenomedullin	1.0 ng/mL
Aldosterone	0.6 ng/mL
Angiotensin I	0.6 ng/mL
Angiotensin II	0.6 ng/mL
Angiotensin III	1.0 ng/mL
Endothelin	20 pg/mL
Renin	50 ng/mL
Urodilatin	3.5 µg/mL
Arg-Vasopressin	1.0 ng/mL

Analytical Sensitivity: ≤ 5 pg/mL

The analytical sensitivity of the NTP method is ≤ 5 pg/mL. Analytical sensitivity is defined as the concentration at two standard deviations (n = 20) above a sample devoid of NT-proBNP such as the LOCI NTP Level 1 Calibrator (0 pg/mL).

Limit of Detection³⁵: 2.3 pg/mL

The Limit of Detection (LoD) for the NTP method is 2.3 pg/mL, determined consistent with CLSI guideline EP17-A and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 30 determinations, with 5 blank and 5 low level samples. LOD is the lowest concentration that can be detected reliably.

Functional Sensitivity: Evaluation of the 20% CV limit

Functional sensitivity was evaluated by determining the total imprecision of natural NT-proBNP serum samples. Total imprecision was measured by a two replicate per day, twenty day ANOVA study. The functional sensitivity/ imprecision profile for the Dimension® NTP method corresponds to a coefficient of variation (CV) of 20% at an NT-proBNP level of ≤ 30 pg/mL.

Dilutional Linearity

A patient plasma sample with a high NT-proBNP concentration was diluted across the assay range. Additional plasma samples were diluted linearly to confirm dilutional accuracy at the low end of the assay range. All samples were diluted with low level patient plasma pools. The following table summarizes the results of these studies:

	Sample A	Sample B
Starting Concentration (pg/mL)	47120	1232
Regression Statistics		
Slope	1.004	0.996
Intercept	12.8	10.9
Correlation Coefficient	0.999	0.999
N	5	5
% Recovery		
Average	100.2	101.9
Range	99 – 101	100 – 106

Manufactured under license from Roche Diagnostics.
US 5786163 and foreign equivalents
US Patent Application: 09/890442 and foreign equivalents

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension Vista®, LOCI®, Dimension®, EXL™ and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía

1. Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001; 22: 1527–1560.
2. Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998; 97: 1921–1929.
3. de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230: 767–770.
4. Valli N, Gobinet A, Bordenave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999; 134: 437–444.
5. de Bold AJ, Boerensstein HB, Varesse AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981; 28: 89–94.
6. Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987; 79: 738–745.
7. Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999; 20: 1374–1375.
8. Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47: 287–296.
9. Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999; 20: 1736–1744.
10. Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996; 78: 284–287.
11. McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998; 351: 9–13.
12. Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89: 879–881.
13. James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108: 275–281.
14. Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotrophin 1 are raised in unstable angina. Heart 2000; 84: 421–424.
15. Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002; 40: 437–445.
16. Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004; 102: 37–40.
17. Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003; 107: 2786–2792.
18. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al. Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method. Clin Chem 42: 9 1996, 1518–1526.
19. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et al. Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl.Acad.Sci. USA, Vol 91, pp. 5426–5430, June 1994 Biochemistry.
20. Ullman EF. Homogeneous Immunoassays, The Immunoassay Handbook, 2nd ed., David Wild editor, 2001 192–194.
21. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 41–58 (specimen collection).
22. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
23. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
24. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
25. Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
26. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. Clin Chem 1992; 38: 1737–1742.
27. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45: 7; 942–956.
28. Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.
29. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
30. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
31. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
32. Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4):247–256.
33. Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 2017;55(6):817–825.
34. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays. Anal Biochem 1982; 127: 308–315.
35. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour “n” tests / Contenuto sufficiente per “n” saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic *in vitro* / Dispositivo medico-diagnostico *in vitro* / Producto sanitario para diagnóstico *in vitro*

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_ERGS