Dimension® EXL™ integrated chemistry system LOCI® Module

Flex® reagent cartridge

TSHL

See shaded sections: Updated information from 2019-05 version.

Issue Date 2020-08

LOCI Thyroid Stimulating Hormone

Intended Use: The TSHL method is an *in vitro* diagnostic test for the quantitative measurement of Thyroid Stimulating Hormone (TSH, thyrotropin) in human serum and plasma on the Dimension® EXL™ integrated chemistry system with LOCI® Module. Measurements of TSH are used in the diagnosis and monitoring of thyroid disease.

Summary: Thyroid stimulating hormone is a glycoprotein secreted by the anterior lobe of the pituitary gland. TSH stimulates the normal thyroid gland to synthesize and secrete thyroxine (T4) and triiodothyronine (T3).¹ Although less sensitive measurements of TSH (or free T4) can be used to diagnose severe, clinically apparent hypo- or hyperthyroidism, only a highly sensitive TSH assay has sufficient clinical sensitivity to detect the minor degrees of thyroxine excess or deficiency associated with early, subclinical phases of hypo- or hyperthyroidism.^{1,2} The TSHL assay meets the criteria of a "third generation" assay, defined as a functional sensitivity of ≤ 0.02 mIU/L with a between-run coefficient of variation (CV) of $\leq 20\%$.²

Principles of Procedure: The TSHL method is a homogeneous, sandwich chemiluminescent immunoassay based on LOCI® technology. The LOCI® reagents include two synthetic bead reagents and a biotinylated anti-TSH monoclonal antibody fragment. The first bead reagent (Sensibeads) is coated with streptavidin and contains a photosensitizer dye. The second bead reagent (Chemibeads) is coated with a second anti-TSH monoclonal antibody and contains chemiluminescent dye. Sample is incubated with biotinylated antibody and Chemibeads to form bead-TSH-biotinylated antibody sandwiches. Sensibeads are added and bind to the biotin to form bead-pair immunocomplexes. Illumination of the complex at 680 nm generates singlet oxygen from Sensibeads which diffuses into the Chemibeads, triggering a chemiluminescent reaction. The resulting signal is measured at 612 nm and is a direct function of the TSH concentration in the sample.^{3,4}

Reagents

Wells ^{a,b}	Form	Ingredient	Concentration ^c	Source
1 – 2	Liquid	Biotinylated TSH antibody	7.5 µg/mL	Mouse monoclonal
3 – 4	Liquid	TSH antibody coated Chemibeads	200 µg/mL	Mouse monoclonal
5 – 6	Liquid	Streptavidin Sensibeads	1400 µg/mL	Recombinant <i>E. coli</i>
7 – 8	Empty			

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Wells 1 – 6 contain buffers, stabilizers and preservatives.

c. Nominal value per well in a cartridge.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used HM reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Recommended specimen types: serum, lithium and sodium heparin plasma, EDTA plasma.

Samples and controls stabilized with sodium azide cannot be used.

Collect serum and plasma specimens using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁵ Separated samples are stable for 1 day at ambient temperature and 7 days at 2 – 8 °C. For longer storage, samples should be frozen at -20 °C for 1 month.^{5,7} Avoid repeated freezing and thawing.

Specimens must be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation. If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of particulate matter.⁸

Follow the instructions provided with your specimen collection device for use and processing.⁹

The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Procedure

TSHL Flex® reagent cartridge, Cat. No. RF612

Materials Provided

Materials Required But Not Provided

LOCI Thyroid Calibrator, Cat. No. RC610/RC610A
MULTI 2 Sample Diluent, Cat. No. KD694
TSH Sample Diluent, Cat. No. KD691
HM reaction vessels, Cat. No. RXV1A
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension® EXL™ with LM System. For details of this processing, refer to your Dimension® EXL™ with LM Operator's Guide.

Test Conditions

Sample Volume (delivered to the HM reaction vessel)	12 µL
Biotinylated Antibody Reagent Volume	20 µL
Chemibeads Reagent Volume	10 µL
Streptavidin Sensibeads Reagent Volume	12 µL
Temperature	37.0 °C
Reaction Time	16 minutes
Wavelength	Illumination 680 nm, Emission 612 nm
Type of Measurement	Chemiluminescence

Calibration

Assay Range 0.007 – 100 µIU/mL [mIU/L]^d
Calibration Material LOCI Thyroid Calibrator, Cat. No. RC610/RC610A
Calibration Scheme 5 levels, n = 3
Units µIU/mL [mIU/L]
(µIU/mL x 1) = [mIU/L]

Typical Calibration Levels

Level 2: 0.0 µIU/mL [mIU/L]
Level 3: 4.0 µIU/mL [mIU/L]
Level 4: 20 µIU/mL [mIU/L]
Level 5: 50 µIU/mL [mIU/L]
Level 6: 105 µIU/mL [mIU/L]

Calibration Frequency

A new calibration is required

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

d. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known Thyroid Stimulating Hormone concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of TSH in µIU/mL [mIU/L] using the calculation scheme described in your Dimension® EXL™ with LM Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.007 – 100 µIU/mL [mIU/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 100 µIU/mL [mIU/L] should be repeated on dilution.

Manual Dilution: Dilute with TSH Sample Diluent, Cat. No. KD691 or MULTI 2 Sample Diluent, Cat. No. KD694, to obtain results within reportable range. The recommended dilution factor is 5.

- Samples with results less than 0.007 µIU/mL [mIU/L] should be reported as "less than 0.007 µIU/mL [mIU/L]".

Limitations of Procedure

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. As with any immuno-recognition measurement of a peptide, extremely rare genetic variants may exhibit varying degrees of detection. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{10,11}

Performance of this assay has not been established with neonatal specimens.

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in TSHL results. Refer to your Dimension® EXL™ with LM Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

United States Federal law restricts this device to sale by or on the order of a physician or other practitioner licensed by the laws of the State in which he practices, to use or order the use of the device.

Do not use in patients taking biotin supplements.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n = 5 replicates at the following TSH concentrations are:

TSH Concentration	Acceptable SD Maximum
4.0 µIU/mL [mIU/L]	0.5 µIU/mL [mIU/L]
20.0 µIU/mL [mIU/L]	2.3 µIU/mL [mIU/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Adult population: 0.358 – 3.74 µIU/mL [mIU/L]

The reference interval was transferred from that previously determined for the TSH method on the Dimension Vista® System. It represents the central 95% of results determined non-parametrically from a population of 297 apparently healthy adults (187 males and 110 females, 18 – 65 years of age). The original determination and transference to the TSHL method was done in accordance with CLSI/NCCLS C28-A2.¹²

Pediatric Population:

Pediatric Age Group	Number of Samples	Reference Intervals	
		µIU/mL	[mIU/L]
Infants (01 – 23 months)	75	0.867 – 6.43	0.867 – 6.43
Children (02 – 12 years)	185	0.704 – 4.01	0.704 – 4.01
Adolescents (13 – 20 years)	147	0.516 – 4.13	0.516 – 4.13

Reference intervals for the pediatric population (infants, children, and adolescents) were established for the Dimension® EXL™ TSHL method in accordance with the CLSI guideline EP28-A3c.¹³ Samples were collected prospectively from apparently healthy (euthyroid) pediatric subjects, using predefined inclusion criteria.

A non-parametric approach was used to establish the reference intervals for children and adolescents where the 2.5th and 97.5th percentiles of the distribution of values were calculated. For the infant population, the reference interval was calculated using an approach to accommodate the smaller sample size.¹⁴ The upper limit (97.5th percentile) of the infant reference interval was 6.43 µIU/mL[mIU/L]. Data from this infant population (n=75) have demonstrated a highly skewed distribution to the right; therefore, the estimate of the upper limit of the reference interval has some uncertainty with a 90% probability that the upper limit of the reference interval can be between 5.40–7.54 µIU/mL[mIU/L].

Each laboratory should establish its own expected values for TSH as performed on the Dimension® EXL™ with LM system.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® EXL™ with LM system.

Material	Precision ^{15,a}		
	Mean	Standard Deviation (% CV)	
	µIU/mL [mIU/L]	Repeatability	Within-Lab
Bio-Rad Lyphochek®			
Level 1	0.457	0.015 (3.30)	0.021 (4.58)
Level 2	6.365	0.212 (3.34)	0.323 (5.07)
Level 3	41.750	1.135 (2.72)	1.616 (3.87)
Serum Pool 1	0.260	0.007 (2.60)	0.012 (4.67)
Serum Pool 2	2.107	0.059 (2.82)	0.103 (4.91)
Serum Pool 3	91.785	4.200 (4.58)	5.486 (5.98)

e. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.
Lyphochek® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Comparative Method	Method Comparison ¹⁶ Regression Statistics ^f			
	Slope	Intercept µIU/mL [mIU/L]	Correlation Coefficient	n
	Dimension Vista® TSH	1.05	-0.04	0.997

f. The method used to fit the linear regression line was ordinary least squares.
g. The range of Dimension Vista® TSH values in the correlation study was 0.009 to 97.2 µIU/mL [mIU/L].

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The TSHL method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁷ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	TSH µIU/mL [mIU/L]	Bias* %
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L]	0.3 [0.3]	< 10
Hemoglobin (hemolysate)	750 mg/dL [0.47 mmol/L]	8 [8]	< 10
Bilirubin (unconjugated)	60 mg/dL [1026 µmol/L]	0.3 [0.3]	< 10
Bilirubin (unconjugated)	60 mg/dL [1026 µmol/L]	8 [8]	< 10
Bilirubin (conjugated)	40 mg/dL [684 µmol/L]	0.3 [0.3]	< 10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	8 [8]	< 10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	0.3 [0.3]	< 10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	8 [8]	< 10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

*Analyte results should not be corrected based on this bias.

Interfering Substances					
The TSHL method was evaluated for interference according to CLSI/NCCLS EP7-A2. ¹⁷ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.					
Analyte Concentration	Biotin Test Level (ng/mL)				
	50	100	250	500	1200
	% Bias				
0.3 µIU/mL	-3	-1.6	-9.2	-49.7	93.6
8 µIU/mL	-3.5	-3.5	-10.6	-42.5	-99.4

- Specimens that contain biotin at a concentration of 100 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely depressed results for patient samples.
- The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.¹⁸ Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.¹⁹

Non-Interfering Substances

The following substances do not interfere with the TSHL method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at Thyroid Stimulating Hormone concentrations of 0.3 µIU/mL [mIU/L] and 8.0 µIU/mL [mIU/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Amikacin	8 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	3000 mg/dL	750 µmol/L
Diazepam	0.51 mg/dL	18.0 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2427.2 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.9 µmol/L
Protein: albumin	6 g/dL	60 g/L
Protein: Total	10.6 g/dL	106 g/L
Rheumatoid Factor	500 IU/mL	500,000 IU/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Hook Effect

One-step sandwich immunometric assays are susceptible to a high dose "hook effect," where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.²⁰ Such samples must be diluted and re-assayed prior to reporting the results. The TSHL method shows no hook effect up to 30000 µIU/mL [mIU/L].

Cross-reactivity

The following substances cross-react with the TSHL method when tested at a TSH level of 0.000 µIU/mL [mIU/L] at the concentrations indicated:

Substance	Concentration Tested	Cross-reactivity
FSH	1000 mIU/mL [IU/L]	≤ 0.005 µIU/mL
LH	1000 mIU/mL [IU/L]	≤ 0.005 µIU/mL
hCG	200,000 mIU/mL [IU/L]	≤ 0.005 µIU/mL

Recovery

Known amounts of TSH (WHO 2nd IRP 80/558) were added to a human serum pool with a baseline TSH value of 0.24 µIU/mL [mIU/L]. The sample concentrations were measured and the calculated percent recovery ranged from 98% to 107% with a mean recovery of 102%.

$$\% \text{ Recovery} = \frac{\text{Value obtained} - \text{Baseline}}{\text{Amount Added}} \times 100$$

Dilution Recovery

Six samples with TSH values from 35.2 to 91.9 µIU/mL [mIU/L] were diluted with TSH Sample Diluent and MULTI 2 Sample Diluent and tested by the Dimension® EXL™ TSHL method. The observed values were plotted against the expected values for each sample and slopes of 0.97 to 1.05 were obtained. The linear relationship was verified by second order polynomial regression (no significant quadratic term at 95% confidence).

Limit of Detection and Limit of Blank²¹

The Limit of Detection (LoD) for the TSHL method is 0.004 µIU/mL [mIU/L], determined consistent with CLSI guideline EP17-A and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 264 determinations, with 120 blank and 144 low level samples.

The Limit of Blank (LoB) is 0.0025 µIU/mL [mIU/L].

- LoD is the lowest concentration of analyte that can be detected reliably.
- LoB is the highest concentration that is likely to be observed for a blank sample.

Functional Sensitivity: 0.007 µIU/mL [mIU/L]

The functional sensitivity represents the lowest concentration with an observed 20% coefficient of variation. The functional sensitivity was determined with multiple serum pools in the range of 0.003 to 0.034 µIU/mL [mIU/L]. All samples were assayed using one replicate per sample, twice a week for 6 weeks, with a recalibration after 3 weeks. Based on this study, the TSHL method meets the definition of a third generation TSH assay.²

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension Vista®, LOCI®, Dimension®, EXL™ and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía

- Burtis CA., Ashwood ER., Bruns DE, eds. Tietz Textbook of Clinical Chemistry and molecular Diagnostics. 4th ed. Philadelphia: Elsevier Inc, 2006 pp. 2053-2087.
- Demers L.M, and Spencer C.A., Eds., Standards of Laboratory Practice, Laboratory support for the diagnosis & monitoring of thyroid disease, National Academy of Clinical Biochemistry, 2002.
- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al. Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad. Sci. USA, Vol 91, pp 5426-5430, June 1994 Biochemistry.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard — Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Tietz NW. Clinical Guide to Laboratory Tests. 4th ed., 2006, Saunders WB Co, Philadelphia, PA., pp 1040.
- Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007: p 867.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection — Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. ClinChem 1999; 45 (7): 942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. Clin Chem 1992; 38:1737-1742.
- Clinical and Laboratory Standards Institute/NCCLS. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline — Second Edition*. CLSI/NCCLS document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Clinical and Laboratory Standards Institute (formerly NCCLS). *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline — Third Edition*. CLSI Document EP28-A3c. CLSI 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
- Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC: AACC Press; 2005
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline — Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline — Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4):247–256.
- Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 2017;55(6):817-825
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays, Anal Biochem 1982;127:308-315.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusqu'à / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contenuto sufficiente per <n> saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EF16S

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare