Dimension® EXL™ integrated chemistry system LOCI® Module

Flex® reagent cartridge

FT3

Issue Date 2020-08

See shaded sections: Updated information from 2019-04 version.

LOCI Free Triiodothyronine

Intended Use: The FT3 method is an *in vitro* diagnostic test for the quantitative measurement of Free Triiodothyronine in human serum and plasma on the Dimension® EXL™ integrated chemistry system with LOCI® Module. Measurements of free triiodothyronine are used in the diagnosis and monitoring of thyroid disease.

Summary: The thyroid gland secretes two hormones, 3,5,3',5'-L-tetraiodothyronine or thyroxine (T4) and 3,5,3'-L-triiodothyronine (T3) in response to the pituitary thyroid stimulating hormone (TSH). At least 85% of normal T3 production comes from peripheral deiodination of T4. Although serum levels of T3 are less than T4, T3 has four to five times higher physiological potency than does T4.¹

In circulation 99.7% of T3 is bound to proteins such as Thyroxine-Binding Globulin (TBG), Thyroxine-Binding Prealbumin (TBPA) and albumin. Changes in binding proteins may affect the interpretation of thyroid hormone levels in the diagnosis of thyroid deficiency states. Measurement of free T3 is useful when total T3 is altered due to changes in thyroid hormone binding proteins. Elevated serum free T3 levels are seen in hyperthyroidism.

Principles of Procedure: The FT3 method is a homogeneous, sequential, chemiluminescent immunoassay based on LOCI® technology. The LOCI® reagents include two synthetic bead reagents and a biotinylated anti-T3 sheep monoclonal antibody. The first bead reagent (Chemibeats) is coated with diiodothyronine (T2), a naturally occurring, weaker binding analog of T3, and contains chemiluminescent dye. The second bead reagent (Sensibeats) is coated with streptavidin and contains a photosensitizer dye. In a first step, sample is incubated with biotinylated antibody which allows T3 from the sample to saturate a fraction of the biotinylated antibody that is directly related to the free T3 concentration. In a second step, T2 chemibeats are added and form bead/biotinylated antibody immunocomplexes with the non-saturated fraction of the biotinylated antibody. Sensibeats are then added and bind to the biotin to form bead pair immunocomplexes. Illumination of the complex at 680 nm generates singlet oxygen from Sensibeats which diffuses into the Chemibeats, triggering a chemiluminescent reaction. The resulting signal is measured at 612 nm and is an inverse function of the FT3 concentration in the sample.^{2,3}

Reagents

Wells ^{a,b}	Form	Ingredient	Concentration ^c	Source
1 – 2	Liquid	Streptavidin Sensibeats	400 µg/mL	Recombinant <i>E. coli</i>
3 – 4	Liquid	T2 Chemibeats	50 µg/mL	
5 – 6	Liquid	Biotinylated FT3 antibody	15 ng/mL	Sheep, monoclonal
7 – 8	Empty			

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Wells 1 – 6 contain buffers, stabilizers and preservatives.

c. Nominal value per well in a cartridge.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!
May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used HM reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Recommended specimen types: serum, lithium and sodium heparin plasma, EDTA plasma. Serum is the preferred specimen.^{4,5}

Samples and controls stabilized with sodium azide cannot be used.

Collect serum and plasma specimens using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁵

Follow the instructions provided with your specimen collection device for use and processing.⁷ Serum separator or plasma separator tubes are acceptable. Variations in collection tubes and separator products may exist. Users should validate the acceptability of these products for patient sample collection.

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation. If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of particulate matter.⁸

Separate serum or plasma from cells as soon as possible. Separated samples are stable for 8 hours at ambient temperature (20 – 25 °C).⁴

Separated samples are stable for 7 days refrigerated at 2 – 8 °C.^{4,5} For longer storage, samples should be frozen at -20 °C for 1 month.^{4,5,9} Avoid repeated freezing and thawing.

The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Procedure

FT3 Flex® reagent cartridge, Cat. No. RF616

Materials Provided

Materials Required But Not Provided

LOCI Thyroid Calibrator, Cat. No. RC610A

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension® EXL™ integrated chemistry system with LOCI® Module. For details of this processing, refer to your Dimension® EXL™ with LM Operator's Guide.

Test Conditions

Sample Volume (delivered to the reaction vessel)	15 µL
T3 Biotinylated Antibody Reagent Volume	25 µL
T2 Chemibeats Reagent Volume	25 µL
Streptavidin Sensibeats Reagent Volume	37 µL
Temperature	37.0 °C
Reaction Time	18.5 minutes
Wavelength	Illumination 680 nm, Emission 612 nm
Type of Measurement	Chemiluminescence

Calibration

LOCI Thyroid Calibrator, Cat. No. RC610A

5 levels, n = 3
pg/mL [pmol/L]^d = [pmol/L]
(pg/mL x 1.54) = [pmol/L]

Typical Calibration Levels

Level 2: 0 pg/mL [0 pmol/L]
Level 3: 2 pg/mL [3 pmol/L]
Level 4: 6 pg/mL [9 pmol/L]
Level 5: 15 pg/mL [23 pmol/L]
Level 6: 32 pg/mL [49 pmol/L]

Calibration Frequency

A new calibration is required:

- Every 15 days for any one lot
- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

d. Système International d'Unités [SI units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze two levels of a Quality Control (QC) material with known free triiodothyronine concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of free triiodothyronine in pg/mL [pmol/L] using the calculation scheme described in your Dimension® EXL™ with LM Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.50 – 30.00 pg/mL [0.8 – 46.2 pmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 30.00 pg/mL [46.2 pmol/L] should be reported as greater than 30.00 pg/mL [46.2 pmol/L]. Samples should NOT be diluted.
- Samples with results less than 0.50 pg/mL [0.8 pmol/L] should be reported as "less than 0.50 pg/mL [0.8 pmol/L]".

Limitations of Procedure

Patient samples may contain human anti-sheep antibodies (HASA) that could give falsely elevated or depressed results with assays that use sheep monoclonal antibodies. This assay has been designed to minimize interference from HASA containing samples. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{10,11} Thyroid hormone autoantibodies in human serum may interfere with the assay.

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in FT3 results. Refer to your Dimension® EXL™ with LM Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported if indicated.

United States Federal law restricts this device to sale by or on the order of a physician or other practitioner licensed by the laws of the State in which he practices, to use or order the use of the device.

Do not use in patients taking biotin supplements.

Interfering Substances

The FT3 method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹² The following substances, tested at concentrations listed, caused elevated FT3 results due to releasing of T3 from serum binding proteins.

Substance	Concentration	SI units	FT3 Concentration pg/mL [pmol/L]	
			2.92 [4.5] Bias (%)	7.79 [12.0] Bias (%)
Diclofenac	50 µg/mL	157 µmol/L	38	34
Furosemide	6 mg/dL	181 µmol/L	57	43
Ibuprofen	50 mg/dL	2425 µmol/L	34	32
Linoleic acid	2.8 mg/mL	10 mmol/L	60	37
Mefenamic acid	0.1 mg/mL	414 µmol/L	49	37
Oleic acid	2.8 mg/mL	10 mmol/L	42	32
Phenylbutazone	15 mg/dL	486 µmol/L	15	11
Phenytoin	5 mg/dL	198 µmol/L	17	13
Salicylic acid	60 mg/dL	4.34 mmol/L	22	25
Valproic acid	50 mg/dL	3467 µmol/L	<10	10

Analyte Concentration	Biotin Test Level (ng/mL)			
	100	250	500	1200
	% Bias			
2.92 pg/mL	2.5	20.9	293.9	Above AMR
7.79 pg/mL	6.8	57.7	Above AMR	Above AMR

* Specimens that contain biotin at a concentration of 100 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely elevated results for patient samples.

* The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.¹³ Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.¹⁴

Expected Values

Adult Population: 2.18 – 3.98 pg/mL [3.4 – 6.1 pmol/L]

The reference interval was transferred from what was previously determined for the FT3 method on the Dimension Vista® System. This reference interval was calculated non-parametrically and represents the central 95% of results determined from a population of 298 apparently healthy adults (188 males and 110 females, 18 – 65 years of age). The original determination and transference of the FT3 method on the Dimension® EXL™ system was done in accordance with CLSI/NCCLS EP28-A2.¹⁵

Pediatric Population:

Pediatric Age Group	Number of Samples	Reference Intervals	
		[pg/mL]	[pmol/L]
Infants (01 – 23 months)	75	3.47 – 5.29	5.3 – 8.2
Children (02 – 12 years)	185	3.35 – 4.82	5.2 – 7.4
Adolescents (13 – 20 years)	147	2.91 – 4.70	4.5 – 7.2

Reference intervals for the pediatric population (infants, children, and adolescents) were established for the Dimension® EXL™ FT3 method in accordance with the CLSI guideline EP28-A3c.¹⁶ Serum samples were collected prospectively from apparently healthy (euthyroid) pediatric subjects, using predefined inclusion criteria.

The reference interval for infants for the FT3 method was calculated by a robust measure of location and spread as developed by Horn and Pesce.¹⁷ A non-parametric approach based on the CLSI guideline was used to establish the reference intervals for children and adolescents. The 2.5th and 97.5th percentiles of the distribution of values were calculated for each age group.

Each laboratory should establish its own expected values for FT3 as performed on the Dimension® EXL™ integrated chemistry system with LOCI® Module.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n = 5 replicates at the following FT3 concentrations are:

FT3 concentration	Acceptable SD Maximum
2.01 pg/mL [3.1 pmol/L]	0.17 pg/mL [0.3 pmol/L]
5.94 pg/mL [9.1 pmol/L]	0.25 pg/mL [0.4 pmol/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® EXL™ integrated chemistry system with LOCI® Module.

Material	Precision ^{18, e}		
	Mean pg/mL [pmol/L]	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Bio-Rad Liquichek™ Immunoassay Plus Control			
Level 1	1.58 [2.4]	0.04 [0.06] (2.5)	0.07 [0.11] (4.4)
Level 2	7.33 [11.3]	0.07 [0.11] (1.0)	0.18 [0.28] (2.5)
Level 3	14.20 [21.9]	0.10 [0.15] (0.7)	0.39 [0.60] (2.7)
Normal Plasma	2.42 [3.7]	0.07 [0.11] (2.9)	0.10 [0.15] (4.1)
Serum Pool 1	7.22 [11.1]	0.08 [0.12] (1.1)	0.17 [0.26] (2.4)
Serum Pool 2	16.48 [25.4]	0.15 [0.23] (0.9)	0.46 [0.71] (2.8)
Serum Pool 3	25.51 [39.3]	0.49 [0.75] (1.9)	0.71 [1.09] (2.8)

e. CLSI/NCCLS EP5-A2¹⁸ was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Method Comparison¹⁹ Regression Statistics¹

Comparative Method	Slope	Intercept pg/mL [pmol/L]	Correlation Coefficient	n

f. CLSI/NCCLS EP9-A2¹⁹ was used. The method used to fit the linear regression line was ordinary least squares.
g. The range of Vista® FT3 values in the correlation study was 0.52 – 28.98 pg/mL [0.8 – 44.6 pmol/L].

Specificity

Albumin, prealbumin, and thyroxine binding globulin (TBG) do not affect the FT3 assay.

One hundred samples were assayed for albumin, prealbumin, and FT3 (albumin range 3.4 – 5.2 g/dL, prealbumin range 17.1 – 48.2 mg/dL, and FT3 range 1.75 – 3.80 pg/mL). Seventy-eight samples were assayed for TBG and FT3 (TBG range 8 – 70 µg/mL and FT3 range 2.09 – 4.72 pg/mL). The FT3 concentration was plotted versus the concentrations of the three binding proteins and a linear regression analysis performed for each. The slope was found to be not significant in all three cases.

Hemolysis, Icterus, Lipemia (HIL) Interference

FT3 was evaluated for interference according to CLSI/NCCLS EP7-A2.¹² Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	pg/mL	FT3 [pmol/L]	Bias* %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L]	2.92	4.5	<10
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L]	7.79	12.0	<10
Bilirubin (unconjugated)	20 mg/dL [342 µmol/L]	2.92	4.5	<10
Bilirubin (unconjugated)	20 mg/dL [342 µmol/L]	7.79	12.0	<10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	2.92	4.5	<10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	7.79	12.0	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	2.69	4.1	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	7.79	12.0	<10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* Analyte results should not be corrected based on this bias.

Non Interfering Substances

The following substances do not interfere with the FT3 method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at an FT3 concentration of approximately 2.92 pg/mL [4.5 pmol/L] and 7.79 pg/mL [12.0 pmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	3.5 mg/dL	232 µmol/L
Amikacin	8 mg/dL	137 µmol/L
Amiodarone	6 µg/mL	8.8 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1.0 mg/dL	33.3 µmol/L
Chlorothiazide	20 µg/mL	67 µmol/L
Chlorpromazine	0.20 mg/dL	6.27 µmol/L
Cholesterol	400 mg/dL	10.3 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2.65 mmol/L
Dextran 40	4500 mg/dL	1125 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Gentamicin	1.0 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Methimazole	4 µg/mL	35 µmol/L
N-acetylcysteine	1.3 mg/mL	8.0 mmol/L
Nicotine	0.10 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Propylthiouracil	50 µg/mL	294 µmol/L
Rheumatoid Factor (RF)	500 IU/mL	500,000 IU/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	1000 mg/dL	11.3 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Vancomycin	10 mg/dL	69 µmol/L

Cross reactivity

The following substances were added into the zero level LOCI Thyroid Calibrator at the concentrations indicated. % cross-reactivity reflects the ratio of the FT3 concentration observed with added substance relative to a reference spike of 10 ng/mL T3.

$$\% \text{ CR} = \frac{\text{FT3 value with substance added}}{\text{FT3 value with T3 added}} \times \frac{\text{Concentration of substance added}}{\text{Concentration of T3 added}} \times 100$$

Substance	Concentration	Cross reactivity
L-Thyroxine	2000 ng/mL	0.09%
D-Thyroxine	2000 ng/mL	0.07%
Reverse T3	2000 ng/mL	0.53%
Monoiodotyrosine	2000 ng/mL	<0.008%
Diiodotyrosine	2000 ng/mL	<0.008%

Limit of Detection and Limit of Blank²⁰

The Limit of Detection (LoD) for the FT3 method is 0.50 pg/mL [0.8 pmol/L], determined consistent with CLSI guideline EP17-A and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 240 determinations, with 120 blank and 120 low level samples.

The Limit of Blank (LoB) is 0.21 pg/mL [0.3 pmol/L].

LoD is the lowest concentration of analyte that can be detected reliably.

LoB is the highest concentration that is likely to be observed for a blank sample.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension Vista®, LOCI®, Dimension®, Flex® and EXL™ are trademarks of Siemens Healthcare Diagnostics.

©2010 Siemens Healthcare Diagnostics
All rights reserved.

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

- Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th edition, Elsevier Inc., St. Louis, Missouri, 2006: p. 2053 – 2054 (The Thyroid: Pathophysiology and Thyroid Function Testing).
- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al. Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method. *Clin Chem* 42:9 1996, 1518 – 1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al. Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. *Proc.Natl.Acad.Sci. USA*, Vol 91, pp. 5426 – 5430, June 1994 Biochemistry.
- Clinical and Laboratory Standards Institute/NCCLS. *Measurement of Free Thyroid Hormones; Approved Guideline*. CLSI/NCCLS document C45-A [ISBN 1-56238-548-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- National Academy of Clinical Biochemistry (NACB). *Laboratory Medicine Practice Guidelines, Laboratory Support for the Diagnosis and Monitoring of Thyroid Disease*, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection—Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Wu A, ed. *Tietz Clinical Guide to Laboratory Tests*, 4th edition. Philadelphia; W.B. Saunders Co., 2006: p. 1076.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. *ClinChem* 1999; 45 (7): 942 – 956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem* 1992; 38:1737–1742.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4): 247–256.
- Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 2017;55(6):817-825
- Clinical and Laboratory Standards Institute/NCCLS. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition*. CLSI/NCCLS document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Clinical and Laboratory Standards Institute (formerly NCCLS). *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. CLSI Document EP28-A3c. CLSI 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
- Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC: AACC Press; 2005.
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_ERGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare