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1. Siekmann L, Ceriotti F, Ferard G, Kanno T, Schumann G, et. al. IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 ° C. Part 2. Reference procedure for the measurement of catalytic concentrations of creatine kinase [ATP: creatine N-phosphotransferase (CK), EC 2.7.3.2] Clin Chem Lab Med. 2002; 40(6):635-642.
2. McPherson RA, Pincus MR, eds. Henry's clinical diagnosis and management by laboratory methods 21st ed. St. Louis: Saunders Elsevier, 2007: 221, 260.
3. Anderson JL, Adams CD, Antman EM, Bridges Cr, Califf RM, Casey DE, et.al. ACC/ACC 2007 Guidelines for the management of patients with unstable angina/non-ST-elevation myocardial infarction-executive summary. JACC.50(7), 2007:221,260.
4. Burtis CA, Ashwood ER, eds. Tietz fundamentals of clinical chemistry, 5th ed. Philadelphia: WB Saunders Co., 2001:358-62.
5. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard- Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard - Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimen; Approved Guideline - Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
8. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition.

CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.

9. Klein G, Berger A, Bertholf R et al. Abstract: Multicenter Evaluation of Liquid Reagents for CK, CK-MB and LDH with Determination of Reference Intervals on Hitachi Systems. Clin Chem 2001; 47:Suppl. A30.

10. Thomas L, Müller M, Schumann G, Weidemann G et al. Consensus of DGKL and VDGH for Interim reference intervals on enzymes in Serum. J Lab Med 2005; 29:301-308.

11. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline- Second Edition. CLSI/NCCLS Document EP5-A2 [ISBN1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.

13. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

13. Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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Siemens Healthcare Diagnostics Inc.

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