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Comparison of PFA-100 and bleeding time testing in pediatric patients with suspected hemorrhagic problems. *J Pediatr Hematol Oncol* 2003;25:474-9.

3. Lippi G, Manzato F, Franchini M, et al. Establishment of reference values for the PFA-100 platelet function analyzer in pediatrics. *Clin Exp Med* 2001; 1:69-70.
4. Koscielny J, Ziemer S, Radtke H, et al. A practical concept for preoperative identification of patients with impaired primary hemostasis. *Clin Appl Thromb Hemo* 2004;10:195-204.
5. Fressinaud E, Veyradier A, Sigaud M, et al. Therapeutic monitoring of von Willebrand disease: interest and limits of a platelet function analyzer at high shear rates. *Brit J Haematol* 1999;106:777-83.
6. Koscielny J, von Tempelhoff G-F, Ziemer S, et al. A practical concept for preoperative management of patients with impaired primary hemostasis. *Clin Appl Thromb Hemo* 2004;10:155-66.
7. Kratzer MAA, Born GVR Simulation of primary haemostasis in vitro. *Haemostasis* 1985;5:357-62.
8. Kratzer MAA, Bellucci S, Caen JP. Detection of abnormal platelet function with an in vitro model of primary hemostasis. *Haemostasis* 1985;15:363-70.
9. Born GVR. Aggregation of blood platelets by adenosine diphosphate and its reversal. *Nature* 1962;194:927-9.
10. Kundu SK, Heilmann EJ, Sio R, et al. Characterization of an in vitro platelet function analyzer, PFA-100. *Clin Appl Thromb Hemo* 1996;2:241-9.
11. MacIntyre DE, Hoover RL, Smith M, et al. Inhibition of platelet function by cis-unsaturated fatty acids. *Blood* 1984;63:848-57.
12. Driss F, Vericel E, Lagarde M, et al. Inhibition of platelet aggregation and thromboxane synthesis after intake of small amount of icosapentaenoic acid. *Thromb Res* 1984;36:389-96.
13. Schick PK. Megakaryocyte and platelet lipids. In: Colman RW, Hirsh J, Marden VJ, Salzman EW, eds. *Hemostasis and Thrombosis; Basic Principles and Clinical Practice*, 3rd ed., Philadelphia: JB Lippincott Co., 1994:574-89.
14. Reinhart W, Felix C. Influence of propofol on erythrocyte morphology, blood viscosity and platelet function. *Clin Hemorheol Microcirc.* 2003;29:33-40.
15. Porter JM, Crowe B, Cahill M, Shorten GD. The effects of ropivacaine hydrochloride on platelet function: an assessment using the platelet function analyser (PFA-100). *Anaesthesia* 2001;56:15-8.
16. Fressinaud E, Veyradier A, Sigaud M, et al. Therapeutic monitoring of von Willebrand disease: interest and limits of a platelet function analyzer at high shear rates. *Brit J Haematol* 1999;106:777-83.
17. Bauer KA, Gerson W, Wright C, et al. Platelet function following administration of a novel formulation of intravenous diclofenac sodium versus active comparators: a randomized, single-dose, crossover study in healthy male volunteers. *J Clin Anesth.* 2010;22:510-8.
18. Van Kraaij D, Hovestad-Witterland A, De Metz M, Vollaard E. A comparison of the effects of nabumetone vs meloxicam on serum

thromboxane B2 and platelet function in healthy volunteers. Br J Clin Pharmacol. 2002;53:644-7.

19. Liu FC, Liao CH, Chang YW, et al. Hydroxyethyl Starch interferes with human blood ex vivo coagulation, Platelet Function and Sedimentation. Acta Anaesthesiol Taiwan 2009;47:71-8.

20. Madan M, Berkowitz SD, Christite DJ, et al. Determination of platelet aggregation inhibition during percutaneous coronary intervention with the platelet function analyzer PFA-100. Am Heart J 2002;144:151-8.

21. Zwerina J, Landsteiner H, Leitgeb U, et al. The influence of VIP and epoprostenol on platelet CD62P expression and primary haemostasis in vitro. Platelets. 2004;15:55-60.

22. Kobsar A, Koessler J, Rajkovic M, et al. Prostacyclin receptor stimulation facilitates detection of human platelet P2Y(12) receptor inhibition by the PFA-100(R) system. Platelets. 2010;19. 112-6.

23. Böck M, de Haan J, Beck KH, et al. Standardisation of the PFA-100® Platelet Function Test in 105 mmol/L buffered citrate: effect of gender, smoking, and oral contraceptives. Br J Haematol 1999;106:898-904.

24. Clinical and Laboratory Standards Institute (CLSI). Defining, establishing, and verifying reference intervals in the clinical laboratory; approved guideline - third edition. CLSI document C28-A3 (ISBN 1-56238-682-4).CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2008.

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Products GmbH