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- 1) Bayer M, et al. Clinical experience with sensitive thyrotropin measurements: diagnostic and therapeutic implications. J Nucl Med 1985;36:1248-56.
- 2) Burger HG, Patel TC. Thyrotrophin releasing hormone—TSH. Clin Endocrinol Metab 1977 March;6(1):83-100.
- 3) Chen I-W, Heminger L. Thyroid-stimulating hormone. In: Kaplan LA, Pesce AJ, editors. Clinical chemistry. St Louis: C.V. Mosby, 1984:1160-64.
- 4) Durham AP. The upper limit of normal for thyrotropin is 3 or 4 milli-int. units/L. Clin Chem 1985;31:296-98.
- 5) Fisher DA, Klein AH. Thyroid development and disorders of thyroid function in the newborn. N Eng J Med 1981;304:702-12.
- 6) Fraser CG, Browning MCK. Measuring serum thyrotropin. Lancet 1985;1:816-7.
- 7) Jackson IMD. Thyrotropin releasing hormone. N Engl J Med 1982;306:145-55.
- 8) Lindstedt G, et al. Thyroid function evaluation in the mid 80s. Scand J Clin Lab Invest 1984;44:465-70.
- 9) Babson, AL. The IMMULITE Automated Immunoassay System. J Clin Immunoassay 1991;14:83-8.
- 10) Ridgway EC, et al. Thyrotropin. In: Rothfeld B, editor. Nuclear medicine in vitro. Philadelphia: J.B. Lippincott, 1974: 205-19.
- 11) Ridgway EC, et al. Peripheral responses to thyroid hormone before and after L-thyroxine therapy in patients with subclinical hypothyroidism. J Clin Endocrinol Metabol 1981;53:1238-42.
- 12) Tsay J-Y, Chen I-W, et al. A statistical method for determining normal ranges from laboratory data including values below the

minimum detectable value. Clin Chem 1979;25:2011-4.

- 13) Walflsch PG. The best way to screen for neonatal hypothyroidism. Diagnostic Medicine 1984 Feb;7(2):67-75.
- 14) Weekie J. The influence of the circadian thyrotropin rhythm on the thyrotropin response to thyrotropin-releasing hormone in normal subjects. Scand J Clin Lab Invest 1974;33:17-20.
- 15) Wehmann RE, et al. Suppression of thyrotropin in the low-thyroxine state of severe nonthyroidal illness. N Engl J Med 1985; 312:546-52.
- 16) Woodhead JS, Weeks I. Circulating thyrotrophin as an index of thyroid function. Ann Clin Biochem 1985;22:455-9.
- 17) Spencer C, et al. Specificity of sensitive assays of thyrotropin (TSH) used to screen for thyroid disease in hospitalized patients. Clin Chem 1987;33:1391-6.
- 18) Nicoloff JT, Spencer CA. The use and misuse of the sensitive thyrotropin assays. J Clin Endocrinol Metab 1990;71:553-8.
- 19) Klee GG, Hay ID. Assessment of sensitive thyrotropin assays for an expanded role in thyroid function testing: proposed criteria for analytic performance and clinical utility. J Clin Endocrinol Metab 1987;64(3):461-71.
- 20) Mandel SJ, Brent GA, Larsen PR. Levothyroxine therapy in patients with thyroid disease. Ann Intern Med 1993;119:492-502.
- 21) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: W.B. Saunders, 1994.
- 22) Spencer CA, Takeuchi M, et al. Interlaboratory/intermethod differences in functional sensitivity of immunometric assays of thyrotropin (TSH) and impact on reliability of measurement of subnormal concentrations of TSH. Clin Chem 1995;41:367-74.
- 23) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
- 24) Baloch Z, Carayon P, Conte-Devolx B, Demers LM, Feldt-Rasmussen U, Henry JF, et al.; Guidelines Committee, National Academy of Clinical Biochemistry (NACB). Laboratory medicine practice guidelines (LMPG). Laboratory support for the diagnosis and monitoring of thyroid disease. Thyroid 2003 Jan;13(1):3-126. Also available at http://www.nacb.org/lmpg/thyroid_LMPG_Word.stm (accessed January 2005).

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Siemens Healthcare Diagnostics Products Limited