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- Bangham DR, Storrin PL. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:390.
- Batzer F. Hormonal evaluation of early pregnancy. *Fertil Steril* 1980;34:1-12.
- Braunstein GD, et al. First-trimester chorionic gonadotropin measurements as an aid in the diagnosis of early pregnancy disorders. *Amer J Obstet Gynecol* 1978;131:25-32.
- Clayton L, et al. Endocrine aspects of trophoblastic neoplasia. *J Reprod Med* 1981;26:192-9.
- Corson SL, Horwitz CA, Lau HL, Soderstrom R. Early diagnosis of pregnancy: a symposium. *J Reprod Med* 1981;26:149-78.
- DeCherney AH, Minkin MJ, Spangler S. Contemporary management of ectopic pregnancy. *J Reprod Med* 1981;26:519-23.
- Goldstein DP, Berkowitz RS, Kosasa TS. Measurement of human chorionic gonadotropin. *Major Problems in Obstet Gynecol* 1982;14:64-97.
- Hager H, Van Weemen BK. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:629.
- Kohorn EI. Hydatidiform mole and gestational trophoblastic disease in southern Connecticut. *Obstet and Gynecol* 1982;59:78-84. Erratum: The first diagram belongs with the caption to Figure 3, and vice versa.
- Kosasa TS. Measurement of human chorionic gonadotropin. *J Reprod Med* 1981;26:201-6.
- Lenton E, Neal L, Sulaiman R. Plasma concentrations of human gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773-8.
- Lindstedt G, Janson PO, Thorburn J. Sensitivity of serum chorionic gonadotropin assay for ectopic pregnancy. *Lancet* 1981;1:781-2.
- McCready J, Braunstein GD, et al. Modification of the chorionic gonadotropin beta-subunit radioimmunoassay for determination of urinary chorionic gonadotropin. *Clin Chem* 1978;24:1958-61.
- Storrin PL, Gaines-Das R, Bangham DR. International Reference Preparation of Human Chorionic Gonadotrophin for Immunoassay: potency estimates in various bioassay and protein binding assay systems. *J Endocrinol* 1980;84:295-310.
- Vaitukaitis J. Human chorionic gonadotropin - a hormone secreted for many reasons. *N Eng J Med* 1979;301:324-6.
- Danzer H, Braunstein GD, et al. Maternal serum human chorionic gonadotropin concentrations and fetal sex predictions. *Fertil Steril* 1980;34:336-40.
- Hussa RO. *The Clinical Marker HCG*. New York: Praeger Publishers, 1987.

18) Greene MF, de Montserrat MF, Tulchinsky D. In: Tietz NW, editor. Textbook of Clinical Chemistry. Philadelphia: Saunders, 1986: 1745.

19) Braunstein GD, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126:678.

20) Babson AL The IMMULITE Automated Immunoassay System. J Clin Immunoassay 1991;14:83-8.

21) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

22) Tietz NW, ed. Clinical Guide to Laboratory Tests. 2nd Ed. Philadelphia: W.B. Saunders, 1990.

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Products Limited