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1. Apple FS. Acute myocardial infarction and coronary reperfusion. Serum cardiac markers for the 1990s. *Am J Clin Pathol*. 1992;97(2):217–226.

2. Adams JE 3rd, Bodor GS, Davila-Romain VG, et al. Cardiac troponin I. A marker for cardiac injury. *Circulation*. 1993;88(1):101–106.

3. Corin SJ, Juhasz O, Zhu L, et al. Structure and expression of the human slow twitch skeletal muscle troponin I gene. *J Biol Chem*. 1994;269(14):10651–10659.

4. Mair J, Wagner I, Puschendorf B, et al. Cardiac troponin I to diagnose myocardial injury. *Lancet*. 1993;341(8848):838–839.

5. Mair J, Larue C, Mair P, et al. Use of cardiac troponin I to diagnose perioperative myocardial infarction in coronary artery bypass grafting. *Clin Chem*. 1994;40(11, pt 1):2066–2070.

6. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. *Clin Chem*. 1992;38(11):2203–2214.

7. Adams JE 3rd, Sicard GA, Allen BT, et al. Diagnosis of perioperative myocardial infarction with measurement of cardiac troponin I. *N Engl J Med*. 1994;330(10):670–674.

8. Adams JE 3rd, Schechtman KB, Landt Y, et al. Comparable detection of acute myocardial infarction by creatine kinase MB isoenzyme and cardiac troponin I. *Clin Chem*. 1994;40(7, pt 1):1291–1295.

9. Cummins B, Auckland ML, Cummins P. Cardiac-specific tro-

ponin-I radioimmunoassay in the diagnosis of acute myocardial infarction. *Am Heart J*. 1987;113(6):1333–1344.

10. Larue C, Calzolari C, Bertinchant JP, et al. Cardiac-specific immunoenzymometric assay of troponin I in the early phase of acute myocardial infarction. *Clin Chem*. 1993;39(6):972–979.

11. The Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. *Circulation*. 1979;59:607–9.

12. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.

13. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med*. 1996;335(18):1342–1349.

14. Galvani M, Ottani F, Ferrini D, et al. Prognostic influence of elevated values of cardiac troponin I in patients with unstable angina. *Circulation*. 1997;95(8):2053–2059.

15. Jaffe AS. Elevations in cardiac troponin measurements: false false-positives. *Cardiovasc Toxicol*. 2001;1(2):87–92.

16. The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology committee for the redefinition of myocardial infarction. *J Am Coll Cardiol*. 2000;36:959–69.

17. Williams EJ, Campbell AK. A homogeneous assay for biotin based on chemiluminescence energy transfer. *Anal Biochem*. 1986;155(2):249–255.

18. Bertinchant JP, Larue G, Pernel I, et al. Release kinetics of serum cardiac troponin I in ischemic myocardial injury. *Clin Biochem*. 1996;29(6):587–594.

19. Galván Ramirez ML, Valdez Alvarado V, Vargas Gutierrez G, et al. Prevalence of IgG and IgM anti-Toxoplasma antibodies in patients with HIV and acquired immunodeficiency syndrome (AIDS). *Rev Soc Bras Med Trop*. 1997;30(6):465–467.

20. Morrow DA, Rifai N, Tanasijevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes: a Thrombolysis in Myocardial Infarction (TIMI) 11B substudy. *Clin Chem*. 2000;46(4):453–460.

21. Keffer JH. The cardiac profile and proposed practice guideline for acute ischemic heart disease. *Am J Clin Pathol*. 1997;107(4):398–409.

22. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.

23. Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.

24. Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.

25. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res.* 1985;45(2):879–885.

26. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem.* 1999;45(7):942–956.

27. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem.* 1992;38(9):1737–1742.

28. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.

29. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c (formerly C28-A3c).

30. Clinical and Laboratory Standards Institute. Assessment of the Clinical Accuracy of Laboratory Tests Using Receiver Operating Characteristic (ROC) Plots; Approved Guideline. Wayne, PA: Clinical and Laboratory Standards Institute; 1995. CLSI Document GP10-A.

31. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.

32. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

33. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.

34. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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住所：511 Benedict Avenue, Tarrytown, New York 10591, USA

生产地址：333 Coney Street, East Walpole, Massachusetts

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Siemens Healthcare Diagnostics Inc.