HDL/LDL Cholesterol
Calibrator

Contents

LOT	538098	 2022-08
------------	--------	---

REF	Contents
00309530 (B03-4763-01)	3 vials of lyophilized calibrator (B13-4763-01) 3 x 0.17 g
HDL/LDL CHOL	CAL

Intended Use

This product is intended for *in vitro* diagnostic use in the calibration of ADVIA® Chemistry systems for the Direct HDL Cholesterol, HDL II and Direct LDL Cholesterol methods.

Calibrator Description

Amount	0.17 g/vial
Ingredients	Lyophilized human serum
Storage	2–8°C
Stability	Unopened – until the expiration date on the vial label. Reconstituted – 3 days

 **CAUTION! POTENTIAL BIOHAZARD:** Contains human source material. While each human serum or plasma donor unit used in the manufacture of this product was tested by FDA-approved methods and found nonreactive for hepatitis B surface antigen (HBsAg), antibody to hepatitis C (HCV), and antibody to HIV-1/2, all products manufactured using human source material should be handled as potentially infectious. Because no test method can offer complete assurance that hepatitis B or C viruses, HIV, or other infectious agents are absent, these products should be handled according to established good laboratory practices.

For *in vitro* diagnostic use.

Preparing the Calibrators

- Pipette exactly 1.0 mL of distilled water into one (1) vial of the lyophilized material.
- Allow to stand for 30 minutes.
- Ensure homogeneity prior to use.

Disposal

Dispose of hazardous and biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with all federal, state, and local requirements.

Standardization

Values are traceable to the standardization of each specific assay. For additional information refer to the method sheet.

Technical Assistance

For customer support, please contact your local technical support provider or distributor.

siemens.com/healthcare

ADVIA is a trademark of Siemens Healthcare Diagnostics.
© 2009 Siemens Healthcare Diagnostics Inc. All rights reserved.

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© 2009 Siemens Healthcare Diagnostics Inc. Med enerett.