ACTH

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 ACTH

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the quantitative measurement of adrenocorticotrophic hormone in EDTA plasma, as an aid in the assessment of adrenal insufficiency and hypersecretion.

Catalog Number: **L2KAC2** (200 tests)

Test Code: **ACT** Color: **Light Blue**

Summary and Explanation

Adrenocorticotrophic hormone (ACTH) is a polypeptide hormone which exists principally as a chain, 39 amino acids long, with a molecular mass of approximately 4,500 daltons. It is produced in the pituitary and serves to stimulate steroid production by the adrenal cortex. ACTH secretion is in turn controlled by the hypothalamic hormone corticotropin releasing factor (CRF) and by negative feedback from cortisol.

ACTH determinations are valuable in the differential diagnosis of adrenal insufficiency and hypersecretion. In Addison's disease (primary adrenal insufficiency), elevated levels are typical, whereas low levels are the rule when adrenal insufficiency is secondary to pituitary dysfunction. ACTH determinations can also help to identify the cause of cortisol hypersecretion in Cushing's syndrome. ACTH levels are typically low when this is due to lesions or hyperplasia of the adrenal cortex, and high when it is due to ectopic ACTH production or hypersecretion of ACTH by the pituitary.

Plasma levels of ACTH exhibit a significant diurnal variation. It is important, therefore, to standardize the time of collection: reference ranges have typically been established for approximately 9 in the morning.

Principle of the Procedure

IMMULITE 2000 ACTH is a solid-phase, two-site sequential chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Collect blood by venipuncture (avoiding hemolysis) into *iced* EDTA tubes, noting the time of collection. The tubes should be immersed in an ice bath following collection. Separate the plasma from the cells by centrifugation in a *refrigerated* centrifuge, then aliquot and freeze immediately in *plastic or silanized glass* tubes. (Do not use nonsilanized glass.)

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 ACTH has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 75 µL EDTA plasma.

Storage: 30 days at –20°C.²⁰ Before assay, samples should be thawed in an ice bath and kept at or below 4°C at all times.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H311, H302,
H411**

**P280, P273,
P301 + P312,
P302 + P312,
P501**

Danger! Toxic in contact with skin. Harmful if swallowed. Toxic to aquatic life with long lasting effects.

Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF

SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: sodium azide; ACTH Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

ACTH Bead Pack (L2AC12)

With barcode. 200 beads, coated with monoclonal murine anti-ACTH. Stable at 2–8°C until expiration date.

L2KAC2: 1 pack.

ACTH Reagent Wedge (L2ACA2)

With barcode. 2 reagents: 11.5 mL protein buffer/serum matrix; 11.5 mL alkaline phosphatase (bovine calf intestine)

conjugated to polyclonal rabbit anti-ACTH in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KAC2: 1 wedge.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

ACTH Adjustors (LACL, LACH)

Two vials (Low and High) of lyophilized ACTH in a bovine protein-based matrix, with preservative. Reconstitute each vial with **4.0 mL** distilled or deionized water. Let stand for 30 minutes, then mix by *gentle* swirling or inversion. Aliquot and freeze. Stable at –20°C for 2 months after reconstitution.

L2KAC2: 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

LACCM: Bi-level ACTH control module (protein-based)

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of ACTH.

Expected Values

Based on its relationship to IMMULITE ACTH (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

A study performed on 59 apparently healthy laboratory volunteers yielded a median of 24 pg/mL and a 95% reference range of

ND to 46 pg/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

ACTH measurement is not recommended for patients who have been administered cosyntropin (ACTH (1-24)) because the results may be falsely elevated or depressed.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on EDTA plasma samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

pg/mL \times 0.222 \rightarrow pmol/L

Calibration Range: Up to 1,250 pg/mL (278 pmol/L).

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 5 pg/mL (1,1 pmol/L).

High-dose Hook Effect: None up to 1,500,000 pg/mL (333,000 pmol/L).

Intraassay Precision (Within-Run):

Statistics were calculated for five samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for five samples assayed in 8 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three ACTH solutions (1,100, 2,400 and 5,550 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: Samples with apparent ACTH concentrations of 0 pg/mL showed some cross-reactivity with ACTH (1-18), ACTH (1-24), ACTH (18-39) and Alpha MSH. These compounds were shown to bind with the antibodies in the assay, leading to an increase in the apparent concentration of ACTH in the sample. The presence of cosyntropin (ACTH (1-24)) in a sample containing ACTH can lead to either an increase or a decrease in concentration, dependent on the amount of ACTH present in the sample. (See "Specificity" table and "Limitations" section.)

The antibody is highly specific for ACTH. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 5,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To determine whether serum may be used in the procedure, blood was collected from 9 laboratory volunteers into EDTA, heparinized and plain vacutainer tubes. All samples were spiked with ACTH then assayed by the IMMULITE 2000 ACTH procedure. Results show that serum and heparinized plasma are not recommended for use with this assay.

(Heparin) = 0.79 (EDTA) + 1.2 pg/mL
 $r = 1.00$
 (Serum) = 0.91 (EDTA) – 3.4 pg/mL
 $r = 0.999$

Means:
 218 pg/mL (EDTA)
 173 pg/mL (Heparin)
 194 pg/mL (Serum)

Method Comparison: The assay was compared to IMMULITE ACTH on 86 samples. (Concentration range: approximately 20 to 1,200 pg/mL. See graph.) By linear regression:

(IML 2000) = 0.95 (IML) + 2.9 pg/mL
 $r = 0.988$

Means:
 241 pg/mL (IMMULITE 2000)
 252 pg/mL (IMMULITE)

References

- 1) Broughton A. Application of adrenocorticotropin assays in a routine clinical laboratory. *Am J Clin Path* 1975;64:618-24.
- 2) Cizza G, Chrousos GP. Adrenocorticotrophic hormone-dependent Cushing's syndrome. *Cancer Treat Res* 1997;89:25-40.
- 3) Demers LM, Whitley RJ. Function of the adrenal cortex. In: Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 3rd ed. Philadelphia: WB Saunders, 1999: 1530-69.
- 4) Donald RA. ACTH and related peptides. *Clin Endocrinol* 1980;12:491-524.
- 5) Findling JW, et al. Selective venous sampling for ACTH in Cushing's syndrome. *Ann Intern Med* 1981;94:647-52.
- 6) Gold EM. The Cushing syndromes: changing views of diagnosis and treatment. *Ann Intern Med* 1979;90:829-44.
- 7) Horrocks PM, London DR. Diagnostic value of 9 AM plasma adrenocorticotrophic hormone concentrations in Cushing's disease. *Br Med J* 1982;285:1302-3.
- 8) Katz FH. Adrenocortical diseases: identification and management. *Postgrad Med* 1979;66(6):52-65.
- 9) Krieger D. Plasma ACTH and corticosteroids. In: DeGroot L, et al, editors. *Endocrinology*, vol 2.

New York: Grune & Stratton, 1979: 1139-56.

- 10) Krieger D, et al. ACTH, β -lipotropin and related peptides in brain, pituitary, and blood. *Recent Prog Horm Res* 1980;36:277-344.
- 11) Nelson DH. The adrenal cortex: physiological function and disease. Philadelphia: W.B. Saunders, 1980.
- 12) Pullan PT, et al. ACTH, LPH and related peptides in the ectopic ACTH syndrome. *Clin Endocrinol* 1980;13:437-45.
- 13) Rees L, Lowry PJ. Adrenocorticotrophin and lipotrophin. In: Gray CH, James V, editors. *Hormones in blood*, vol 3. London: Academic Press, 1979: 129-78.
- 14) Schwyzer R. ACTH: a short introductory review. *Ann N Y Acad Sci* 1977;297:3-26.
- 15) Wolfsen A. Aging and the adrenals. In: Korenman SG, editor. *Endocrine aspects of aging*. New York: Elsevier Biomedical, 1982: 55-79.
- 16) Wolfsen A, Odell WD. The dose-response relationship of ACTH and cortisol in Cushing's disease. *Clin Endocrinol* 1980;12:557-68.
- 17) Woodard BH, et al. Adrenocorticotropin production by a mammary carcinoma. *Cancer* 1981;47:1823-7.
- 18) Yalow R. Ectopic ACTH in carcinoma of the lung. In: Muggia FM, Rozenzweig M, editors. *Lung cancer: progress in therapeutic research*. New York: Raven Press, 1978: 209-16.
- 19) Yalow R, et al. Plasma and tumor ACTH in carcinoma of the lung. *Cancer* 1979;44:1789-92.
- 20) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs


Intraassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	23	2.0	8.7%
2	30	2.0	6.7%
3	40	2.7	6.8%
4	208	14.3	6.9%
5	421	39.8	9.5%

Interassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	24	2.4	10.0%
2	44	3.6	8.2%
3	89	7.7	8.7%
4	229	21.3	9.3%
5	496	30.0	6.1%

Method Comparison

$$(IML\ 2000) = 0.95 (IML) + 2.9\ \text{pg/mL}$$

$$r = 0.988$$

Specificity

Compound/ Amount pg/mL ¹	Initial Concentration pg/mL ²	Change in Concentration pg/mL ³	% Cross- reactivity ⁴
ACTH (1-18)			
500	0	ND	ND
5,000	0	ND	ND
50,000	0	7	0.01%
500,000	0	5	0.001%
ACTH (1-24)			
500	0	5	1%
5,000	0	26	0.5%
50,000	0	86	0.2%
50,000	43	48	0.1%
50,000	263	-148	-0.3%
500,000	0	190	0.04%
ACTH (18-39)			
500	0	87	17.3%
5,000	0	829	16.6%
Alpha MSH			
500	0	ND	ND
5,000	0	ND	ND
50,000	0	7	0.01%
500,000	0	ND	ND

ND: Not detectable.⁵

Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	16	—	—
	A	66	70	94%
	B	130	135	96%
	C	259	290	89%
2	—	27	—	—
	A	86	81	106%
	B	147	146	101%
	C	294	301	98%
3	—	99	—	—
	A	159	149	107%
	B	237	214	111%
	C	364	369	99%
4	—	233	—	—
	A	302	276	109%
	B	355	341	104%
	C	513	496	103%

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	63	—	—
	4 in 8	31	32	97%
	2 in 8	17	16	106%
	1 in 8	7.0	7.9	89%
2	8 in 8	225	—	—
	4 in 8	118	113	104%
	2 in 8	58	56	104%
	1 in 8	26	28	93%
3	8 in 8	321	—	—
	4 in 8	167	161	104%
	2 in 8	75	80	94%
	1 in 8	36	40	90%
4	8 in 8	326	—	—
	4 in 8	184	163	113%
	2 in 8	89	82	109%
	1 in 8	45	41	110%
5	8 in 8	526	—	—
	4 in 8	272	263	103%
	2 in 8	131	132	99%
	1 in 8	65	66	99%
6	8 in 8	658	—	—
	4 in 8	374	329	114%
	2 in 8	180	165	109%
	1 in 8	82	82	100%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2008-2020 Siemens Healthcare Diagnostics. All rights reserved.