IMMULITE®

ACTH

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE®/IMMULITE® 1000 ACTH

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of adrenocorticotrophic hormone (ACTH) in EDTA plasma, as an aid in the assessment of adrenal insufficiency and hypersecretion.

Catalog Number: **LKAC1** (100 tests)

Test Code: **ACT** Color: **Light Blue**

Summary and Explanation

Adrenocorticotrophic hormone (ACTH) is a polypeptide hormone which exists principally as a chain, 39 amino acids long, with a molecular mass of approximately 4,500 daltons. It is produced in the pituitary and serves to stimulate steroid production by the adrenal cortex. ACTH secretion is in turn controlled by the hypothalamic hormone corticotropin releasing factor (CRF) and by negative feedback from cortisol.

ACTH determinations are valuable in the differential diagnosis of adrenal insufficiency and hypersecretion. In Addison's disease (primary adrenal insufficiency), elevated levels are typical, whereas low levels are the rule when adrenal insufficiency is secondary to pituitary dysfunction. ACTH determinations can also help to identify the cause of cortisol hypersecretion in Cushing's syndrome. ACTH levels are typically low when this is due to lesions or hyperplasia of the adrenal cortex, and high when it is due to ectopic ACTH production or hypersecretion of ACTH by the pituitary.

Plasma levels of ACTH exhibit a significant diurnal variation. It is important, therefore, to standardize the time of collection: reference ranges have typically been established for approximately 9 in the morning.

Principle of the Procedure

IMMULITE/IMMULITE 1000 ACTH is a solid-phase, two-site sequential chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

Collect blood by venipuncture²⁰ (avoiding hemolysis) into *iced* EDTA tubes, noting the time of collection. The tubes should be immersed in an ice bath following collection. Separate the plasma from the cells by centrifugation in a *refrigerated* centrifuge, then aliquot and freeze immediately in *plastic or silanized glass* tubes. (Do not use nonsilanized glass.)

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 ACTH has not been tested with all possible variations of tube types.

Volume Required: 75 µL EDTA plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 30 days at –20°C.²¹ Before assay, samples should be thawed in an ice bath and kept at or below 4°C at all times.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H311, H302,
H411**

**P280, P273,
P301 + P312,
P302 + P312,
P501**

Danger! Toxic in contact with skin. Harmful if swallowed. Toxic to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: sodium azide; ACTH Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

ACTH Test Units (LAC1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-ACTH. Stable at 2–8°C until expiration date.

LKAC1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

ACTH Reagent Wedges (LACA, LACB)

LACA: 7.5 mL of a protein buffer/serum matrix. **LACB:** 7.5 mL of an alkaline phosphatase (bovine calf intestine) conjugated to polyclonal rabbit anti-ACTH antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKAC1: 1 set

ACTH Adjustors (LACL, LACH)

Two vials (Low and High) of lyophilized ACTH in a bovine protein-based matrix, with preservative. Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Aliquot and freeze. Stable at –20°C for 2 months after reconstitution.

LKAC1: 2 sets

Kit Components Supplied Separately

ACTH Sample Diluent (LACZ)

For the manual dilution of patient samples. One vial containing 25 mL of an ACTH-free bovine protein-based matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LACCM: Bi-level ACTH control module (protein-based)

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
4 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of ACTH.

Expected Values

A study performed on 59 apparently healthy laboratory volunteers yielded a median of 24 pg/mL and a 95% reference range of

ND to 46 pg/mL

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

ACTH measurement is not recommended for patients who have been administered cosyntropin (ACTH (1-24)) because the results may be falsely elevated or depressed.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1998;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test

components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all results were generated on EDTA plasma samples.)

Calibration Range: Up to 1,250 pg/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 9 pg/mL

High-dose Hook Effect:
None up to 908,000 pg/mL.

Intraassay Precision (Within-Run):
Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):
Statistics were calculated for samples assayed in *n* different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with two ACTH solutions (2,160 and 7,420 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: Samples with apparent ACTH concentrations of 0 pg/mL showed no cross-reactivity with ACTH (22-39) or Alpha MSH; however, samples with apparent ACTH concentrations of 0 pg/mL showed some cross-reactivity with ACTH (1-18), ACTH (1-24), and ACTH (18-39). ACTH (1-18), ACTH (1-24), and ACTH (18-39) were shown to bind with the antibodies in the assay, leading to an increase in the apparent concentration of ACTH in the sample. The presence of cosyntropin (ACTH (1-24)) in a sample containing ACTH can lead to either an increase or a decrease in concentration, dependent on the amount of ACTH present in the sample. (See "Specificity" table and "Limitations" section.)

The antibody is highly specific for ACTH. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 μ L/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to a commercially available immunoradiometric assay from Nichols on 75 plasma samples from apparently healthy volunteers (28 unspiked samples and 47 spiked samples). (Concentration range: approximately 12 to 1,200 pg/mL. See graph.) By linear regression:

(IML) = 0.99 (Nichols IRMA) + 9.5 pg/mL
r = 0.989

Means:

164 pg/mL (IMMULITE)

156 pg/mL (Nichols IRMA)

References

- 1) Broughton A. Application of adrenocorticotropin assays in a routine clinical laboratory. *Am J Clin Path* 1975;64:618-24.
- 2) Cizza G, Chrousos GP. Adrenocorticotrophic hormone-dependent Cushing's syndrome. *Cancer Treat Res* 1997;89:25-40.
- 3) Demers LM, Whitley RJ. Function of the adrenal cortex. In: Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 3rd ed. Philadelphia: WB Saunders, 1999: 1530-69.
- 4) Donald RA. ACTH and related peptides. *Clin Endocrinol* 1980;12:491-524.
- 5) Findling JW, et al. Selective venous sampling for ACTH in Cushing's syndrome. *Ann Intern Med* 1981;94:647-52.
- 6) Gold EM. The Cushing syndromes: changing views of diagnosis and treatment. *Ann Intern Med* 1979;90:829-44.
- 7) Horrocks PM, London DR. Diagnostic value of 9 AM plasma adrenocorticotrophic hormone concentrations in Cushing's disease. *Br Med J* 1982;285:1302-3.
- 8) Katz FH. Adrenocortical diseases: identification and management. *Postgrad Med* 1979;66(6):52-65.
- 9) Krieger D.

Plasma ACTH and corticosteroids. In: DeGroot L, et al, editors. *Endocrinology*, vol 2. New York: Grune & Stratton, 1979: 1139-56.

- 10) Krieger D, et al. ACTH, β -lipotropin and related peptides in brain, pituitary, and blood. *Recent Prog Horm Res* 1980;36:277-344.
- 11) Nelson DH. The adrenal cortex: physiological function and disease. Philadelphia: W.B. Saunders, 1980.
- 12) Pullan PT, et al. ACTH, LPH and related peptides in the ectopic ACTH syndrome. *Clin Endocrinol* 1980;13:437-45.
- 13) Rees L, Lowry PJ. Adrenocorticotrophin and lipotropin. In: Gray CH, James V, editors. *Hormones in blood*, vol 3. London: Academic Press, 1979: 129-78.
- 14) Schwyzer R. ACTH: a short introductory review. *Ann N Y Acad Sci* 1977;297:3-26.
- 15) Wolfsen A. Aging and the adrenals. In: Korenman SG, editor. *Endocrine aspects of aging*. New York: Elsevier Biomedical, 1982: 55-79.
- 16) Wolfsen A, Odell WD. The dose-response relationship of ACTH and cortisol in Cushing's disease. *Clin Endocrinol* 1980;12:557-68.
- 17) Woodard BH, et al. Adrenocorticotropin production by a mammary carcinoma. *Cancer* 1981;47:1823-7.
- 18) Yalow R. Ectopic ACTH in carcinoma of the lung. In: Muggia FM, Rozenzweig M. editors. *Lung cancer: progress in therapeutic research*. New York: Raven Press, 1978: 209-16.
- 19) Yalow R, et al. Plasma and tumor ACTH in carcinoma of the lung. *Cancer* 1979;44:1789-92.
- 20) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
- 21) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Intraassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	24	2.3	9.6%
2	51	3.1	6.1%
3	101	6.2	6.1%
4	247	9.6	3.9%
5	504	15.6	3.1%
6	999	49	4.9%

Interassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³	<i>n</i>
1	26	2.3	8.8%	16
2	52	4.9	9.4%	16
3	103	8.6	8.3%	16
4	435	39	9.0%	16
5	629	37	5.9%	20
6	1,121	57	5.1%	19

Specificity

Compound/ Amount pg/mL ¹	Initial Concentration pg/mL ²	Change in Concentration pg/mL ³	% Cross- reactivity ⁴
ACTH (1-18)			
500	0	ND	ND
5,000	0	ND	ND
50,000	0	ND	ND
500,000	0	15	0.003%
ACTH (1-24)			
500	0	ND	ND
5,000	0	ND	ND
50,000	0	27	0.05%
50,000	53	38	0.1%
50,000	272	-159	-0.3%
500,000	0	70	0.01%
ACTH (18-39)			
500	0	66	13%
5,000	0	752	15%
ACTH (22-39)			
500	0	ND	ND
5,000	0	ND	ND
50,000	0	ND	ND
500,000	0	ND	ND
Alpha MSH			
500	0	ND	ND
5,000	0	ND	ND
50,000	0	ND	ND
500,000	0	ND	ND

ND: not detectable.⁵


Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	407	—	—
	4 in 8	207	204	101%
	2 in 8	105	102	103%
	1 in 8	53	51	104%
2	8 in 8	627	—	—
	4 in 8	350	314	111%
	2 in 8	167	157	106%
	1 in 8	77	78	99%
3	8 in 8	855	—	—
	4 in 8	393	428	92%
	2 in 8	197	214	92%
	1 in 8	106	107	99%
4	8 in 8	1,136	—	—
	4 in 8	573	568	101%
	2 in 8	284	284	100%
	1 in 8	141	142	99%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2008-2020 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB