

Haptoglobin (Hapt)

Current Revision and Date ^a	Rev. 03, 2020-03
Product Name	Atellica CH Haptoglobin (Hapt) REF 11097643 (300 tests)
Abbreviated Product Name	Atellica CH Hapt
Test Name/ID	Hapt
Systems	Atellica CH Analyzer
Materials Required but Not Provided	Atellica CH LSP CAL REF 11099434
Specimen Types	Serum
Sample Volume	3.4 µL
Measuring Interval	1–340 mg/dL (0.01–3.40 g/L)

^a A vertical bar in the page margin indicates technical content that differs from the previous version.



Intended Use

The Atellica® CH Haptoglobin (Hapt) assay is for *in vitro* diagnostic use in the quantitative measurement of haptoglobin in human serum using the Atellica® CH Analyzer. Such measurements are used as an aid in the diagnosis of hemolytic disorders.¹

Summary and Explanation

The use of polyethylene glycol (PEG) to accelerate the antigen-antibody interaction is described in the work of Hellsing. The enhancing effect of polymers on the precipitin reaction was first reported in 1964 and is dependent upon such factors as the concentration and molecular weight of the added polymer.^{1,2}

Principles of the Procedure

The Atellica CH Haptoglobin (Hapt) assay is based upon the reaction between antibody and haptoglobin in a serum sample. When the antibody reacts with human haptoglobin in dilute solution, insoluble complexes are formed, which can be quantitated by the turbidity they produce. The amount of complex formed varies in proportion to the concentration of haptoglobin in the sample. Turbidity is measured at 340/694 nm. A calibration curve is generated from a set of 6 calibrators of known haptoglobin concentration. Values of unknowns are calculated by interpolating their absorbance readings versus that obtained during calibration.

Reagents

Material Description	Storage	Stability ^a
Atellica CH Hapt	Unopened at 2–8°C	Until expiration date on product
Pack 1 (P1)	Onboard per pack	30 days
Well 1 (W1) Reagent 1 (R1) 18.0 mL Polyethylene glycol (4%); sodium chloride (150 mmol/L); preservative; surfactant		
Well 2 (W2) Empty		
Pack 2 (P2)		
Well 1 (W1) Reagent 2 (R2) 5.3 mL Buffer		
Well 2 (W2) Empty		
Vial (Hapt R2) 1.04 mL Anti-human haptoglobin antibody (goat); sodium azide ($\leq 0.1\%$)	Unopened at 2–8°C	Until expiration date on product

^a Refer to *Storage and Stability*.

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Safety data sheets (SDS) available on siemens.com/healthineers.

Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. (R1 and R2)

CAUTION

This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Contains sodium azide as a preservative. Sodium azide can react with copper or lead plumbing to form explosive metal azides. On disposal, flush reagents with a large volume of water to prevent buildup of azides. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

Note For information about reagent preparation, refer to *Preparing the Reagents* in the *Procedure* section.

Storage and Stability

Unopened reagents are stable until the expiration date on the product when stored at 2–8°C. Do not use products beyond the expiration date printed on the product labeling.

Onboard Stability

Reagents are stable onboard the system for 30 days per pack. Discard reagents at the end of the onboard stability interval. Do not use products beyond the expiration date printed on the product labeling.

Specimen Collection and Handling

Serum is the recommended sample type for this assay.

Collecting the Specimen

- Observe universal precautions when collecting specimens. Handle all specimens as if they are capable of transmitting disease.³
- Follow recommended procedures for collection of diagnostic blood specimens by venipuncture.⁴
- Follow the instructions provided with your specimen collection device for use and processing.⁵
- Allow blood specimens to clot completely before centrifugation.⁶
- Keep tubes capped at all times.⁶

Storing the Specimen

Specimens may be stored frozen for up to 14 days at -20°C.⁷

The handling and storage information provided here is based on data or references maintained by the manufacturer. It is the responsibility of the individual laboratory to use all available references and/or its own studies when establishing alternate stability criteria to meet specific needs.

Transporting the Specimen

Package and label specimens for shipment in compliance with applicable federal and international regulations covering the transport of clinical specimens and etiological agents.

Preparing the Samples

This assay requires 3.4 µL of sample for a single determination. This volume does not include the unusable volume in the sample container or the additional volume required when performing duplicates or other tests on the same sample. For information about determining the minimum required volume, refer to the online help.

Note Do not use specimens with apparent contamination.

Before placing samples on the system, ensure that samples are free of:

- Bubbles or foam.
- Fibrin or other particulate matter.

Note Remove particulates by centrifugation according to CLSI guidance and the collection device manufacturer's recommendations.⁶

Note For a complete list of appropriate sample containers, refer to the online help.

Procedure

Materials Provided

The following materials are provided:

REF	Contents	Number of Tests
11097643	<p>Pack 1 (P1) Well 1 (W1) 18.0 mL of Atellica CH Hapt Reagent 1 Well 2 (W2) Empty</p> <p>Pack 2 (P2) Well 1 (W1) 5.3 mL of Atellica CH Hapt Reagent 2 Well 2 (W2) Empty</p> <p>Vial (Hapt R2) 2 x 1.04 mL of Atellica CH Hapt R2</p>	2 x 150

Materials Required but Not Provided

The following materials are required to perform this assay, but are not provided:

REF	Description
	Atellica CH Analyzer ^a
11099434	Atellica CH LSP CAL (calibrator) <ul style="list-style-type: none"> 1 x 1.0 mL calibrator level 1 CAL 1 1 x 1.0 mL calibrator level 2 CAL 2 1 x 1.0 mL calibrator level 3 CAL 3 1 x 1.0 mL calibrator level 4 CAL 4 1 x 1.0 mL calibrator level 5 CAL 5 1 x 1.0 mL calibrator level 6 CAL 6 Calibrator lot-specific value sheet CAL LOT VAL
	Commercially available quality control materials

^a Additional system fluids are required to operate the system: Atellica CH Diluent, Atellica CH Wash, Atellica CH Conditioner, Atellica CH Cleaner, Atellica CH Reagent Probe Cleaner 1, Atellica CH Reagent Probe Cleaner 2, Atellica CH Reagent Probe Cleaner 4, Atellica CH Lamp Coolant, and Atellica CH Water Bath Additive. For system fluid instructions for use, refer to the Document Library.

Assay Procedure

The system automatically performs the following steps:

1. For serum, dispenses 50 µL of primary sample and 200 µL of Atellica CH Diluent into a dilution cuvette.
2. Dispenses 84 µL of Reagent 1 into a reaction cuvette.
3. Dispenses 3.4 µL of pre-diluted sample into a reaction cuvette.
4. Measures the absorbance after sample addition.
5. Dispenses 16.8 µL of Reagent 2 into a reaction cuvette.
6. Mixes and incubates the mixture at 37°C.
7. Measures the absorbance after Reagent 2 addition.
8. Reports results.

Test Duration: 10 minutes

Preparing the Reagents

Prepare Hapt Reagent 2:

1. Dilute the contents of one Hapt R2 vial of Reagent 2 Concentrate with the buffer contained in well 1 of the P2 pack.
2. Mix well to ensure homogeneity.
3. Pour the solution back into well 1 of the P2 pack and mix well.
4. Carefully rinse the Hapt R2 vial several times with the contents of well 1 and empty the contents back into well 1 of the P2 pack.

Preparing the System

Ensure that the system has sufficient reagent packs loaded in the reagent compartment. For information about loading reagent packs, refer to the online help.

Performing Calibration

For calibration of the Atellica CH Hapt assay, use Atellica CH LSP CAL. Use the calibrators in accordance with the calibrator instructions for use.

Calibration Frequency

Perform a calibration if one or more of the following conditions exist:

- When changing lot numbers of primary reagent packs.
- At the end of the lot calibration interval, for a specified lot of calibrated reagent on the system.
- At the end of the pack calibration interval, for calibrated reagent packs on the system.
- When indicated by quality control results.
- After major maintenance or service, if indicated by quality control results.

At the end of the onboard stability interval, replace the reagent pack on the system with a new reagent pack. Recalibration is required.

Stability Interval	Days
Lot Calibration	30
Pack Calibration	30
Reagent Onboard Stability	30

For information about lot calibration and pack calibration intervals, refer to the online help.

Follow government regulations or accreditation requirements for calibration frequency. Individual laboratory quality control programs and procedures may require more frequent calibration.

Performing Quality Control

For quality control of the Atellica CH Hapt assay, use at least two levels (low and high) of the appropriate quality control material of known analyte concentration. Use the quality control material in accordance with the quality control instructions for use.

For the assigned values, refer to the lot-specific value sheet provided. A satisfactory level of performance is achieved when the analyte values obtained are within the expected control range for the system or within your range, as determined by an appropriate internal laboratory quality control scheme. Follow your laboratory's quality control procedures if the results obtained do not fall within the acceptable limits. For information about entering quality control definitions, refer to the online help.

Follow government regulations or accreditation requirements for quality control frequency. Individual laboratory quality control programs and procedures may require more frequent quality control testing.

Taking Corrective Action

If the quality control results do not fall within the assigned values, do not report results. Perform corrective actions in accordance with established laboratory protocol. For suggested protocol, refer to the online help.

Results

Calculation of Results

The system determines the result using the calculation scheme described in the online help. The system reports results in mg/dL (common units) or g/L (SI units), depending on the units defined when setting up the assay.

Conversion formula: $\text{mg/dL} \times 0.01 = \text{g/L}$

For information about results outside the specified measuring interval, refer to *Measuring Interval*.

Interpretation of Results

Results of this assay should always be interpreted in conjunction with the patient's medical history, clinical presentation, and other findings.

Limitations

The Atellica CH Hapt assay is limited to the detection of haptoglobin in human serum.

Hemoglobin concentrations above 15.6 mg/dL may result in a negative bias in serum specimens. Hemoglobin at 31.3 mg/dL decreases the haptoglobin result in serum at 59 mg/dL (0.59 g/L) by 12%. Hemoglobin at 250 mg/dL decreases the haptoglobin result in serum at 191 mg/dL (1.91 g/L) by 14%.

Do not use hemolyzed samples.

A number of substances cause physiological changes in serum analyte concentrations. A comprehensive discussion of possible interfering substances, their serum concentrations, and their possible physiological involvements is beyond the scope of this document. Consult the listed reference for specific details on known potential interfering substances.⁸

As with any chemical reaction, you must be alert to the possible effect on results of unknown interferences from medications or endogenous substances. The laboratory and physician must evaluate all patient results in light of the total clinical status of the patient.

Expected Values

Reference Interval

A reference interval for healthy adults was established in accordance with CLSI Document EP28-A3c and verified on the Atellica CH Analyzer.⁹

The reference interval for haptoglobin is 40–280 mg/dL (0.40–2.80 g/L) for adults.¹⁰

As with all *in vitro* diagnostic assays, each laboratory should determine its own reference interval for the diagnostic evaluation of patient results. Consider these values as guidance only.⁹

Performance Characteristics

Measuring Interval

The Atellica CH Hapt assay provides results from 1 mg/dL (0.01 g/L) to 340 mg/dL (3.40 g/L). The system flags all values that are outside the specified measuring interval.

Detection Capability

Detection capability was determined in accordance with CLSI Document EP17-A2.¹¹ The assay is designed to have a limit of blank (LoB) \leq limit of detection (LoD) and LoD \leq 1 mg/dL (\leq 0.01 g/L).

The LoD corresponds to the lowest concentration of haptoglobin that can be detected with a probability of 95%. The LoD for the Atellica CH Hapt assay is 1 mg/dL (0.01 g/L), and was determined using 480 determinations, with 240 blank and 240 low level replicates, and a LoB of 0 mg/dL (0.00 g/L).

Assay results obtained at individual laboratories may vary from the data presented.

Precision

Precision was determined in accordance with CLSI Document EP05-A3.¹² Samples were assayed on an Atellica CH Analyzer in duplicate in 2 runs per day for 20 days (N \geq 80 for each sample).

The assay is designed to have the following precision:

- Repeatability: CV \leq 6.0% at 50–110 mg/dL (0.50–1.10 g/L) and CV \leq 5.0% at 111–300 mg/dL (1.11–3.00 g/L)
- Within-Laboratory: CV \leq 8.0% at 50–110 mg/dL (0.50–1.10 g/L) and CV \leq 7.0% at 111–300 mg/dL (1.11–3.00 g/L)

The following results were obtained:

Sample Type	N	Mean mg/dL (g/L)	Repeatability		Within-Lab Precision	
			SD ^a mg/dL (g/L)	CV ^b (%)	SD mg/dL (g/L)	CV (%)
QC	80	80 (0.80)	0.6 (0.006)	0.7	1.8 (0.018)	2.2
QC	80	125 (1.25)	0.8 (0.008)	0.7	2.7 (0.027)	2.2
Serum	80	214 (2.14)	2.5 (0.025)	1.2	6.2 (0.062)	2.9

^a Standard deviation.

^b Coefficient of variation.

Assay results obtained at individual laboratories may vary from the data presented.

Assay Comparison

The Atellica CH Hapt assay is designed to have a correlation coefficient of > 0.950 and a slope of 1.00 ± 0.10 compared to ADVIA® Chemistry HAPT. Assay comparison was determined using the Deming linear regression model in accordance with CLSI Document EP09-A3.¹³ The following results were obtained:

Specimen	Comparative Assay (x)	Regression Equation	Sample Interval	N ^a	r ^b
Serum	ADVIA Chemistry 1800 HAPT	$y = 1.06x + 0$ mg/dL ($y = 1.06x + 0.00$ g/L)	8–310 mg/dL (0.08–3.10 g/L)	105	0.997

^a Number of samples tested.

^b Correlation coefficient.

The agreement of the assay may vary depending on the study design, comparative assay, and sample population. Assay results obtained at individual laboratories may vary from the data presented.

Interferences

Hemolysis, Icterus, and Lipemia (HIL)

The Atellica CH Hapt assay is designed to have $\leq 10\%$ interference from hemoglobin, bilirubin, and lipemia. Interfering substances at the levels indicated in the table below were tested in accordance with CLSI Document EP07-A2 using the Atellica CH Hapt assay.¹⁴

Bias is the difference in the results between the control sample (does not contain the interferent) and the test sample (contains the interferent) expressed in percent. Bias $> 10\%$ is considered interference. Analyte results should not be corrected based on this bias.

Substance	Substance Test Concentration Common Units (SI Units)	Analyte Concentration mg/dL (g/L)	Percent Bias
Hemoglobin	15.6 mg/dL (0.010 mmol/L)	59 (0.59)	0
	125 mg/dL (0.078 mmol/L)	191 (1.91)	1
Bilirubin, conjugated	25 mg/dL (428 μ mol/L)	56 (0.56)	2
	25 mg/dL (428 μ mol/L)	194 (1.94)	0
Bilirubin, unconjugated	25 mg/dL (428 μ mol/L)	57 (0.57)	0
	25 mg/dL (428 μ mol/L)	194 (1.94)	0
Lipemia (Intralipid®)	250 mg/dL (2.8 mmol/L)	58 (0.58)	0
	250 mg/dL (2.8 mmol/L)	195 (1.95)	0
Lipemia (Triglycerides)	1000 mg/dL (11.3 mmol/L)	55 (0.55)	0
	1000 mg/dL (11.3 mmol/L)	186 (1.86)	2

Assay results obtained at individual laboratories may vary from the data presented.

High-Dose Hook Effect

No prozone effect (high-dose hook effect) was observed for haptoglobin concentrations up to at least 600 mg/dL (6.00 g/L).

Standardization

The Atellica CH Hapt assay is traceable to the Institute for Reference Materials and Measurements (IRMM) CRM-470 from IFCC.

Assigned values for calibrators are traceable to this standardization.¹⁵

Technical Assistance


















For customer support, contact your local technical support provider or distributor.
siemens.com/healthineers



















References









1. Hellsing K. Influence of polymers on the antigen-antibody reaction in a continuous flow system. In: Hamm JD, ed. *Automated Immunoprecipitin Reactions*. Tarrytown, NY: Technicon Instruments Corp; 1972:17-20.
2. Hellsing K, Laurent TC. The influence of dextran on the precipitin reaction. *Acta Chem Scand*. 1964;18(5):1303-1304.
3. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
4. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
5. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
6. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
7. Wu AHB, ed. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Saunders; 2006:512.
8. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Washington, DC: AACC Press; 2000.
9. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
10. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*. 3rd ed. Washington, DC: AACC Press; 2007:469.
11. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
12. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
13. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
14. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.
15. Data on file at Siemens Healthcare Diagnostics.

Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Symbol Title and Description
	Consult instructions for use
 Rev. 01	Version of instructions for use
 siemens.com/healthcare	Internet URL address to access the electronic instructions for use
 siemens.com/document-library	
Rev. 	Revision
	Caution Consult instructions for use or accompanying documents for cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device.
	Biological risks Potential biological risks are associated with the medical device.
	Corrosive
	Dangerous to environment
	Irritant Oral, dermal, or inhalation hazard
	Inhalation hazard Respiratory or internal health
	Flammable Flammable to extremely flammable
	Oxidizing
	Explosive
	Toxic
	Compressed gas
	Keep away from sunlight Prevent exposure to sunlight and heat.

Symbol	Symbol Title and Description
	Up Store in an upright position.
	Do not freeze
	Temperature limit Upper and lower limits of temperature indicators are adjacent to the upper and lower horizontal lines.
	Handheld barcode scanner
	<i>In vitro</i> diagnostic medical device
	Contains sufficient for <n> tests Total number of IVD tests the system can perform with the IVD kit reagents appears adjacent to the symbol.
RxOnly	Prescription device (US only) Applies only to United States-registered IVD assays. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.
	Mixing of substances Mix product before use.
	Reconstitute and mix lyophilized product before use.
	Target
	Interval
	Legal Manufacturer
	Authorized Representative in the European Community
	Use-by date Use by the designated date.
	Batch code
	Catalog number
	Recycle
	Printed with soy ink
	CE Mark


Symbol	Symbol Title and Description
	CE Mark with notified body ID number Notified body ID number can vary.
YYYY-MM-DD	Date format (year-month-day)
	Variable hexadecimal number that ensures the Master Curve and Calibrator definition values entered are valid.
	Common Units
	International System of Units
	Material
	Unique material identification number
	Name of control
	Type of control

Legal Information

Atellica and ADVIA are trademarks of Siemens Healthcare Diagnostics.

All other trademarks are the property of their respective owners.

© 2017–2020 Siemens Healthcare Diagnostics. All rights reserved.

 Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthineers

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthineers