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1. Remme WJ, Swedberg K; Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Guidelines for the diagnosis and treatment of chronic heart failure. Eur Heart J. 2001;22(17):1527-1560.
2. Richards AM, Nicholls MG, Yandle TG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: new neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation. 1998;97(19):1921-1929.
3. de Bold AJ. Atrial natriuretic factor: a hormone produced by the heart. Science. 1985;230(4727):767-770.
4. Valli N, Gobinet A, Bordernave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med. 1999(5);134:437-444.
5. de Bold AJ, Boerenstein HB, Varesse AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial myocardial extract in rats. J Am Soc Nephrol. 2001;12(2):403-409. Reprinted in: Life Sci. 1981;28:89-94.
6. Epstein M, Loutzenhiser R, Friedland E, et al. Relationship of increased plasma atrial natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest. 1987;79(3):738-745.
7. Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J.

1999;20(19):1374-1375.

8. Hunt PJ, Richards AM, Nicholls MG, et al. Immunoreactive amino-terminal pro-brain natriuretic peptide (NT-proBNP): A new marker of cardiac impairment. *Clin Endocrinol*. 1997;47(3):287-296.

9. Talwar S, Squire IB, Davies JE, et al. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. *Eur Heart J*. 1999;20(23):1736-1744.

10. Darbar D, Davidson NC, Gillespie N, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. *Am J Cardiol*. 1996;78(3):284-287.

11. McDonagh TA, Robb SD, Murdoch DR, et al. Biochemical detection of left-ventricular systolic dysfunction. *Lancet*. 1998;351(9095):9-13.

12. Fisher C, Berry C, Blue L, et al. N-terminal pro B type natriuretic peptide, but not the new putative cardiac hormone relaxin, predicts prognosis in patients with chronic heart failure. *Heart*. 2003;89(8):879-881.

13. James SK, Lindahl B, Siegbahn A, et al. N-terminal pro-brain natriuretic peptide and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease: a Global Utilization of Strategies To Open occluded arteries (GUSTO)-IV substudy. *GUSTO IV Substudy*. *Circulation*. 2003;108(3):275-281.

14. Talwar S, Squire IB, Downie PF, Davies JE, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotrophin 1 are raised in unstable angina. *Heart*. 2000;84(4):421-424.

15. Jernberg T, Stridsberg M, Venge P, Lindahl B. N-terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. *J Am Coll Cardiol*. 2002;40(3):437-445.

16. Zeller M, Collin Y, Laurent Y, et al. N-terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. *Cardiology*. 2004;102(1):37-40.

17. Richards AM, Nicholls MG, Espiner E, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. *Circulation* 2003;107(22):2786-2792.

18. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections*; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.

19. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture*; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.

20. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection*; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.

21. Clinical and Laboratory Standards Institute. *Procedures for*

the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.

22. Clinical and Laboratory Standards Institute. *Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves*; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Institute; 2011. CLSI Document EP24-A2.

23. Kondratovich M. Matched receiver operating characteristic (ROC) analysis and propensity scores. Paper presented at: US Food and Drug Administration Joint Statistical Meetings, Biopharmaceutical Section; 2003.

24. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45(7):942-956.

25. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38(9):1737-1742.

26. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory*; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.

27. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry*; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.

28. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures*; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

29. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures*; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.

30. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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住所：Glyn Rhonwy, Llanberis, Gwynedd, LL55 4EL, United Kingdom

生产地址：333 Coney Street, East Walpole, MA, 02032, United States

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Siemens Healthcare Diagnostics Products Limited