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1) Becker W, Reimers C, Brner W. Immunological criteria for

the differential diagnosis of thyroid disease. *Therapiewoche* 1985;35:1167-76.

2) Czarnocka B, Ruf J, Ferrand M, et al. Purification of the human thyroid peroxidase and its identification as the micro-somal antigen involved in autoimmune thyroid diseases. *FEBS Lett* 1985;190:147-52.

3) Czarnocka B, Ruf J, Ferrand M, et al. Interaction of highly purified thyroid peroxidase with anti-microsomal antibodies in autoimmune thyroid diseases. *J Endocrinol Invest* 1986;9:135-8.

4) Falk U, et al. Radioimmunoassay for serum thyroglobulin designed for early detection of metastases and recurrences in the follow-up of patients with differentiated thyroid carcinoma. *J Clin Chem Clin Biochem* 1984;22:661-70.

5) Feldt-Rasmussen U. Serum thyroglobulin and thyroglobulin autoantibodies in thyroid diseases. *Allergy* 1983;38:369-87.

6) Gons MH, et al. Concentration of plasma thyroglobulin and urinary excretion of iodinated material in the diagnosis of thyroid disorders in congenital hypothyroidism. *Acta Endocrinol* 1983;104:27-34.

7) Hashimoto H. Understanding the lymphatic changes of the thyroid. *Arch Klin Chir* 1912;97:219-48.

8) Horster FA. The importance of MAK, TAK, TRAK and thyroglobulin in the diagnosis of thyroid disease. *Internist* 1988;29:538-40.

9) Ket JL, et al. Serum thyroglobulin levels: physiological decrease in infancy and the absence in athyroidism. *J Clin Endocrinol Metab* 1981;58:1301-3.

10) Larsen PR, Ingbar SH. The thyroid gland. In: Wilson JD, Foster DW, editors. *Williams textbook of endocrinology*. 8th ed. Philadelphia: WB Saunders, 1992: 357-487.

11) Libert F, Ruel J, Ludgate M, et al. Thyroperoxidase, an autoantigen with a mosaic structure made of nuclear and mitochondrial gene modules. *EMBO J* 1987;6:4193-6.

12) McKenzie JM, Zakarija M, Sato A. Humoral immunity in Graves' disease. *Clin Endocrinol Metab* 1978;7:1-45.

13) Ruf J, Czarnocka B, Ferrand M, et al. Thyroid peroxidase is the organ-specific 'microsomal' antigen involved in thyroid autoimmunity. *Acta Endocrinol (Copenh)* 1987; Suppl. 281:49-55.

14) Ruf J, Czarnocka B, Ferrand M, et al. Relationship between immunological structure and biochemical properties of human thyroid peroxidase. *Endocrinol* 1988;125(3):1211-8.

15) Schatz H, Löbig The diagnostic and prognostic importance of antibodies against thyroid microsomes and thyroglobulin. *Akt Endokr Stoffw* 1989;10:146-53.

16) Scherbaum WA, Berg PA. The importance of autoantibodies in the diagnosis of endocrine disease. *Dtsch Med Wschr* 1981;106:308-13.

17) Scherbaum WA, Berg PA. Thyroid autoimmune disease – new aspects in development and diagnosis. *Dtsch Med Wschr* 1984;109:1574-81.

18) Spencer CA, Wang C-C. Thyroglobulin measurement: techniques, clinical benefits, and pitfalls. *Endocrinol Metab Clin North Am* 1995;24:841-63.

19) Spencer CA, Takeuchi M, Kazarosyan M. Current status and performance goals for serum thyroglobulin assays. *Clin Chem* 1996;42:164-73.

20) Spencer CA. Recoveries cannot be used to authenticate thyroglobulin (Tg) measurements when sera contain Tg autoantibodies. *Clin Chem* 1996;42:661-3.

21) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994.

22) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

23) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8.

24) Clinical and Laboratory Standard Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline – Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.

25) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

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Products Limited