Dimension Vista® System

Flex® reagent cartridge

TRIG

See shaded sections: Updated information from 2019-05 version.


Issue Date 2020-11

Triglycerides

Intended Use: The TRIG method is an *in vitro* diagnostic test for the quantitative measurement of triglycerides in human serum and plasma on the Dimension Vista® System. Measurements obtained are used in the diagnosis and treatment of patients with diabetes mellitus, nephrosis, liver obstruction, other diseases involving lipid metabolism, or various endocrine disorders.

Summary: Triglycerides are water-insoluble lipids consisting of three fatty acids linked to one glycerol molecule. Triglycerides are transported in the blood as core constituents of all lipoproteins, but the greatest concentration of these molecules is carried in the triglycerides-rich chylomicrons and very low density lipoproteins (VLDL).¹ Through the action of lipases and bile acids, triglycerides are hydrolyzed into glycerol and fatty acids which are absorbed by adipose tissue for storage or by other tissues requiring a source of energy. A peak concentration of chylomicron-associated triglycerides occurs within 3–6 hours after ingestion of a fat-rich meal; however, the rate of absorption of fats is highly variable, depending on the individual and dietary composition of the fat. After absorption, triglycerides are resynthesized in the epithelial cells and combined with cholesterol and a number of apolipoproteins to form chylomicrons.²

Principles of Procedure: The triglycerides method is based on an enzymatic procedure in which combinations of enzymes are employed for the measurement of serum or plasma triglycerides. The sample is incubated with lipoprotein lipase (LPL) enzyme reagent that converts triglycerides into free glycerol and fatty acids. Glycerol kinase (GK) catalyzes the phosphorylation of glycerol by adenosine-5-triphosphate (ATP) to glycerol-3-phosphate. Glycerol-3-phosphate-oxidase (GPO) oxidizes glycerol-3-phosphate to dihydroxyacetone phosphate and hydrogen peroxide (H₂O₂). The catalytic action of peroxidase (POD) forms quinoneimine from H₂O₂, aminoantipyrine and 4-chlorophenol. The change in absorbance due to the formation of quinoneimine is directly proportional to the total amount of glycerol and its precursors in the sample and is measured using a bichromatic (510, 700 nm) endpoint technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^{b,c}
1–12	Liquid	Lipoprotein Lipase	7.5 KU/L
		ATP	3 mmol/L
		Glycerol kinase	0.5 KU/L
		Glycerol-3-phosphate-oxidase	2.2 KU/L
		4-aminoantipyrine	0.75 mmol/L
		4-chlorophenol	6 mmol/L
		Peroxidase	5 KU/L
		Mg ²⁺	22.5 mmol/L
		Buffer pH 7.2	50 mmol/L

- a. Wells are numbered consecutively from the wide end of the cartridge.
 b. This represents the nominal value in final reaction mixture.
 c. Contains bovine serum albumin.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Contains sodium azide (< 0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2–8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 7 days for wells 1–12

Specimen Collection and Handling: Recommended specimen types: serum, lithium heparin plasma.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.^{3,4} Blood collection tubes containing glycerol lubricated stoppers should be avoided since they will cause erroneously elevated results.

Follow the instructions provided with your specimen collection device for use and processing.⁷

Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁵ Samples should be kept at 4 °C and analyzed within 24 hours. Samples can be stored refrigerated for up to 7 days. For longer storage, samples may be frozen at -20 °C for up to 3 months or many years at -70 °C.⁶

Procedure

TRIG Flex® reagent cartridge, Cat. No. K2069

CHEM 2 CAL, Cat. No. KC120
 Quality Control Materials

Materials Provided

Materials Required But Not Provided

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

Test Conditions

Sample Volume (delivered to the cuvette)	1.6 µL
Reagent Volume	55 µL
Temperature	37.0 °C
Reaction time	5.6 minutes
Wavelength	510 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Calibration Material	CHEM 2 CAL, Cat. No. KC120
Calibration Scheme	2 levels, n=5
Units	mg/dL [mmol/L] ^a (mg/dL x 0.0113) = [mmol/L]
Typical Calibration Levels	Level 1 (Calibrator A): 0 mg/dL [0.00 mmol/L] Level 2 (Calibrator B): 1064 mg/dL [12.02 mmol/L]
Calibration Frequency	Every 90 days for any one lot Calibration interval may be extended based on acceptable verification of calibration.
A new calibration is required:	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations

d. Système International d'Unités [SI units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known triglyceride concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of triglyceride in mg/dL [mmol/L] using the calculation scheme described in your Dimension Vista® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 2–1000 mg/dL [0.02–11.30 mmol/L]^a

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 1000 mg/dL [11.30 mmol/L] should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water to obtain results within reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 13 µL (dilution factor = 4) for serum/plasma. Refer to your Dimension Vista® Operator's Guide.

- Samples with results less than 2 mg/dL [0.02 mmol/L] will be reported as "less than 2 mg/dL [0.02 mmol/L]" by the instrument.

e. Triglycerides-based samples (not glycerol samples) should be used to verify the assay range.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in triglyceride results. Refer to your Dimension Vista® Operator's Guide for the meaning of report slip flags and comments. Any report slip containing flags and/or comments should be held for follow-up and addressed according to your laboratory's procedure manual.

Venipuncture should occur prior to N-Acetyl Cysteine or Metamizole (Sulpyrine) administration due to the potential for falsely depressed results.

In the presence of etamsylate at 2 mg/dL [76 µmol/L], falsely depressed results > 10% for triglycerides may be observed.

Interfering Substances:

- Ascorbic acid at a concentration of 5 mg/dL [227 µmol/L] decreases triglyceride results by 11.8% at a triglyceride concentration of 180 mg/dL [2.03 mmol/L].
- Small amounts of free glycerol may be found in blood samples from healthy individuals due to natural lipolysis. The concentration of free glycerol may be increased by stress, disease states or administration of intravenous infusates.⁸ Free glycerol or other polyols may cause a positive interference.
- ‡ Glycerol-based quality control products should not be used with this method.

Expected Values:

The National Cholesterol Education Program Adult Treatment Panel III (NCEP- ATP III)⁹ provides the following categories of triglycerides concentrations:

Category	Serum Triglycerides	
	mg/dL	[mmol/L]
Normal	< 150	< 1.70
Borderline high	150 – 199	1.70 – 2.25
High	200 – 499	2.26 – 5.64
Very High	≥ 500	≥ 5.65

Each laboratory should establish its own expected values for TRIG as performed on the Dimension Vista® System.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n=5 replicates at the following analyte concentrations are:

TRIG concentration	Acceptable SD Maximum
70 mg/dL [0.79 mmol/L]	9 mg/dL [0.10 mmol/L]
375 mg/dL [4.24 mmol/L]	26 mg/dL [0.29 mmol/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Material	Precision ^{10,f}		
	Mean mg/dL [mmol/L]	Standard Deviation (% CV)	
		Repeatability	Within-Lab
Multiquel® Control			
Level 1	68 [0.77]	2 [0.02] (3)	3 [0.03] (4)
Level 2	384 [4.34]	6 [0.07] (2)	9 [0.11] (2)

- f. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.
Multiquel® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Comparative Method	Method Comparison ¹¹			
	Regression Statistics ^g			n
	Slope	Intercept mg/dL [mmol/L]	Correlation Coefficient	
Dimension® Rxl System Serum/plasma	1.04	7.8 [0.09]	0.991	113 ^h

- g. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.
h. The range of serum/plasma triglyceride values in the correlation study was 50–968 mg/dL [0.57–10.94 mmol/L].

Specificity

HIL Interference
The TRIG method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹² Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Triglyceride mg/dL [mmol/L]	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer)		
	1000 mg/dL [0.62 mmol/L]	180 [2.03]	<10
Bilirubin (unconjugated)	5 mg/dL [86 µmol/L]	180 [2.03]	<10
	10 mg/dL [171 µmol/L]		11
	20 mg/dL [342 µmol/L]		20
	60 mg/dL [1026 µmol/L]		24
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	180 [2.03]	<10

* Analyte results should not be corrected based on this bias.

Non Interfering Substances

The following substances do not interfere with the TRIG method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at triglycerides concentration of 180 mg/dL [2.03 mmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	5 ng/mL	6.15 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25,000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

‡ The following substances do not interfere with the TRIG method when present in serum and plasma at the concentrations indicated.

Substance	Test Concentration	SI Units
Dextran 75	2500 mg/dL	333 mmol/L
Sodium Heparin	8000 U/L	8000 U/L
Lithium (lithium chloride)	3.5 mg/dL	5.07 mmol/L
Rheumatoid factor	500 IU/mL	500 IU/mL

Analytical Sensitivity: 2 mg/dL [0.02 mmol/L]

The analytical sensitivity represents the lowest concentration of TRIG that can be distinguished from zero. This sensitivity is defined as the mean value (n=20) plus two standard deviations of the Level 1 (0 mg/dL [0 mmol/L]) CHEM 2 Calibrator.

‡ The Dimension Vista® TRIG method (REF K2069) and the Dimension® TGL (REF DF69A) method utilize the same reagents under equivalent reaction conditions. Interference testing was performed using the Dimension® TGL (REF DF69A) method and the results are representative of both methods.

Symbols Key: See Adjacent Panel.

Bibliography: See Adjacent Panel.

Dimension Vista®, Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

1. Rifai N, Russell Warnick G, Dominiczak MH. Handbook of Lipoprotein Testing, AACC Press, Washington, DC 1997: p. 115.
2. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders Co., Philadelphia, PA 1994: p 1017.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Tietz NW. Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1986, pp. 52-53 (techniques and procedures to minimize laboratory infections) and pp 478-497 (specimen collection and storage recommendations).
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Tietz NM. Clinical Guide to Laboratory Tests, 2nd Edition, W.B. Saunders Co., Philadelphia, PA 1990: pp 554.
7. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
8. Kaplan AK, Amadeo JP. Clinical Chemistry: Theory, Analysis, and Correlation, 3rd Edition, Mosby, St. Louis Missouri, 1996, pp. 680.
9. National Cholesterol Education Program: Third Report of the Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). Journal of American Medical Association, May 16, 2001: p 18.
10. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
11. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
12. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFSS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare