Dimension® clinical chemistry system

Flex® reagent cartridge

TGL

See shaded sections: Updated information from 2019-04 version.


Issue Date 2020-11

Triglycerides

Intended Use: The TGL method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of triglycerides in human serum and plasma. Measurements obtained are used in the diagnosis and treatment of patients with diabetes mellitus, nephrosis, liver obstruction, other diseases involving lipid metabolism, or various endocrine disorders.

Summary: Triglycerides are water-insoluble lipids consisting of three fatty acids linked to one glycerol molecule. Triglycerides are transported in the blood as core constituents of all lipoproteins, but the greatest concentration of these molecules is carried in the triglycerides-rich chylomicrons and very low density lipoproteins (VLDL).¹ Through the action of lipases and bile acids, triglycerides are hydrolyzed into glycerol and fatty acids which are absorbed by adipose tissue for storage or by other tissues requiring a source of energy. A peak concentration of chylomicron-associated triglycerides occurs within 3 – 6 hours after ingestion of a fat-rich meal; however, the rate of absorption of fats is highly variable, depending on the individual and dietary composition of the fat. After absorption, triglycerides are resynthesized in the epithelial cells and combined with cholesterol and a number of apolipoproteins to form chylomicrons.²

Principles of Procedure: The triglycerides method is based on an enzymatic procedure in which a combination of enzymes are employed for the measurement of serum or plasma triglycerides. The sample is incubated with lipoprotein lipase (LPL) enzyme reagent that converts triglycerides into free glycerol and fatty acids. Glycerol kinase (GK) catalyzes the phosphorylation of glycerol by adenosine-5-triphosphate (ATP) to glycerol-3-phosphate. Glycerol-3-phosphate-oxidase oxidizes glycerol-3-phosphate to dihydroxyacetone phosphate and hydrogen peroxide (H₂O₂). The catalytic action of peroxidase (POD) forms quinoneimine from H₂O₂, aminoantipyrine and 4-chlorophenol. The change in absorbance due to the formation of quinoneimine is directly proportional to the total amount of glycerol and its precursors in the sample and is measured using a bichromatic (510, 700 nm) endpoint technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^{b,c}
1 – 6	Liquid	Lipoprotein Lipase	7.5 KU/L
		ATP	3 mmol/L
		Glycerol kinase	0.5 KU/L
		Glycerol-3-Phosphate-oxidase	2.2 KU/L
		4-aminoantipyrine	0.75 mmol/L
		4-chlorophenol	6 mmol/L
		Peroxidase	5 KU/L
		Mg ²⁺	22.5 mmol/L
		Buffer pH 7.2	50 mmol/L

a. Wells are numbered consecutively from the wide end of the cartridge.

b. This represents the nominal value in final reaction mixture.

c. Contains bovine serum albumin.

Risk and Safety

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Contains sodium azide (< 0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 10 days for wells 1 – 6

Specimen Collection and Handling: Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.³

Follow the instructions provided with your specimen collection device for use and processing.⁴

Blood collection tubes containing glycerol lubricated stoppers should be avoided since they will cause erroneously elevated result.

Corvac®, and SST®, collection tubes do not affect the TGL method.

Separated specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.⁵

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.

SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

TGL Flex® reagent cartridge, Cat. No. DF69A

Materials Required But Not Provided

CHEM II Calibrator, Cat. No. DC20

Quality Control Materials

Test Steps

Sampling,⁴ reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume.

Test Conditions

Sample Size	4 µL
Reagent 1 Volume	133 µL
Temperature	37 °C ± 0.1 °C
Wavelengths	510 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range ^a	15 – 1000 mg/dL [0.17 – 11.3 mmol/L] ^f
Calibration Material	CHEM II Calibrator, Cat. No. DC20
Calibration Scheme	3 Levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.0113) = [mmol/L]
Typical Calibration Levels	120, 240, 485 mg/dL [1.37, 2.74, 5.54 mmol/L]
Calibration Frequency	Every 90 days for any one lot.
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -2.6 C ₁ 1.5

e. Triglycerides-based samples (not glycerol samples) should be used to verify the assay range.

f. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze at least two levels of a quality control material with known triglycerides concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of triglycerides in mg/dL [mmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 15 – 1000 mg/dL [0.17 – 11.3 mmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 1000 mg/dL [11.3 mmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water to obtain a result within the assay range. Enter the dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): If using the auto-dilution feature, results above 1000 mg/dL [11.3 mmol/L] will automatically be repeated (for serum, plasma).

Results less than 15 mg/dL [0.17 mmol/L] should be reported as "less than 15 mg/dL [0.17 mmol/L]".

Limitations of Procedure

The instrument reporting system contains error messages to warn the user of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

Venipuncture should occur prior to N-Acetyl Cysteine or Metamizole (Sulpyrine) administration due to the potential for falsely depressed results.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
100 mg/dL [1.13 mmol/L]	>5 mg/dL [0.06 mmol/L]
400 mg/dL [4.53 mmol/L]	>16 mg/dL [0.18 mmol/L]

In the presence of etamsylate at 2 mg/dL [76 µmol/L], falsely depressed results >10% for triglycerides may be observed.

Interfering Substances

Small amounts of free glycerol may be found in blood samples from healthy individuals due to natural lipolysis. The concentration of free glycerol may be increased by stress, disease states or administration of intravenous infusates.⁶ Free glycerol or other polyols may cause a positive interference.

Glycerol-based quality control products should not be used with this method.

Hemoglobin (hemolysate) of 500 mg/dL [0.31 mmol/L] (monomer) will increase a triglycerides result of 155 mg/dL [1.75 mmol/L] by 12%.

Bilirubin (unconjugated) of 20 mg/dL [342 µmol/L] will increase a triglycerides result of 156 mg/dL [1.76 mmol/L] by 11%.

Expected Values

The National Cholesterol Education Program Adult Treatment Panel III (NCEP- ATP III)⁷ provides the following categories of triglycerides concentrations:

Category	Serum mg/dL	Triglycerides [mmol/L]
Normal	< 150	[< 1.70]
Borderline high	150 – 199	[1.70 – 2.25]
High	200 – 499	[2.26 – 5.64]
Very High	≥ 500	[≥ 5.65]

Each laboratory should establish its own reference interval for triglycerides as performed on the Dimension® system.

Specific Performance Characteristics^a

Material	Precision ^{b,i}		
	Mean mg/dL [mmol/L]	Standard Deviation (% CV)	
		Within-run	Between-day
Multiqua [®]			
Level 2	132 [1.49]	0.7 [0.01] (0.5)	1.1 [0.01] (0.8)
Level 3	216 [2.44]	0.9 [0.01] (0.4)	1.5 [0.02] (0.7)
Plasma Pool	364 [4.11]	1.4 [0.02] (0.4)	3.6 [0.04] (1.0)
Serum Pool	425 [4.80]	1.5 [0.02] (0.4)	5.5 [0.06] (1.3)

- g. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed on the Dimension® Rxl system. (Refer to your Dimension® Operator's Guide).
- h. Reproducibility testing was done in accordance with the NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999).
- i. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Multiqua[®] is a registered trademark of Bio-Rad Laboratories, Irvine, CA, 92618-2017, USA.

Method Comparison

Regression Statistics^j

Comparative Method	Slope	Intercept	Correlation Coefficient	n ^k
Dimension® TRIG	1.01	-4.2 [-0.5]	0.999	230

- j. Model equation for regression statistics is: Result of Dimension® analyzer = [Slope x Comparative method result] + Intercept.
- k. The range of triglycerides in the correlation study was 30 – 906 mg/dL [0.34 – 10.24 mmol/L]. Triglycerides results greater than 500 mg/dL were obtained from diluted samples when analyzed by the Dimension® TRIG method (Cat. No. DF69). All triglycerides results analyzed by the Dimension® TGL method (Cat. No. DF69A) were obtained from undiluted samples.

Specificity

HIL Interference

The TGL method was evaluated for interference from hemolysis and icterus according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	TGL Concentration mg/dL [mmol/L]	Bias % ^l
Hemoglobin (hemolysate)	300 mg/dL [0.19 mmol/L] (monomer)	155 [1.75]	<10
Bilirubin (unconjugated)	5 mg/dL [86 µmol/L]	156 [1.76]	<10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	200 [2.26]	<10

l. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the TGL method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at triglycerides concentrations of 200 mg/dL [2.26 mmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Albumin	6.8 g/dL	68 g/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic Acid	3 mg/dL	170.3 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol	500 mg/dL	13.0 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Sodium Heparin	8000 U/L	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium (lithium chloride)	3.5 mg/dL	5.07 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein: Total	3.6 g/dL	36 g/L
Protein: Total	11.8 g/dL	118 g/L
Rheumatoid factors	500 IU/mL	500 IU/mL
Salicylic Acid	50 mg/dL	3.62 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 15 mg/dL [0.17 mmol/L]

The analytical sensitivity represents the lowest activity of TGL that can be reported.

Symbols Key: See Adjacent Panel

Bibliography: See Adjacent Panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Rifai N, Russell Warnick G, Dominiczak MH. Handbook of Lipoprotein Testing. AACC Press, Washington, DC 1997: p. 115.
2. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders Co., Philadelphia, PA 1994: p 1017.
3. Tietz NW. Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1986, pp 478-497 (specimen collection and storage recommendations).
4. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Kaplan AK, Amadeo JP. Clinical Chemistry: Theory, Analysis, and Correlation, 3rd Edition, Mosby, St. Louis Missouri, 1996, pp. 680.
7. National Cholesterol Education Program: Third Report of the Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). Journal of American Medical Association, May 16, 2001: p 18.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFSS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare