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1. Burtis CA., Ashwood ER., Bruns DE, eds. Tietz Textbook of Clinical Chemistry and molecular Diagnostics.4th ed. Philadelphia: Elsevier Inc, 2006 pp. 2053-2087.
2. Demers LM, and Spencer C.A., Eds., Standards of Laboratory Practice, Laboratory support for the diagnosis& monitoring of thyroid disease, National Academy of Clinical Biochemistry, 2002.
3. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al. Luminescent oxygen channeling assay (LOCI™):sensitive, broadly applicable homogeneous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
4. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al Luminescent oxygen channeling immunoassay:Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad. Sci. USA, Vol 91,pp 5426-5430, June 1994 Biochemistry.
5. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard — Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Tietz NW. Clinical Guide to Laboratory Tests. 4th ed., 2006, Saunders WB Co, Philadelphia, PA., pp 1040
7. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007: p 867.
8. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
9. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection — Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
10. Kricka LJ. Human anti-animal antibody interference in immunological assays. ClinChem 1999; 45 (7):942-956.
11. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine

kinase MB by using F(ab'), conjugate and polyclonal mouse IgG.
Clin Chem 1992; 38:1737-1742.

12. Clinical and Laboratory Standards Institute/NCCLS. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline — Second Edition. CLSI/NCCLS document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.

13. Clinical and Laboratory Standards Institute (formerly NCCLS). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. CLSI 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.

14. Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC: AACC Press; 2005

15. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline — Second Edition. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

16. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

17. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline — Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.

18. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays, Anal Biochem 1982;127:308-315.

19. Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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Siemens Healthcare Diagnostics Inc.

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©2008 Siemens Healthcare Diagnostics