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安全数据 (MSDS/SDS) 见 www.siemens.com/diagnostics

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1. Tucker JF, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. Annals of Emerg Med 1994;4:704-708.
2. Woo J, et al. Is myoglobin useful in the diagnosis of acute myocardial infarction in the emergency department setting? Am J Clin. Pathol. 1995;103:No.6;725-729.
3. Brogan GX, et al. Evaluation of a new rapid quantitative immunoassay for serum myoglobin versus ck-mb for ruling out acute myocardial infarction in the emergency department. Anals of Emerg. Med. 1994;4:665-671.
4. Carraro P, et al. A new immunoassay for the measurement of myoglobin in serum. Journal of Clinical Laboratory Analysis 1994;8:70-75.
5. de Winter RJ, et al. Value of Myoglobin, Troponin T, CK-MB mass in ruling out an acute myocardial infarction in the emergency room. Circulation 1995; 92(12): 3401-3407.
6. Newby LK, Storow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain units: The chest pain evaluation by Creatine Kinase-MB, myoglobin, and Troponin I (CHECKMATE) Study. Circulation 2001;103:1832-1837.
7. Montague C, Kircher T. Myoglobin in the early evaluation of acute chest pain. A.J.C.P. 1995;4:472-476.
8. Birkmeyer RC., Diaco, Hutson DK., et al. Applications of novel chromium dioxide magnetic particles to immunoassay development. Clin Chem. 1987;33:1543-1547.
9. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
10. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
11. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
12. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition, AACC Press, Washington, DC 1997; PP. 4-501.
13. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for Creatine Kinase MB by using F(ab') 2 conjugate and polyclonal mouse IgG. Clin Chem 1992;38:1737-1742.
14. Kricka LJ. Human Anti-animal antibody interferences in immunological assays. Clin Chem 1999;45(7):942-956.
15. Kontos MC, et al. Use of the combination of myoglobin of acute

myocardial infarction. Amer Journal of Emer Med 1997;1:14-19.

16. Ryan TJ, Anderson JL, Antman EM, Braniff BA, Brooks NH, Califf RM, Hillis LD, Hiratzka LF, Rappaport E, Riegel BJ, Russell RO, Smith EE III, Weaver WD. ACC/AHA guidelines for management of patients with acute myocardial infarction: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Acute Myocardial Infarction). J Am Coll Cardiol 1996; 28:1356-1358.

17. Braunwald E, Antman EM, Beasley JW, Califf RM, Chaitlin MD, Hochman JS, Jones RH, Kereiakes D, Kupersmith J, Levin TN, Pepine CJ, Schaeffer JW, Smith EE III, Steward DE, Theroux P. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction: executive summary and recommendations: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Patients with Unstable Angina). Circulation 2000;102:1193-1209.

18. World Health Organization criteria for the diagnosis of acute myocardial infarction. Proposal for the multinational monitoring of trends and determinants in cardiovascular disease. Geneva Cardiovascular Disease Unit of World Health Organization;1981.

19. Wu HB, Apple FS, Warshaw MM. Recommendations for the use of cardiac markers in coronary artery diseases. National Academy of Clinical Biochemistry Standard of Laboratory Practice 1999;(5):17-18.

20. Thygesen K, Alpert JS, et al. Myocardial Infarction Redefined – A Consensus Document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology 2000;36(3):959-969.

21. Ryall RG, Story CJ, DR. Reappraisal of the causes of “hook effect” in two-site immunoradiometric assays. Anal Biochem 1982;127:308.

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Siemens Healthcare Diagnostics Inc.