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1. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 547-548, (clinical significance).
2. Whicher J, Spence C. When is serum albumin worth measuring? Ann Clin Biochem 1987; 24:572–80.
3. Carter P. Ultramicroestimation of human serum albumin: binding cationic dye 5,5'-dibromo-o cresolsulphonphthalein, Microchem J 1970; 15:531-539.
4. Louderback A, Measley EH, Taylor NA. A new dye-binder technic using bromocresol purple for Determination of albumin in serum, Clin Chem 1968; 14:793-794.
5. Laskey FD, Li ZMC, Shaver DD, Savory MG, Willey DG, Mikolak BJ, Lantry CL. Evaluation of a Bromocresol purple method for the determination of albumin adapted to the Du Pont aca® discrete-Clinical analyzer, Clin Biochem 1985; 18:290-296.
6. Tietz NW. Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1986, pp 478-497 (specimen collection and storage recommendations).
7. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 1898 USA, 2003.
8. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection: Approved Standard – Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
9. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 (ISBN 1-56238 555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
10. Mabuchi H, Nakahashi H. Underestimation of serum albumin by the bromocresol purple method and a major endogenous ligand in uremia. Clin Chim Acta 1987; 167: 89-96.
11. Maguire GA, Price CP. Bromocresol purple method for serum albumin gives falsely low values in patients with renal insufficiency. Clin Chim Acta 1986; 155: 83-88.
12. Willey DA, Savory J, Lasky F. An Evaluation of a Revised Albumin Method, for the aca® discrete clinical analyzer, Du Pont Company, Wilmington, DE, August 1982.
13. Humphrey KR, Gruemer HD, Lott JA, Impact of posture on the reference range for serum proteins and calcium, Clin Chem 1977; 23:1343–1345.

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Siemens Healthcare Diagnostics Inc.

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