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[siemens.com/healthineers](https://www.siemens.com/healthineers)

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2. Libby RL. New and rapid quantitation technique for the determination of potency of types I and II antipneumococcal serum. J Immunol. 1938;34:269.

3. Gitlin D, Edelhock H. A study of the reaction between human serum and its homologous equine antibody through the medium of light scattering. J Immunol. 1951;66:67-77.

4. Goldberg RJ, Campbell DH. The light scattering properties of an antigen-antibody reaction. J Immunol. 1951;66:70.

5. Johnson P, Ottewill RH. A light scattering study of diphtheria toxin-antitoxin interaction. Disc Faraday Soc. 1954;18:327.

6. Marrack JR, Grant RA. The interaction of antigen and antibody in low concentrations of salt. Brit J Exp Path. 1953;34:263.

7. Marrack JR, Richards CB. Light scattering studies of the formation of aggregates in mixtures of antigen and antibody. J Immunol 1971;20:1019-1038.

8. Ritchie RF. A simple, direct and sensitive technique for measurement of specific protein in dilute solution. J Lab Clin Med. 1967;70:512-517.

9. Hellsing K. Influence of polymers on the antigen-antibody reaction in a continuous flow system. In: Automated Immunoprecipitin Reaction. Colloquium on AIP. Tarrytown, NY: Technicon Inst. Corp.; 1972:17.

10. Hellsing K, Laurent TC. The influence of dextran on the precipitin reaction. Acta Chem Scand. 1964;18:1303.

11. Hellsing K. The effect of dextran on the reaction between iodine-125 labeled human serum albumin and gamma G-globulin from rabbit anti-albumin sera. Acta Chem Scand. 1966;20:1251.

12. Hellsing K. The effect of hyaluronate, chondroitin sulfate and chondroitin sulphate-protein complex on the precipitin reaction. Biochem J. 1969;112:475.

13. Hellsing K. Polysaccharide - enhanced precipitin reactions with antigens of various sizes. Biochem J. 1969;114:145.

14. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.

15. Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.

16. Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.

17. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document

GP44-A4.

18. Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. Clin Chim Acta. 2005;355:5407

19. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. Washington, DC: AACC Press, 2007:537

20. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.

21. Wu AHB, ed. Tietz Clinical Guide to Laboratory Tests. 4th ed. St. Louis, MO: Saunders; 2006:608..

22. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

23. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.

24. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

25. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.

26. Data on file at Siemens Healthcare Diagnostics.

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住所: 511 Benedict Avenue, Tarrytown, New York 10591, USA

生产地址: Ardmore, Diamond Road, Crumlin, Co Antrim, BT29
4QY, United Kingdom

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© 2014 Siemens Healthcare Diagnostics. 保留所有权利。
Siemens Healthcare Diagnostics Inc.