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安全数据表 (MSDS/SDS) 见 siemens.com/Healthineers

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- Shimojo N, Naka K, Nakajima C, Yoshikawa C, Okuda K, Okada K. Test-strip method for measuring lactate in whole blood. Clin Chem. 1989;35(9):1992-1994.
- Shimojo N, Fujino K, Kitahashi S, Nakao M, Naka K, Okuda K. Lactate analyzer with continuous blood sampling for monitoring blood lactate during physical exercise. Clin Chem. 1991;37(11):1978-1980
- Lehninger AL. Principles of Biochemistry. New York, NY: Worth Publishers; 1993:416.
- Mascini M. Lactate and pyruvate electrochemical biosensors for whole blood in extracorporeal experiments with an endocrine artificial pancreas. Clin Chem. 1987;33(4): 591-593.
- Jacobs DS, Kasten BL, DeMott WR, eds. Laboratory Test Handbook. 2nd ed. Stow, OH: LexiComp Inc; 1990:245.
- Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Tietz NW. Clinical Guide to Laboratory Tests. WB Saunders Co., Philadelphia, PA.1990;348.
- Henry RS, Cannon DC, Winkelman JW. Clinical Chemistry Principles and Techniques. New York: Harper and Row; 1974;1330-1333.
- Data on file at Siemens Healthcare Diagnostics.
- Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Tietz NW. Clinical Guide to Laboratory Tests, 3rd Ed. Philadelphia: WB Saunders Company; 1995: 382.
- Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved

Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.

17. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

18. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.

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住所：511 Benedict Avenue, Tarrytown, New York 10591, USA

生产地址：Ardmore, 55 Diamond Road, Crumlin, Co Antrim,
BT29 4QY, UK

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© 2017 Siemens Healthcare Diagnostics. 版权所有。
Siemens Healthcare Diagnostics Inc.