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- Mackey SE, Creasman WT. Ovarian cancer screening. J Clin Oncol. 1995;13(3):783–793.
- Schwartz MK. Cancer markers. In: DeVita VT, Hellman S, Rosenberg SA, eds. Cancer: Principles and Practice of Oncology. Philadelphia: JB Lippincott Company, 1993;534.
- Lavin PT, Knapp RC, Malkasian G, et al. CA 125 for the monitoring of ovarian carcinoma during primary therapy. Obstet Gynecol. 1987;69(2):223–227.
- O'Brien TJ, Hardin JW, Bannon GA, et al. CA 125 antigen in human amniotic fluid and fetal membranes. Am J Obstet Gynecol. 1986;155(1):50–55.
- Kenemans P, Verstraeten AA, van Kamp GJ, von Mendonck-Pouilly S. The second generation CA 125 assays. Ann Med. 1995;27(1):107–113.
- Bonfrer JM, Baan AW, Jansen E, Lentfer D, Kenemans P. Technical evaluation of three second generation CA 125 assays. Eur J Clin Chem Clin Biochem. 1994;32(3):201–207.
- Kenemans P, van Kamp GJ, Oehr P, Verstraeten RA. Heterologous double-determinant immunoradiometric assay CA 125 II: reliable second-generation immunoassay for determining CA 125 in serum. Clin Chem. 1993;39(12):2509–2513.

8. Parker D, Bradley C, Bogle SM, et al. Serum albumin and CA 125 are powerful predictors of survival in epithelial cancer. *Br J Obstet Gynaecol.* 1994;101(10):888-893.
9. Niloff JM, Knapp RC, Lavin PT, et al. The CA 125 assay as a predictor of clinical recurrence in epithelial ovarian cancer. *Am J Obstet Gynecol.* 1986;155(1):56-60.
10. Sevela P, Schemper M, Spona J. CA 125 as an independent prognostic factor for survival in patients with epithelial ovarian cancer. *Am J Obstet Gynecol.* 1989;161(5):1213-1216.
11. Berek JS, Knapp RC, Malkasian GD, et al. CA 125 serum levels correlated with second-look operations among ovarian cancer patients. *Obstet Gynecol.* 1986;67(5):685-689.
12. Niloff JM, Bast RC, Schetzl EM, et al. Predictive value of CA 125 antigen levels in second-look procedures for ovarian cancer. *Am J Obstet Gynecol.* 1985;151(7):981-986.
13. Helzlsouer KJ, Bush TL, Alberg AJ, et al. Prospective study of serum CA-125 levels as markers of ovarian cancer. *JAMA.* 1993;269(9):1123-1126.
14. Mogensen O. Prognostic value of CA 125 in advanced ovarian cancer. *Gynecol Oncol.* 1992;44(3):207-212.
15. Malkasian GD, Knapp RC, Lavin PT, et al. Preoperative evaluation of serum CA 125 levels in premenopausal and postmenopausal patients with pelvic masses: discrimination of benign from malignant disease. *Am J Obstet Gynecol.* 1988;159(2):341-346.
16. Chou CY, Chang CH, Yao BL, et al. Color Doppler ultrasonography and serum CA 125 in the differentiation of benign and malignant ovarian tumors. *J Clin Ultrasound.* 1994;22(8): 491-496.
17. Cannistra SA. Cancer of the ovary. *N Engl J Med.* 1993;329(21):1550-1559.
18. Bast RC, Klug TL, St. John E, et al. A radioimmunoassay using a monoclonal antibody to monitor the course of epithelial ovarian cancer. *N Engl J Med.* 1983;309(15):883-887.
19. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
20. Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
21. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
22. Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
23. Inloes R, Clark D, Drobenies A. Interference of fluorescein, used in retinal angiography, with certain clinical laboratory tests. *Clin Chem.* 1987;33(11):2126-2127.
24. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem.* 1999;45(7):942-956.
25. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem.* 1992;38(9):1737-1742.
26. Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28-A2.
27. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
28. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
29. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
30. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
31. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07-ed3.
32. Wild D, ed. *The Immunoassay Handbook*. 4th ed. Oxford, UK: Elsevier Science; 2013:842.

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Siemens Healthcare Diagnostics Inc.