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- Holick MF. Vitamin D deficiency. *N Engl J Med.* 2007;357(3):266–281.
- Blutt SE, Weigel NL. Vitamin D and prostate cancer. *Proc Soc Exp Biol Med.* 1999; 221(2): 89–98.
- Feskanich D, Ma J, Fuchs CS, et al. Plasma vitamin D metabolites and risk of colorectal cancer in women. *Cancer Epidemiol Biomarkers Prev.* 2004; 13(9):1502–1508.
- Shin MH, Holmes MD, Hankinson SE, et al. Intake of dairy products, calcium, and vitamin D and risk of breast cancer. *J Natl Cancer Inst.* 2002; 94(17):1301–1311.
- Krause R, Bühring M, Hopfenmüller W, et al. Ultraviolet B and blood pressure. *Lancet.* 1998; 352(9129):709–710.
- Lips P. Vitamin D deficiency and secondary hyperparathyroidism in the elderly: consequences for bone loss and fractures and therapeutic implications. *Endocr Rev.* 2001; 22(4):477–501.
- Lips P, Hosking D, Lippuner K, et al. The prevalence of vitamin D inadequacy amongst women with osteoporosis: an international epidemiological investigation. *J Intern Med.* 2006; 260(3):245–254.
- Hyppönen E, Läärä E, Reunanen A, et al. Intake of vitamin D and risk of type 1 diabetes: a birth-cohort study. *Lancet.* 2001; 358(9292):1500–1503.
- Munger KL, Zhang SM, O'Reilly E, et al. Vitamin D intake and incidence of multiple sclerosis. *Neurology.* 2004; 62(1):60–65.
- Merlino LA, Curtis J, Mikuls TR, et al. Vitamin D intake is inversely associated with rheumatoid arthritis: results from the Iowa Women's Health Study. *Arthritis Rheum.* 2004; 50(1):72–77.
- Hollis BW, Horst RL. The assessment of circulating 25(OH)D and 1,25(OH)2D: where we are and where we are going. *J Steroid Biochem Mol Biol.* 2007;103(3-5):473–476.
- Hart GR. Medical conditions associated with vitamin D deficiency and the clinical consequences. *Immunodiagnostic Systems (IDS).* 2004;1:1-10.
- Centers for Disease Control. Perspectives in disease prevention and health promotion update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR.* 1988;37(24):377–382, 387–388.
- Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Data on file at Siemens Healthcare Diagnostics.
- Inloes R, Clark D, Drobenies A. Interference of fluorescein, used in retinal angiography, with certain clinical laboratory tests. *Clin Chem.* 1987;33(11):2126–2127.
- Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem.* 1999;45(7):942–956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polyclonal mouse IgG. *Clin Chem.* 1992;38(9):1737–1742.
- Holick MF. MrOs is D-ficient. *J Clin Endocrinol Metab.* 2009;94(4):1092–1093.
- Rollins G. Vitamin D testing—What's the right answer? Labs grapple with confusing analytics, evidence. *Clin Lab News.* 2009;35(7):1–9.
- Freeman R. Vitamin D: The sunshine hormone. How and when to treat deficiencies. *Menopausal Medicine.* 2009;17(2):S8–S11.
- Misra M, Pacaud D, Petryk A, et al. Vitamin D deficiency in children and its management: review of current knowledge and

recommendations. Pediatrics. 2008;122(2):398–417.

26. Holick MF. The role of vitamin D for bone health and fracture prevention. Curr Osteoporosis Rep. 2006;4(3):96–102.

27. Özkan B, Hatun Ş, Bereket A. Vitamin D intoxication. Turk J Pediatr. 2012;54(2):93–98.

28. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.

29. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

30. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.

31. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

32. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.

33. Sempos CT, Vesper HW, Phinney KW, et al; Vitamin D Standardization Program (VDSP). Vitamin D status as an international issue: national surveys and the problem of standardization. Scand J Clin Lab Invest. 2012;72(suppl 243):32–40.

34. Thienpont LM, Stepman HC, Vesper HW. Standardization of measurements of 25-hydroxyvitamin D3 and D2. Scand J Clin Lab Invest. 2012;72(suppl 243):41–49.

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Siemens Healthcare Diagnostics Inc.