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1. Holick MF. Vitamin D Deficiency. N Engl J Med. 2007;357:266–81.
2. Blutt SE, Weigel NL. Vitamin D and prostate cancer. Proc Soc Exp Biol Med. 1999; 221(2):89–98.
3. Feskanich D, Ma J, Fuchs CS, et al. Plasma vitamin D metabolites and risk of colorectal cancer in women. Cancer Epidemiol Biomarkers Prev.2004;13(9):1502–1508.
4. Shin MH, Holmes MD, Hankinson SE, Wu K, Colditz GA, Willett WC. Intake of dairy products,calcium, and vitamin D and risk of breast cancer. J Natl Cancer Inst. 2002; 94(17):1301–1311.
5. Krause R, Buhring M, Hopfenmuller W, Holick MF, Sharma AM. Ultraviolet B and blood pressure. Lancet. 1998; 352(9129):709–710.

6. Lips P. Vitamin D deficiency and secondary hyperparathyroidism in the elderly: consequences for bone loss and fractures and therapeutic implications. *Endocr Rev.* 2001; 22(4):477–501.
7. Lips P, Hosking D, Lippuner K, et al. The prevalence of vitamin D inadequacy amongst women with osteoporosis: an international epidemiological investigation. *J Unrestricted Med.* 2006;260(3):245–254.
8. Hypponen E, Laara E, Reunanen A, Jarvelin MR, Virtanen SM. Intake of vitamin D and risk of type 1 diabetes: a birth-cohort study. *Lancet.* 2001; 358(9292):1500–1503.
9. Munger KL, Zhang SM, O'Reilly E, et al. Vitamin D intake and incidence of multiple sclerosis. *Neurology.* 2004; 62(1):60–65.
10. Merlino LA, Curtis J, Mikuls TR, Cerhan JR, Criswell LA, Saag KG. Vitamin D intake is inversely associated with rheumatoid arthritis: results from the Iowa Women's Health Study. *Arthritis Rheum.* 2004; 50(1):72–77.
11. Hollis BW, Horst RL. The Assessment of Circulating 25(OH)D and 1,25 (OH)2D: Where We Are and Where We Are Going. *J Steroid Biochem Mol Biol.* March 2007;103(3-5):473–76.
12. Hart GR. Medical Conditions Associated with Vitamin D Deficiency and the Clinical Consequences. Review Series Volume 1, Sept 2004;1–8.
13. Centers for Disease Control. 1988. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR.* 37:377–382, 387, 388.
14. Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
15. Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.
16. Clinical and Laboratory Standards Institute (formerly NCCLS). Procedures for the Handling and Processing of Blood Specimens; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document H18-A3.
17. Data on file at Siemens Healthcare Diagnostics.
18. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem.* 1988;34:27–33.
19. Inloes R, Clark, D and Drobnies A: Interference of fluorescein, used in retinal angiography with certain clinical laboratory tests. *Clin Chem* 1987, 33:2126-2127.
20. Holick MF. MrOs is D-ficient. *J Clin Endocrinol Metab.* 2009;94(4):1092–3.
21. Rollins G. Vitamin D Testing—What's the Right Answer? Labs Grapple with Confusing Analytics, Evidence. *Clinical Laboratory News.* July 2009;35(7): 1,6.
22. Freeman R. Vitamin D: The sunshine hormone. How and when to treat deficiencies. *Menopausal Medicine.* May 2009; S8–11.
23. Misra M, Pacaud D, Petryk A, Collett-Solberg PF, Kappy M. Vitamin D Deficiency in Children and Its Management: Review of Current Knowledge and Recommendations. *Pediatrics.* 2008;122:398–417.
24. Holick MF. The role of vitamin D for bone health and fracture prevention. *Curr Osteoporosis Rep.* 2006;4(3):96–102.
25. Özkan B, Hatun, Ş, Bereket A. Vitamin D intoxication. *Turk J Pediatr.* 2012;54:93–98.
26. Clinical and Laboratory Standards Institute (formerly NCCLS). How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. NCCLS Document C28-A2.
27. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. NCCLS Document EP17-A2.
28. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP5-A2.
29. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
30. Clinical and Laboratory Standards Institute (formerly NCCLS). Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document EP7-A2.
31. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. Wayne, PA: Clinical and Laboratory Standards Institute; 2003. NCCLS Document EP6-A.
32. Sempas CT, Vesper HW, Phinney KW, Thienpont LM, Coates PM, The Vitamin D Standardization Program (VDSP). Vitamin D Status as an International Issue: National Surveys and the Problem of Standardization. *Scandinavian Journal of Clinical & Laboratory Investigation*, 2012; 72(Suppl 243): 32–40.
33. Thienpont LM, Stepman HCM, Vesper HW. Standardization of Measurements of 25-Hydroxyvitamin D3 and D2. *Scandinavian Journal of Clinical & Laboratory Investigation*, 2012; 72(Suppl 243): 41–49.

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Siemens Healthcare Diagnostics Inc.