Anti-CCP IgG

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Anti-CCP IgG

English

Intended Use: The IMMULITE® 2000 Anti-CCP IgG assay is an *in vitro* diagnostic immunoassay for the semi-quantitative determination of the IgG class of autoantibodies specific to cyclic citrullinated peptide (CCP) in human serum (including Serum Separator Tubes) or plasma (EDTA or lithium heparin) on the IMMULITE 2000 Systems Analyzers. Detection of anti-CCP antibodies is used as an aid in the diagnosis of Rheumatoid Arthritis (RA) and should be used in conjunction with other clinical information. Autoantibody levels represent one parameter in a multi-criteria diagnostic process, encompassing both clinical and laboratory-based assessments.

Catalog Number: **L2KPI2** (200 tests)

Test Code: **CCP** Color: **Red**

Summary and Explanation

Anti-cyclic citrullinated peptide (anti-CCP) IgG is a specific marker for rheumatoid arthritis, a disease that affects one percent of the population worldwide.¹ Studies show that the presence of anti-CCP IgG can be detected in the early stages of rheumatoid arthritis²⁻⁴ and that it is indicative of a more progressive form of the disease.⁴

The epitope for CCP is derived from a precursor protein belonging to filaggrin, a group of filament aggregating proteins. During the final stage of epidermal cell differentiation, arginine residues in filaggrin are deiminated to citrulline.⁵ Autoantibodies targeting citrulline-containing peptide sequences in the major antigenic epitopes are markers for rheumatoid arthritis in affected patients.^{5,6}

Various studies have shown between 41–89% sensitivity and 89–99% specificity of anti-CCP for RA diagnosis.⁷ Given that anti-cyclic citrullinated protein antibodies (ACPA) is more specific than Rheumatoid Factor (RF), it is used to distinguish various causes of arthritis.⁸

Principle of the Procedure

IMMULITE 2000 Anti-CCP IgG is a solid-phase, two-cycle sequential chemiluminescent immunometric assay. The solid phase, a polystyrene bead, is coated with a cyclic citrullinated peptide (CCP) antigen. The liquid phase consists of two reagents: 1) protein-based buffer and 2) alkaline phosphatase (bovine calf intestine) conjugated to a monoclonal murine anti-human IgG antibody in buffer.

In the first cycle, the patient sample and the buffer are incubated together with the coated bead for 30 minutes. During this time, human IgG in the sample binds to CCP on the bead. Unbound sample is then removed by centrifugal washes.

In the second cycle, the enzyme-conjugated monoclonal murine anti-human IgG is added to the original reaction tube for additional 30 minutes incubation. The enzyme conjugated monoclonal murine anti-human IgG antibody binds to immobilized anti-CCP IgG on the bead. The unbound enzyme conjugate is removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 2 × 30 minutes

Time to First Result: 65 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Anti-CCP IgG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 5 µL serum or plasma (EDTA and Lithium Heparin). After automatic predilution, 10 µL of diluted sample is used for the immunoassay.

Automatic Predilution Factor: 40

Storage: Serum for 2 days at room temperature, 7 days at 2–8°C; or at or below –20°C up to 12 months.⁹

Plasma for 2 days at room temperature, 7 days at 2–8°C; or at or below –20°C up to 8 months.⁹

When both serum and plasma samples were subject to 4 freeze/thaw cycles, no clinically significant differences were observed.⁹

Thoroughly mix thawed specimens and centrifuge.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹³⁻¹⁵

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H302 + H312,
H411**

**P280, P273,
P301 + P312,
P302 + P312,
P501**

Warning! Harmful if swallowed or in contact with skin. Toxic to aquatic life with long lasting effects.

Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF

SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: sodium azide, 2-methyl-2H-isothiazol-3-one; Anti-CCP IgG Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Caution: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Safety data sheets (MSDS/SDS) available on siemens-healthineers.com.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Anti-CCP IgG Bead Pack (L2PI12)

With barcode. 200 beads, coated with cyclic citrullinated peptide (CCP) antigen. Stable at 2–8°C until expiration date.

L2KPI2: 1 pack

Anti-CCP IgG Reagent Wedge (L2PIA2)

With barcode. 11.5 mL protein-buffer matrix and 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to a monoclonal murine anti-human IgG antibody in buffer. Stable at 2–8°C until expiration date.

L2KPI2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Anti-CCP IgG Adjustors (LPIL, LPIH)

Two vials (Low and High) of lyophilized human serum with IgG reactive to CCP, in buffer with preservative. At least 30 minutes before use reconstitute each vial with 2.0 mL distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. (No further dilution is required.) Stable at 2–8°C for 30 days after reconstitution or for 12 months (aliquotted) at –20°C.

L2KPI2: 1 set

Anti-CCP IgG Controls (LPIC1, LPIC2)

Two vials containing 2 mL each. **LPIC1 (Negative Control):** human serum nonreactive to anti-CCP IgG, with preservative. **LPIC2 (Positive Control):** human serum with IgG reactive to CCP, with preservative. Stable at 2–8°C for 30 days after opening or for 12 months (aliquotted) at –20°C.

L2KPI2: 1 set

For the current control ranges in U/mL, please refer to the control insert.

Aliquot Labels with barcodes are supplied with the kit, for use with the Adjustors and Controls. Before use, place the appropriate Aliquot Labels on test tubes, so the barcodes can be read by the barcode reader on the IMMULITE 2000.

Autoantibody Sample Diluent (L2A1Z)

For the on-board dilution of patient samples and controls. 50 mL concentrated (ready-to-use) protein/buffer matrix. Stable

at 2–8°C for 30 days after opening or for 6 months (aliquotted) at –20°C.

L2KPI2: 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KPI2: 3 labels

**Kit Components
Supplied Separately**

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: The controls supplied with the kit should be used as quality control material to monitor assay performance.

Follow government regulations or accreditation requirements for quality control frequency.

Interpretation of Results

The cutoff of the IMMULITE 2000 Anti-CCP IgG assay was determined with positive and negative patient samples by a ROC analysis, with a balanced consideration of sensitivity and specificity.

Reactive: A result greater than or equal to 4.00 U/mL is considered to be “reactive”, and indicates that anti-CCP IgG antibodies were detected in the sample.

Non-reactive: A result less than 4.00 U/mL is considered to be “nonreactive”, and indicates that anti-CCP IgG antibodies were not detected in the sample.

The Anti-CCP IgG results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Anti-CCP IgG assay. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Expected Values

A total of 197 serum samples from presumed healthy male (n = 99) and female (n = 98) donors were analyzed using the IMMULITE 2000 Anti-CCP IgG assay. The results from this study suggest a median of < 1.50 U/mL and a nonparametric 2.5th–97.5th percentile range of < 1.50–1.93 U/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Anti-CCP IgG test results alone are not sufficient to diagnose a patient with rheumatoid arthritis. Diagnosis should be made by the physician only after all clinical and laboratory findings have been evaluated.

Lipemic, hemolyzed or grossly contaminated samples may give erroneous results.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from neonates, cord blood, or pretransplant patients.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic

antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Reportable Range: 1.50–200 U/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: Limit of Blank (highest value expected for a sample with no analyte; determined in accordance with CLSI EP17-A¹⁰): 0.26 U/mL

Limit of Detection (lowest detectable concentration; determined in accordance with CLSI EP17-A¹⁰): 1.50 U/mL

Functional Sensitivity: (concentration with 20% coefficient of variation (CV) determined):¹⁰ 2.34 U/mL

Precision: Precision was evaluated according to CLSI guideline EP5-A2.¹¹ Serum samples were assayed in duplicate over the course of 20 days, two runs per day, using 3 lots for a total of 40 runs and 80 replicates. (See Precision table.)

Specificity: The specificity assessment of the IMMULITE 2000 Anti-CCP IgG assay was conducted on 8 subgroups of patients (n = 40) with potentially cross-reacting conditions and disease states where other autoantibodies may be present in the patient samples. The following disease states were evaluated: Systemic Lupus Erythematosus (SLE), Sjogren's Syndrome, Lyme Disease, Osteoarthritis, Scleroderma, Polymyositis, Thyroiditis and Mononucleosis.

Disease State	Total Tested (n)	Positive (n)
SLE	5	0
Sjogren's Syndrome	5	1
Lyme Disease	5	1
Osteoarthritis	5	0
Scleroderma	5	0
Polymyositis	5	1
Thyroiditis	5	0
Mononucleosis	5	0

Three out of 40 samples (1 Sjogren's syndrome; 1 Lyme Disease, and 1 Polymyositis) yielded positive results. All samples were also evaluated on a commercially available FDA-cleared predicate testing device for Anti-CCP. The three positive samples were confirmed positive by the predicate testing device.

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 20 mg/dL has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 3500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 500 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 1500 mg/dL has no effect on results, within the precision of the assay.

Rheumatoid Factor: Evaluation of high levels of endogenous rheumatoid factor (RF) was achieved by screening RF positive samples and diluting them to 200 IU/mL. No clinically significant changes were observed.

Alternate Sample Type: To assess the effect of alternate sample types, 39 matched sets of apparently healthy donors were drawn into, heparinized (Li Heparin), EDTA, and BD SST Vacutainer® tubes. Twenty percent of these samples were spiked to achieve anti-CCP levels across the reportable range. All samples were analyzed by the IMMULITE 2000 Anti-CCP IgG assay, with the following results:

(Li Heparin) = 1.02 (Serum) – 0.46 U/mL
r = 1.000

(SST) = 1.02 (Serum) – 0.52 U/mL
r = 1.000

(EDTA) = 1.01 (Serum) – 0.31 U/mL
r = 1.000

Means:

52.3 U/mL (Serum)
52.9 U/mL (Li Heparin)
53.0 U/mL (SST)
52.2 U/mL (EDTA)

This study showed 100% positive and negative agreements between all samples collected in Serum, SST, EDTA and Lithium Heparin.

Method Comparison : The IMMULITE 2000 Anti-CCP IgG assay was compared to a commercially available Axis-Shield DIASTAT Anti-CCP IgG assay on 457 patient samples.

IMMULITE 2000 Anti-CCP IgG

	Non-Reac	% Positive Agreement	% Negative Agreement
DIAGNOSTIC	Reac		
Positive	187	13	93.5%
Negative	14	243	94.6%

Total Agreement: 94.1%

Clinical Performance: The assay was performed on the IMMULITE 2000 to assess the clinical sensitivity and specificity of well-characterized samples.

A total of 457 patient serum samples were collected for the study. 257 samples were rheumatoid arthritis (RA) positive¹² and the remaining 200 samples were from apparently healthy donors (non-RA). RA positive patients were classified according to the ACR criteria.¹²

IMMULITE 2000 Anti-CCP IgG

Clinical Status	Reac	Non-Reac	Sensitivity	Specificity
Positive	201	56	78.2%	99.5%
Negative	1	199		

References

1) Yelin E, Callahan LF. The economic cost and social and psychological impact of musculoskeletal conditions. *Arthritis Rheum.* 1995;38:1351-62. 2) Kroot EJJA, de Jong BAW, van Leeuwen MA, et al. The prognostic value of anti-cyclic citrullinated peptide antibody in patients with recent-onset rheumatoid arthritis. *Arthritis Rheum.* 2000;43:1831-5. 3) Bizzaro N, Mazzanti G, Tonutti E, et al. Diagnostic accuracy of the anti-citrulline antibody assay for rheumatoid arthritis. *Clin Chem.* 2001;47:1089-93. 4) Visser H, Cessie S, Vols K, Breedveld F, Hazes J. How to diagnose rheumatoid arthritis early. *Arthritis Rheum.* 2002;46:357-65. 5) Girbal-Neuhauser E, Durieux JJ, Arnaud M, et al. The epitopes targeted by the rheumatoid arthritis-associated anti-filaggrin autoantibodies are postrationally generated on various sites of (pro)filaggrin by deimination of arginine residues. *J Immunol.* 1999;162:585. 6) Vincet C, Simon M, Sebbag M, et al. Immunoblotting detection of autoantibodies to human epidermis filaggrin: a new diagnostic test for rheumatoid arthritis. *J Rheumatol.* 1998;25:838-46. 7) Silveira I.G., Burlingame R.W., von Muhlen C.A., Bender A.L. and Staub, H.L. (Feb 2007). "Anti-CCP antibodies have more diagnostic impact than rheumatoid factor (RF) in a population tested RF." *Clin. Rheumatol.* DOI 10.1007/s10067-007-0601-6. 8) Avouac J, Gossec L, Dougados M (July 2006). "Diagnostic and predictive value of anti-cyclic citrullinated protein antibodies in rheumatoid arthritis: a systematic literature review." *Ann. Rheum. Dis.* 65 (7): 845-51. 9) Data on file at Siemens Healthcare Diagnostics. 10) CLSI. Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI document EP17-A Vol. 24 (No. 34). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004. 11) Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Clinical Devices; Approved Guideline. Wayne PA: Clinical and Laboratory Standards Institute; 1999. NCCLS Document EP5-A. 12) Arnett FC, Edworthy SM, Bloch DA, et al. The American Rheumatism Association 1987 revised criteria for the classification of rheumatoid arthritis. *Arthritis Rheum* 1988;31(3):315-24. 13) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8.

14) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 15) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

siemens-healthineers.com

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE is a trademark of Siemens Healthineers.
All other trademarks are the property of their respective owners.
© 2018-2022 Siemens Healthineers. All rights reserved.
CONTROL AB